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RECURSOS

1. Technical Manuel, 15th Edition, American Association of Blood Banks, Bethesda, 2005

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INDICATION

Across Identi Cell Extended, Across Identi Cell Extended P products complement Across Identi Cell and Across Identi Cell P to identify unexpected antibodies in immunohematological tests. During antibody identification one very often requires extra cells of certain phenotypes to confirm possible antibodies. Across Identi Cell Extended, Across Identi Cell Extended P offer such additional cells.

INTRODUCTION

Identification of unexpected antibodies helps to ensure detection of clinically important antibodies in patient/donor sample. The antibody detected as a result of the test is an indication of the presence of an autoantibody, an alloantibody or both (1). The cold agglutinin identification test are used to help detection of antibodies that react at 4 °C or at 18 – 25 °C (room temperature), such as anti-M, anti-N, anti-P1, anti-Lea, anti-Leb, anti-I which are freely available in patient or donor serum. They are also used to define hemolytic anemia associated with cold antibodies.

PRINCIPLE

The principle of this test is based on the gel centrifugation technique to detect agglutination reactions that result from binding of antigens in red blood cells with the corresponding antibodies present in reagent, serum or plasma samples.

CONTENT

Each vial contains human erythrocytes of O blood group with or without papain, in a 0,8 % suspension in a buffer solution containing preservatives (sulfamethoxazole and trimethoprim). Across Identi Cell Extended (4 test cells), Across Identi Cell Extended P (4 test cells P) are reagent red blood cells with polyvalent antigens of four single blood donors in separate vials for the identification of red blood cell antibodies. Reagent red blood cells is given in the supplement table (specific for each product lot). Panel of red blood cell reagents show the most important antigenic determinants of most blood group systems. These red blood cells are different in terms of antigenic configuration and have been selected to help detection of most clinically important antibodies.

STORAGE CONDITIONS

- Keep at 2-8°C. Prevent contamination after opening vials.
- Do not freeze cards.
- Stable until the expiry date on the label if stored in a suitable conditions.

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- Recently extracted plasma/serum samples collected in tubes with or usual blood bank anticoagulants (containing plasma) or without anticoagulants (containing serum).
- Do not use hemolyzed, cloudy, contaminated blood samples or with clot presence.
- Blood samples stored at 2-8°C and taken 48 hours ago can be used.
- Plasma/serum samples stored between -20 / -80°C can be used.
- Reagent red blood cells are ready to use.
- Perform necessary centrifuging procedure for clean plasma or serum.

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be controlled/checked.

TEST PROCEDURE

Across Identi Cell Extended and Across Identi Cell Extended P are used with manual system and Across Auto System octo-m. For Across Auto System octo-m apply the procedure in the user manual.

Allow the samples and reagents to reach room temperature (18-25°C). Inspect the state of red blood cell reagent before use (do not use cloudy hemolyzed/hemolysed reagents.) In manual equipment and Across Auto System octo-m devices, Working principles of Across Identi Cell Extended, Across Identi Cell Extended P have equivalent properties and example of working method with manual equipments is given below.

Pool Cell cells are used only for donor antibody identification (Not used for the patient).

MANUAL METHOD

1. Paste sample identification onto the card by barcode or writing manually.
2. Open the aluminium folio carefully keeping the card upright to avoid reagent contamination.

Cold Agglutinin Test

- Homogenize the papain free Across red blood cell (820114) reagents
- Pipette 50 µl Across red blood cell reagent to corresponding microtubes (Neutral) to be used.
- Pipette 25 µl of patient serum or plasma to correspond - ing microtubes.
- Incubate for 15 minutes at 4°C. Incubate for 15 minutes at 4°C.

Antibody Identification

- Homogenize the papainized and/or papain free Across red blood cell reagents (820114 and/or 820115).
- Add 50 µl Across red blood cell reagent to correspond - ing microtubes (Neural and/or AHG)
- Add 25 µl of patient serum or plasma to corresponding microtubes
- Incubate for 15 minutes at 37 °C
- 3. Centrifuge card or cards for 9 min at 990 rpm.
- 4. Read the results

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Note: Test cell reagents for antibody identification containing pooled cells can show a double cell population appearance depending on the antibody present. This is considered to be a positive result.

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www.diapro.com.tr

Revision No:.....

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ACROSS CELL SCREEN 2 / ACROSS CELL SCREEN 2P / ACROSS CELL SCREEN 3 / ACROSS CELL SCREEN 3P / ACROSS CELL SCREEN 4 / ACROSS CELL SCREEN 4P / ACROSS IDENTI CELL / ACROSS IDENTI CELL P / ACROSS POOL CELL

For 'in vitro' diagnostic use only. Read carefully before using the product.

INDICATION

Across Cell Screen 2, Across Cell Screen 2P, Across Cell Screen 3, Across Cell Screen 3P, Across Cell Screen 4, Across Cell Screen 4P, Across Pool Cell are used for screening of unexpected antibodies and Across Identi Cell, Across Identi Cell P are used for identification of unexpected antibodies in immunohematologic test.

INTRODUCTION

Screening of unexpected antibodies ensures detection of clinically important antibodies in patient/donor sample. Positive test result is an indication of presence of an autoantibody, an alloantibody or both⁽¹⁾.

The cold agglutinin test are used to detect antibodies that react at 4 °C or at 18–25 °C (room temperature), such as anti-M, anti-N, anti-P1, anti-Lea, anti-Leb, anti-I which are freely available in patient or donor serum. They are also used to define haemolytic anemia associated with cold antibodies.

PRINCIPLE

The principle of this test is based on the gel centrifugation technique to detect agglutination reactions that result from binding of antigens on Red Blood Cells Reagent with the corresponding antibodies present in serum or plasma samples.

CONTENT

Each vial contains human erythrocytes of O blood group with or without Papain, in a 0,8 % suspension in a buffer solution containing preservatives.

Across Cell Screen 2 (R1R1/ R1wR1 + R2R2) Across Cell Screen 2P (R1R1P/ R1wR1P + R2R2P), Across Cell Screen 3 (R1wR1 + R2R2 + rr), Across Cell Screen 3P (R1wR1P + R2R2P + rrP), Across Cell Screen 4 (R1R1 + R2R2 + R1wR1 + rr), Across Cell Screen 4P (R1R1P + R2R2P + R1wR1P + rrP), Across Pool Cell (R1R1/R1wR1 + R2R2), Across Identi Cell (11 test cells), Across Identi Cell P (11 test cells P)

Across Cell Screen 2, Across Cell Screen 2P, Across Cell Screen 3, Across Cell Screen 3P, Across Cell Screen 4, Across Cell Screen 4P are Red Blood Cells reagent with polyvalent antigens of two, three or four single blood donors in separate vials for the detection of red blood cell antibodies. Across Pool Cell is a single vial containing reagent Red Blood Cells with polyvalent antigens of two single blood donors. Across Identi Cell (11 test cells), Across Identi Cell P (11 test cells P) are reagent red blood cells with polyvalent antigens of eleven single blood donors in separate vials for the identification of red blood cell antibodies. Reagent Red Blood Cells is given in the supplement table (specific for each product lot).

Panel of Red Blood Cell reagents show the most important antigenic determinants of most blood group systems. These red blood cells are different in terms of antigenic configuration and have been selected to help detection of most clinically important antibodies.

STORAGE CONDITIONS

- Keep at 2-8 °C. Prevent contamination after opening vials.
- Do not freeze cards.
- Stable until the expiry date on the label if stored in a suitable conditions.

Notes:

- Do not use the product if it is detected that the temperature conditions are wrong during the storage and transportation of the product.
- The used vial should be stored at a re-storage temperature.
- The caps of vials should be kept closed to prevent contamination during the period when the product is not used.

MATERIAL REQUIRED FOR TEST BUT NOT PROVIDED IN THE KIT

- Across Auto System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Incubator
- Across System Card Reader
- Across System Pipette
- Across System Workstation
- Disposable Pipette Tips
- Across Gel Cards

Please refer to the Instructions for Use of Across Gel Cards.

COLLECTION AND PREPARATION OF SAMPLES

Samples

For detection of unexpected antibodies in the serum or plasma;

- Serum or plasma samples collected in tubes with or without anticoagulant can be used. Anticoagulants used in the blood bank should be used as anticoagulants.
- Do not use serum or plasma samples that are hemolyzed, cloudy, contaminated or clot presence.
- Serum or plasma samples stored at 2-8 °C up to 48 hours after extraction can be used.
- Plasma/serum samples stored between -20 / -80 °C can be used.
- Reagent Red Blood Cells are ready to use.
- Perform necessary centrifuging procedure for clean plasma or serum.

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be controlled.

TEST PROCEDURE

Across Cell Screen 2, Across Cell Screen 2P, Across Cell Screen 3, Across Cell Screen 3P, Across Cell Screen 4, Across Cell Screen 4P, Across Pool Cell, Across Identi Cell and Across Identi Cell P are used with manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments apply the procedure in the user manual.

Allow the samples and reagents to reach room temperature (18-25 °C). Check the condition of Red Blood Cell Reagent before use (do not use cloudy, change of color or haemolysed reagents.) Follow the instructions for use of the Across Gel Cards used. Do not use the card if there are bubbles in the microtubes on the relevant cards, if there is no supernatant, the amount of gel decreases, the gel cracks or changes in the microtube color. In manual equipment and Across Auto System semi-automatic or automatic instruments, working principles of the Across Cell Screen 2, Across Cell Screen 2P, Across Cell Screen 3, Across Cell Screen 3P, Across Cell Screen 4, Across Cell Screen 4P, Across Pool Cell, Across Identi Cell, Across Identi Cell P have equivalent properties and example of working method with manual equipments is given below.

Pool Cell cards are used only for donor antibody screening (Not used for the patient).

MANUAL METHOD

1. Paste sample identification onto the card by barcode or writing manually.
2. Open the aluminium foil carefully keeping the card upright to avoid reagent contamination.

Cold Agglutinin Test

- Homogenize the Papainized and Papain free Across Red Blood Cell Reagents.
- Pipette 50 µl Across Red Blood Cell Reagent to corresponding microtubes to be used.
- Pipette 25 µl of patient serum or plasma to corresponding microtubes.
- Incubate for 15 minutes at 4 °C.

Indirect Coombs / Antibody Screening

- Homogenize the Papainized and Papain free Across Red Blood Cell Reagents (820103, 820104, 820105, 820106, 820111, 820112, 820113)
- Add 50 µl Across Red Blood Cell Reagent to corresponding microtubes.
- Add 25 µl of patient/donor serum or plasma to corresponding microtubes (Only donor samples are used for Pool Cell).
- Incubate for 15 minutes at 37 °C.

Indirect Coombs / Antibody Identification

- Homogenize the Papainized and Papain free Across Red Blood Cell Reagents (820107, 820108).
- Add 50 µl Across Red Blood Cell Reagent to corresponding microtubes.
- Add 25 µl of patient serum or plasma to corresponding microtubes.
- Incubate for 15 minutes at 37 °C.
- 3. Centrifuge card or cards for 9 mins at 990 rpm.
- 4. Read the results.

RESULTS

Result reading:

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

Note: Test cell reagents for antibody screening containing pooled cells can show a double cell population appearance depending on the antibody present. This is considered to be a positive result.

It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position. Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8 °C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results

Red Blood Cell Reagent	Interpretation
0	Absence of alloantibodies and/or autoantibodies
+	Absence of alloantibodies and/or autoantibodies

= positive

0= negative

Positive or negative reactions should be evaluated according to the antigen present in the antigen table.

PERFORMANCE CHARACTERISTICS

The studies were carried out based on antibody screening results by selecting 200 samples from random population using Across Screening and Identification Cell products and 202 donor samples for Across Pool Cell product. Although there is not any defined procedure or technique can precisely determine all of the unexpected antibodies in a sample population performance evaluation studies realized have shown that Across Cell products are suitable for intended use. These studies contain tests for screening and identification of unexpected antibodies and gave results that are equivalent to or better than the results of other products with the same intended use.

WARNINGS AND PRECAUTIONS

- Results cannot be used alone for diagnostic use. They should be evaluated with other clinical and laboratory findings of the patient.
- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Do not use cloudy or haemolysed reagents.
- If poor anti-coagulated plasma or partially coagulated serum is used, the fibrin residue may trap non-agglutinated Red Blood Cells at the top of the gel, so that a pinkish or reddish layer can be observed. As a result, although the results can be interpreted correctly, the false appearance of the double population could lead to misinterpretation when the result is negative. If serum samples that are not fully clotted are used, it is recommended to centrifuge the serum after 10 minutes of incubation at 37 °C and repeat the test.
- The reagents are of human origin. Red Blood Cells used in the preparation of this product obtained from HBsAg, anti-HIV, anti-HCV negative donors. However, all products and samples should be handled as if they were potentially capable of transmitting infectious diseases.
- Due to the low-incidence antigens are not included in the Across Red Blood Cells, negative reactions with them do not always indicate the absence of antibodies in the sample.
- It is in the nature of a Pool cell that the antigen dose of antigens may be reduced. The reduced antigen density leads to reduced sensitivity for the respective antigen. Therefore, the sensitivity of the pool cell may be slightly lower compared to the sensitivity of a screening cell of a two cell.
- Due to the pool of two cells of two distinct donors, it is possible that the Pool Cell is positive and negative for a specific antigen, which leads in the case of a positive result in antibody screening to a double population. Such a result is to be interpreted as positive.
- The presence of high concentration of proteins, abnormal proteins, medications, plasma volume intensifier, etc. In the serum causes "rouleaux formation", leading to false positive results.
- The reagent Red Blood Cells do not detect the presence of anti-A or anti-B antibodies.
- In case of presence of multiple antibodies or the antibodies corresponding to frequently encountered antigens in the sample, all Red Blood Cell reagents can be agglutinated.
- If the plasma is used, complement dependent haemolytic reactions may not be detected.
- When interpreting the reactions, ensure that the lot number of the antigen table corresponds with the lot number of the reagent vials.
- Do not use any product beyond the expiry date.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent.

REFERENCES

1. Technical Manuel, 15th edition, American Association of Blood Banks, Bethesda, 2005

PACKAGING


Across Cell Screen 2, Across Cell Screen 2P, Across Cell Screen 3, Across Cell Screen 3P, Across Cell Screen 4, Across Cell Screen 4P, Across Identi Cell, Across Identi Cell P, Across Pool Cell

Ref No:

820103 Across Cell Screen 2	2x10 ml
820104 Across Cell Screen 2P	2x10 ml
820111 Across Cell Screen 3	3x10 ml
820112 Across Cell Screen 3P	3x10 ml
820105 Across Cell Screen 4	4x10 ml
820106 Across Cell Screen 4P	4x10 ml
810107 Across Identi Cell	11x5 ml
820108 Across Identi Cell P	11x5 ml
820113 Across Pool Cell	1x10 ml

Document No: TK3435

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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