STANDARD F TnI/CK-MB Combo FIA

STANDARD™ F TnI/CK-MB Combo FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer

SPECIMEN COLLECTION AND STORAGE

Whole blood

- Collect the venous whole blood into the commercially available EDTA tube by venipuncture.
- It is recommended that collected venous whole blood specimens are used immediately. If venous whole blood in an anticoagulant tube is stored at room temperature and refrigerated condition (2-8°C/36-46°F), the specimen can be used for testing within 8 hours after collection.
- Do not use hemolyzed blood specimens.



Serum

- Collect the whole blood into the commercially available plain tube, NOT containing anti-coagulants such as heparin, EDTA or sodium citrate by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
- Serum specimen may be stored at room temperature and refrigerator conditions (2-8°C/36-46°F) for up to 24 hours prior to testing.
- It should be brought to room temperature prior to use.

PREPARATION AND TEST PROCEDURE


Preparation

- Allow kit components and collected specimen to room temperature at least 30 minutes before starting the test.
- Carefully read instructions for the STANDARD F TnI/CK-MB Combo FIA.
- Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
- Open the foil pouch and check the test device in the foil pouch.

Analysis of specimen

'STANDARD TEST' mode	STANDARD F200 and F2400 Analyzer
1. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.	
2. Insert the test device to the test slot of the analyzer. When inserting the analyzer read the barcode data, and check the test device is valid.	
3. Using a disposable dropper, collect 100µl of the specimen (whole blood, serum) in the black line of the disposable dropper. Add the collected specimen into the extraction buffer and mix thoroughly.	
4. Apply 100µl of specimen mixture to the specimen well of the test device.	

System Accuracy	Serum	Whole Blood		
TnI	CK-MB	TnI	CK-MB	
Below -25%	2 / 210 (1%)	0 / 210 (0%)	3 / 210 (1.4%)	0 / 210 (0%)
Within ±25%	206 / 210 (98%)	210 / 210 (100%)	207 / 210 (98.6%)	210 / 210 (100%)
Over +25%	2 / 210 (1%)	0 / 210 (0%)	0 / 210 (0%)	0 / 210 (0%)

PERFORMANCE CHARACTERISTICS

Analysator performance

1) Accuracy (Method Comparison)
Results comparing the STANDARD F TnI/CK-MB Combo FIA with an automated chemiluminent method are shown below.

Regression Analysis				
	Serum		Whole Blood	
	TnI	CK-MB	TnI	CK-MB
Slope	0.9558	1.0016	1.0039	0.9914
Y-intercept	105.9	-0.0267	0.3879	0.2689
R	0.9954	0.9965	0.9961	0.9969
R ²	0.9909	0.9931	0.9922	0.9939
n	210	210	210	210

System Accuracy				
	Serum		Whole Blood	
	TnI	CK-MB	TnI	CK-MB
Below -25%	2 / 210 (1%)	0 / 210 (0%)	3 / 210 (1.4%)	0 / 210 (0%)
Within ±25%	206 / 210 (98%)	210 / 210 (100%)	207 / 210 (98.6%)	210 / 210 (100%)
Over +25%	2 / 210 (1%)	0 / 210 (0%)	0 / 210 (0%)	0 / 210 (0%)

INTERPRETATION OF TEST RESULTS

Troponin I

Measuring Range	10 - 20,000 ng/L (0.01 - 20 ng/mL)
Out of measuring range	- If below 10 ng/L, "1 < 10 ng/mL" message will be displayed. - If above 20,000 ng/L, "1 > 20,000" message will be displayed.
CV ≤ 10% at 99th Percentile URL	70.0 ng/L (0.07 ng/mL)

TnI Concentration (ng/L)				
Population	N	99th Percentile URL	90% Confidence Interval	Maximum Observed
	732	70.0	54.0 - 81.0	100.0

Troponin I concentrations were determined using serum specimens obtained from 732 apparently healthy individuals. The upper reference limit (99th percentile) for TnI was 70 ng/L, and the lowest TnI concentration with CV 10% was found to be 70 ng/L.

Elevated troponin is specific to myocardial injury and may be due to conditions other than MI. Acute MI can be established if there is a rise and/or fall in cardiac troponin with at least one value above the 99th percentile URL, together with a clinical likelihood or ECG evidence of myocardial ischemia.

The default results are measured in ng/L. Should an alternate result unit be desired, the user can go to the "Settings" menu of the analyzer or manually calculate using the following conversion factors:

(Concentration in ng/L) x (0.001) = ng/mL = µg/L
(Concentration in ng/L) x (1) = pg/mL

CK-MB

Measuring Range	1-200 ng/mL
Out of measuring range	- If below 1 ng/mL, "1 < 1" message will be displayed. - If above 200 ng/mL, "200 > 1" message will be displayed.
Cut-off value for MI	< 5.0 ng/mL
CK Index calculation	People with skeletal muscle damage may have elevated CK-MB levels. To help differentiate heart attack from skeletal muscle damage, a CK index can be calculated using CK-MB and total CK as follows: CK index = (CK-MB, ng/mL) x100 / (total CK activity, IU/L). If CK-MB is elevated and the CK index is higher than 2.5 to 3.0, heart damage is likely.

The Troponin I and CK-MB reference ranges are provided for orientational purpose only. Clinicians should use the test results in conjunction with the patient's other disease findings and clinical signs, and interpret the concrete values in the context of the patient's clinical situation.

EXPLANATION AND SUMMARY

Introduction

The Troponin (Tn), one of the subunits of the troponin regulatory complex, binds to actin and inhibits interactions between actin and myosin. Troponin complex plays an important role in the regulation of skeletal and cardiac muscle contraction. Troponin complex consists of three subunits: Troponin T (TnT), Troponin I (TnI) and Troponin C (TnC). The subunits are held together by non-covalent interactions. TnI has two important functions: it regulates interactions between the troponin subunits and the thin filament. The TnI subunit is responsible for inhibiting actomyosin formation at low intracellular Ca2+ concentrations. The third subunit, TnC, binds Ca2+ ions during the excitation of the muscle and changes the conformation of troponin complex, thus enabling the formation of actomyosin complex and the subsequent muscle contraction. Cardiac Troponin I (cTnI) has been known as reliable markers of cardiac cell death. It is now widely used for the diagnosis of acute myocardial infarction (AMI), postmyocardium trauma, chemotherapy cardiotoxicity, as well as some other diseases related to cardiac muscle injury. Being expressed cTnI only in cardiac tissue, cardiac TnI have been the preferred biomarkers for myocardial infarction for a long time. Creatine kinase (CK), also known as Creatine Phosphokinase or Phospho-creatine Kinase is an enzyme expressed by various tissues and cell types. Disruption of cell membranes due to hypoxia or other injury releases CK from the cellular cytosol into the systemic circulation. CK enzyme is a dimer formed by the association of two subunits from muscle (Muscle type) or nerve cells (Brain type). These subunits associate to form three isoenzymic forms: CK-BB, CK-MM and CK-MB. These isoenzymes are expressed at different levels in various human tissues. Though CK-MM is the most abundant CK isoenzyme in the cardiac muscles, CK-MB constitutes about 20% of the total CK in the cardiac muscle tissue. Elevated levels of total CK is not specific to the myocardial tissue and may be observed in patients with skeletal muscle injury and certain other disorders but as CK-MB is more specific to myocardial tissue, CK-MB levels along with total CK can be considered as an important diagnostic indicator of myocardial infarction. The concentration of CK-MB in the healthy adult maintains low concentration but it shows great increases in several malignant diseases, mostly primary coronary syndrome, myocardial injury and infarction. CK-MB has been found to be more sensitive and early indicator of myocardial injury because it has a lower basal level and a much narrower normal range. Medical literature commonly reveals that following an acute myocardial infarction, CK-MB levels become elevated in 4 to 9 hours after the onset of chest pain, attain peak at 10 to 24 hours, and return to normal within 2 to 3 days. Use of CK-MB level as a percentage of total CK in the diagnosis of myocardial infarction is the most important clinical application of CK measurements in clinical chemistry. STANDARD F TnI/CK-MB Combo FIA, employing immunofluorescent detection system with STANDARD F Analyzer, provides significantly fast, easy and accurate system to identify cTnI and CK-MB levels in human serum and whole blood as an aid in the diagnosis of myocardial infarction and enables supportive treatment decisions.

Intended use

STANDARD F TnI/CK-MB Combo FIA is a fluorescent immunoassay for the quantitative determination of cardiac troponin I and total creatine kinase isoenzyme-MB(CK-MB) levels in human serum and whole blood using STANDARD F Analyzers manufactured by SD BIOSENSOR. This test is an *in vitro* professional diagnostic use and intended for use as an aid in the screening and monitoring of myocardial infarction(MI).

Test Principle

STANDARD F TnI/CK-MB Combo FIA uses antibodies that specifically bind to the TnI and CK-MB proteins. When a human specimen is applied into the specimen well of the test device, the specimen migrates through the membrane from the specimen well. If Troponin I (TnI) and/or CK-MB is present, it will be bound to the detector coupled to europium microparticle that migrates through the membrane. The europium microparticle complex will be captured by capture antibody on the test line where it is detected by STANDARD F Analyzer. If is not present, the europium microparticle will not be trapped by the capture antibody nor detected by STANDARD F Analyzer. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzer can analyze the presence of the TnI and CK-MB simultaneously in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

KIT STORAGE AND STABILITY

For 18 months when stored at 2-30°C/36-86°F and 10-90% RH.

WARNINGS AND PRECAUTIONS

- STANDARD F TnI/CK-MB Combo FIA is for *in vitro* diagnostics use only.
- Carefully follow the instructions and procedures described in this Instructions for use before testing.
- STANDARD F TnI/CK-MB Combo FIA should be used with STANDARD F Analyzer.
- STANDARD F TnI/CK-MB Combo FIA should remain in its original sealed pouch until ready to use. Do not use the test device if pouch is damaged or the seal is broken.
- STANDARD F TnI/CK-MB Combo FIA is single use only. Do not re-use it.
- Do not use hemolyzed specimens or frozen specimens.
- Do not use any artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm and soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Desiccant in foil pouch is absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.

- Check the expiration date printed at the top or package.
- Check the volume (100µl) of extraction buffer.
- Use the STANDARD F TnI/CK-MB Combo FIA at room temperature (15-32°C/59-89.6°F).
- All kit components are must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

LIMITATION OF TEST

- The test should be used for the detection of Troponin I and CK-MB levels in human whole blood and serum specimens.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid result.
- Invalid results may occur if a poor-quality specimen is obtained.
- The test result must always be evaluated with other clinical data available to the physician.

QUALITY CONTROL

Calibration

The calibration set test of STANDARD F Analyzer should be conducted according to the analyzer's manual.

When to use calibration set

- Before using the analyzer for the first time.
- When you drop the analyzer.
- Whenever you do not agree with the final result.
- When you want to check the performance of the analyzer and test device.

How to use calibration set

- Calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.
- Select the "Calibration" in main menu
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 for white calibration, CAL-2 for UV LED calibration, and CAL-3 for RGB LED calibration in sequence.

STANDARD F Analyzer automatically calibrate and identify the optical performance through measuring the membrane of the test device whenever the test is conducted in STANDARD TEST™ mode. If "EEF" message displays on the screen, it means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the "EEF" message still appears.

Internal quality control

- The internal procedural control zone is on the membrane of the test device. STANDARD F Analyzers read the fluorescence signal of the internal procedural control zone and decide whether the result is valid or invalid.
- The invalid result denotes that the fluorescence signal is not within the pre-set range. If the screen of STANDARD F Analyzer shows "Invalid", turn off and turn on the analyzer again and re-test with a new test device.

External quality control

Quality control testing should be run to check the performance of STANDARD F TnI/CK-MB Combo FIA and STANDARD F Analyzer. STANDARD F TnI Control and STANDARD F CK-MB Control manufactured by SD BIOSENSOR can be used for the external quality control. The test should be conducted in accordance with the Instructions for use of STANDARD F TnI Control and STANDARD F CK-MB Control. External quality control test should be run:

- once for each new lot.
- once for each untrained operator.
- as required by test procedures in the Instructions for use of STANDARD F TnI Control and STANDARD F CK-MB Control in accordance with local, state and federal regulations or accreditation requirements.

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

Within Run												
TnI (ng/L)	Level 1			Level 2			Level 3					
	Site 1	Site 2	Site 3	Site 1	Site 2	Site 3	Site 1	Site 2	Site 3	Site 1	Site 2	Site 3
N	100	100	100	100	100	100	100	100	100	100	100	100
Ref.	70	70	70	6500	6100	5800	13420	12920	13130			
AVG.	62.67	62.67	62.14	6271.2	6659.3	6479.7	14086.7	13000.6	14854.4			
CV (%)	7.04%	7.62%	6.79%	6.34%	7.41%	11.25%	7.44%	5.39%	9.84%			
DIFF (%)	-10.47%	-10.31%	-9.80%	3.40%	9.17%	11.72%	4.97%	6.02%	13.10%			
N	100	100	100	100	100	100	100	100	100	100	100	100
Ref.	70	70	70	6500	6100	5600	13240	15120	14310			
AVG.	62.38	63.02	62.22	6984.9	6413.6	6063.5	14097.2	15651.1	14632.2			
CV (%)	7.68%	7.64%	7.49%	6.04%	5.73%	6.69%	7.87%	5.80%	8.85%			
DIFF (%)	-10.89%	-9.97%	-11.11%	7.46%	5.14%	8.28%	6.47%	3.51%	3.06%			

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