CA50

Carbohydrate Antigen 50 (CLIA)

Order Information

Catalog No.	Package Size
1001446	2x50 tests
1001447	2x100 tests

Intended Use

The CLIA series CA50 assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of carbohydrate antigen 50 (CA50) in human serum or plasma. This assay is to be used as an aid in monitoring disease of patients with pancreatic cancer.

Summary

Carbohydrate antigen 50 (CA50) is a tumor related antigen discovered by Sweden scientist Lindholm in 1983. The antigenic determinant of CA50 consists of monosialoganglioside and sialoglycoprotein.¹⁻³ The CA50 is not expressed in healthy tissues, however, it is highly expressed in many malignant tumor especially gastrointestinal tumor and non-gastrointestinal epithelial cell tumor.⁴⁻⁶ The increase of serum CA50 level is related with pancreas cancer, gallbladder cancer, colorectal cancer, hepatobiliary cancer, gastric cancer and ovarian cancer. The level of CA50 will increase in the occurrence of pancreatitis, colitis and pulmonitis. It will decrease once the inflammation is eliminated.^{7,8}

Assay Principle

The CLIA series CA50 assay is a two-site sandwich assay.

In the first step, sample, paramagnetic microparticle coated with monoclonal anti-CA50 antibody (mouse) and monoclonal anti-CA50 antibody (mouse)-alkaline phosphatase conjugate are added into a reaction cuvette. After incubation, CA50 present in the sample binds to both anti-CA50 antibody coated microparticle and anti-CA50 antibody alkaline phosphatase-labeled conjugate to form a sandwich complex. Microparticle is magnetically captured while other unbound substances are removed by washing.

In the second step, the substrate solution is added to the reaction cuvette. It is catalyzed by anti-CA50 antibody (mouse)-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of CA50 present in the sample is proportional to the light units (RLUs) generated during the reaction. The CA50 concentration can be determined via a calibration curve.

Reagent Components

Components in different lots of reagent kit cannot be interchanged:

Ra	Paramagnetic microparticles coated with monoclonal anti-CA50 antibody (mouse). Minimum concentration: 0.5 g/L solids. TRIS ^{a)} buffer: 50 mmol/L. Preservatives: 0.02% BND and 0.05% ProClin 300.
Rb	Monoclonal anti-CA50 antibody (mouse)-alkaline phosphatase conjugate. Minimum concentration: 0.05 mg/L. MES ^{b)} buffer: 50 mmol/L. Preservatives: 0.02% BND and 0.05% ProClin 300.

- a) TRIS=Tris (hydroxymethyl)-aminomethane
b) MES=2-(N-morpholino) ethanesulfonic acid

Storage and Stability

The unopened CA50 (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The CA50 (CLIA) reagent kit can be stored on-board and used for a maximum of 56 days after opening at 2-8°C.

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Specimen Types

- Human serum or plasma collected in EDTA-K₂, EDTA-K₃, sodium heparin and lithium heparin are recommended for this assay.
- Blood collection tubes from various manufacturers may contain different raw materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Wiener lab. Each laboratory should determine the acceptability of the tubes and serum/plasma separation products.

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Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter have been removed prior to analysis.
- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

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Assay Procedure

For optimal performance of this assay, operators should read the related system operator's manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as

well.

Before loading the CA50 (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Wiener lab. Customer Service for help. Do not invert opened reagent bottle.

This assay requires 20 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operator's manual and specific requirement of the assay to determine the minimum sample volume.

Calibration

CLIA series CA50 (CLIA) has been standardized against a commercial CA50 test (CLIA).

The specific information of master calibration curve of CA50 (CLIA) reagent kit is stored in the two-dimensional barcode on the reagent pack, which is used in combination with CA50 calibrators for calibration of the specific reagent lot. Before initiating calibration on each new lot of reagent, load the assay master curve by scanning the two-dimensional barcode on the reagent pack. When performing the calibration, scan the two-dimensional barcode on the calibrator pack, and then test the CA50 calibrators at three levels. Valid calibration curve is required before any CA50 test. Recalibration is recommended every 28 days, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

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Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three-level calibrators of defined concentration values. The results are shown in the unit of U/mL.

Dilution

Samples with CA50 concentrations above the upper limit can be diluted with Wiener lab. Sample Diluent. The recommended dilution is 1:40 (either automatically by the analyzer or manually). The maximum manual dilution is 1:100. The concentration of the diluted sample must be > 2 U/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzers, the system automatically multiplies the result by the dilution factor when calculating the sample concentration.

Expected values

An extensive study on a cohort of 130 healthy males and 120 healthy females has determined the reference range of CLIA series CA50 assay.

Category	N	95 th percentile
Healthy	250	25 U/mL

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Limitation

The upper limit of this assay is 500 U/mL. A specimen with a CA50 concentration lower than the upper limit can be quantitatively determined, while the specimen with a concentration higher than the upper limit will be reported as > 500 U/mL.

The concentration of CA50 in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

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Performance Characteristics

Lower limits of measurement

The CA50(CLIA) reagent kit has a Limit of Blank (LoB) of ≤ 0.3 U/mL, a Limit of Detection (LoD) of ≤ 1.0 U/mL and a Limit of Quantitation (LoQ) of ≤ 1.0 U/mL (with a total allowable relative error of ≤ 25%).

Measuring range

1.0-500 U/mL (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 1.00 U/mL. Values above the measuring range are reported as > 500 U/mL (or up to 50,000 U/mL for 100-fold diluted samples).

Analytical Specificity

Hemoglobin up to 3,500 mg/dL, bilirubin up to 70 mg/dL, triglycerides up to 3,500 mg/dL, total protein up to 15 g/dL biotin up to 5,000 ng/mL, and antinuclear antibody will not interfere with the CLIA series CA50 assay. Criterion: Recovery within ± 10% of initial value.

Rheumatoid factors up to 1,500 IU/mL will not interfere with the CLIA series CA50 assay. Criterion: Recovery within ± 15% of initial value.

In vitro tests were performed on 27 commonly used pharmaceuticals. No interference was observed from these substances at the levels indicated below. Criterion: Recovery within ± 10% of initial value.

Substance	Concentration
Acetylcysteine	1,660 µg/mL
Ampicillin sodium	1,000 µg/mL
Ascorbic acid	1,000 µg/ml
Cyclosporine	5 µg/mL
Cefoxitin	2,500 µg/mL
Heparin	5,000 IU/L
Levodopa	20 µg/mL
Methyldopa	20 µg/mL
Metronidazole	200 µg/mL
Phenylbutazone	400 µg/mL
Doxycycline	50 µg/mL
Acetylsalicylic acid	1,000 µg/mL
Rifampicin	60 µg/mL

Substance	Concentration
Acetaminophen	200 µg/mL
Ibuprofen	500 µg/mL
Theophylline	100 µg/mL
Doxorubicin	75 mg/L
Cyclophosphamide	1,000 mg/L
Cisplatin	225 mg/L
5-fluorouracil	500 mg/L
Methotrexate	1,000 mg/L
Tamoxifen	50 mg/L
Mitomycin	25 mg/L
Carboplatin	1,000 mg/L
Etoposide	400 mg/L
Vincristine	70 mg/dl
Taxol	5.5 mg/L

No cross reactivity was observed when Wiener lab. CA50 Calibrator C0 was spiked with other tumor markers at specific levels indicated in the table below. Criterion: Reported CA50 ≤ 10.0 U/mL. The results are summarized in the table below*.

Tumor Marker	Cross-reactant Concentration	Reported CA50 (U/mL)
CA125	1,000 U/mL	0.17
CEA	1,000 ng/mL	0.36
CA72-4	300 U/mL	0.73

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High Dose Hook

For the CLIA series CA50 assay, no high dose hook effect was observed when samples contain up to approximately 100,000 U/mL of CA50 were assayed.

Trueness

Spiking recovery test was performed to verify the trueness of this assay. Two different samples with high levels of CA50 were spiked into a base sample that has low level of CA50 respectively, at a volume ratio of 1:9. The results showed that this assay had a mean recovery of 100 ± 15%. The results are listed in the following table*.

Undiluted CA50 (U/mL)	Spiked CA50(U/mL)	Mean Recovery
Low level: 7.08	18.25	107.79%
High level: 110.21		

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Precision

Precision was determined by following Clinical and Laboratory Standards Institute (CLSI) Protocol EP5-A2.¹² Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below*.

Sample	Mean CA50 (U/mL)	Within-run CV	Between-run CV	Within-Device CV
Level 1	19.89	1.44%	1.26%	2.54%
Level 2	161.36	1.73%	5.66%	5.92%

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Linearity

A high concentration CA50 sample (approximately 500 U/mL) was mixed with a low concentration sample (< 1 U/mL) at different ratios, generating a series of dilutions. The CA50 of each dilution was determined using the Wiener lab. CLIA series CA50 Assay. Linearity was demonstrated in the range of 1 U/mL to 500 U/mL, the correlation coefficient r ≥ 0.9900. The linearity data are summarized in the table below*.

Sample	Expected CA50(U/mL)	Measured CA50 (U/mL)
1	0.00	0.00
2	92.04	96.75
3	184.26	204.77
4	276.40	309.11
5	368.16	389.00
6	460.20	474.72
7	552.24	571.99
8	644.28	644.28

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Method Comparison

The Wiener lab. CLIA series CA50 Assay was compared to a commercially available diagnostic kit in a correlation study with about 1161 specimens. The statistical data obtained by linear regression are shown in the table below.

Concentration Range (U/mL)	Slope	Intercept	Correlation Coefficient
2.5-500	1.0010	-0.0516	0.9978

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- The concentration of CA50 in a given specimen determined with different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the CA50 assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining CA50 values is changed, additional sequential testing should be carried out to confirm baseline values.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 56 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimens and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages.
- If the reagents are opened unintentionally before use, they shall be used as soon as possible.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
- Do not freeze. The results can't be assured when the reagents are stored at inappropriate condition.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

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Graphical Symbols

References

- Bunworasate U, Voravud N. CA 50: a tumor marker for gastrointestinal malignancies. J. Med Assoc Thai. 1995; 78: 255-270.
- Mansson JE, Fredman P, Nilsson O, et al. Chemical structure of carcinoma ganglioside antigens defined by monoclonal antibody C-50 and some allied gangliosides of human pancreatic adenocarcinoma. Biochim Biophys Acta. 1985; 834: 110-117.
- Nilsson O, Mansson JE, Lindholm L, et al. Sialosylactotetraosylceramide, a novel ganglioside antigen detected in human carcinomas by a monoclonal antibody. FEBS Lett. 1985; 182: 398-402.
- Kuusela P, Haglund C, Roberts PJ, et al. Comparison of CA50, a new tumor marker, with carcinoembryonic antigen (CEA) and alpha-fetoprotein (AFP) in patients with gastrointestinal diseases. Br J Cancer. 1987; 55: 673-676.
- Genolla J, Moragas G, Mane S, et al. CA50 serum levels in patients with liver diseases. Int J Biol Markers. 1987; 2:207-208.
- Paganuzzi M, Onetto M, Marroni P, et al. CA19-9 and CA50 in benign and malignant pancreatic and biliary diseases. Cancer. 1988; 61: 2100-2108.
- Plsson B, Masson P, Andrén-Sandberg A. Tumour marker CA 50 levels compared to signs and symptoms in the diagnosis of pancreatic cancer. EJSO. 1997; 23: 151-156.
- Lucarotti ME, Rothnie ND, Kelly SB, et al. Clinical evaluation of tumour marker combinations in the differential diagnosis of benign and malignant liver disease. HPB Surgery. 1991; 5: 23-28.

- Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline -Second Edition.

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- Kuusela P, Haglund C, Roberts PJ, et al. Comparison of CA50, a new tumor marker, with carcinoembryonic antigen (CEA) and alpha-fetoprotein (AFP) in patients with gastrointestinal diseases. Br J Cancer. 1987; 55: 673-676.
- Genolla J, Moragas G, Mane S, et al. CA50 serum levels in patients with liver diseases. Int J Biol Markers. 1987; 2:207-208.
- Paganuzzi M, Onetto M, Marroni P, et al. CA19-9 and CA50 in benign and malignant pancreatic and biliary diseases. Cancer. 1988; 61: 2100-2108.
- Plsson B, Masson P, Andrén-Sandberg A. Tumour marker CA 50 levels compared to signs and symptoms in the diagnosis of pancreatic cancer. EJSO. 1997; 23: 151-156.
- Lucarotti ME, Rothnie ND, Kelly SB, et al. Clinical evaluation of tumour marker combinations in the differential diagnosis of benign and malignant liver disease. HPB Surgery. 1991; 5: 23-28.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Endereço: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Wiener Laboratórios S.A.I.C.
Riobamba 2944
2000 Rosario - Argentina
<http://www.wiener-lab.com>

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Contenido

Referencias

- Bunworasate U, Voravud N. CA 50: a tumor marker for gastrointestinal malignancies. J. Med Assoc Thai. 1995; 78: 255-270.
- Mansson JE, Fredman P, Nilsson O, et al. Chemical structure of carcinoembryonic antigens defined by monoclonal antibody C-50 and some allied gangliosides of human pancreatic adenocarcinoma. Biochim Biophys Acta. 1985; 834: 110-117.
- Nilsson O, Mansson JE, Lindholm L, et al. Sialosylactotetraosylceramide, a novel ganglioside antigen detected in human carcinomas by a monoclonal antibody. FEBS Lett. 1985; 182: 398-402.
- Kuusela P, Haglund C, Roberts PJ, et al. Comparison of CA50, a new tumor marker, with carcinoembryonic antigen (CEA) and alpha-fetoprotein (AFP) in patients with gastrointestinal diseases. Br J Cancer. 1987; 55: 673-676.
- Genolla J, Moragas G, Mane S, et al. CA50 serum levels in patients with liver diseases. Int J Biol Markers. 1987; 2:207-208.
- Paganuzzi M, Onetto M, Marroni P, et al. CA19-9 and CA50 in benign and malignant pancreatic and biliary diseases. Cancer. 1988; 61: 2100-2108.
- Plisson B, Masson P, Andrén-Sandberg A. Tumour marker CA 50 levels compared to signs and symptoms in the diagnosis of pancreatic cancer. EJSO. 1997; 23: 151-156.
- Lucarotti ME, Rothnie ND, Kelly SB, et al. Clinical evaluation of tumour marker combinations in the differential diagnosis of benign and malignant liver disease. HPB Surgery. 1991; 5: 23-28.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline -Second Edition.

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047-046068-00 (1.0)

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.

Mindray Building, Keji 12th Road South, High-Tech Industrial Park,
Nanshan, Shenzhen, 518057, P.R. China

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
CA50 Calibrators

Order Information

Catalog No.	Package Size
1001449	C0:1×1.2 mL/vial C1:1×1.0 mL/vial C2:1×1.0 mL/vial

Intended Use

Wiener lab. CA50 Calibrators (CA50 CAL) are intended to calibrate the quantitative Carbohydrate Antigen 50 (CA50) assay on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Summary

Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, Wiener lab. reagent kits, calibrators, and controls. CA50 Calibrators contain Carbohydrate Antigen 50 (CA50) analyte. The mathematical relationship between the measured responses and the known analyte concentration establishes the calibration curve. It is used to convert Relative Light Units (RLUs) measured from samples to quantitative analyte concentration.

Components

C0, C1, C2	Carbohydrate antigen 50 in buffered bovine serum albumin (BSA) at three concentration levels with 0.05% ProClin 300 and 0.02% 5-Bromo-5-Nitro-1,3-Dioxane (BND) as preservatives.
------------	---

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Laboratory professionals must be trained to take all necessary precautions while handling laboratory reagents. Do not use reagents beyond the expiration date.
- The calibrators' values are lot-specific with the matched set listed in the Calibration Card.
- Calibration is recommended after each reagent lot change, special maintenance, or troubleshooting procedure.
- Do not mix calibrators from different lots.
- Please perform maintenance and standard operation as scheduled, including calibration and analysis to assure the performance of the measurement system.
- Recalibrate the system when controls are out of specified range.
- The calibrators were tested with CE-approved kits for HIV-1/2 Antibody, HCV Antibody, Treponema pallidum Antibody and HBsAg. All results are negative. However, as no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazards.¹

- If there is any evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- Material safety data sheet is available upon request.
- Disposal of all waste material should be in accordance with local guidelines.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
- If the reagents are opened unintentionally before use, they shall be tightly capped and stored at 2-8°C, and then used within in-use stability time.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

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Prevention:

P261 Avoid breathing dust/fume/ gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P273 Avoid release to the environment.

Response:

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
P302 + P352 IF ON SKIN: Wash with plenty of water.
P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/ container in accordance with local regulation.

Testing Method

- Refer to Wiener lab. System Operator's Manual and / or customer service to acquire information of calibration principle, use of calibrators, and validity of data of calibrators.
- Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation.
- Three levels of CA50 calibrators are provided, which are 0, ~20 and ~200 U/mL. Test C0 three times and C1, C2 twice for calibration.

- Refer to calibration card for accurate information about lot-specific CA50 concentration.

Preparation

- Take the vial out of refrigerator and let it stand for no less than 30 minutes at room temperature.
- According to the calibration procedure described in the instrument manual, use the Calibrators' lot-specific value to set calibration parameters.

Storage and Stability

The calibrators are stable up to the expiration date as indicated on the label when stored in unopened vial at 2-8°C. Once opened, it is stable for 30 days at 2-8°C or 90 days at -20°C.

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- Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, reagent kits and controls.
- General laboratory equipment.

Traceability of Calibrator

The traceability process is based on EN ISO 17511:2021,² the measurand (analyte) in the calibrators is traceable to a commercial CA50 test (CLIA). The concentrations of calibrators are lot-specific on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, and cannot be used on the other system or method, as bias may exist among different systems and methods. Details about traceability are available upon request.

Quality Control

Use Wiener lab. controls as the samples to monitor the performance of system. The control result should be within the defined range as shown in the control value sheet. If the controls are out of specified range, the measurement system should be checked. For instance, expiration date or storage condition of the calibrators, reagent and control; contamination of reagent; position of reagent or sample placed on the analyzer; the parameter setting in software; or performance of the analyzer. The assigned values of controls are fitted to Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer and are lot-specific.

Graphical Symbols

In vitro diagnostic medical device

Authorized representative in the European Community

Consult instructions for use

European Conformity

Catalogue number

Temperature limit

Manufacturer

Use-by date

Biological risks

Batch code

Unique device identifier

Caution

Contents

References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
- EN ISO 17511:2021 In vitro diagnostic medical devices – Requirements for establishing metrological traceability of values assigned to calibrators, trueness control materials and human samples.

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Riobamba 2944, 2000 Rosario - Argentina

<http://www.wiener-lab.com>

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.
Mindray Building, Keji 12th Road South, High-Tech Industrial Park,
Nanshan, Shenzhen, 518057, P. R. China

Shanghai International Holding Corp. GmbH(Europe)
Eiffestraße 80, 20537 Hamburg, Germany

0177917780430564102023010111

UDI

(01)06903053058582 (97)105-039942-A0

(10)2023010111

(99)1001449

(17)241019

(98)PBQ202301011100047

LOT

2023010111


2024-10-19

IVD

CE 0123

Wiener lab.

047-046071-00(1.0)

CA50 Calibrators

1.2 mL

IVD

C0

CE 0123

LOT 2023010111

2024-10-19

2°C - 8°C

Wiener lab.

CA50 Calibrators

1.0 mL

IVD

C1

CE 0123

LOT 2023010111

2024-10-19

2°C - 8°C

Wiener lab.

CA50 Calibrators

1.0 mL

IVD

C2

CE 0123

LOT 2023010111

2024-10-19

2°C - 8°C

Wiener lab.

CA242

Carbohydrate Antigen 242 (CLIA)

Order Information

Catalog No.	Package Size
1001444	2x50 tests
1001445	2x100 tests

Intended Use

The CLIA series CA242 assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of carbohydrate antigen 242 (CA242) in human serum or plasma. This assay is to be used as an aid in monitoring disease of cancer patients.

Summary

Carbohydrate antigen 242 (CA242) is a tumor related antigen discovered by Sweden scientist Lindholm in 1985.¹ In serum, CA242 is found on the same mucin-complex as CA19-9; however, their expression levels are different in benign and malignant diseases or different tumorigenesis organs. The chemical structure of the antigen is not exactly known, but the determinant has been shown to be sialylated carboxydrate structure, similar but different from sialyl-lewsia.^{2,3} The elevated levels of CA242 are commonly found in pancreas cancer and gastrointestinal cancer.⁴

Assay Principle

The CLIA series CA242 assay is a two-site sandwich assay.

In the first step, sample, paramagnetic microparticle coated with monoclonal anti-CA242 antibody (mouse) and monoclonal anti-CA242 antibody (mouse)-alkaline phosphatase conjugate are added into a reaction cuvette. After incubation, CA242 present in the sample binds to both anti-CA242 antibody coated microparticle and anti-CA242 antibody alkaline phosphatase-labeled conjugate to form a sandwich complex. Microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, the substrate solution is added to the reaction cuvette. It is catalyzed by anti-CA242 antibody (mouse)-alkaline phosphatase conjugate in the immunocomplex retained on the microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of CA242 present in the sample is proportional to the light units (RLUs) generated during the reaction. The CA242 concentration can be determined via a calibration curve.

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Ra	Paramagnetic microparticles coated with monoclonal anti-CA242 antibody (mouse). Minimum concentration: 0.5 g/L solids. TRIS ^{a)} buffer: 50 mmol/L. Preservatives: 0.02% BND and 0.05% ProClin 300.
Rb	Monoclonal anti-CA242 antibody (mouse)-alkaline phosphatase conjugate. Minimum concentration: 0.14 mg/L. MES ^{b)} buffer: 50 mmol/L. Preservatives: 0.02% BND and 0.05% ProClin 300.

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b) MES=2-(N-morpholino) ethanesulfonic acid

Storage and Stability

The unopened CA242 (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The CA242 (CLIA) reagent kit can be stored onboard and used

for a maximum of 56 days after opening at 2-8°C.

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Materials Required but Not Provided

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Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.
- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample is covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

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Assay Procedure

For optimal performance of this assay, operators should read the related system Operator's Manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the CA242 (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Wiener lab. Customer Service for help. Do not invert opened reagent bottle.

This assay requires 20 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system Operator's Manual and specific requirement of the assay to determine the minimum sample volume.

Calibration

CLIA series CA242 (CLIA) has been standardized against a commercial CA242 test (CLIA).

The specific information of master calibration curve of CA242 (CLIA) reagent kit is stored in the two-dimensional barcode on the reagent pack, which is used in combination with CA242 calibrators for calibration of the specific reagent lot. Before initiating calibration on each new lot of reagent, load the assay master curve by scanning the two-dimensional barcode on the reagent pack. When performing the calibration, scan the two-dimensional barcode on the calibrator pack, and then test the CA242 calibrators at three levels. Valid calibration curve is required before any CA242 test. Recalibration is recommended every 28 days, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system Operator's Manual.

Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Wiener lab. Tumor Marker Multi Control II (L) and Wiener lab. Tumor Marker Multi Control II (H). In addition, other suitable control materials can be used.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system Operator's Manual. If the quality control results are still out of the specified range, please contact Wiener lab. Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three level calibrators of defined concentration values. The results are shown in the unit of U/mL.

Dilution

Samples with CA242 concentrations above the upper limit can be diluted with Wiener lab. Sample Diluent. The recommended dilution is 1:10 (either automatically by the analyzer or manually). The maximum manual dilution is 1:40. The concentration of the diluted sample must be > 5 U/mL. After

manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzer, the system automatically multiplies the result by the dilution factor when calculating the sample concentration.

Expected values

An extensive study on a cohort of 266 males and 239 females has determined the reference range of CLIA series CA242 assay.

Category	N	95 th percentile	97.5 th percentile	99 th percentile
Healthy	505	10 U/mL	20 U/mL	30 U/mL

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Limitation

The upper limit of this assay is 200 U/mL. A specimen with a CA242 concentration lower than the upper limit can be quantitatively determined, while the specimen with a concentration higher than the upper limit will be reported as > 200 U/mL.

The concentration of CA242 in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

The heterophilic antibodies in serum or plasma may react with the immunoglobulin in reagent components to interfere the test.⁵ Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies.^{6,7} However, no obvious interference of HAMA was observed in this assay.

Performance Characteristics

Lower limits of measurement

The CA242(CLIA) reagent kit has a Limit of Blank (LoB) of ≤ 0.3 U/mL, a Limit of Detection (LoD) of ≤ 0.5 U/mL and a Limit of Quantitation (LoQ) of ≤ 0.5 U/mL (with a total allowable relative error of ≤ 25%).

Measuring range

0.5-200 U/mL (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 0.50 U/mL. Values above the measuring range are reported as > 200 U/mL (or up to 8,000 U/mL for 40-fold diluted samples).

Analytical Specificity

Hemoglobin up to 3,500 mg/dL, bilirubin up to 70 mg/dL, triglycerides up to 3,500 mg/dL, total protein up to 15 g/dL, biotin up to 5,000 ng/mL and antinuclear antibody will not interfere with the CLIA series CA242 assay. Criterion: Recovery within ± 10% of initial value.

Rheumatoid factors up to 1,500 IU/mL will not interfere with the CLIA series CA242 assay.

Criterion: Recovery within ± 15% of initial value.

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Substance	Concentration
Acetylcysteine	1,660 µg/mL
Ampicillin sodium	1,000 µg/mL
Ascorbic acid	1,000 µg/ml
Cyclosporine	5 µg/mL

Substance	Concentration
Cefoxitin	2,500 µg/mL
Heparin	5,000 IU/L
Levodopa	20 µg/mL
Methyldopa	20 µg/mL
Metronidazole	200 µg/mL
Phenylbutazone	400 µg/mL
Doxycycline	50 µg/mL
Acetylsalicylic acid	1,000 µg/mL
Rifampicin	60 µg/mL
Acetaminophen	200 µg/mL
Ibuprofen	500 µg/mL
Theophylline	100 µg/mL
Doxorubicin	75 mg/L
Cyclophosphamide	1,000 mg/L
Cisplatin	225 mg/L
5-fluorouracil	500 mg/L
Methotrexate	1,000 mg/L
Tamoxifen	50 mg/L
Mitomycin	25 mg/L
Carboplatin	1,000 mg/L
Etoposide	400 mg/L
Vincristine	70 mg/dl
Taxol	5.5 mg/L

No cross reactivity was observed when Wiener lab. CA242 Calibrator C0 was spiked with other tumor markers at specific levels indicated in the table below. Criterion: Reported CA242 ≤ 2.0 U/mL. The results are summarized in the table below*.

Tumor Marker	Cross-reactant Concentration	Reported CA242 (U/mL)
CA125	1,000 U/mL	0.09
CEA	1,000 ng/mL	0.23
CA72-4	300 U/mL	1.00

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High Dose Hook

For the CLIA series CA242 assay, no high dose hook effect was observed when samples containing up to approximately 200,000 U/mL of CA242 were assayed.

Trueness

Spiking recovery test was performed to verify the trueness of this assay. Two different samples with high levels of CA242 were spiked into a base sample that has low level of CA242 respectively, at a volume ratio of 1:9. The results showed that this assay had a mean recovery of 100 ± 15%. The results are listed in the following table*.

Undiluted CA242 (U/mL)	Spiked CA242 (U/mL)	Mean Recovery
Low level: 3.98	20.04	100.05%
High level: 164.51		

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Precision

Precision was determined by following Clinical and Laboratory

Standards Institute (CLSI) Protocol EP5-A2.⁸ Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below*.

Sample	Mean CA242 (U/mL)	Within-run CV	Between-run CV	Within-Device CV
Level 1	7.67	1.25%	2.52%	3.59%
Level 2	53.97	1.46%	1.87%	2.65%

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Linearity

A high concentration CA242 sample (approximately 200 U/mL) was mixed with a low concentration sample (< 0.5 U/mL) at different ratios, generating a series of dilutions. The CA242 of each dilution was determined using the Wiener lab. CLIA series CA242 Assay. Linearity was demonstrated in the range of 0.5 U/mL to 200 U/mL, the correlation coefficient $r \geq 0.9900$. The linearity data are summarized in the table below*.

Sample	Expected CA242 (U/mL)	Measured CA242 (U/mL)
1	0.00	0.00
2	33.11	34.40
3	66.30	65.61
4	99.44	98.03
5	132.46	128.14
6	165.57	162.00
7	198.69	197.09
8	231.80	231.80

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Method Comparison

The Wiener lab. CLIA series CA242 Assay was compared to a commercially available diagnostic kit in a correlation study with about 955 specimens. The statistical data obtained by linear regression are shown in the table below.

Concentration Range (U/mL)	Slope	Intercept	Correlation Coefficient
0.5-200	0.9992	0.02364	0.998

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- The concentration of CA242 in a given specimen determined with different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the CA242 assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining CA242 values is changed, additional sequential testing should be carried out to confirm baseline values.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 56 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be

considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.

- The Material Safety Data Sheet (MSDS) is available upon request.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages.
- If the reagents are opened unintentionally before use, they shall be used as soon as possible.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
- Do not freeze. The results can't be assured when the reagents are stored at inappropriate condition.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

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Graphical Symbols

References

- Nilsson O, Jansson EL, Johansson C, et al. CA242, a novel tumour-associated carbohydrate antigen with increased tumour specificity and sensitivity. J Tumour Marker Oncol 1988; 3: 314.
- Haglund C, Lindgren J, Roberts PJ, Tissue expression of the tumour associated antigen CA242 in benign and

- malignant pancreatic lesions. A comparison with CA 50 and CA 19-9. Br. J. Cancer 1989; 60: 845-851.
- Rye PD, Bovin NV, Vlasova EV, et al. Summary report on the ISOBM TD-6 Workshop: Analysis of 20 monoclonal antibodies against sialyl Lewis x and related antigens. Tumour Biol 1998; 19: 390-420.
- Kuusela P, Haglund C, Robert PJ, Comparison of a new tumour marker CA242 with CA19-9, CA50 and carcinoembryonic antigen (CEA) in digestive tract diseases. Br J Cancer, 1991, 63:636-640.
- Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

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Conteúdo

Referências

1. Nilsson O, Jansson EL, Johansson C, et al. CA242, a novel tumour-associated carbohydrate antigen with increased tumour specificity and sensitivity. J Tumour Marker Oncol 1988; 3: 314.
2. Haglund C, Lindgren J, Roberts PJ, Tissue expression of the tumour associated antigen CA242 in benign and malignant pancreatic lesions. A comparison with CA 50 and CA 19-9. Br. J. Cancer 1989; 60: 845-851.
3. Rye PD, Bovin NV, Vlasova EV, et al. Summary report on the ISOBM TD-6 Workshop: Analysis of 20 monoclonal antibodies against sialyl Lewisia and related antigens. Tumour Biol 1998; 19: 390-420.
4. Kuusela P, Haglund C, Robert P.J. Comparison of a new tumour marker CA242 with CA19-9, CA50 and carcinoembryonic antigen (CEA) in digestive tract diseases. Br J Cancer, 1991; 63:636-640.
5. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
6. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
7. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
8. CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Endereço: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Wiener Laboratórios S.A.I.C.
Riobamba 2944
2000 Rosario - Argentina
<http://www.wiener-lab.com>

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Símbolos gráficos

Referencias

- Nilsson O, Jansson EL, Johansson C, et al. CA242, a novel tumour-associated carbohydrate antigen with increased tumour specificity and sensitivity. J Tumour Marker Oncol 1988; 3: 314.
- Haglund C, Lindgren J, Roberts PJ, Tissue expression of the tumour associated antigen CA242 in benign and malignant pancreatic lesions. A comparison with CA 50 and CA 19-9. Br. J. Cancer 1989; 60: 845-851.
- Rye PD, Bovin NV, Vlasova EV, et al. Summary report on the ISOBM TD-6 Workshop: Analysis of 20 monoclonal antibodies against sialyl Lewis x and related antigens. Tumour Biol 1998; 19: 390-420.
- Kuusela P, Haglund C, Robert PJ, Comparison of a new tumour marker CA242 with CA19-9, CA50 and carcinoembryonic antigen (CEA) in digestive tract diseases. Br J Cancer, 1991, 63:636-640.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

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CA242 Calibrators

Order Information

Catalog No.	Package Size
1001448	C0:1×1.2 mL/vial C1:1×1.0 mL/vial C2:1×1.0 mL/vial

Intended Use

Wiener lab. CA242 Calibrators (CA242 CAL) are intended to calibrate the quantitative Carbohydrate Antigen 242 (CA242) assay on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Summary

Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, Wiener lab. reagent kits, calibrators, and controls. CA242 Calibrators contain Carbohydrate Antigen 242 (CA242) analyte. The mathematical relationship between the measured responses and the known analyte concentration establishes the calibration curve. It is used to convert Relative Light Units (RLUs) measured from samples to quantitative analyte concentration.

Components

C0, C1, C2	Carbohydrate antigen 242 in buffered bovine serum albumin (BSA) at three concentration levels with 0.05% ProClin 300 and 0.02% 5-Bromo-5-Nitro-1,3-Dioxane (BND) as preservatives.
------------	--

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Testing Method

- Refer to Wiener lab. System Operator's Manual and / or customer service to acquire information of calibration principle, use of calibrators, and validity of data of calibrators.
- Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation.
- Three levels of CA242 calibrators are provided, which are 0, ~10 and ~70 U/mL. Test C0 three times and C1, C2 twice for calibration.
- Refer to calibration card for accurate

information about lot-specific CA242 concentration.

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Traceability of Calibrator

The traceability process is based on EN ISO 17511:2021,² the measurand (analyte) in the calibrators is traceable to a commercial CA242 test (CLIA). The concentrations of calibrators are lot-specific on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, and cannot be used on the other system or method, as bias may exist among different systems and methods. Details about traceability are available upon request.

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
Intended Use

The CLIA series HE4 assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of human epididymal protein 4 (HE4) in human serum or plasma. This assay can be used as an aid in monitoring disease of patients with ovarian cancer.

Summary

The human epididymal protein 4 (HE4) encoded by WFDC2 is an acidic cysteine rich secretion protein, belongs to the family of whey acidic four-disulfide core (WFDC) proteins. The WFDC2 codes for a 13 kD protein, however the mature glycosylated form protein has a molecular weight of approximately 25 kD.¹⁻⁴ HE4 was first detected in the epithelium of the distal epididymal and vas deferens, supposed to play an important role in the process of sperm maturation. Deep research suggests HE4 shows high expression in epithelial of female reproductive system and proximal trachea, parotid gland, proximal and distal convoluted tubules, anterior pituitary, mammary duct as well as colonic epithelium. It is also found in thyroid, prostate, tear gland, exocrine gland and pancreas. However, it doesn't express in ovarian tissue, gastrointestinal tissue, liver, spleen, lymph gland, musculature, brain tissue, lymph system, hematopoietic system and nervous system.⁵⁻⁷ High secreted levels can be found in ovarian cancer and endometrial carcinoma.^{4,8}

Assay Principle

The CLIA series HE4 assay is a two-site sandwich assay.

In the first step, sample, paramagnetic microparticle coated with monoclonal anti-HE4 antibody (mouse) and monoclonal anti-HE4 antibody (mouse)-alkaline phosphatase conjugate are added into a reaction cuvette. After incubation, HE4 present in the sample binds to both anti-HE4 antibody coated microparticle and anti-HE4 antibody alkaline phosphatase-labeled conjugate to form a sandwich complex. Microparticle is magnetically captured while other unbound substances are removed by washing.

In the second step, the substrate solution is added to the reaction cuvette. It is catalyzed by anti-HE4 antibody (mouse)-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of HE4 present in the sample is proportional to the light units (RLUs) generated during the reaction. The HE4 concentration can be determined via a calibration curve.

Reagent Components

Ra	Paramagnetic microparticles coated with monoclonal anti-HE4 antibody (mouse). Minimum concentration: 0.5 g/L solids. TRIS ^{a)} buffer: 50 mmol/L. Preservatives: 0.02% BND and 0.05% ProClin 300.
Rb	Monoclonal anti-HE4 antibody (mouse)-alkaline phosphatase conjugate. Minimum concentration: 0.08 mg/L. MES ^{b)} buffer: 50 mmol/L. Preservatives: 0.02% BND and 0.05% ProClin 300.

a) TRIS=Tris (hydroxymethyl)-aminomethane

b) MES=2-(N-morpholino) ethanesulfonic acid

Storage and Stability

The unopened HE4 (CLIA) reagent kit is stable up to the stated

expiration date when stored at 2-8°C.

The HE4 (CLIA) reagent kit can be stored onboard and used for a maximum of 56 days after opening at 2-8°C.

Reagent Preparation

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

Materials Required but not Provided

Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Cat.No. 1001452: HE4 Calibrators, 1x1.2 mL for calibrator C0, 1x1.0 mL for each of calibrator C1 and C2.

Cat.No. 1001453: Tumor Marker Multi Control II (L), 3x2.0 mL.

Cat.No. 1001455: Tumor Marker Multi Control II (H), 3x2.0 mL.

Cat.No. 1001454: Tumor Marker Multi Control II (L), 6x2.0 mL.

Cat.No. 1001456: Tumor Marker Multi Control II (H), 6x2.0 mL.

Cat.No.1001185: Wash Buffer, 1x10 L.

Cat.No.1001139: Substrate Solution, 4x75 mL.

Cat.No.1001138: Substrate Solution, 4x115mL.

Reaction Cuvette.

Applicable Instrument

Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer

Specimen Collection and Preparation

Specimen Types

- Human serum or plasma collected in k₂EDTA, k₃EDTA, sodium heparin and lithium heparin are recommended for this assay.
- Blood collection tubes from various manufacturers may contain different raw materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Wiener lab. Each laboratory should determine the acceptability of the tubes and serum/plasma separation products.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.
- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

- Specimens should be tested as soon as possible after sample collection. If testing is not completed within 24 hours, specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 7 days, specimens should be frozen at -20°C or colder. The specimens can be stored at -20°C for up to 90 days.
- The number of repeated freezing and thawing shall not exceed 3 times.

Assay Procedure

For optimal performance of this assay, operators should read the related system operator's manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the HE4 (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Wiener lab. Customer Service for help. Do not invert opened reagent bottle.

This assay requires 5 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operator's manual and specific requirement of the assay to determine the minimum sample volume.

Calibration

CLIA series HE4 (CLIA) has been standardized against a commercial HE4 test (CLIA).

The specific information of master calibration curve of HE4 (CLIA) reagent kit is stored in the two-dimensional barcode on the reagent pack, which is used in combination with HE4 calibrators for calibration of the specific reagent lot. Before initiating calibration on each new lot of reagent, load the assay master curve by scanning the two-dimensional barcode on the reagent pack. When performing the calibration, scan the two-dimensional barcode on the calibrator pack, and then test the HE4 calibrators at three levels. Valid calibration curve is required before any HE4 test. Recalibration is recommended every 28 days, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operator's manual.

Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Wiener lab. Tumor Marker Multi Control II (L) and Wiener lab. Tumor Marker Multi Control II (H). In addition, other suitable control material can be used.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operator's manual. If the quality control results are still out of the specified range, please contact Wiener lab. Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three level calibrators of defined concentration values. The results are shown in the unit of pmol/L.

Dilution

Samples with HE4 concentrations above the upper limit can be diluted with Wiener lab. Sample Diluent. The recommended dilution is 1:20 (either automatically by the analyzer or manually). The concentration of the diluted sample must be > 30 pmol/L. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzers, the

system automatically multiplies the result by the dilution factor when calculating the sample concentration.

Expected values

An extensive study on a cohort of 696 healthy females has determined the reference range of CLIA series HE4 assay.

Age (years)	N	HE4 (pmol/L)	
		Median	95 th percentile
<40	205	42.11	60.53
40-49	123	44.47	76.33
50-59	125	47.65	74.31
60-69	122	56.12	83.26
>70	121	61.52	105.69

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

Limitation

The upper limit of this assay is 1,500 pmol/L. A specimen with a HE4 concentration lower than the upper limit can be quantitatively determined, while the specimen with a concentration higher than the upper limit will be reported as > 1,500 pmol/L.

The concentration of HE4 in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

The heterophilic antibodies in serum or plasma may react with the immunoglobulin in reagent components to interfere the test.⁹ Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies.^{10,11} However, no obvious interference of HAMA was observed in this assay.

Performance Characteristics

Measuring range

15-1,500 pmol/L (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 15.00 pmol/L. Values above the measuring range are reported as > 1,500 pmol/L (or up to 30,000 pmol/L for 20-fold diluted samples).

Lower limits of measurement

The HE4 (CLIA) reagent kit has a Limit of Blank (LoB) of ≤ 5 pmol/L, a Limit of Detection (LoD) of ≤ 15 pmol/L and a Limit of Quantitation (LoQ) of ≤ 20 pmol/L (with a total allowable relative error of ≤ 30%).

Analytical Specificity

Hemoglobin up to 1,000 mg/dL, bilirubin up to 70 mg/dL, triglycerides up to 3,000 mg/dL, total protein up to 12 g/dL and antinuclear antibody will not interfere with the CLIA series HE4 assay. Criterion: Recovery within ± 10% of initial value.

Rheumatoid factors up to 1,500 IU/mL will not interfere with the CLIA series HE4 assay. Criterion: Recovery within ± 15% of initial value.

In vitro tests were performed on 10 commonly used pharmaceuticals. No interference was observed from these substances at the levels indicated below. Criterion: Recovery within ± 10% of initial value.

Substance	Concentration
Carboplatin	500 µg/mL
Cisplatin	165 µg/mL
Clotrimazole	0.3 µg/mL
Cyclophosphamide	500 µg/mL

Substance	Concentration
Doxorubicin	1.16 µg/mL
Methotrexate	45 µg/mL
Paclitaxel	0.0035 µg/mL
Dexamethasone	10 µg/mL
Folinate	2.68 µg/mL
Melphalan	2.8 µg/mL

No cross reactivity was observed when Wiener lab. HE4 Calibrator C0 was spiked with other tumor markers at specific levels indicated in the table below. Criterion: Reported HE4 ≤ 50.0 pmol/L. The results are summarized in the table below*.

Tumor Marker	Cross-reactant Concentration	Reported HE4 (pmol/L)
AFP	1,000 ng/mL	8.87
CA125	1,000 U/mL	2.39
CA15-3	100 U/mL	40.34
CA19-9	1,000 U/mL	0.24
CEA	1,000 ng/mL	5.10
CA72-4	300 U/mL	18.50
NSE	350 ng/mL	2.21
CYFRA 21-1	500 ng/mL	11.67
FERR	1,000 ng/mL	11.67
CA50	500 U/mL	29.84
CA242	200 U/mL	4.88
ProGRP	5,000 pg/mL	2.93
SCCA	70 ng/mL	2.02

*Representative data; results in individual laboratories may vary.

High Dose Hook

For the CLIA series HE4 assay, no high dose hook effect was observed when samples containing up to approximately 40,000 pmol/L of HE4 were assayed.

Accuracy

Spiking recovery test was performed to verify the accuracy of this assay. Two different samples with high levels of HE4 were spiked into a base sample that has low level of HE4 respectively, at a volume ratio of 1:9. The results showed that this assay had a mean recovery of 100 ± 15%. The results are listed in the following table*.

Undiluted HE4 (pmol/L)	Spiked HE4 (pmol/L)	Mean Recovery
Low level: 38.52	145.51	101.86%
High level: 1,088.18		

*Representative data; results in individual laboratories may vary.

Precision

Precision was determined by following Clinical and Laboratory Standards Institute (CLSI) EP5-A2.12 Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below*.

Sample	Mean HE4 (pmol/L)	Within-run CV	Between-run CV	Within-Device CV
Level 1	79.54	0.91%	0.96%	2.26%
Level 2	367.65	0.87%	0.90%	2.12%

*Representative data; results in individual laboratories may vary.

Linearity

A high concentration HE4 sample (approximately 1,500 pmol/L) was mixed with a low concentration sample (< 15 pmol/L) at different ratios, generating a series of dilutions. The HE4 of each dilution was determined using the Wiener lab. CLIA series HE4 Assay. Linearity was demonstrated in the range of 15 pmol/L to 1,500 pmol/L, the correlation coefficient r ≥ 0.9900.

The linearity data are summarized in the table below*.

Sample	Expected HE4 (pmol/L)	Measured HE4 (pmol/L)
1	0.00	0.00
2	281.89	304.05
3	563.77	611.03
4	845.66	914.40
5	1127.54	1223.72
6	1409.43	1504.04
7	1691.31	1691.31

*Representative data; results in individual laboratories may vary.

Method Comparison

The Wiener lab. CLIA series HE4 Assay was compared to a commercially available diagnostic kit in a correlation study with about 1183 specimens. The statistical data obtained by linear regression are shown in the table below.

Concentration Range (pmol/L)	Slope	Intercept	Correlation Coefficient
15-1,500	1.0007	1.3750	0.9976

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- The concentration of HE4 in a given specimen determined with different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the HE4 assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining HE4 values is changed, additional sequential testing should be carried out to confirm baseline values.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 56 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages.
- If the reagents are opened unintentionally before use, they shall be used as soon as possible.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
- Do not freeze. The results can't be assured when the reagents are stored at inappropriate condition.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Warning

H317 May cause an allergic skin reaction.
H412 Harmful to aquatic life with long lasting effects.

Prevention:

P261 Avoid breathing dust /fume/gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P273 Avoid release to the environment.


Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/container in accordance with local regulation.

Graphical Symbols

References

- Anstasi, E., et al., HE4 in the differential diagnosis of a pelvi mass: a case report. International journal of molecular sciences 2011; 12:627-632.
- Anstasi, E., et al., Ovarian tumor marker HE4 is differently expressed during the phases of the menstrual cycle in healthy young women. Tumor Biology 2010; 31:411-415.
- Anstasi, E., et al., HE4: a new potential early biomarker for the recurrence of ovarian cancer. Tumor Biology 2010; 31:113-119.
- Hellström, I., et al., The HE4(WFDC2) protein is a biomarker for ovarian carcinoma. Cancer research 2003; 63:3695-3700.
- Bingle, L., V. Singleton, and C.D. Bingle, The putative ovarian tumor marker gene HE4 (WFDC2), is expressed in normal tissues and undergoes alternative splicing to yield multiple protein isoforms. Oncogene 2002; 21:2768-2773.
- Hellström, I., et al., Detection of the HE4 protein in urine as a biomarker for ovarian neoplasms. Cancer letters 2010; 296: 43-48.
- Hellström, I., and K.E. Hellstrom, Two new biomarkers, mesothelin and HE4, for diagnosis of ovarian carcinoma. Expert opinion on medical diagnostics 2011; 5:227-240.
- Moore RG, Brown AK, Miller MC, et al. Utility of a novel serum tumor biomarker HE4 in patients with endometrioid adenocarcinoma of the uterus. Gynecol Oncol 2008; 110:196-201.
- Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays- still a threat. Clin Chem 2000;46:1037-1038.
- Bjerner J, et al. Immunometric assay interference:incidence and prevention. Clin Chem 2002; 48:613-621.

- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Second Edition.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Conteúdo

Referências

- Anstasi, E., et al., HE4 in the differential diagnosis of a pelvi mass: a case report. International journal of molecular sciences 2011; 12:627-632.
- Anstasi, E., et al., Ovarian tumor marker HE4 is differently expressed during the phases of the menstrual cycle in healthy young women. Tumor Biology 2010; 31:411-415.
- Anstasi, E., et al., HE4: a new potential early biomarker for the recurrence of ovarian cancer. Tumor Biology 2010; 31:113-119.

- Hellström, I., et al., The HE4(WFDC2) protein is a biomarker for ovarian carcinoma. Cancer research 2003; 63:3695-3700.
- Bingle, L., V. Singleton, and C.D. Bingle, The putative ovarian tumor marker gene HE4 (WFDC2), is expressed in normal tissues and undergoes alternative splicing to yield multiple protein isoforms. Oncogene 2002; 21:2768-2773.
- Hellström, I., et al., Detection of the HE4 protein in urine as a biomarker for ovarian neoplasms. Cancer letters 2010; 296: 43-48.
- Hellström, I., and K.E. Hellstrom, Two new biomarkers, mesothelin and HE4, for diagnosis of ovarian carcinoma. Expert opinion on medical diagnostics 2011; 5:227-240.
- Moore RG, Brown AK, Miller MC, et al. Utility of a novel serum tumor biomarker HE4 in patients with endometrioid adenocarcinoma of the uterus. Gynecol Oncol 2008; 110:196-201.
- Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays- still a threat. Clin Chem 2000;46:1037-1038.
- Bjerner J., et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48:613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Second Edition.

Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Contenido

Referencias

- Anstasi, E., et al., HE4 in the differential diagnosis of a pelvi mass: a case report. International journal of molecular sciences 2011; 12:627-632.
- Anstasi, E., et al., Ovarian tumor marker HE4 is differently expressed during the phases of the menstrual cycle in healthy young women. Tumor Biology 2010; 31:411-415.
- Anstasi, E., et al., HE4: a new potential early biomarker for the recurrence of ovarian cancer. Tumor Biology 2010; 31:113-119.
- Hellström, I., et al., The HE4(WFDC2) protein is a biomarker for ovarian carcinoma. Cancer research 2003; 63:3695-3700.
- Bingle, L., V. Singleton, and C.D. Bingle, The putative ovarian tumor marker gene HE4 (WFDC2), is expressed in normal tissues and undergoes alternative splicing to yield multiple protein isoforms. Oncogene 2002; 21:2768-2773.
- Hellström, I., et al., Detection of the HE4 protein in urine as a biomarker for ovarian neoplasms. Cancer letters 2010; 296: 43-48.
- Hellström, I., and K.E. Hellstrom, Two new biomarkers, mesothelin and HE4, for diagnosis of ovarian carcinoma. Expert opinion on medical diagnostics 2011; 5:227-240.
- Moore RG, Brown AK, Miller MC, et al. Utility of a novel serum tumor biomarker HE4 in patients with endometrioid adenocarcinoma of the uterus. Gynecol Oncol 2008; 110:196-201.
- Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays- still a threat. Clin Chem 2000;46:1037-1038.
- Bjerner J., et al. Immunometric assay interference:incidence and prevention. Clin Chem 2002; 48:613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Second Edition.

Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Dirección: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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047-046072-00 (1.0)

SHENZHEN MINDRAY BIO-MEDICAL
ELECTRONICS CO., LTD.
Mindray Building, Keji 12th Road South, High-Tech
Industrial Park,
Nanshan, Shenzhen, 518057, P.R. China

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
HE4 Calibrators

Order Information

Catalog No.	Package Size
1001452	C0:1×1.2 mL/vial C1:1×1.0 mL/vial C2:1×1.0 mL/vial

Intended Use

Wiener lab. HE4 Calibrators (HE4 CAL) are intended to calibrate the quantitative Human Epididymal Protein 4 (HE4) assay on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Summary

Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, Wiener lab. reagent kits, calibrators, and controls. HE4 Calibrators contain Human Epididymal protein 4 (HE4) analyte. The mathematical relationship between the measured responses and the known analyte concentration establishes the calibration curve. It is used to convert Relative Light Units (RLUs) measured from samples to quantitative analyte concentration.

Components

C0, C1, C2	Human Epididymal protein 4 in buffered bovine serum albumin (BSA) at three concentration levels with 0.05% ProClin 300 and 0.02% 5-Bromo-5-Nitro-1,3-Dioxane (BND) as preservatives.
------------	--

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Laboratory professionals must be trained to take all necessary precautions while handling laboratory reagents. Do not use reagents beyond the expiration date.
- The calibrators' values are lot-specific with the matched set listed in the Calibration Card.
- Calibration is recommended after each reagent lot change, special maintenance, or troubleshooting procedure.
- Do not mix calibrators from different lots.
- Please perform maintenance and standard operation as scheduled, including calibration and analysis to assure the performance of the measurement system.
- Recalibrate the system when controls are out of specified range.
- The calibrators were tested with CE-approved kits for HIV-1/2 Antibody, HCV Antibody, Treponema pallidum Antibody and HBsAg. All results are negative. However, as no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazards.¹
- If there is any evidence of microbial contamination or excessive turbidity in the

product, discard the vial.

- Material safety data sheet is available upon request.
- Disposal of all waste material should be in accordance with local guidelines.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
- If the reagents are opened unintentionally before use, they shall be tightly capped and stored at 2-8°C, and then used within in-use stability time.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

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Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P273 Avoid release to the environment.

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Testing Method

- Refer to Wiener lab. System Operator's Manual and / or customer service to acquire information of calibration principle, use of calibrators, and validity of data of calibrators.
- Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation.
- Three levels of HE4 calibrators are provided, which are 0, ~70 and ~500 pmol/L. Test C0 three times and C1, C2 twice for calibration.
- Refer to calibration card for accurate information about lot-specific HE4 concentration.

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Traceability of Calibrator

The traceability process is based on EN ISO 17511:2021,² the measurand (analyte) in the calibrators is traceable to a commercial HE4 test (CLIA). The concentrations of calibrators are lot-specific on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, and cannot be used on the other system or method, as bias may exist among different systems and methods. Details about traceability are available upon request.

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<http://www.wiener-lab.com>

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.

Mindray Building, Keji 12th Road South, High-Tech Industrial Park,

Nanshan, Shenzhen, 518057, P. R. China

Shanghai International Holding Corp. GmbH(Europe)

Eiffestraße 80, 20537 Hamburg, Germany

0177917780430878102023050100

UDI

(01)06903053058643 (97)105-039945-A0

(10)2023050100 (99)1001452

(17)241214

(98)PA5202305010000045

LOT

2023050100


2024-12-14

IVD

CE 0123

Wiener lab.

047-046074-00(1.0)

HE4 Calibrators

1.2 mL

C0

CE

0123

LOT

2023050100

2024-12-14

IVD

Wiener lab.

HE4 Calibrators

1.0 mL

C1

CE

0123

LOT

2023050100

2024-12-14

IVD

Wiener lab.

HE4 Calibrators

1.0 mL

C2

CE

0123

LOT

2023050100

2024-12-14

IVD

Wiener lab.

Tumor Marker Multi Control II

Order Information

Product	Catalog No.	Package Size
Tumor Marker Multi Control II (L)	1001453	3×2.0 mL/vial
	1001454	6×2.0 mL/vial
Tumor Marker Multi Control II (H)	1001455	3×2.0 mL/vial
	1001456	6×2.0 mL/vial

Intended Use

Wiener lab. Tumor Marker Multi Control II is used for quality control by monitoring the accuracy and precision of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer and test ability of clinical laboratory in the quantitative measurement of tumor marker analytes.

Summary

Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, Wiener lab. reagent kits, calibrators, and controls. Wiener lab. Tumor Marker Multi Control II can be used to check whether QC result is within specified range or not. It is used with Wiener lab. chemiluminescent immunoassays, including Human Epididymal Protein 4 (CLIA), Carbohydrate Antigen 50 (CLIA), Carbohydrate Antigen 242 (CLIA) and Squamous Cell Carcinoma Antigen (CLIA).

Components

Tumor Marker Multi Control II contains Human Epididymal Protein 4 (HE4), Carbohydrate Antigen 50 (CA50), Carbohydrate Antigen 242 (CA242), Squamous Cell Carcinoma Antigen (SCCA) analytes, 0.05% ProClin 300 and 0.02% 5-Bromo-5-Nitro-1,3-Dioxane (BND) in buffered bovine serum albumin (BSA).

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Laboratory professionals must be trained to take all necessary precautions while handling laboratory reagents. Do not use reagents beyond the expiration date.
- The control values are lot-specific with the matched set listed in the value sheet.
- The controls should be run daily in parallel with the patient samples, after each calibration, or after the reagent lot has been changed.
- Please perform maintenance and standard operation as scheduled, including calibration and system analysis to assure the performance of measurement.
- The controls were tested with CE-approved kits for HIV-1/2 Antibody, HCV Antibody, Treponema pallidum Antibody and HBsAg. The results are all negative. However, as no testing method can rule out the potential risk

of infection with absolute certainty, this material should be handled as potential biohazard.1

- If there is any evidence of microbial contamination or excessive turbidity in the product, contact Wiener lab. Customer Service for help.
- Material safety data sheet is available upon request.
- Disposal of all waste material should be handled in accordance with local guidelines.
- The target values and ranges of Wiener lab. Tumor Marker Multi Control II are for reference only. Each laboratory can establish its own internal quality control procedure and set up its own target values and ranges.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
- If the reagents are opened unintentionally before use, they shall be tightly capped and stored at 2-8°C, and then used within in-use stability time.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

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Test Method

Tumor Marker Multi Control II has similar components as test samples. Therefore, the controls can be tested by following the same procedure as tested samples, the target values and ranges are then calculated and compared to see if the controls meet specifications.

Preparation

- Refer to Wiener lab. System Operator's Manual and / or Customer Help Desk to acquire information about control principle, control use, and control validity evaluation.
- Take the vial out of fridge, and let it stand for no less than 30 minutes at room temperature.
- Mix contents by gently inverting the vial at least 10 times before use. Avoid bubble formation.
- Dispense the required volume into a sample cup, and then run QC test you have requested, or have it analyzed in the same way as patient samples.

Storage and Stability

Tumor Marker Multi Control II is stable up to the expiration date as indicated on the label when stored in unopened vial at 2-8°C. Once opened, it is stable for 30 days at 2-8°C or 60 days at -20°C.

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- Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, reagent kits, and matched calibrators.
- General laboratory equipment.

Assay Values

The control values (the target value and range) determined by Wiener lab. routine measurement procedure are listed in the target value sheet. The target value was obtained from Wiener lab. measurement system and the range was calculated as the target values ± 3 standard deviations. The control values are lot-specific. Please check the lot number of each control before use.

Quality Control

Use Wiener lab. controls as the samples to monitor the performance of system. The values of controls are lot-specific on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer. The control result should be within the defined range as shown in the control target value sheet. If the controls are out of specified range, the measurement system and reagents should be checked, for instance, expiration date or storage condition of the calibrators, reagent and controls; possible contamination of reagent; position of reagent, sample placement on the analyzer; the parameter setting in software; or operation status of the analyzer. Each laboratory should establish its own internal quality control scheme and procedure for corrective and preventive actions.

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Contents

References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.

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WIENER Laboratorios S.A.I.C.
Referencia: WIENER Laboratorios S.A.I.C.