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

Roche Diabetes Care India Pvt. Ltd.  
146-B, 166A, Unit no. 7, 8, 9,  
8th Floor, R City Office, R City Mall,  
Lal Bahadur Shastri Marg, Ghatkopar,  
Mumbai Suburban - 400086,  
Maharashtra, India.  
Customer Service: Accu-Chek Extra Care:  
Toll-Free: 1800 120 6020  
SMS - HELP <Message> to 5757535  
Email: india.accu-chek@roche.com  
Monday to Saturday 7.30am to 9.30pm

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Roche

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MEDIDOR DE GLUCEMIA

This file may not print or view at 100%.  
Die lines and color breaks do not print.  
“Holding Area” text, box and rules do  
not print.

Roche USA – 83856-  
V1/1 – 07791208061A /  
07791208061(01) – Black

## PMS-MANUAL FOLD OUT

Approved according to  
SHR\_SOP 04.07.08\_Att A

**Category 2**

The Signature will not be printed!

LAN

IPM

PM/LAB

GD/LAY

Farm. ROBERTA MILEAZZA  
PRODUCED BY ROCHÉ S.A. del.  
Division Diagnostica  
DT & APOTECARIA LEGAL

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TESTS

REF 07453701 / 07453710 / 07453736

Device for self-testing  
Device for near-patient testing

## Intended Use

The Accu-Chek Guide test strips with the Accu-Chek Guide family of meters are indicated to quantitatively measure glucose in fresh capillary whole blood from the finger, palm, forearm, and upper arm as an aid in monitoring the effectiveness of glucose control.

The Accu-Chek Guide test strips with the Accu-Chek Guide family of meters are intended for in vitro diagnostic self-testing by people with diabetes.

The Accu-Chek Guide test strips with the Accu-Chek Guide family of meters are intended for in vitro diagnostic near-patient testing by healthcare professionals in clinical settings. Venous, arterial, and neonatal blood testing is limited to healthcare professional use.

The systems are not for diagnosing diabetes mellitus and not for testing neonate cord blood samples.

## Consumer Information

Read this package insert and the User's Manual before performing a blood glucose test. The package insert contains warnings and precautions:  
A **WARNING** indicates a foreseeable serious hazard.

A **PRECAUTION** describes a measure you should take to use the product safely and effectively or to prevent damage to the product.

### ⚠ WARNING

#### Risk of suffocation

This product contains small parts that can be swallowed. Keep the small parts away from small children and people who might swallow small parts.

#### Risk of a serious health incident

Failure to follow testing instructions or test strip storage and handling instructions can lead to an incorrect test result that may lead to improper therapy.

Carefully read and follow the instructions in the User's Manual and package inserts for the test strips and control solutions.

#### Risk of infection

Human blood is a potential source for the transmission of infection. Avoid exposing other people to contaminated components.

Discard a used test strip as infectious material according to the regulations applicable in your country.

#### Contents of the pack

Pack containing test strips and package inserts.

All components of the pack can be discarded in domestic waste. Discard used test strips according to local regulations. If you have any questions, contact Roche.

#### Test strip storage and handling

- If the container is open or damaged before using the test strips for the first time, if the cap is not fully closed, if you see any damage to the cap or container, or if anything prevents the cap from closing properly, do not use the test strips. Contact Roche.
- Store the test strips at temperatures between 4 and 30 °C. Do not freeze the test strips.
- Use the test strips at temperatures between 4 and 45 °C.
- Use the test strips between 10 and 90 % humidity. Do not store the test strips in high heat and moisture areas such as the bathroom or kitchen.
- Store the unused test strips in their original test strip container with the cap closed.
- Close the test strip container tightly immediately after removing a test strip to protect the test strips from humidity.
- Use the test strip immediately after removing it from the test strip container.
- Discard the test strips if they are past the use by date. Expired test strips can produce incorrect results. The use by date is printed on the test strip box and on the label of the test strip container next to . The test strips can be used until the printed use by date when they are stored and used correctly. This applies for test strips from a new, unopened test strip container and for test strips from a test strip container already opened by the user.
- Use a test strip only once. Test strips are for single use only.

## Performing a Blood Glucose Test

Refer to the meter User's Manual for instructions on obtaining a blood sample and performing a blood glucose test.

## Understanding Test Results

The normal fasting glucose level for a non-diabetic adult is below 100 mg/dL (5.6 mmol/L).<sup>1</sup> The normal glucose level for a non-diabetic adult 2 hours post-meal, e.g. simulated by 75 g Oral Glucose Tolerance Test (OGTT), is less than 140 mg/dL (7.8 mmol/L).<sup>2</sup> A criterion for the diagnosis of diabetes in adults is a fasting glucose level of 126 mg/dL or higher (7.0 mmol/L or higher) confirmed in two tests.<sup>1,2,4</sup> Adults with a fasting glucose level between 100 and 125 mg/dL (5.6 and 6.9 mmol/L) are defined as having impaired fasting glucose (prediabetes).<sup>1</sup> Other diagnostic criteria for diabetes exist. Consult your healthcare professional to determine if you have diabetes or not. For people with diabetes: Consult your healthcare professional for the blood glucose range appropriate for you. You should treat your low or high blood glucose as recommended by your healthcare professional.

#### Unusual test results

If **LO** is displayed on the meter, your blood glucose may be below 10 mg/dL (0.6 mmol/L). If **HI** is displayed on the meter, your blood glucose may be over 600 mg/dL (33.3 mmol/L). For detailed information on error messages, refer to the User's Manual.

### ⚠ PRECAUTION

#### Risk of a serious health incident

Never ignore symptoms or make significant changes to your diabetes therapy without talking to your healthcare professional. If your blood glucose result does not match how you feel, follow these steps:

- Repeat the blood glucose test with a new test strip.
- Perform a control test as described in the User's Manual.
- Refer to the User's Manual for other causes.

If your symptoms still do not match your blood glucose results, contact your healthcare professional.

## Accu-Chek Guide Control Ranges

Control 1: 30–60 mg/dL (1.7–3.3 mmol/L)

Control 2: 252–342 mg/dL (14.0–19.0 mmol/L)

## Healthcare Professional Information

The system can be used in professional healthcare environments such as doctors' offices, general wards, in suspected cases of diabetes, and in emergency cases.

#### Sample collection and preparation by healthcare professionals

- When using the Accu-Chek Guide family of meters, always follow the recognized procedures for handling objects that are potentially contaminated with human material. Practice the hygiene and safety policy of your laboratory or institution.
- A blood drop is required to perform a blood glucose test. Capillary blood can be used. Venous, arterial, or neonatal blood may be used, but must be obtained by healthcare professionals.
- Take caution to clear arterial lines before the blood sample is obtained and applied to the test strip.
- The system has been tested with neonatal blood. As a matter of good clinical practice, caution is advised in the interpretation of neonate blood glucose values below 50 mg/dL (2.8 mmol/L). Follow the recommendations for follow-up care that have been set by your institution for critical blood glucose values in neonates.
- To minimize the effect of glycolysis, venous or arterial blood glucose tests need to be performed within 30 minutes of obtaining the blood samples.
- Avoid air bubbles when using pipettes.
- Capillary, venous, and arterial blood samples containing these anticoagulants or preservatives are acceptable: EDTA, lithium heparin, or sodium heparin. Anticoagulants containing iodacetate or fluoride are not recommended.
- Refrigerated samples should be brought to room temperature slowly prior to testing.

#### Additional information for healthcare professionals

If the blood glucose result does not reflect the patient's clinical symptoms, or seems unusually high or low, perform a control test. If the control test confirms that the system is working properly, repeat the blood glucose test. If the second blood glucose result still seems unusual, follow facility guidelines for further action.

Discard components of the pack per facility guidelines. Consult local ordinances as they may vary by country.

## Limitations

Certain health conditions can lead to an incorrect test result. If you know that one or more of the following health conditions apply to you, do not use the test strip. If you are unsure whether any of the health conditions apply to you, contact your healthcare professional.

- Lipemic samples (triglycerides) >1,800 mg/dL (>20.3 mmol/L) may produce elevated blood glucose results.
- Do not use this system during xylose absorption test.
- Do not use this system if you are undergoing intravenous administration of ascorbic acid.
- If peripheral circulation is impaired, collection of capillary blood from the approved sample sites is not advised as the results might not be a true reflection of the physiological blood glucose level. This may apply in the following circumstances: Severe dehydration as a result of diabetic ketoacidosis or due to hyperglycemic hyperosmolar non-ketotic syndrome, hypotension, shock, decompensated heart failure NYHA Class IV, or peripheral arterial occlusive disease.
- Hematocrit should be between 10 and 65 %. Ask your healthcare professional if you do not know your hematocrit.
- This system has been tested at altitudes up to 3,094 meters. Do not use this system at altitudes above 3,094 meters.

## Performance Characteristics

The Accu-Chek Guide family of meters complies with the requirements of ISO 15197:2013 (In vitro diagnostic test systems – Requirements for blood glucose monitoring systems for self-testing in managing diabetes mellitus).

**Calibration and traceability:** The system (meter and test strips) is calibrated with venous blood containing various glucose concentrations as a calibrator. The reference values are obtained using the hexokinase method which is calibrated using the ID-GCMS method. The ID-GCMS method as the method of highest metrological quality (order) is traceable to a primary NIST standard. Using this traceability chain, the results obtained with these test strips for control solutions can also be traced back to the NIST standard.

**Detection limit (lowest value displayed):** 10 mg/dL (0.6 mmol/L) for the test strip

**System measuring range:** 10–600 mg/dL (0.6–33.3 mmol/L)

**Sample size:** 0.6 µL

**Test time:** <4 seconds

#### System accuracy:

System accuracy results for glucose concentrations less than 100 mg/dL (less than 5.55 mmol/L)

within ±5 mg/dL (within ±0.28 mmol/L)	within ±10 mg/dL (within ±0.56 mmol/L)	within ±15 mg/dL (within ±0.83 mmol/L)
175/186 (94.1 %)	186/186 (100 %)	186/186 (100 %)

System accuracy results for glucose concentrations equal to or greater than 100 mg/dL (equal to or greater than 5.55 mmol/L)

within ±5 %	within ±10 %	within ±15 %
296/414 (71.5 %)	404/414 (97.6 %)	413/414 (99.8 %)

System accuracy results for glucose concentrations between 14 mg/dL (0.78 mmol/L) and 578 mg/dL (32.1 mmol/L)

within ±15 mg/dL or within ±15 % (within ±0.83 mmol/L or within ±15 %)
599/600 (99.8 %)

#### Repeatability:

Mean value	[mg/dL]	40.5	81.7	132.1	206.7	330.2
	[mmol/L]	2.3	4.5	7.3	11.5	18.3
Standard deviation	[mg/dL]	1.4	2.0	2.8	5.4	8.6
	[mmol/L]	0.1	0.1	0.1	0.3	0.5
Coefficient of variation [%]	—	—	2.1	2.6	2.6	2.6

#### Intermediate precision:

Mean value	[mg/dL]	44.9	116.6	297.4
	[mmol/L]	2.5	6.5	16.5
Standard deviation	[mg/dL]	1.4	2.8	6.8
	[mmol/L]	0.1	0.2	0.4
Coefficient of variation [%]	—	2.4	2.3	

**Performance assessment by the user:** A representative study was performed using the Accu-Chek Guide blood glucose monitoring system evaluating glucose values from fingertip capillary blood samples obtained by 100 lay persons and showed the following results:

- For glucose concentrations less than 100 mg/dL (less than 5.55 mmol/L), 100 % of the test results were within ±15 mg/dL (within ±0.83 mmol/L) of the results obtained through laboratory testing.
- For glucose concentrations equal to or greater than 100 mg/dL (equal to or greater than 5.55 mmol/L), 98.9 % of the test results were within ±15 % of the results obtained through laboratory testing.

**Test principle:** The enzyme on the test strip, a FAD-dependent glucose dehydrogenase (GDH) expressed in *A. oryzae*, converts the glucose in the blood sample to gluconolactone. This reaction creates a harmless DC electrical current that the meter interprets for the blood glucose result. The sample and environmental conditions are evaluated using AC and DC signals.

These test strips deliver results that correspond to blood glucose concentrations in plasma as per the recommendation of the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).<sup>3</sup> Therefore, the meter displays blood glucose concentrations that refer to plasma although whole blood is always applied to the test strip.

#### Reagent composition<sup>3</sup>

Mediator	6.6 %
FAD-GDH enzyme	21.3 %
Buffer	22.6 %
Stabilizer	2.3 %
Non-reactive ingredients	47.2 %

©Minimum at time of manufacture

**Note:** For an explanation of symbols used and a list of references, refer to the end of this package insert.

#### Control and linearity test kits (if available)

Accu-Chek Guide control solution – Refer to the control solution package insert for details.  
Accu-Chek Guide linearity test kit – Refer to the linearity test kit package insert for details.  
Visit our website at [www.accu-chek.com](http://www.accu-chek.com) or contact the local Roche representative for more information.

## Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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- American Diabetes Association. 2. Classification and Diagnosis of Diabetes: Standards of Medical Care in Diabetes-2021. *Diabetes Care*. 2021 Jan; 44(Supplement 1): S15-S33.
- American Diabetes Association website: Diagnosing Diabetes and Learning about Prediabetes. <http://www.diabetes.org/diabetes-basics/diagnosis/>. Accessed April 22, 2019.
- IDF Clinical Guidelines Task Force. Global guideline for Type 2 diabetes. Brussels: International Diabetes Federation, 2012.
- Definition and diagnosis of diabetes mellitus and intermediate hyperglycemia: report of a WHO/IDF consultation. WHO, Geneva 2006 (ISBN 92 4 159493 4, ISBN 978 92 4 159493 6).
- D’Orazio et al.: “Approved IFCC Recommendation on Reporting Results for Blood Glucose (Abbreviated),” *Clinical Chemistry* 51:9 1573-1576 (2005).

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
KEEP OUT AREA  
does not print,  
for reference only.

989039576680

08394563068A

KEEP OUT AREA  
does not print,  
for reference only.

# Accu-Chek® Guide

<b>CONTROLS</b>	
	REF 07748906

**EN**
**Device for self-testing**  
**Device for near-patient testing**  
**For in vitro diagnostic use**

**Intended Use**

The Accu-Chek Guide control solution is intended for performing control tests on the Accu-Chek Guide family of blood glucose meters and the Accu-Chek Guide test strips. It is intended for self-testing by people with diabetes and for near-patient testing by healthcare professionals. Testing control solutions with known glucose levels establishes that the operator and the system are performing acceptably. Control results must be within the defined acceptable ranges before valid patient testing is allowed.

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Do not ingest or inject the control solution or use the control solution as eye drops.

EUH208 - Contains DIAZOLIDINYL UREA. May produce an allergic reaction.\*

\*Regulation (EC) No 1272/2008 (CLP)

##### Contents of the pack

Pack containing 2 × 2.5 mL control solution and package inserts. The Accu-Chek Guide control solution pack contains two control solutions, one for the hypoglycemic range (control solution 1, gray cap) and one for the hyperglycemic range (control solution 2, white cap).

##### Disposal

All components of the pack can be discarded in domestic waste. Consult local ordinances as they may vary by country. If you have any questions, contact Roche.

##### Control solution storage and handling

- Refer to the test strip package insert for system operating conditions.
- The printed use by date is valid if the unopened control solutions are stored between 2 and 32 °C
- The control solution does not need to be kept in the refrigerator. Do not freeze the control solution.
- Write the date the control solution bottle was opened on the bottle label. The control solution must be discarded 3 months from the date the control solution bottle was opened (discard date) or on the use by date on the bottle label, whichever comes first.
- Control solution taken directly from the refrigerator must be allowed to adjust to room temperature (without opening the control solution bottle).
- The control solution can stain fabric. Remove stains by washing with soap and water.

#### Performing a Control Test

Refer to the meter User's Manual for instructions on performing a control test.

#### Sources of Error

If the control results are outside the acceptable range, do not use the meter until you solve the problem. Check this list to help solve the problem.

- Were the test strips or control solutions expired?
- Did you wipe the tip of the control solution bottle before use?
- Were the caps on the test strip container and control solution bottle always closed tightly?
- Was the test strip used immediately after it was removed from the test strip container?
- Were the test strips and control solutions stored in a cool, dry place?
- Was the correct control solution level selected when the test was performed?
- Did you follow the directions?

#### Control Testing Intervals

Follow your facility's policy for control testing intervals.

Control tests should be performed:

- the first time before using the meter for patient blood glucose testing.
- at intervals established by the facility.
- when a new test strip box is opened.
- if the test strip container was left open.
- if test strips were improperly stored.
- if there is a question about a blood glucose result.
- to check the performance of the system.
- if the meter was dropped.

Your facility may require that control tests be successfully performed **after** any of the following occur and **before** patient testing resumes:

- Previous control results were outside the acceptable range.
- Control tests were not performed at the proper interval.

Control results must be within the designated range on the test strip container label, or as defined by your facility, before being considered acceptable. Patients can be tested after control tests have been acceptably performed at the proper testing interval.

Ingredients		
Component	Control Solution Level 1 (Low)	Control Solution Level 2 (High)
Glucose	0.06 <span> </span> %	0.30 <span> </span> %
Buffer	4.83 <span> </span> %	4.82 <span> </span> %
Biological salt	3.39 <span> </span> %	3.38 <span> </span> %
Preservative	0.30 <span> </span> %	0.30 <span> </span> %
Non-reactive ingredients	9.99 <span> </span> %	9.97 <span> </span> %
FD&C Blue #1	0.08 <span> </span> %	0.08 <span> </span> %
Percentages in weight/weight		

**Note:** For an explanation of symbols used, refer to the end of this package insert. Visit our website at [www.accu-chek.com](http://www.accu-chek.com) or contact the local Roche representative for more information.

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are trademarks of Roche,
ACCU-CHEK, ACCU-CHEK GUIDE, ACCU-CHEK GUIDE ME, 罗氏, and แอดคิว-เช็ก are trademarks of Roche.
**Referencia:** PRODUCTOS ROCHE S.A.