Seccion 1 – Sobrerótulos

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Fondo: Blanco

Tamaño: 4,5cm x 3cm

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Bioline™ Zika IgM

Importador: Abbott Rapid Diagnostics
Argentina S.A.

Dom.: Ing. Eiffel 4180, El Triángulo, Bs. As.

DT. Farm. Carolina Tchicourel

Autorizado por ANMAT: PM-1275-274

NO APTO PARA BANCO DE SANGRE

[Código de barras]

12FK21

[Revisión]

Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

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INTENDED USE

The Bioline™ Zika IgM test is an *in vitro* immunochromatographic assay for the qualitative detection of IgM antibodies to Zika virus in human serum, plasma or whole blood. The test is intended for professional use, to aid in the presumptive diagnosis of Zika infections.

INTRODUCTION

Zika virus infection has recently received great scrutiny and clinical attention since the linkages to fetal birth defects such as microcephaly and the linkage to Guillain-Barré syndrome have been observed. Zika virus is a member of the *Flaviviridae* family and transmitted to humans primarily through the bite of an infected *Aedes* species mosquito (*A. aegypti* and *A. albopictus*); the same vectors of dengue and chikungunya viruses. Most Zika infections are asymptomatic. However, the symptoms when present are similar to those of dengue and chikungunya, and may include fever, rash, joint pain, conjunctivitis (red eyes), muscle pain and headache. Thus it is difficult to symptomatically differentiate Flavivirus infections including dengue, chikungunya, West Nile, Japanese Encephalitis and yellow fever virus. Diagnosis of Zika virus infections can be made by direct detection of viral RNA or antibodies (IgM). The detection of viral RNA from Zika virus is often detectable in serum during the acute phase of infection (<14 days after illness onset). There is limited evidence to suggest that some patients may have RNA detectable for a longer period of time. Detectable levels of anti-Zika IgM antibodies in infected individuals typically appear within the first week after symptoms onset and persist for approximately 8-12 weeks.

TEST PRINCIPLE

The Bioline™ Zika IgM Test is a qualitative assay for the detection of IgM antibodies for the Zika virus. When patient specimen is added to the specimen well, Zika specific-IgM antibodies present in the specimen bind to anti-human IgM antibodies on the membrane. Colloidal gold complexes containing inactivated Zika antigen (Zika virus PRVABC59 strain, classified as *Flaviviridae*, *Flavivirus* and identified from human serum specimen, Puerto Rico, December 2015) migrate along the membrane chromatographically and are captured by the bound patient's IgM to produce colored Test Lines adjacent to "T" (Test Line). A second line, the Control Line ("C") is used for procedural control and should always appear adjacent to "C" (Control Line) if the test procedure is performed correctly and the active ingredients of the main components on the strip are working. When Zika antibodies are present in the specimen, both the Control ("C") and Test ("T") Lines appear in the result window, indicating a reactive result. Only the Control Line ("C") appears for specimens that do not contain Zika IgM, or that contain antibodies below the detection limit of the assay, indicating a non-reactive result.

MATERIALS PROVIDED AND ACTIVE INGREDIENTS OF MAIN COMPONENTS

1. The Bioline™ Zika IgM test kit contains the following items to perform the assay.

1-1 Bioline™ Zika IgM : Cat. No. 12FK20

- 25 Test devices with desiccant in individual foil pouches
- Assay diluent (1 x 5 ml/vial) ⚠
- 25 Capillary pipettes (10 µl)
- 1 Instructions for Use

1-2 Bioline™ Zika IgM : Cat. No. 12FK21

- 25 Test devices with desiccant in individual foil pouches
- Assay diluent (1 x 5 ml/vial) ⚠
- 25 Capillary pipettes (10 µl)
- 25 Safety lancets
- 25 Alcohol swabs
- 1 Instructions for Use

1-3 Bioline™ Zika IgM : Cat. No. 12FK26

- 25 Test devices with desiccant in individual foil pouches
- Assay diluent (1 x 5 ml/vial) ⚠
- 25 Capillary pipettes (10 µl)
- 25 Sterile lancets
- 25 Alcohol swabs
- 1 Instructions for Use

Farm. Carolina A. Tchicourei
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

2. Active ingredients of main components

1 test strip includes :

Gold conjugate : Mouse monoclonal anti-Zika virus antibody conjugated gold colloid ($0.080 \pm 0.024 \mu\text{g}$)

Test Line : Mouse monoclonal anti-human IgM antibody ($0.640 \pm 0.192 \mu\text{g}$),

Control Line : Goat anti-mouse IgG ($0.580 \pm 0.174 \mu\text{g}$)

Assay diluent includes :

Tris buffer (q.s.), Tween 20 (q.s.), Sodium chloride (q.s.), Proclin™ 300 (q.s.), Sodium azide (q.s.), Normal human serum (q.s.)

※ Classification according to Regulation (EC) No 1272/2008:

Product identifier	Trade name	Bioline™ Zika IgM Assay diluent
	Hazardous substance	Proclin™ 300 (5-chloro-2-methyl-4- isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)

Classification	Skin sensitization Category 1
Hazard Pictogram	
Signal Word	Warning
Hazard Statement	H317: May cause an allergic skin reaction
Precautionary statements	
Prevention	P261: Avoid breathing dust/fume/gas/mist/vapours/spray P272: Contaminated work clothing should not be allowed out of the workplace P280: Wear protective gloves/protective clothing/eye protection/face protection
Response	P302+P352: IF ON SKIN: Wash with plenty of soap and water P321: Specific treatment P333+P313: If skin irritation or rash occurs: Get medical advice/attention P363: Wash contaminated clothing before reuse
Disposal	P501: Dispose of contents/container in accordance with local/regional/national/international regulation

MATERIALS REQUIRED BUT NOT PROVIDED

- Protective gloves, Timer, Biohazard container, Micropipette

KIT STORAGE AND STABILITY

The test device is sensitive to both heat and humidity and must be used and stored properly. Store the test kit at temperatures between 2 °C and 30 °C. Do not expose the kit or its components to freezing or high temperatures. The Bioline™ Zika IgM kit and assay diluent are stable until the expiration dates marked on their outer packaging and containers. Assay diluent may be opened and resealed for each assay. Cap should be firmly sealed between each use. Assay diluent is stable until expiration date if kept at 2 - 30 °C.

Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

WARNINGS

1. For *in vitro* diagnostic use only.
2. Do not reuse test device and kit components.
3. The instructions must be followed exactly to achieve accurate results. Any individual using this product must be trained in its use.
4. Do not eat or smoke while handling specimens and kit.
5. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
6. Avoid splashing or aerosol formation of specimen and assay diluent.
7. Clean up spills thoroughly using an appropriate disinfectant.
8. Decontaminate and dispose of all specimens, tested devices and potentially contaminated materials in a biohazard container, the same as if they were infectious waste.
9. Use one test per specimen. Do not combine specimens; pooled specimens may give erroneous results.
10. Do not use reagents from different lots or those for other products.
11. Do not drink assay diluent.
12. The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.
13. Perform the test immediately after opening the foil pouch, to avoid exposure to humidity or moisture.
14. Do not use kit past its expiration date.
15. Do not use the test if the aluminum pouch is damaged.

QUALITY CONTROL

The Bioline™ Zika IgM test has built-in procedural controls (internal quality controls). The Bioline™ Zika IgM test has “T” (Test Line) and “C” (Control Line) marked on the surface of the device. Neither the Test Line nor the Control Line is visible in the result window before applying a specimen. The Control Line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a reactive specimen control.

Farm. Carolina A. Tchicourei
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

SPECIMEN COLLECTION AND HANDLING

Use freshly collected specimens for optimal test performance. Inadequate specimen collection or improper specimen handling/transport may yield incorrect results.

1. Whole blood

[Collection by venipuncture]

- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
- Whole blood specimens should be tested as soon as possible after collection. If whole blood specimens cannot be tested immediately, they must be refrigerated at 2 - 8 °C.
- If stored at 2 - 8 °C, the whole blood specimen must be tested within 3 days of collection.
- Allow all whole blood specimens to warm to room temperature between 15 °C and 30 °C prior to testing in the Bioline™ Zika IgM test.

[Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.
- Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.
- Immerse the open end of a new capillary pipette (10 µl) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

2. Plasma or Serum

- **[Plasma]** Using venipuncture, draw whole blood into the collection tube (containing heparin, EDTA or sodium citrate anticoagulants) and then centrifuge the tube to get a plasma specimen.
- **[Serum]** Using venipuncture, draw whole blood into the collection tube (NOT containing anticoagulant). Leave blood to coagulate for 30 minutes and then centrifuge the tube to get a serum supernatant.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Use standard laboratory techniques to clarify such specimens prior to testing.
- If plasma or serum specimens cannot be tested immediately, they must be refrigerated at 2 - 8 °C. If stored at 2 - 8 °C, the serum or plasma specimen must be tested within 14 days of collection. For storage periods longer than 2 weeks, freezing (below -20 °C) is required.
- Allow all plasma or serum specimens to warm to room temperature between 15 °C and 30 °C prior to testing in the Bioline™ Zika IgM test.

3. Precautions

- Repeated frozen-thawed cycle for specimen should be avoided.
- Use of heparin, EDTA, and citrate anticoagulants do not affect the test result. Use of other anticoagulants has not been validated. Their use may alter the test result and should not be used.
- Use a new disposable capillary pipette or pipette tip for each specimen in order to avoid cross-contamination of specimens, which could produce erroneous results.

TEST PROCEDURE (REFER TO FIGURE)

1. Bring all kit components and specimens to room temperature between 15 °C and 30 °C prior to opening and testing.
2. Invert the specimen collection tube several times to mix.
3. Open the foil pouch and place the test device on a flat, dry surface. Label the test device with a patient identifier. Perform the test immediately to avoid exposing the test device to humidity and moisture.
4. **[Using a capillary pipette]** Draw serum, plasma or whole blood specimen to the black line (10 µl), then carefully place the pipette tip into the square specimen well and dispense the specimen.
OR,
[Using a micropipette] Draw 10 µl of serum, plasma or whole blood specimen, then carefully place the pipette tip into the square specimen well and dispense the specimen.
5. Holding the assay diluent bottle vertically, dispense 4 drops into the round assay diluent well. Do not let the bottle tip touch the device in order to avoid cross-contamination.
6. Interpret test results 15 minutes after adding assay diluent.

Caution: Do not read the test results after 20 minutes; reading after 20 minutes can yield erroneous results.

RESULT INTERPRETATION (REFER TO FIGURE)

1. **Non-reactive:** The presence in the result window of one colored line at "C" (Control Line) indicates that IgM antibodies to Zika are not present in the specimen or are present below the minimum level of detection.
2. **Zika IgM reactive:** The presence in the result window of two colored lines; one each at "C" (Control Line) and "T" (Test Line) indicates that the specimen is reactive for IgM antibodies to Zika virus.

Caution : If any "T" (Test Line) appears, no matter how faint, the result is considered reactive.

3. **Invalid:** If no colored line is visible in the result window at "C" (Control Line), the result is invalid. The instructions may not have been followed correctly or the test may have been compromised. When an invalid result is obtained, it is recommended that the specimen is checked for integrity and retested using a new device.

TEST LIMITATIONS

1. A non-reactive result can occur if the Zika antibody is not present in the specimen or is below the detection limit of this assay. A non-reactive test result does not exclude the possibility of Zika infection.
2. As with all diagnostic tests, results must be considered with other clinical information available to the physician.
3. Serological cross reactivity across the Flavivirus genus (e.g. Dengue virus, West Nile virus, St. Louis Encephalitis virus, Japanese Encephalitis virus, Yellow Fever virus,) is common. Limited evaluation of cross-reactivity with flaviviruses has been conducted. Cross reactivity was only observed with dengue virus. Cross-reactivity with the Flavivirus genus may be possible.
4. The Bioline™ Zika IgM test has been evaluated using only the specimen types indicated. Performance with other specimen types has not been evaluated.
5. Test results can vary according to the time at which the specimen was collected following onset of symptoms, the specimen type, serotype present in the population tested, reference method utilized and other factors.
6. The Bioline™ Zika IgM test provides a qualitative result. Therefore, correlation cannot be drawn between the line intensity of a reactive assay signal and the concentration of IgM in the patient specimen.

PERFORMANCE CHARACTERISTICS

1. CLINICAL PERFORMANCE

- A total of 156 plasma specimens collected from Brazil were tested on the Bioline™ Zika IgM in Brazil. The specimens were comprised of 98 positive and 58 negative generated by CDC MAC ELISA. The Bioline™ Zika IgM showed a sensitivity of 90.8 % (89/98) and a specificity of 98.3 % (57/58). Also, no cross-reactivity (0/68) for Dengue positive plasma specimens on the Bioline™ Zika IgM.

Farm. Carolina A. Tchicourei
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

		CDC MAC ELISA	
		Positive	Negative
Bioline™ Zika IgM	Positive	89	1
	Negative	9	57
Total		98	58
Sensitivity		90.8 % (89/98) (95% CI: 83.3 – 95.7 %)	
Specificity		98.3 % (57/58) (95% CI: 90.8 – 100 %)	

• The clinical performance of the Bioline™ Zika IgM test was evaluated with 228 serum or plasma specimens obtained from five commercial resources. 130 specimens suspected of acute Zika infection and 98 normal human serum specimens were tested with Bioline™ Zika IgM and the InBios ZIKV Detect™ IgM Capture ELISA.

The results from the Bioline™ Zika IgM test were compared to the InBios ZIKV Detect™ IgM Capture ELISA. For the purposes of this evaluation, a specimen was classified as “Zika Positive” when the InBios ZIKV Detect™ IgM Capture ELISA was positive. A specimen was classified as “Negative” when the InBios ZIKV Detect™ IgM Capture ELISA was negative. The summary of results is shown in the table below.

		InBios ZIKV Detect™ IgM Capture ELISA	
		Reactive	Non-reactive
Bioline™ Zika IgM	Reactive	65	3
	Non-reactive	3	157
Total		68	160
Positive percent agreement (PPA) (%)		95.6 % (65/68) (95 % CI: 87.8 – 98.5 %)	
Negative percent agreement (NPA) (%)		98.1 % (157/160) (95 % CI: 94.6 – 99.4 %)	
Overall percent agreement (%)		97.4 % (222/228)	

• Of the 228 specimens tested, 55 were confirmed positive by Zika Real Time PCR (test results provided by the vendor). Bioline™ Zika IgM demonstrated 95.2 % PPA and 97.1 % NPA compared to the Zika Real Time PCR assay.

		Zika Real Time PCR	
		Reactive	Non-reactive
Bioline™ Zika IgM	Reactive	20	1
	Non-reactive	1	33
Total		21	34
Positive percent agreement (PPA) (%)		95.2 % (20/21) (95 % CI: 77.3 – 99.2 %)	
Negative percent agreement (NPA) (%)		97.1 % (33/34) (95 % CI: 85.1 - 99.5 %)	

• Of the 228 specimens tested, 22 were from pregnant women. Bioline™ Zika IgM demonstrated 100 % PPA and 84.6 % NPA compared to the InBios ZIKV Detect™ IgM Capture ELISA.

		InBios ZIKV Detect™ IgM Capture ELISA	
		Reactive	Non-reactive
Bioline™ Zika IgM	Reactive	9	2
	Non-reactive	0	11
Total		9	13
Positive percent agreement (PPA) (%)		100 % (9/9) (95 % CI: 70.1 – 100 %)	
Negative percent agreement (NPA) (%)		84.6 % (11/13) (95 % CI: 57.8 - 95.7 %)	

The stated positive and negative agreement performance of the assay is derived from the above studies. Users should note that performance of individual laboratories may vary because of characteristics including flavivirus cross-reactivity, sampling time after onset of illness, patient ethnicity, dietary, environmental and other factors of the served population.

2. ANALYTICAL SPECIFICITY

• CROSS-REACTIVITY

To determine the analytical specificity of the Bioline™ Zika IgM test, a panel of clinical serum or plasma specimens confirmed positive for viruses and organisms known to cause similar signs and symptoms of the Zika virus infection were tested. All of the following positive clinical specimens were non-reactive by the Bioline™ Zika IgM test, unless otherwise noted.

No.	Disease/ Infectious agent	Number of panel	Bioline™ Zika IgM reactive	Cross Reactivity (%)
1	Anti- West Nile Virus	18	0	0 %
2	Anti-tick borne encephalitis (TBE)	3	0	0 %
3	Anti-Japanese Encephalitis IgM	10	0	0 %
4	Yellow fever virus post-immunization	12	0	0 %
5	Anti- Chikungunya virus IgM	40	0	0 %
6	Anti-Cytomegalovirus (CMV)	20	0	0 %
7	Anti-Epstein Barr virus(EBV)	20	0	0 %
8	Anti- Parvovirus B19	18	0	0 %
9	Anti-varicella zoster virus	10	0	0 %
10	Anti-Malaria	10	0	0 %
11	Adenovirus IgM	3	0	0 %
12	Enterovirus IgM	3	0	0 %
13	Anti- Leptospira	7	0	0 %
14	Anti-Hepatitis C virus(HCV)	15	0	0 %
15	Anti-Nuclear antibody (ANA)	10	0	0 %
16	Rheumatoid Factor (RF)	10	0	0 %
17	Human Anti-Mouse Antibody (HAMA)	9	0	0 %
18	Zika IgG	19	0	0 %

9/89 (10.1 %) anti-Dengue specimens generated reactive results with the Bioline™ Zika IgM test. High cross-reactivity of antibodies to Dengue and Zika viruses is well known.^{3,4}

Although cross reactivity in the results of internal studies was only with anti-dengue virus, the possibility of other cross reactivity cannot be excluded.

• INTERFERING SUBSTANCES

The impact of potentially interfering substances on Bioline™ Zika IgM was evaluated. No interference with potential interfering substances was observed at the concentrations tested, Hemoglobin (2 mg/ml), Bilirubin (conjugate) (0.3 mg/ml), Albumin from human serum (60 mg/ml) and Human leukocyte antigen (1.7 mg/ml) and Human Anti-Mouse Antibody (34 ng/ml).

• REPRODUCIBILITY


Reproducibility of the Bioline™ Zika IgM test was established using in-house reference panels containing strong positive, medium positive, weak positive and negative specimens. There were no differences observed within-run (replicates tested by one operator), between days, between operators (three different test operators), between sites (three sites), or lot-to-lot (3 different lots).

Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

TEST PROCEDURE

I. Blood (by venipuncture), Plasma or Serum specimen

10 µl : 10 µl of Serum, Plasma or Whole blood

4 drops : 4 drops of Assay diluent

1

[Using a micropipette]

Take 10 µl of serum, plasma or whole blood specimen using a micropipette.


2

Carefully place the pipette tip into the square specimen well and dispense the specimen.

Farm. Carolina A. Tchicounei
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.


10 µl

3

Holding the assay diluent bottle vertically, dispense **4 drops** into the round assay diluent well.

Do not let the bottle tip touch the device in order to avoid cross-contamination.

4

Interpret test results 15 minutes after adding assay diluent.

Caution: Do not read the test results after 20 minutes; reading after 20 minutes can yield erroneous results.


Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

TEST PROCEDURE


II. Blood specimen (with a lancet)

1 Clean the area to be lanced with an alcohol swab.


2 Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.


[12FK21]

[12FK26]


3 **[Using a capillary pipette provided]**
Immerse the open end of a new capillary pipette (10 µl) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

Fermi Carolina A. Tchikourou
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

4


Carefully place the pipette tip into the square specimen well and dispense the specimen.

5

Holding the assay diluent bottle vertically, dispense **4 drops** into the round assay diluent well.

Do not let the bottle tip touch the device in order to avoid cross-contamination.

6

Interpret test results 15 minutes after adding assay diluent.

Caution: Do not read the test results after 20 minutes; reading after 20 minutes can yield erroneous results.


Farm. Carolina A. Tchicouret
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

RESULT INTERPRETATION

Non-reactive

The presence in the result window of one colored line at "C" (Control Line) indicates that IgM antibodies to Zika are not present in the specimen or are present below the minimum level of detection.


Zika IgM reactive

The presence in the result window of two colored lines; one each at "C" (Control Line) and "T" (Test Line) indicates that the specimen is reactive for IgM antibodies to Zika virus.

Caution: If any "T" (Test Line) appears, no matter how faint, the result is considered reactive.


Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

Invalid

If no colored line is visible in the result window at "C" (Control Line), the result is invalid. The instructions may not have been followed correctly or the test may have been compromised.

When an invalid result is obtained, it is recommended that the specimen is checked for integrity and retested using a new device.

Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

Product. Disclaimer :

While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the analysis should consult a doctor for further confirmation of the test result.

Warning :

The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect reactive or non-reactive result using this product.

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
1. CDC - National Center for Emerging and Zoonotic Infectious Diseases (NCEZID); Division of Vector-Borne Diseases (DVBD); Zika virus
2. Cécile Beck et al. Flaviviruses in Europe: Complex Circulation Patterns and Their Consequences for the Diagnosis and Control of West Nile Disease. Int. J. Environ. Res. Public Health 2013, 10, 6049-6083
3. Dejnirattisai W. et al. Dengue virus sero-cross-reactivity drives antibody-dependent enhancement of infection with Zika virus, 2016, Nature Immunology. 17, 9: 1102-1108
4. Priyamvada L et al. Human antibody responses after dengue virus infection are highly cross-reactive to Zika virus, 2016, PNAS, 113,28:7852-7857

Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

Authorized Representative
MT Promedt Consulting GmbH
Altenhofstrasse 80, 66386 St. Ingbert, Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

Abbott Diagnostics Korea Inc.
65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, 17099,
Republic of Korea
Tel : 82-31-899-2800 | Fax : 82-31-899-2840
abbott.com/poct

© 2020 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Farm. Carolina A. Tchicourel
M. N. 17617
Directora Técnica
Abbott Rapid Diagnostics Argentina S.A.

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12FK21-40-A0 © 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.
Referencia: ABBOTT RAPID DIAGNOSTICS ARGENTINA S.A. rótulos e instrucciones de uso