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VIASURE West Nile Virus Real Time PCR Detection Kit

I. Sin control de extracción

I. 1. Perfil bajo

I. 1. a. 6 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus REAL TIME PCR DETECTION KIT 6x8-well strip - low profile REF VS-WNV106L LOT WNV106L-xxx yyyy-mm CE IVD 40°C 48 010843544020370417aammd10WNV106Lxxx <table border="0"><tr><td>VS-WNV1SL</td><td>West Nile Virus 8- well strips</td><td>6x8-well strip</td></tr><tr><td>VS-RB02</td><td>Rehydration Buffer</td><td>1 blue vial</td></tr><tr><td>VS-WNV1C</td><td>West Nile Virus Positive control</td><td>1 red vial</td></tr><tr><td>VS-NC1</td><td>West Nile Virus Negative control</td><td>1 violet vial</td></tr><tr><td>VS-H2O</td><td>Water RNase/DNase free</td><td>1 white vial</td></tr><tr><td>VS-OCS</td><td>Tear-off 8-cap strips</td><td>6x8-cap strip</td></tr></table> West Nile Virus REAL TIME PCR DETECTION KIT CERTEST BIO TEC S.L Pol. Industrial Rio Gallego II, Calle J, Nº 1, 50840 San Mateo de Gallego, Zaragoza (SPAIN) VS-WNV106L v.02	VS-WNV1SL	West Nile Virus 8- well strips	6x8-well strip	VS-RB02	Rehydration Buffer	1 blue vial	VS-WNV1C	West Nile Virus Positive control	1 red vial	VS-NC1	West Nile Virus Negative control	1 violet vial	VS-H2O	Water RNase/DNase free	1 white vial	VS-OCS	Tear-off 8-cap strips	6x8-cap strip	Importador: BIOARS S.A. Estomba 961 - C.A.B.A. C.P.: C1427COU Tel.: (011) 4555 4601 www.bioars.com.ar Directora Técnica: Dra. C. Etchevés - Bioquímica M.N.: 7028 Venta exclusiva a laboratorios de análisis clínicos USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496 VS-WNV106L VIASURE West Nile Virus Real Time PCR Detection Kit VS-WNV106L WNV106L-xxx DD-MM-AAAA
VS-WNV1SL	West Nile Virus 8- well strips	6x8-well strip																	
VS-RB02	Rehydration Buffer	1 blue vial																	
VS-WNV1C	West Nile Virus Positive control	1 red vial																	
VS-NC1	West Nile Virus Negative control	1 violet vial																	
VS-H2O	Water RNase/DNase free	1 white vial																	
VS-OCS	Tear-off 8-cap strips	6x8-cap strip																	

I. 1. b. 12 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus REAL TIME PCR DETECTION KIT 12x8-well strip - low profile REF VS-WNV112L LOT WNV112L-xxx yyyy-mm CE IVD 40°C 96 010843544020371117aammd10WNV112Lxxx <table border="0"><tr><td>VS-WNV1SL</td><td>West Nile Virus 8- well strips</td><td>12x8-well strip</td></tr><tr><td>VS-RB02</td><td>Rehydration Buffer</td><td>1 blue vial</td></tr><tr><td>VS-WNV1C</td><td>West Nile Virus Positive control</td><td>1 red vial</td></tr><tr><td>VS-NC1</td><td>West Nile Virus Negative control</td><td>1 violet vial</td></tr><tr><td>VS-H2O</td><td>Water RNase/DNase free</td><td>1 white vial</td></tr><tr><td>VS-OCS</td><td>Tear-off 8-cap strips</td><td>12x8-cap strip</td></tr></table> West Nile Virus REAL TIME PCR DETECTION KIT CERTEST BIO TEC S.L Pol. Industrial Rio Gallego II, Calle J, Nº 1, 50840 San Mateo de Gallego, Zaragoza (SPAIN) VS-WNV112L v.02	VS-WNV1SL	West Nile Virus 8- well strips	12x8-well strip	VS-RB02	Rehydration Buffer	1 blue vial	VS-WNV1C	West Nile Virus Positive control	1 red vial	VS-NC1	West Nile Virus Negative control	1 violet vial	VS-H2O	Water RNase/DNase free	1 white vial	VS-OCS	Tear-off 8-cap strips	12x8-cap strip	Importador: BIOARS S.A. Estomba 961 - C.A.B.A. C.P.: C1427COU Tel.: (011) 4555 4601 www.bioars.com.ar Directora Técnica: Dra. C. Etchevés - Bioquímica M.N.: 7028 Venta exclusiva a laboratorios de análisis clínicos USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496 VS-WNV112L VIASURE West Nile Virus Real Time PCR Detection Kit VS-WNV112L WNV112L-xxx DD-MM-AAAA
VS-WNV1SL	West Nile Virus 8- well strips	12x8-well strip																	
VS-RB02	Rehydration Buffer	1 blue vial																	
VS-WNV1C	West Nile Virus Positive control	1 red vial																	
VS-NC1	West Nile Virus Negative control	1 violet vial																	
VS-H2O	Water RNase/DNase free	1 white vial																	
VS-OCS	Tear-off 8-cap strips	12x8-cap strip																	

I. 1. c. Placa de 96 pocillos

Sobre-Rótulo

West Nile Virus REAL TIME PCR DETECTION KIT 96-well plate - low profile REF VS-WNV113L LOT WNV113L-xxx yyyy-mm CE IVD 40°C 96 010843544020372817aammd10WNV113Lxxx <table border="0"><tr><td>VS-WNV1PL</td><td>West Nile Virus 96-well plate</td><td>1 plate</td></tr><tr><td>VS-RB02</td><td>Rehydration Buffer</td><td>1 blue vial</td></tr><tr><td>VS-WNV1C</td><td>West Nile Virus Positive control</td><td>1 red vial</td></tr><tr><td>VS-NC1</td><td>West Nile Virus Negative control</td><td>1 violet vial</td></tr><tr><td>VS-H2O</td><td>Water RNase/DNase free</td><td>1 white vial</td></tr><tr><td>VS-OCS</td><td>Tear-off 8-cap strips</td><td>12x8-cap strip</td></tr></table> West Nile Virus REAL TIME PCR DETECTION KIT CERTEST BIO TEC S.L Pol. Industrial Rio Gallego II, Calle J, Nº 1, 50840 San Mateo de Gallego, Zaragoza (SPAIN) VS-WNV113L v.02	VS-WNV1PL	West Nile Virus 96-well plate	1 plate	VS-RB02	Rehydration Buffer	1 blue vial	VS-WNV1C	West Nile Virus Positive control	1 red vial	VS-NC1	West Nile Virus Negative control	1 violet vial	VS-H2O	Water RNase/DNase free	1 white vial	VS-OCS	Tear-off 8-cap strips	12x8-cap strip	Importador: BIOARS S.A. Estomba 961 - C.A.B.A. C.P.: C1427COU Tel.: (011) 4555 4601 www.bioars.com.ar Directora Técnica: Dra. C. Etchevés - Bioquímica M.N.: 7028 Venta exclusiva a laboratorios de análisis clínicos USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496 VS-WNV113L VIASURE West Nile Virus Real Time PCR Detection Kit VS-WNV113L WNV113L-xxx DD-MM-AAAA
VS-WNV1PL	West Nile Virus 96-well plate	1 plate																	
VS-RB02	Rehydration Buffer	1 blue vial																	
VS-WNV1C	West Nile Virus Positive control	1 red vial																	
VS-NC1	West Nile Virus Negative control	1 violet vial																	
VS-H2O	Water RNase/DNase free	1 white vial																	
VS-OCS	Tear-off 8-cap strips	12x8-cap strip																	

BIOARS S.A.
Patricia del Carmen Etchevés
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVÉS
DIRECTOR TÉCNICO

I. 2. Perfil alto

I. 2. a. 6 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
6x8-well strip - high profile

REF VS-WNV106H LOT WNV106H-xxx yyyy-mm

CE IVD 40°C 48

010843544020373517aamdd10WNV106Hxxx

VS-WNV1SH	West Nile Virus 8- well strips	6x8-well strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	6x8-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

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VS-WNV106H VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV106H WNV106H-xxx DD.MM.AAAA

I. 2. b. 12 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
12x8-well strip - high profile

REF VS-WNV112H LOT WNV112H-xxx yyyy-mm

CE IVD 40°C 96

010843544020374217aamdd10WNV112Hxxx

VS-WNV1SH	West Nile Virus 8- well strips	12x8-well strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	12x8-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

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VS-WNV112H VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV112H WNV112H-xxx DD.MM.AAAA

I. 2. c. Placa de 96 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
96-well plate - high profile

REF VS-WNV113H LOT WNV113H-xxx yyyy-mm

CE IVD 40°C 96

010843544020375917aamdd10WNV113Hxxx

VS-WNV1PH	West Nile Virus 96-well plate	1 plate
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	12x8-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

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VS-WNV113H VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV113H WNV113H-xxx DD.MM.AAAA

Patricia del Carmen Etchevés
BIOARS S.A.
Patricia del Carmen Etchevés
Presidente

Claudia Etchevés
BIOARS S.A.
BIOQ CLAUDIA ETCHÉVES
DIRECTOR TÉCNICO

I. 3. Perfil para instrumentos Rotor-Gene

I. 3. a. 9 tiras de 4 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
9x4-well strip

REF VS-WNV136 LOT WNV136-xxx yyyy-mm

IVD CE   2°C 40°C  36

010843544021410617yyymmdd10WNV136xxx

VS-WNV1RG	West Nile Virus 4-well strips	9x4 well-strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNAse free	1 white vial
VS-OCS-RG	4-cap strips	9x4-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gallego II, Calle J, Nº 1,
50840 San Mateo de Gallego, Zaragoza (SPAIN)

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VS-WNV136 VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV136 WNV136-xxx DD-MM-AAAA

I. 2. b. 18 tiras de 4 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
18x4-well strip

REF VS-WNV172 LOT WNV172-xxx yyyy-mm

IVD CE   2°C 40°C  72

010843544021411317yyymmdd10WNV172xxx

VS-WNV1RG	West Nile Virus 4-well strips	18x4 well-strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNAse free	1 white vial
VS-OCS-RG	4-cap strips	18x4-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gallego II, Calle J, Nº 1,
50840 San Mateo de Gallego, Zaragoza (SPAIN)

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VS-WNV172 VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV172 WNV172-xxx DD-MM-AAAA

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I. 4. a. 4 tubos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
4 tubes x 24 reactions

REF VS-WNV196T LOT WNV196T-xxx yyyy-mm

CE IVD   2°C 40°C  96

010843544022639017yyymmdd10WNV196Txxx

VS-WNV1TB	West Nile Virus Reaction-Mix tube	4 white vials
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNAse free	1 white vial

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gallego II, Calle J, Nº 1,
50840 San Mateo de Gallego, Zaragoza (SPAIN)

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

VS-WNV196T VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV196T WNV196T-xxx DD-MM-AAAA

BIOARS S.A.
Patricia del Carmen Etchevés
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVES
DIRECTOR TÉCNICO

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II. 1. a. 6 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
6x8-well strip - low profile

REF VS-WNV106LE LOT WNV106LE-xxx yyyy-mm

CE IVD    

010843544022388717aamdd10WNV106LExxx

VS-1WNV1SL	West Nile Virus 8- well strips	6x8-well strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	6x8-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT
CERTEST BIOTEC S.L. Pol. Industrial Rio Gállego II, Calle J, Nº 1,
50840 San Mateo de Gallego, Zaragoza (SPAIN) VS-WNV106LE v.00

Importador: BIOARS S.A. | Estomba 961 - C.A.B.A.
C.P.: C1427COU | Tel.: (011) 4555 4601 | www.bioars.com.ar
Directora Técnica: Dra. C. Etchevés - Bioquímica M.N.:7028

Venta exclusiva a laboratorios de análisis clínicos
USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496
VS-WNV106LE VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV106LE WNV106LE-xxx DD.MM.AAAA

II. 1. b. 12 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
12x8-well strip - low profile

REF VS-WNV112LE LOT WNV112LE-xxx yyyy-mm

CE IVD    

010843544022402017yyymmdd10WNV112LE000

VS-1WNV1SL	West Nile Virus 8- well strips	12x8-well strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	12x8-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT
CERTEST BIOTEC S.L. Pol. Industrial Rio Gállego II, Calle J, Nº 1,
50840 San Mateo de Gallego, Zaragoza (SPAIN) VS-WNV112LE v.00

Importador: BIOARS S.A. | Estomba 961 - C.A.B.A.
C.P.: C1427COU | Tel.: (011) 4555 4601 | www.bioars.com.ar
Directora Técnica: Dra. C. Etchevés - Bioquímica M.N.:7028

Venta exclusiva a laboratorios de análisis clínicos
USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496
VS-WNV112LE VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV112LE WNV112LE-xxx DD.MM.AAAA

II. 1. c. Placa de 96 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
96-well plate - low profile

REF VS-WNV113LE LOT WNV113LE-xxx yyyy-mm

CE IVD    

010843544022404417aamdd10WNV113LExxx

VS-1WNV1PL	West Nile Virus 96-well plate	1 plate
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	12x8-cap strip


West Nile Virus REAL TIME PCR DETECTION KIT
CERTEST BIOTEC S.L. Pol. Industrial Rio Gállego II, Calle J, Nº 1,
50840 San Mateo de Gallego, Zaragoza (SPAIN) VS-WNV113LE v.00

Importador: BIOARS S.A. | Estomba 961 - C.A.B.A.
C.P.: C1427COU | Tel.: (011) 4555 4601 | www.bioars.com.ar
Directora Técnica: Dra. C. Etchevés - Bioquímica M.N.:7028

Venta exclusiva a laboratorios de análisis clínicos
USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496
VS-WNV113LE VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV113LE WNV113LE-xxx DD.MM.AAAA

BIOARS S.A.
Patricia del Carmen Etchevés
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVES
DIRECTOR TÉCNICO

II. 2. Perfil alto

II. 2. a. 6 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
6x8-well strip - high profile

REF VS-WNV106HE LOT WNV106HE-xxx yyyy-mm

CE IVD   40°C  48

010843544022401317aammdd10WNV106HExxx

VS-1WNV1SH	West Nile Virus 8- well strips	6x8-well strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	6x8-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT
CERTEST BIOTEC S.L Pol. Industrial Rio Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN) VS-WNV106HE v.00

Importador: BIOARS S.A. | Estomba 961 - C.A.B.A.
C.P.: C1427COU | Tel.: (011) 4555 4601 | www.bioars.com.ar
Directora Técnica: Dra. C. Etchevés - Bioquímica M.N. :7028

Venta exclusiva a laboratorios de análisis clínicos
USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496
VS-WNV106HE VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV106HE WNV106HE-xxx DD.MM.AAAA

II. 2. b. 12 tiras de 8 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
12x8-well strip - high profile

REF VS-WNV112HE LOT WNV112HE-xxx yyyy-mm

CE IVD   40°C  96

010843544022403717yyymmdd10WNV112HE000

VS-1WNV1SH	West Nile Virus 8- well strips	12x8-well strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	12x8-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT
CERTEST BIOTEC S.L Pol. Industrial Rio Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN) VS-WNV112HE v.00

Importador: BIOARS S.A. | Estomba 961 - C.A.B.A.
C.P.: C1427COU | Tel.: (011) 4555 4601 | www.bioars.com.ar
Directora Técnica: Dra. C. Etchevés - Bioquímica M.N. :7028

Venta exclusiva a laboratorios de análisis clínicos
USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496
VS-WNV112HE VIASURE West Nile Virus Real Time PCR Detection Kit

VS-WNV112HE WNV112HE-xxx DD.MM.AAAA

II. 2. c. Placa de 96 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
96-well plate - high profile

REF VS-WNV113HE LOT WNV113HE-xxx yyyy-mm

CE IVD   40°C  96

010843544022405117aammdd10WNV113HExxx

VS-1WNV1PH	West Nile Virus 96-well plate	1 plate
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS	Tear-off 8-cap strips	12x8-cap strip


West Nile Virus REAL TIME PCR DETECTION KIT
CERTEST BIOTEC S.L Pol. Industrial Rio Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN) VS-WNV113HE v.00

Importador: BIOARS S.A. | Estomba 961 - C.A.B.A.
C.P.: C1427COU | Tel.: (011) 4555 4601 | www.bioars.com.ar
Directora Técnica: Dra. C. Etchevés - Bioquímica M.N. :7028

Venta exclusiva a laboratorios de análisis clínicos
USO PROFESIONAL EXCLUSIVO - Cert. ANMAT PM-1127-496
VS-WNV113HE VIASURE West Nile Virus Real Time PCR Detection Kit

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II. 3. a. 9 tiras de 4 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
9x4-well strip

REF VS-WNV136E LOT WNV136E-xxx yyyy-mm

IVD CE    2°C 40°C  36

(01)08435440226437(17)yymmdd(10)WNV136Exxx

VS-1WNV1RG	West Nile Virus 4-well strips	9x4 well-strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS-RG	4-cap strips	9x4-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

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VS-WNV136E VIASURE West Nile Virus Real Time PCR Detection Kit

(01)08435440226437(17)yymmdd(10)WNV136Exxx

VS-WNV136E WNV136E-xxx DD-MM-AAAA

II. 2. b. 18 tiras de 4 pocillos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
18x4-well strip

REF VS-WNV172E LOT WNV172E-xxx yyyy-mm

IVD CE    2°C 40°C  72

(01)08435440226444(17)yymmdd(10)WNV172Exxx

VS-1WNV1RG	West Nile Virus 4-well strips	18x4 well-strip
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial
VS-OCS-RG	4-cap strips	18x4-cap strip

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

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VS-WNV172E VIASURE West Nile Virus Real Time PCR Detection Kit

(01)08435440226444(17)yymmdd(10)WNV172Exxx

VS-WNV172E WNV172E-xxx DD-MM-AAAA

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II. 4. a. 4 tubos

Sobre-Rótulo

West Nile Virus
REAL TIME PCR DETECTION KIT
4 tubes x 24 reactions

REF VS-WNV196TE LOT WNV196TE-xxx yyyy-mm

CE IVD    2°C 40°C  96

(01)08435440226413(17)yymmdd(10)WNV196TExxx

VS-1WNV1TB	West Nile Virus Reaction-Mix tube	4 white vials
VS-RB02	Rehydration Buffer	1 blue vial
VS-WNV1C	West Nile Virus Positive control	1 red vial
VS-EC01	Extraction Control	1 green vial
VS-NC1	West Nile Virus Negative control	1 violet vial
VS-H2O	Water RNase/DNase free	1 white vial

West Nile Virus REAL TIME PCR DETECTION KIT

CERTEST BIOTEC S.L. Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840 San Mateo de Gállego, Zaragoza (SPAIN)

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VS-WNV196TE VIASURE West Nile Virus Real Time PCR Detection Kit

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








III. Control positivo

Control positivo – Sobre de aluminio	Control positivo - Tubo
West Nile Virus Positive Control REF VS-WNV1C LOT WNV1C-xxx yyy-mm IVD   2°C 40°C 20  CAUTION Do not open in pre-PCR environment VIASURE	West Nile Virus Positive Control REF VS-WNV1C 20  LOT WNV1C-xxx yyy-mm IVD  2°C 40°C

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


VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-WNV106L
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-WNV106H
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-WNV112L
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-WNV112H
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-WNV113L
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-WNV113H
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-WNV136
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-WNV172

ENGLISH

1. Intended use

VIASURE *West Nile Virus* Real Time PCR Detection Kit is designed for specific identification of West Nile virus in clinical samples from patients with signs and symptoms of West Nile virus infection. This test is intended for use as an aid in the diagnosis of the West Nile virus in combination with clinical and epidemiological risk factors. RNA is extracted from specimens, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for West Nile virus.

2. Summary and Explanation

Since its initial isolation in Uganda in 1937 through the present, West Nile virus (WNV) has propagated to a vast region of the globe and is now considered the most important causative zoonotic pathogen of viral encephalitis worldwide. WNV belongs to the genus *Flavivirus* and is a member of the Japanese encephalitis serogroup. Phylogenetic analysis has revealed two main lineages of WNV, lineage 1, which is widely distributed and further subdivided in three clades (A, B and C) and lineage 2. While several additional lineages have been proposed during the last years.

In nature, there is an enzootic transmission cycle of WNV between birds (hosts) and mosquitoes (bridge vectors) of the genus *Culex* (mainly *Cx. pipiens*, *Cx. perexiguus* and *Cx. modestus*) and *Aedes*. But this virus can infect other vertebrates and humans, between which it has been reported additional methods of transmission: blood transfusion, organ transplantation, transplacental transmission, and via breastmilk. Most human infections with WNV (~80%) are asymptomatic or characterized by a mild-febrile syndrome and flu-like malaise. Nevertheless, fewer than 1% infection may progress to more severe neuroinvasive disease (WNND) and even death, depending on the host immune response and the viral strain. WNND usually encompasses three different syndromes: encephalitis, meningitis, and acute flaccid paralysis.

WNV has a single stranded positive-polarity RNA genome which encodes a single polyprotein that is co- and posttranslationally cleaved into 3 structural proteins and 7 nonstructural proteins. Laboratory diagnosis is generally accomplished by testing of serum or cerebrospinal fluid (CSF) to detect WNV-specific IgM antibodies and confirmed by plaque-reduction neutralization tests (PRNTs). However, the cross-reactivity with infections caused by other flaviviruses remains the major problem with most serological diagnostic tests. In this regard, different PCR-based protocols have been developed to detect viral RNA and can be performed on serum and CSF, both collected early in the course of illness, and additional tissue specimens as urine, where it can be detected much longer and at higher concentrations.

3. Principle of the procedure

VIASURE *West Nile Virus* Real Time PCR Detection Kit is designed for the diagnosis of the West Nile virus in clinical samples. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed

generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved sequence of genomic region 5'UTR using specific primers and a fluorescent-labeled probe.

VIASURE *West Nile Virus* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the RNA target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *West Nile Virus* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPs, buffer, polymerase, retrotranscriptase) in an stabilized format, as well as an internal control to monitor PCR inhibition. West Nile virus RNA targets are amplified and detected in FAM channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).


4. Reagents provided

VIASURE *West Nile Virus* Real Time PCR Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>West Nile Virus</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>West Nile Virus</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *West Nile Virus* Real Time PCR Detection Kit with Ref. VS-WNV106L, VS-WNV106H, VS-WNV112L and VS-WNV112H.

BIOARS S.A.
Patricia del Carmen Etchev s
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCH VES
DIRECTOR TECNICO

Reagent/Material	Description	Color	Amount
West Nile Virus 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
West Nile Virus Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE West Nile Virus Real Time PCR Detection Kit with Ref VS-WNV113L and VS-WNV113H.

Reagent/Material	Description	Colour	Amount
West Nile Virus 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
West Nile Virus Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table 3. Reagents and materials provided in VIASURE West Nile Virus Real Time PCR Detection Kit with Ref. VS-WNV136 and VS-WNV172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).


5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE West Nile Virus Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

BIOARS S.A.
Patricia del Carmen Etchevès
Presidente

VIASURE West Nile Virus Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen)

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVES
DIRECTOR TÉCNICO

SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommended to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-WNV113L, VS-WNV113H, VS-WNV136 and VS-WNV172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-WNV136 and VS-WNV172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. RNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For RNA extraction from clinical samples (serum, CSF, urine, and others) you can use your manually or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- QIAamp Viral RNA Mini Kit, using QIAcube instrument (Qiagen).

8.2. Lyophilized positive control

West Nile Virus Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized West Nile Virus Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.


Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of RNA sample, reconstituted West Nile Virus Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

BIOARS S.A.
Patricia del Carmen Etchevès
Presidente

BIOARS S.A.
BIOCLAUDIA ETCHEVES
DIRECTOR TECNICO

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (West Nile Virus) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for West Nile Virus in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:


West Nile Virus (FAM)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+/-	-	+	West Nile Virus Positive
-	+	-	+	West Nile Virus Negative
-	-	-	+	Experiment fail
+	+	+	-	Experiment fail

Table 5. Sample interpretation

+: Amplification curve

-: No amplification curve


BIOARS S.A.
Patricia del Carmen Etchevès
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVÉS
DIRECTOR TÉCNICO

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated with RNA extracted from serum, and urine samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.

[Signature]
BIOARS S.A.
Patricia del Carmen Etchev s
Presidente

[Signature]
BIOARS S.A.
BIOQ CLAUDIA ETCHEVES
DIRECTOR TECNICO

- There is a possibility of false positive results due to cross-contamination by West Nile Virus either by samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *West Nile Virus* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *West Nile Virus* Real Time PCR Detection Kit was tested using 28 samples (26 serum and 2 urine samples) from symptomatic patients. The results were compared with those obtained by commercial Real Time RT-PCR kit (RealStar® WNV RT-PCR Kit (Altona Diagnostics)). None of 28 samples were positive for West Nile Virus by Real Time PCR assays.

VIASURE *West Nile Virus* Real Time PCR Detection Kit was evaluated with different EQA panels (QCMD and INSTAND). These panels consist of 60 clinical specimens dissolved in transport medium. The results were compared with those presented by EQA programme final reports. All West Nile Virus (44/60) samples could be detected. Samples contained different dilutions of West Nile Virus strains NY99 (lineage 1), Ug37 and Heja (lineage 2). In addition, WNV Negative and Non-WNV flaviviruses samples (which include Japanese Encephalitis virus/Dengue virus types 1-4/ Tick borne Encephalitis virus/ Yellow Fever virus) could be confirmed as a negative.

The results show a high sensitivity and specificity to detect West Nile virus using VIASURE *West Nile Virus* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *West Nile Virus* Real Time PCR Detection Kit has a detection limit of ≥ 10 RNA copies per reaction. (Figure 2).

BIOARS S.A.
Patricia del Carmen Etchev s
Presidente



BIOARS S.A.
BIOQ CLAUDIA ETCHEV S
DIRECTOR TECNICO

Figure 2. Dilution series of West Nile Virus (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).

12.3. Analytical specificity

The specificity of the West Nile virus assay was confirmed by testing a panel consisting of different microorganisms representing the most common arbovirus. No cross-reactivity was detected against any of the following microorganisms tested.

Cross-reactivity testing			
Zika virus strain MR 766	-	Dengue 3 virus strain H87	-
Chikungunya virus strain S27 Petersfield	-	Dengue 4 virus strain H241	-
Dengue 1 virus strain Hawaii	-	St Louis Encephalitis virus strain 17D	-
Dengue 2 virus strain New Guinea C	-	Yellow Fever virus strain 17D	-

Table 6. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE West Nile Virus Real Time PCR Detection Kit was evaluated against West Nile virus H160/99, Heja,Ug37 and NY99 showing positive results.

BIOARS S.A.
Patricia del Carmen Etchev s
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHEV S
DIRECTOR TECNICO

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".

(2) See Annex 3 to check optical measurement exposure setting.

(3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

(4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5) No detection in Cy5 channel.

(6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3


OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTLite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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13. Bibliography/Bibliografía


1. C. Chancey *et al.* The global ecology and epidemiology of West Nile virus. *BioMed Research International* 2015; 2015:376230.
2. A. Rizzoli *et al.* The challenge of West Nile virus in Europe: knowledge gaps and research priorities. *Eurosurveillance* 2015; 20(20): 21135.
3. S. Linke *et al.* Detection of West Nile virus lineages 1 and 2 by real-time PCR. *Journal of Virological Methods* 2007; 146(1-2): 355-358.
4. Centers for Disease Control and Prevention. West Nile Virus (<https://www.cdc.gov/westnile/>).

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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
BIOARS S.A.
Patricia del Carmen Etchevés
Presidente


BIOARS S.A.
BIOQ CLAUDIA ETCHÉVES
DIRECTOR TÉCNICO

- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

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VIASURE *West Nile Virus* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPs, buffer, polymerase, retrotranscriptase) in an stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. *West Nile virus* RNA targets are amplified and detected in FAM channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

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Reagent/Material	Description	Colour	Amount
<i>West Nile Virus</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>West Nile Virus</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table 1. Reagents and materials provided in VIASURE *West Nile Virus* Real Time PCR Detection Kit with Ref. VS-WNV196T.

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- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

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
6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Store the kit at 2-40°C until expiration date stated in the label, except the rehydrated Reaction-Mix tube. Once the West Nile Virus Reaction-Mix tube has been reconstituted, it may be kept at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, we recommend store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).
- The rehydrated Reaction-Mix should be recapped and stored inside the pouch with desiccant material in order to protect from light and humidity.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed.

BIOARS S.A.
Patricia del Carmen Etchevès
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVÈS
DIRECTOR TÉCNICO

- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

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Recommendation is to open and manipulate the West Nile Virus Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of

Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

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- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *West Nile Virus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tube in the thermocycler.

Set up the thermocycler (to check compatibility see Annex 1).

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Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*West Nile Virus*) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *West Nile Virus* in the positive control well. Check Internal

Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

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

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated with RNA extracted from serum, and urine samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by West Nile Virus either by samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

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VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-WNV106LE
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-WNV106HE
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-WNV112LE
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-WNV112HE
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-WNV113LE
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-WNV113HE
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-WNV136E
VIASURE <i>West Nile Virus</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-WNV172E

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VIASURE *West Nile Virus* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in an stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity. West Nile virus RNA targets are amplified and detected in FAM channel and the Extraction control (EC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

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Reagent/Material	Description	Colour	Amount
<i>West Nile Virus</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>West Nile Virus</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *West Nile Virus* Real Time PCR Detection Kit with Ref. VS-WNV106LE, VS-WNV106HE, VS-WNV112LE and VS-WNV112HE.

Reagent/Material	Description	Color	Amount
West Nile Virus 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
West Nile Virus Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE West Nile Virus Real Time PCR Detection Kit with Ref VS-WNV113LE and VS-WNV113HE.

Reagent/Material	Description	Colour	Amount
West Nile Virus 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
West Nile Virus Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip


Table 3. Reagents and materials provided in VIASURE West Nile Virus Real Time PCR Detection Kit with Ref. VS-WNV136E and VS-WNV172E. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

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BIOARS S.A.
Patricia del Carmen Etchevès
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVÉS
DIRECTOR TÉCNICO

VIASURE *West Nile Virus* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

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7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-WNV113LE, VS-WNV113HE, VS-WNV136E and VS-WNV172E). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-WNV136E and VS-WNV172E (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.

- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *West Nile Virus* Positive Control in an area away from the other components.

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If the Extraction Control is used to monitor nucleic acid isolation and as PCR inhibition control, add 5 µL of the EC to the specimen and/or lysis buffer mixture. (clinical specimen, as well as, positive control and/or negative control). Close each tube and vortex for 10 seconds.

If the Extraction Control is used only as a PCR inhibition control, 1 µL of the EC should be added to the reconstituted Reaction-Mix.

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- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) in the reserved wells for negative control.

Add 5 µL of RNA sample in different wells.

Add 5 µL of reconstituted West Nile Virus Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, add 1 µL of the EC (green vial) to the negative control/sample/positive control wells.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

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Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (West Nile Virus) and HEX, JOE or VIC channels (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for West Nile Virus in the positive control well. Check Extraction Control (EC) signal to verify the extraction procedure and correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions. Using the following table read and analyze the results:

West Nile Virus (FAM)	Extraction control (HEX)	Negative Control	Positive Control	Interpretation
+	+/-	-	+	West Nile Virus Positive
-	+	-	+	West Nile Virus Negative
-	-	-	+	Experiment fail
+	+	+	-	Experiment fail

Table 5. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the Extraction Control shows or not an amplification signal. Sometimes, the detection of Extraction Control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the Extraction Control is positive. A failure in the extraction procedure and/or an inhibition of the PCR reaction can be excluded by the amplification of Extraction Control.

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In case of absence of signal in both, extraction control and target in sample wells, we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of nucleic acids purification and/or inhibition.

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Differences can be observed in the values of Ct in the extraction controls between the control samples (negative control and positive control) and the clinical samples, due to the extraction process.

For samples at the limit of detection (LOD) or negative with a value of Ct for the extraction control ≥ 35 , it is recommended to filter the sample, or dilute the extracted (1:10 and 1:100) to avoid possible interference and inhibitors in the amplification reaction. It is recommended to review the instructions for using the extraction process used by the user.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated with RNA extracted from serum, and urine samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected but results may not be reproducible.

- There is a possibility of false positive results due to cross-contamination by West Nile Virus either by samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and West Nile Virus Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.

11. Quality control

VIASURE West Nile Virus Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the Extraction control (EC) in each well confirms the correct performance of the technique.

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The clinical performance of VIASURE West Nile Virus Real Time PCR Detection Kit was tested using 28 samples (26 serum and 2 urine samples) from symptomatic patients. The results were compared with those obtained by commercial Real Time RT-PCR kit (RealStar® WNV RT-PCR Kit (Altona Diagnostics)). None of 28 samples were positive for West Nile Virus by Real Time PCR assays.

VIASURE West Nile Virus Real Time PCR Detection Kit was evaluated with different EQA panels (QCMD and INSTAND). These panels consist of 60 clinical specimens dissolved in transport medium. The results were compared with those presented by EQA programme final reports. All West Nile Virus (44/60) samples could be detected. Samples contained different dilutions of West Nile Virus strains NY99 (lineage 1), Ug37 and Heja (lineage 2). In addition, WNV Negative and Non-WNV flaviviruses samples (which include Japanese Encephalitis virus/Dengue virus types 1-4/ Tick borne Encephalitis virus/ Yellow Fever virus) could be confirmed as a negative.


The results show a high sensitivity and specificity to detect West Nile virus using VIASURE West Nile Virus Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE West Nile Virus Real Time PCR Detection Kit has a detection limit of ≥ 10 RNA copies per reaction. (Figure 2).

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12.4. Analytical reactivity

The reactivity of VIASURE West Nile Virus Real Time PCR Detection Kit was evaluated against West Nile virus H160/99, Heja, Ug37 and NY99 showing positive results.

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
VIASURE *West Nile Virus* Real Time PCR Detection Kit

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generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved sequence of genomic region 5'UTR using specific primers and a fluorescent-labelled probe.

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VIASURE *West Nile Virus* Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in an stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity. West Nile virus RNA targets are amplified and detected in FAM channel and the Extraction control (EC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *West Nile Virus* Real Time PCR Kit includes the following materials and reagents detailed in Table 1.

Reagent/Material	Description	Colour	Amount
<i>West Nile Virus</i> Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>West Nile Virus</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL

Table 1. Reagents and materials provided in VIASURE *West Nile Virus* Real Time PCR Detection Kit with Ref. VS-WNV196TE.

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The following list includes the materials that are required for use but not included in the VIASURE *West Nile Virus* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).

- Filter tips.
- Powder-free disposable gloves.

VIASURE West Nile Virus Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

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

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Store the kit at 2-40°C until expiration date stated in the label, except the rehydrated Reaction-Mix tube. Once the West Nile Virus Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, we recommend store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).
- The rehydrated Reaction-Mix should be recapped and stored inside the pouch with desiccant material in order to protect from light and humidity.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.

BIOARS S.A.
Patricia del Carmen Etchevès
Presidente

BIOARS S.A.
BIOQ CLAUDIA ETCHÉVÉS
DIRECTOR TÉCNICO

- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

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If the Extraction Control is used to monitor nucleic acid isolation and as PCR inhibition control, add 5 µL of the EC to the specimen and/or lysis buffer mixture (clinical specimen, as well as, positive control and/or negative control). The EC must not be added directly to the sample. Close each tube and vortex for 10 seconds.

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For RNA extraction from clinical samples (serum, CSF, urine, and others) you can use your manually or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions for use. We have recommended the following extraction kit:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- QIAamp Viral RNA Mini Kit, using QIAcube instrument (Qiagen).

8.3. Lyophilized positive control

West Nile Virus Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *West Nile Virus* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

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8.4. Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the *West Nile Virus* Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

8.5. PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated *West Nile Virus* Reaction-Mix (white vial) into each tube.

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Add 5 µL of reconstituted *West Nile Virus* Positive Control (red vial) in the reserved wells for positive control.

If the EC is only used as PCR inhibition control, we recommend to add 1 µl of the EC (green vial) to the negative control/sample/positive control wells.

Close the wells with the caps provided. Centrifuge briefly.

Load the plate, the strips or tube in the thermocycler.

3) Set up your thermocycler.

Program your thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

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Table A1/A2. Compatible low and high profile real Time PCR systems.

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- DTprime Real-time Detection Thermal Cyclers (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTLite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

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Referencia: BIOARS S.A. ROTULOS E INSTRUCCIONES DE USO