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For 'in vitro' diagnostic use only. Read carefully before using the product.

INDICATION

It is used for determination of Reverse ABO blood grouping.

INTRODUCTION

Blood grouping is vital before blood transfusion. Before blood transfusion, blood group compatibility between the patient and the donor is one of the things that should be done before blood transfusion. These cells are used to determine the reverse ABO blood group. They are not used for diagnostic purposes and are not specific for the diagnosis of a specific disease. Reverse (reverse) ABO grouping is based on the presence of anti-A and/or anti-B isoagglutinins in serum/plasma, known as A₁ and B isoagglutinins. or A₁, A₂, B, and O Red Blood Cell Reagents.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in Red Blood Cells Reagents are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples.

CONTENT

Across RBC A₁/B and Across RBC A₁/A₂/B/O Cell reagents contain 0.8% suspension of A₁ and B or A₁, A₂, B and O group human red blood cells.

STORAGE CONDITIONS

- Keep at 2-8°C. Prevent contamination after opening vials. Do not freeze cards.
- Stable until the expiry date on the label if stored in a suitable conditions.

Notes:

- Do not use the product if it is detected that the temperature conditions are wrong during the storage and transportation of the product.
- The used vial should be stored at a re-storage temperature.
- The caps of vials should be kept closed to prevent contamination during the period when the product is not used.

MATERIAL REQUIRED FOR TEST BUT NOT PROVIDED IN THE KIT

- Across Auto System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Incubator
- Across System Card Reader
- Across System Card Reader
- Across System Pipette
- Across System Workstation
- Disposable Pipette Tips
- Across Gel Cards

Please refer to the Instructions for Use of Across Gel Cards.

COLLECTION AND PREPARATION OF SAMPLES

Samples

For detection of unexpected antibodies in the serum or plasma;

- Serum or plasma samples collected in tubes with or without anticoagulant can be used. Anticoagulants used in the blood bank should be used as anticoagulants.
- Do not use serum or plasma samples that are hemolyzed, cloudy, contaminated or clot presence.
- Serum or plasma samples stored at 2-8°C up to 48 hours after extraction can be used.
- Plasma/serum samples stored between -20 / -80°C can be used.
- Reagent Red Blood Cells are ready to use.
- Perform necessary centrifuging procedure for clean plasma or serum.

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and materials were used must be controlled.

TEST PROCEDURE

Across RBC A₁/B or Across RBC A₁/A₂/B/O are used with manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments apply the procedure in the user manual.

Allow the samples and reagents to reach room temperature (18-25°C). Check the condition of Red Blood Cell Reagents before use (Do not use cloudy, change of color or haemolysed reagents.) Follow the instructions for use of the Across Gel Cards used. Do not use the card if there are bubbles in the microtubes on the relevant cards, if there is no supernatant, the amount of gel decreases, the gel cracks or changes in the microtube color. Across RBC A₁/B and Across RBC A₁/A₂/B/O test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments manual method is located below.

MANUAL METHOD

1. Paste sample identification onto the card by barcode or writing manually.
2. Open the aluminium foil carefully keeping the card upright to avoid reagent contamination.
3. Homogenize the Red Blood Cell Reagents by gently upside down.
4. Determine needed Across Gel Card microtubes for each serum/plasma to be tested.
5. Add 50 µl of reverse grouping Red Blood Cell Reagent to each microtubes.
6. Add 50 µl of serum or plasma to be tested to each microtubes.
7. Centrifuge card or cards.
8. Read the results.

RESULTS

Result reading:

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

It is recommended that the results are read immediately after centrifuging the cards.

Do not leave processed cards in horizontal position.

Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8 °C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Neutral Mikrotube Across RBC A ₁	Neutral Mikrotube Across RBC A ₂	Neutral Mikrotube Across RBC B	Neutral Mikrotube Across RBC O	ABO Group
0	0	+	0	A
+	+	0	0	B
0	0	0	0	AB
+	+	+	0	0
+	0	+	0	A ₂ (with A ₁ antibody)
+	0	0	0	A ₂ B(with A ₁ antibody)

+ = positive

0= negative

PERFORMANCE CHARACTERISTICS

The studies were carried out by selecting 1030 samples from the random population using Across RBC A₁/B, RBC A₁/A₂/B/O cell products. The results are as follows:

Specificity %	A ₁	A ₂	B	O
Positive Result	100%	98.09%	100%	NA
Negative Result	100%	100%	100%	100%

NA : Not Applicable

WARNINGS AND PRECAUTIONS

- Results cannot be used alone for diagnostic use. They should be evaluated with other clinical and laboratory findings of the patient.
- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Do not use any product with past expiry date.
- If plasma is used, complement dependent haemolytic reactions may not be detected.
- If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residue may trap non-agglutinated Red Blood Cells Reagents at the top of the gel, appearing as a pinkish or reddish layer. In that case, it is recommended to incubate the serum for 10 minutes at 37 °C, centrifuge and repeat the test.
- The reagents are of human origin. Red blood cells used in the preparation of this product obtained from HBsAg, anti-HIV, anti-HCV negative donors. However, all products and samples should be handled as if they were potentially capable of transmitting infectious diseases.
- The presence of high concentrations of proteins, abnormal proteins, medications, plasma volume intensifier, etc. In the serum causes "rouleaux formation", leading to false positive results.
- Since the optimal temperature for isoagglutinins is 4°C. If weak or suspicious reactions are observed, the test is repeated by adding incubation for 15 minutes at 2-8 °C to the test procedure.
- Some of A₂ or A₂B individuals may carry anti-A₁ antibody, leading to discrepancy between direct and reverse groups. While a determination for reverse group is made, it is recommended to also add Across RBC A₂ reagent.
- There can be unexpected antibodies which react with other antigens of Across RBC A₁/B or A₁/A₂/B/O Red Blood Cells. In order to ensure whether they are autoantibody or not, add an autocontrol (the patient's serum tested with the patient's own Red Blood Cells) if autocontrol is negative, check for alloantibodies presence in serum using Across Cell Screen Red Blood Cell reagent.
- In newborns, elderly persons, patients with immune deficiency or very dilute antibody because of plasma changing procedures, reverse grouping reactions maybe weak or missing⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent.

REFERENCES

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 10th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Technical Manuel, 15th edition, American Association of Blood Banks, Bethesda, 2005

PACKAGING

Across RBC A₁/B - Across RBC A₁/A₂/B/O

Ref No:

820101 Across RBC A₁/B 2 X 10 ml


820102 Across RBC A₁/A₂/B/O 4 X 10 ml

Document No: TK3434

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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CE XXXX

INTRODUCTION

To work with samples with known content is a recognized quality control method.

When expected test results are obtained, this means test procedures are properly applied and reagents and equipment give correct results. If unexpected results are observed, giving the patient's test results must be discontinued and the source of erroneous result must be determined. After the problem is found and necessary corrections are made, giving the patient's test results can be proceeded.

INDICATION

Quality Control cells are not used for diagnostic purposes and are not specific to the diagnosis of a specific disease. It is used to verify the performance of materials (cards, devices, procedures applied) in routine tests.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in Red Blood Cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples.

CONTENT

Concentration of red blood cells in each of Across Quality Control whole blood tubes was adjusted to 25%-30%. The kit includes 4 blood samples of 4 ml with details given below.

PRODUCT	BLOOD GROUP	RH, KELL PHENOTYPE				ANTIBODIES	
						Regular	Irregular
QC 1	A	R1R1	D+C+E-c-e+	KEL1 (+)	ANTI B	-	
QC 2	B	R1R2	D+C+E+c-e+	KEL1 (-)	ANTI-A	ANTI-K	
QC 3	AB	rr	D-C-E-c-e+	KEL1 (-)	-	ANTI-D	
QC 4	0	R2R2	D+C-E+c-e-	KEL1 (-)	ANTI-A ANTI-B	-	

* DAT (Direct Coombs) results of all tubes are negative.

** Any discrepant results must be investigated.

STORAGE CONDITIONS

- Keep at 2-8°C. Prevent contamination after opening vials.
- Do not freeze cards.
- Stable until the expiry date on the label if stored in a suitable conditions.

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- Across Auto System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Incubator
- Across System Card Reader
- Across System Pipette
- Across System Workstation
- Across System Dispenser
- Tube for preparing of suspension
- Across Gel Cards
- Disposable pipette tips

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QUALITY CONTROL

Instructions for use must be followed carefully and the equipment used for tests must be checked regularly.

VERONICA E. GIAMPERI
FARMACUTICA
MIL 1615

Valgioni R. Nicolas
Socio Generale
NO MED S.R.L.
30-11697909-3

TEST PROCEDURE

Across Quality Control is used with manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments apply the procedure in the user manual.

Allow the samples and reagents to reach room temperature before use (18-25°C). Check the condition of Red Blood Cell reagents before use (Do not use reagents that display high amounts of haemolysis or cloudy). Across Quality Control test method is the same principle for manual system and Across Auto System automatic or semi-automatic instruments manual method is located below.

MANUAL METHOD

1. Paste sample identification onto the card by barcode or write manually.

2. Open the aluminium foil carefully keeping the card upright to avoid reagent contamination.

Reverse ABO Grouping

- Ensure the resuspension of the Q1, Q2, Q3, Q4 blood samples before use.
- Homogenize the Across RBC A1/B or Across RBC A1/A2/B/O reagents.
- Mark the microtubes of Across Gel Card to be used for testing serum or plasma.
- Pipette 50 µl Across Red Blood Cell reagent into the corresponding microtubes.
- Pipette 50 µl of serum from Across Quality Control to microtubes to be used.

Forward ABO Grouping

- Ensure the resuspension of the Q1, Q2, Q3, Q4 blood samples before use.
- Pipette 10 µl, 5% Red Blood Cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).

Rh / Kell (K)

- Ensure the resuspension of the Q1, Q2, Q3, Q4 blood samples before use.
- Pipette 10 µl, 5% Red Blood Cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).

Indirect Coombs / Antibody Identification

- Ensure the resuspension of the Q1, Q2, Q3, Q4 blood samples before use.
- Homogenize the Across RBC reagents (with or without Papaine).
- Pipette 50 µl of Across RBC reagents to microtubes to be used. (Papaine containing RBC should be added to neutral microtubes, while Papaine free RBC should be added to AHG microtubes).
- Pipette 25 µl of serum from the Across QC product you will be using.
- 3. Incubate for 15 minutes at 37°C.
- 4. Centrifuge card(s).
- 5. Read the results.

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It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position.

Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8°C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results:

ABO System

Forward ABO Group				Reverse ABO Group		
A Microtube	B Microtube	AB Microtube	Ctl. Microtube	Microtube N+ A1 Red Blood Cells Reagent	Microtube N+ B Red Blood Cells Reagent	ABO Group
0	0	0	0	+	+	0
+	0	+	0	0	+	A
0	+	+	0	+	0	B
+	+	+	0	0	0	AB

Rh System (D antigen)

D ⁺ microtube	D ⁺ microtube	Ctl. microtube	Interpretation
+	+	0	D positive
0	0	0	D negative
0	+	0	Weak or partial D
+	0	0	

+ = positive

0 = negative

Across Rh Phenotyping with Kell (K) Card (D⁺, C, E, c, e, C^W, Kell antigens): D⁺, C, E, c, e, C^W, Kell Kell microtubes indicates presence of D⁺, C, E, c, e, C^W, Kell antigens.

Antibody Screening - Indirect Coombs tests: If the results obtained in Neutral and AHG microtubes are positive, presence of antibodies in free-state in serum/plasma is established.

PERFORMANCE CHARACTERISTICS

The total number of tests (n) for determination of repeatability and reproducibility is as follows. Sensitivity and specificity were found to be 100%.

Antigen Parameters: Control for Antibody Reactives

Method	Detection n									
	A	B	D	C	E	c	e	C ^W	K	
Gel	112	112	112	112	112	112	112	112	112	
Tube	56	56	56	56	56	56	56	56	56	

Antibody Parameters: Control for Reverse Grouping Cells

Method	Detection n							
	Anti-A	Anti-B	Anti-AB	No isoagglutinin				
Gel	112	112	112	112				
Tube	56	56	56	56				

Antibody Parameters: Control for Antibody Screening Cell

Method	Detection n	
	Anti-Kell	Anti-D
Gel	112	112
Tube	56	56

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Do not use any product with past expiry date.
- Across Quality Control products are of human origin. Red Blood Cells used in the preparation of this product obtained from HBsAg, anti-HIV, anti-HCV negative donors.
- All blood products and samples should be handled carefully as if they were potentially transmitting.
- Do not use for donor or patients tests except for quality control purposes.
- False positive or false negative results could occur due to contamination of test material, improper reaction temperature and/or time, improper centrifugation, improper storage conditions of material and omission of test reagents.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent.

PACKAGING

Across Quality Control
Ref No: 820110 --- 4x4ml
Document No: TK3436

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IVD	"In Vitro" Diagnostic Medical Device
LOT	Batch Code
Use by date	
Temperature Limit	
Manufacturer	
Consult Instructions for Use	
REF	Catalogue Number
Do not use if package is damaged and consult	
Contains biological material of animal origin	

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