hs-cTnI

High Sensitivity Cardiac Troponin I (CLIA)

Order Information

Catalog No.	Package Size
1001468	2×50 tests
1001469	2×100 tests

Intended Use

The CLIA series high sensitivity cardiac troponin I (hs-cTnI) assay is an automated Chemiluminescent Immunoassay (CLIA) for the quantitative determination of cardiac troponin I (cTnI) in human serum and plasma. This assay can be used as an aid in the diagnosis of myocardial infarction (MI).

Summary

Troponin I (TnI) is a regulatory subunit of the troponin complex which plays a crucial role in regulation of muscle contraction. Troponin I is expressed in three isoforms: fast and slow skeletal isoforms and specific cardiac isoform. The cardiac-specific TnI (cardiac troponin I, cTnI) isoform has a molecular weight of approximately 24,000 Da and contains a post-translational tail of 31 amino acids on the N-terminus of the molecule that is not found in the skeletal isoforms. As a result of its high tissue specificity, cTnI is a cardiac-specific, highly sensitive marker for myocardial injury^{1,2}.

cTnI and cTnT in the blood of patients with AMI present as full-sized form and proteolytic fragments with various lengths, and the major portion of cTnI presents as a binary complex with TnC (IC)³. Several studies showed that at the early stages (approximately first 10 hours after the onset of symptoms) of AMI, cTn is present in blood mainly as a ternary ITC or binary IC complex and later (approximately 20–30 hours after AMI) as a binary complex IC and free fragments of cTnT³.

The International Federation of Clinical Chemistry (IFCC) has issued guidance on high sensitivity cardiac troponin assays. In order to be classified as a high sensitivity assay, two performance criteria should be met: (1) the assay should have analytical imprecision $\leq 10\%$ CV at the 99th

percentile of a healthy population; (2) the assay must be able to measure cTn at or above the limit of detection (LOD) in $\geq 50\%$ of a healthy population⁴. Several studies have demonstrated that high sensitivity cardiac troponin assay can detect cTn release earlier than contemporary cardiac troponin assay and improve sensitivity in the diagnosis of AMI, which may help to "rule in" and "rule out" myocardial infarction earlier⁵.

According to the Fourth Universal Definition of Myocardial Infarction document, myocardial injury is defined by cTn values above 99th percentile, and cTn dynamic change is the key to distinguish acute myocardial injury from chronic myocardial injury. Patients may be diagnosed as myocardial infarction if they are accompanied by clinical evidences of ischemia⁶. In addition, numerous pathologies can potentially cause cTn elevations without overt ischemic heart disease, such as heart failure, diabetes, hypertension, arrhythmias, pulmonary embolism, sepsis, critical illness, myocarditis, stroke, non-cardiac surgery, extreme exercise and drug toxicity. Therefore, it is crucial to combine clinical symptoms, cTn dynamic changes and ischemic evidences to diagnose AMI⁷.

Assay Principles

The CLIA series hs-cTnI assay is a multi-site sandwich assay to quantitatively detect cardiac troponin I (cTnI).

In the first step, sample, pretreatment solution, paramagnetic microparticles coated with anti-cTnI antibodies and anti-cTn-complex antibody-alkaline phosphatase conjugates are added into a reaction cuvette. After incubation, cTn-complex present in the sample binds to both anti-cTnI antibodies coated microparticles and anti-cTn-complex antibody alkaline phosphatase (ALP)-labeled conjugates to form sandwich complexes. Microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, chemiluminescent substrate (sodium 3-((1R,3S,5r,7r)-4'-methoxyspiro[adamantane-2,3'-[1,2]dioxetan]-4'-yl)phenyl phosphate, AMPPD) is added to the reaction cuvette. A phosphate group of AMPPD was removed by ALP to form an unstable intermediate structure, which produces methyl methoxybenzoate anion through intramolecular electron transfer. Chemiluminescence is generated when the methyl

methoxybenzoate anion in the excited state returns to the ground state. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier in the instrument.

The amount of cTnI present in the sample is proportional to the RLUs generated during the reaction. The cTnI concentration can be determined via a calibration curve, which is established on an encoded Calibration Master Curve and three-level product calibrators.

Reagent Components

Components in different lots of reagent kit cannot be interchanged:

Ra	Paramagnetic microparticles coated with monoclonal anti-cTnI antibodies in 50mM TRIS buffer. Minimum concentration: 0.7 g/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.
Rb	Anti-cTn-complex antibody-alkaline phosphatase conjugates in 50mM MES buffer. Minimum concentration: 1.6 mg/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.
Rc	Pretreatment solution in 50mM TRIS buffer. Preservative: 0.05% ProClin 300 and 0.09% sodium azide.

Storage and Stability

The unopened hs-cTnI (CLIA) reagent is stable up to the stated expiration date indicated on the label when stored at 2-8°C.

The hs-cTnI (CLIA) reagent can be stored onboard at 2-8°C and used for a maximum of 28 days after opening for use.

Reagent Preparation

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

Materials Required but not Provided

Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Cat.No.1001470: hs-cTnI Calibrators, C0: 4×0.50 mL, C1: 4×0.40 mL, C2: 4×0.40 mL.

Cat.No.1001471: hs-cTnI Control (L1), 3×1.00 mL.

Cat.No.1001472: hs-cTnI Control (L2), 3×1.00 mL.

Cat.No.1001473: hs-cTnI Control (L3), 3×1.00 mL.

Cat.No.1001474: hs-cTnI Control (L1), 6×1.00 mL.
Cat.No.1001475: hs-cTnI Control (L2), 6×1.00 mL.
Cat.No.1001476: hs-cTnI Control (L3), 6×1.00 mL.
Cat.No.1001139: Substrate Solution, 4×75 mL.
Cat.No.1001138: Substrate Solution, 4×115mL.
Cat.No.1001185: Wash Buffer, 1×10 L.
Cat.No.1001136: Sample Diluent, 6×8 mL.
Cat.No.1001137: Sample Diluent, 2×30 mL.
Reaction Cuvette.

Instrument System

Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Specimen Collection and Preparation

Specimen Types

- Human serum and plasma (sodium heparin, lithium heparin) are recommended for this assay.
- Blood collection tubes from various manufacturers may contain different raw materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Wiener lab. Each laboratory needs to judge the applicability of the tubes and plasma separation products. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Specimen Conditions

Do not use specimen with the following conditions:

- heat inactivated
- grossly hemolyzed
- obvious microbial contamination
- visible fibrin or other cell debris

For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure

residual fibrin and cellular matters have been removed prior to analysis.

- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

- Specimens should be tested as soon as possible after sample collection. If testing is not completed within 8 hours, specimens should be tightly capped and refrigerated at 2-8°C. The sample is stable for 2 days at 2-8°C and 30 days at -20°C.
- Avoid repeated freeze and thaw cycles. The samples may be frozen just once.
- Frozen or refrigerated samples, controls and calibrators must be restored to room temperature before analysis.

Assay Procedure

For optimal performance of this assay, operators should read the related system operator's manual carefully to get enough information such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well.

Before loading the hs-cTnI reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to resuspend the microparticles, which have settled during shipment or storage. Avoid foam formation. Visually inspect the bottle to ensure the microparticles have been well mixed. If the microparticles remain adhered to the bottle, continue inverting until suspended completely. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Wiener lab. Customer Service for help. Avoid inverting opened reagent bottle. Wear clean gloves to avoid contaminating the reagents.

This assay requires 80 µL of sample volume for a single test not including dead volume of the sample container. Additional volume is required when performing additional

tests from the same sample. Operators should refer to the system operator's manual and follow the specific requirements of the assay to determine the minimum sample volume.

Calibration

CLIA series hs-cTnI (CLIA) has been standardized against a commercial hs-cTnI test (CLIA).

The specific information of calibration master curve of hs-cTnI (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used in combination with hs-cTnI calibrators for the calibration of the specific reagent lot. When performing the calibration, please scan the information of calibration master curve from the barcode into the system, then test the three-level calibrators.

Recalibration is required as follows:

1. The previous calibration was performed 28 days before.
2. Quality controls lie outside the expected ranges.
3. A reagent kit with a new lot number is used.

For detailed instruction of calibration, refer to the operator's manual.

Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use and following each calibration. The quality control frequency should be adapted to each laboratory's individual requirements.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system by referring to the operator's manual. If the quality control results are still out of the specified range, please contact Wiener lab. Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the basis of calibration master curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three-level calibrators of defined concentration values. The results are shown in the unit of ng/L.

Conversion factors:	$\text{ng/l} \cdot 0.001 = \text{ng/ml}$
	$\text{ng/l} \cdot 0.001 = \mu\text{g/l}$
	$\text{ng/l} \cdot 1 = \text{pg/ml}$

Measuring Interval

Measuring interval is defined as the range of values in ng/L which meets the limits of acceptable performance for linearity, imprecision, and bias. The analytical measuring interval of the CLIA series hs-cTnI assay is 2.4 to 25,000.0 ng/L. When sample results exceed the measuring interval, refer to dilution.

Specimen Dilution

Samples with concentrations above the upper limit (25,000.0 ng/L) can be diluted with Wiener lab. Sample Diluent. The recommended dilution is 10-fold (either automatically by the analyzer or manually). After automated dilution by the analyzer, the system automatically multiplies the result by the dilution factor when calculating the sample concentration. After manual dilution, multiply the result by the dilution factor. The concentration must be >100.0 ng/L before the dilution factor is applied.

Expected values

A reference range study was conducted according to the guidance from Clinical and Laboratory Standards Institute (CLSI) document EP28-A3⁸. Serum specimens were collected from 832 apparently healthy adults with normal levels of NT-proBNP, HbA1c and estimated GFR values. The 99th percentile upper reference limits (URLs) of CLIA series hs-cTnI assay were determined using the non-parametric statistical method described in CLSI document EP28-A3⁸.

Population	N	Age Range (years)	99th URL (ng/L)	95% CI (ng/L)
Females	407 ^a	21-87	15.3	8.6-29.2
Males	424	21-83	31.3	17.8-47.3
Overall	832	21-87	24.2	16.6-43.5

^aDuring the gender-specific analysis, one result from female subject was identified as outliers using the Dixon method, but it was included in the overall analysis.

The 99th percentile URLs also apply to sodium heparin and lithium heparin plasma samples. The CLIA series hs-cTnI assay has an analytical precision of $\leq 10\%$ CV at the 99th

percentile URL and can measure cTnI concentrations above the limit of detection (LoD) in $\geq 50\%$ of healthy population, which met the criteria of high sensitivity assays stated by IFCC guidance.

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

Limitation

The measuring range of reagent kit is 2.4-25,000.0 ng/L. A specimen with hs-cTnI concentration lower than the lower limit will be reported as < 2.4 ng/L. A specimen with hs-cTnI concentration higher than the upper limit will be reported as $> 25,000.0$ ng/L or repeat the test after dilution.

Heterophilic antibodies in samples can react with immunoglobulin in the kits and interfere with the immunoassay results⁹⁻¹⁰. Patients who have been regularly exposed to animals or have been treated with immunoglobulins or immunoglobulin fragments may produce such antibodies, resulting in abnormal results. For example, specimens from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values of hs-cTnI. Additional clinical or diagnostic information may be required to determine patient status.

Autoantibodies may react with cardiac troponin I or T and the presence of autoantibodies may result in the observation of depressed or elevated troponin values¹¹⁻¹². Cardiac troponin results should always be used in conjunction with other clinical information.

Although the CLIA series hs-cTnI assay is specifically designed to minimize the effects of heterophilic antibodies, autoantibodies and rheumatoid factor (RF), there may still be interference problems in a few samples. Assay results that are not consistent with other clinical observations may require additional clinical examinations or other relevant information for diagnosis.

The hs-cTnI result in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity.

Performance Characteristics

Lower limits of measurement

The Limit of Blank (LoB), Limit of Detection (LoD) and Limit of Quantitation (LoQ) were determined in accordance with the CLSI document EP17-A2 requirements¹³. Testing was conducted using 2 lots of the hs-cTnI reagent kit on several instruments over 4 days.

The Limit of Blank, Limit of Detection and Limit of Quantitation values are summarized below.

	Range (ng/L)	Maximum (ng/L)
LoB ^a	0.1-0.5	0.5
LoD ^b	0.5-0.7	0.7
LoQ ^c	1.1-2.4	2.4

^aThe LoB is the 95th percentile value from $n \geq 60$ replicates of analyte-free samples.

^bThe LoD is the lowest concentration which the analyte can be detected with 95% probability based on $n \geq 60$ replicates of low-analyte level samples.

^cThe LoQ was determined from $n \geq 60$ replicates of low-analyte level samples and is defined as the lowest analyte concentration that can be reproducibly measured with a total error of $\leq 20\%$ CV. In addition, the lowest analyte concentration that met a total error of $\leq 10\%$ CV ranged from 1.1 to 3.9 ng/L.

Precision

Precision was determined by following CLSI document EP05-A3¹⁴. Three controls and three serum samples were tested 2 runs each day in duplicate for totally 20 days ($n = 80$). Results are shown in the table below* (criterion: CV $\leq 10\%$).

Sample	Mean	Within-Run	Within-Laboratory
	(ng/L)	CV	CV
QC-L1	10.9	2.01%	2.97%
QC-L2	36.4	1.84%	2.49%
QC-L3	10087.0	1.61%	1.97%
Sample 1	15.8	2.92%	2.90%
Sample 2	111.9	2.52%	2.77%

Sample	Mean	Within-Run	Within-Laboratory
	(ng/L)	CV	CV
Sample 3	10398.5	1.23%	1.97%

**Representative data; results in individual laboratories may vary.*

Linearity

Linearity study was performed based on guidance from CLSI document EP06-Ed2¹⁵. A high concentration hs-cTnI sample (> 25,000.0 ng/L) was mixed with a low concentration sample (< 2.4 ng/L) at different ratios, generating a series of dilutions. The Linearity of hs-cTnI was demonstrated in the range of 2.4 to 25,000.0 ng/L, the correlation coefficient $r \geq 0.9900$.

Trueness

Two samples defined to the reference reagent were tested to verify the accuracy of hs-cTnI assay. Results showed relative deviation within $\pm 10\%$ of defined value.

Sample	Measured Value (ng/L)	Target Value (ng/L)	Relative Deviation
Trueness control L	97.9	100.6	-2.68%
Trueness control H	9919.0	10354.6	-4.21%

High Dose Hook

For the CLIA series hs-cTnI assay, no high dose hook effect was observed when samples contain up to approximately 1,000,000.0 ng/L of cTnI.

Analytical Specificity

The CLIA series hs-cTnI assay has an analytical specificity of < 0.1% cross-reactivity with the following substances.

Substance	Concentration
Skeletal troponin I	1,000 ng/mL
Skeletal troponin T	1,000 ng/mL
Cardiac troponin T	1,000 ng/mL
Troponin C	1,000 ng/mL
Actin	1,000 ng/mL
Tropomyosin	1,000 ng/mL

Substance	Concentration
Myosin light chain	1,000 ng/mL
Myoglobin	1,000 ng/mL
CK-MB	1,000 ng/mL

Interference

The test results of CLIA series hs-cTnI assay are not observed to be interfered with the potential endogenous substances up to the following concentrations (criterion: recovery within $\pm 10\%$ of initial value)*.

Substance	Concentration
Triglycerides	3,000 mg/dL
Cholesterol	500 mg/dL
Bilirubin (conjugated)	40 mg/dL
Bilirubin (unconjugated)	60 mg/dL
Fibrinogen	490 mg/dL
Hemoglobin	608 mg/dL
Protein gamma globulin	2.5 g/dL
Total protein	14.5 g/dL
Protein albumin	7 g/dL
Rheumatoid factors	1,200 IU/mL

*Representative data; results in individual laboratories may vary.

In vitro tests were performed on 73 commonly used pharmaceuticals. No interference with the assay was found.

Method comparison

The CLIA series hs-cTnI Assay was compared to a commercially available diagnostic kit in a correlation study with about 560 specimens. The statistical data obtained by linear regression are shown in the table below*.

Concentration Range (ng/L)	Slope	Intercept	Correlation Coefficient
2.4-25,000.0	0.9079	59.84	0.973

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.

3. Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Wiener lab. system, or using Wiener lab. reagent kits on other systems.
4. Do not use reagent kits beyond the expiration date.
5. Do not use reagents mixed from different reagent lots.
6. Avoid bubble formation in all reagents. Wear clean gloves to avoid contaminating the reagents.
7. Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
8. Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. Inappropriate storage will lead to incorrect results.
9. If the reagents are opened unintentionally before use, they shall be used as soon as possible.
10. Reagent pack opened for more than 28 days is not recommended for use.
11. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
12. All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
13. The Material Safety Data Sheet (MSDS) is available upon request.
14. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
15. Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
16. Do not freeze. The results can't be assured when the reagents are stored at inappropriate condition.
17. Swallowing or getting the reagents on yourself must be avoided. Wash hands thoroughly after handling. The device must be used in a room with adequate ventilation. The contaminated clothing must be removed immediately.
18. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

**Warning**

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Prevention:

P280 Wear protective gloves and protective clothing.

P261 Avoid breathing mist/vapours/spray.

P273 Avoid release to the environment.

P272 Contaminated work clothing should not be allowed out of the workplace.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.

Graphical Symbols

In vitro
diagnostic
medical
device

Authorized
representative
in the
European
Community

Consult
instructions
for use

European
Conformity

Catalogue
number

Temperature
limit

Manufacturer

Use-by
date

Caution

Batch code

This way up

Unique
device
identifier

Contents

References

1. Katrukha IA. Human cardiac troponin complex. Structure and functions. *Biochemistry (Mosc)* 2013; 78: 1447-1465.
2. Perry SV. The regulation of contractile activity in muscle. *Biochem Soc Trans* 1979; 7: 593-617.
3. Katrukha IA, Katrukha AG. Myocardial Injury and the Release of Troponins I and T in the Blood of Patients. *Clin Chem.* 2021; 67:124-130.
4. Apple FS, Jaffe AS, Collinson P, *et al.* IFCC educational materials on selected analytical and clinical applications of high sensitivity cardiac troponin assays. *Clin Biochem* 2015; 48: 201-203.
5. Roffi M, Patrono C, Collet JP, *et al.* 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation. *European Heart Journal* 2016; 37: 267-315.
6. Thygesen K, Alpert JS, Jaffe AS, *et al.* Fourth universal definition of myocardial infarction (2018). *Circulation* 2018; 138: e618-e651.
7. Members of the cardiovascular specialist branch of the Chinese Medical Association. The consensus of Chinese experts on laboratory detection and clinical application of cardiac troponin. *Chinese Medical Journal* 2021; 101: 2947-2961.
8. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved

- Guideline – Third Edition. CLSI Document EP28-A3. Wayne, PA: CLSI; 2008.
9. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-885.
 10. Heterophilic antibodies; A problem for all immunoassays. Clin Chem 1988;34(1):27-33.
 11. Prevalence of autoantibodies to cardiac troponin T in healthy blood donors. Clin Chem 2009; 55(8):1592-1593.
 12. Troponin autoantibodies: from assay interferent to mediator of cardiotoxicity. Clin Chem 2017;63(1):30-32.
 13. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
 14. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline - Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
 15. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Linearity of Quantitative Measurement Procedures. Second Edition. CLSI guideline EP06. CLSI; 2020.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Contenido

Referencias

1. Katrukha IA. Human cardiac troponin complex. Structure and functions. *Biochemistry (Mosc)* 2013; 78: 1447-1465.
2. Perry SV. The regulation of contractile activity in muscle. *Biochem Soc Trans* 1979; 7: 593-617.
3. Katrukha IA, Katrukha AG. Myocardial Injury and the Release of Troponins I and T in the Blood of Patients. *Clin Chem*. 2021; 67:124-130.
4. Apple FS, Jaffe AS, Collinson P, et al. IFCC educational materials on selected analytical and clinical applications of high sensitivity cardiac troponin assays. *Clin Biochem* 2015; 48: 201-203.
5. Roffi M, Patrono C, Collet JP, et al. 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation. *European Heart Journal* 2016; 37: 267-315.
6. Thygesen K, Alpert JS, Jaffe AS, et al. Fourth universal definition of myocardial infarction (2018). *Circulation* 2018; 138: e618-e651.
7. Members of the cardiovascular specialist branch of the Chinese Medical Association. The consensus of Chinese experts on laboratory detection and clinical application of cardiac troponin. *Chinese Medical Journal* 2021; 101: 2947-2961.
8. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI Document EP28-A3. Wayne, PA: CLSI; 2008.
9. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885.
10. Heterophilic antibodies; A problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.

11. Prevalence of autoantibodies to cardiac troponin T in healthy blood donors. Clin Chem 2009; 55(8):1592-1593.
12. Troponin autoantibodies: from assay interferent to mediator of cardiotoxicity. Clin Chem 2017;63(1):30-32.
13. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
14. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline - Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
15. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Linearity of Quantitative Measurement Procedures. Second Edition. CLSI guideline EP06. CLSI; 2020.

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
Intended Use

Wiener lab. hs-cTnI Calibrators (hs-cTnI Cal) are intended to calibrate the quantitative CLIA series High Sensitivity Cardiac Troponin I (hs-cTnI) assay on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Summary

Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, Wiener lab. reagent kits, calibrators and controls. hs-cTnI Calibrators contain cardiac troponin I-C. The mathematical relationship between the measured responses and the known analyte concentration establishes the calibration curve. It is used to convert Relative Light Units (RLUs) measured from samples to quantitative analyte concentration.

Components

C0, C1, C2	Cardiac Troponin I-C in buffered bovine serum albumin (BSA) at three concentration levels with 0.05% ProClin 300 and 0.02% BND as preservatives.
------------	--

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Laboratory professionals must be trained to take all necessary precautions while handling laboratory reagents. Do not use reagents beyond the expiration date.
- The calibrators' values are lot-specific with the matched set listed in the Calibration Card.
- Calibration is recommended after each reagent lot change, special maintenance, or troubleshooting procedure.
- Do not mix calibrators from different lots.
- Please perform maintenance and standard operation as scheduled, including calibration and analysis to assure the performance of the measurement system.
- Recalibrate the system when controls are out of specified range.
- The calibrators were tested with CE-approved kits for HIV-1/2 Antibody, HCV Antibody, Treponema pallidum Antibody and HBsAg. All results are negative. However, as no testing method can rule out the potential

risk of infection with absolute certainty, this material should be handled as potential biohazards.¹

- If there is any evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- Material safety data sheet is available upon request.
- Disposal of all waste material should be in accordance with local guidelines.
- Please confirm the integrity of the package before use. Do not use the calibrators with damaged packages. Inappropriate storage will lead to incorrect results.
- If the calibrators are opened unintentionally before use, they shall be tightly capped and stored at 2-8°C, and then used within in-use stability time.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Warning
H317 May cause an allergic skin reaction.

Prevention:
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response:
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
P302 + P352 IF ON SKIN: Wash with plenty of water.
P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:
P501 Dispose of contents/ container in accordance with local regulation.

Testing Method

- Refer to Wiener lab. Operator's Manual and / or customer service to acquire information of calibration principle, use of calibrators, and validity of data of calibrators.
- Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation.
- Because of possible evaporation effects, the calibrators on the analyzers should be measured as soon as possible. One vial of

- calibrator can be used only once.
- Three levels of hs-cTnI calibrators are provided, which are 0, ~100 and ~15000 ng/L. Test C0 three times and C1, C2 twice for calibration.
 - Refer to calibration card for accurate information about lot-specific hs-cTnI concentration.

Preparation

- Take the vial out of refrigerator and let it stand for no less than 30 minutes at room temperature.
- According to the calibration procedure described in the instrument manual, use the Calibrators' lot-specific values to set calibration parameters.

Storage and Stability

The calibrators are stable up to the expiration date as indicated on the label when stored in unopened vial at 2-8°C. Once opened, it is stable for 30 days at 2-8°C.

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- Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, reagent kits and controls.
- General laboratory equipment.

Applicable Instrument

Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Traceability of Calibrator

The traceability process is based on EN ISO 17511:2021², the measurand (analyte) in the calibrators is traceable to a commercial hs-cTnI test. The concentrations of calibrators are lot-specific on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, and cannot be used on the other system or method, as bias may exist among different systems and methods. Details about traceability are available upon request.

Quality Control

Use Wiener lab. controls as the samples to monitor the performance of system. The control result should be within the defined range as shown in the control value sheet. If a control is out of specified range, the measurement system and reagents should be checked. For instance, expiration date or storage condition of the calibrators, reagents and controls; contamination of reagents; position of reagents or samples placed on the analyzer; the parameter setting of software; or performance of the analyzer. The assigned values of controls are lot-specific on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer.

Graphical Symbols

In vitro diagnostic medical device

Authorized representative in the European Community

Consult instructions for use

Catalogue number

Temperature limit

Manufacturer

Biological risks

Batch code

This way up

Caution

European Conformity

Use-by date

Unique device identifier

Contents

References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
- EN ISO 17511: 2021 In vitro diagnostic medical devices - Requirements for establishing metrological traceability of values assigned to calibrators, trueness control materials and human samples.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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









Intended Use

Wiener lab. hs-cTnI Control is used for quality control by monitoring the precision and the detection of systematic analytical deviations of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer and test ability of clinical laboratory internally in the quantitative measurement of cardiac troponin I analyte.

Summary

Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, Wiener lab. reagent kits, calibrators, and controls. Wiener lab. hs-cTnI Control can be used to check whether quality control (QC) result is within specified range or not. It is used with CLIA series High Sensitivity Cardiac Troponin I (CLIA).

Components

hs-cTnI Control contains Cardiac Troponin I-C analyte in buffered bovine serum albumin (BSA) with 0.02% BND and 0.05% ProClin 300 as the preservative.

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Laboratory professionals must be trained to take all necessary precautions while handling laboratory reagents. Do not use reagents beyond the expiration date.
- The control values are lot-specific with the matched set listed in the value sheet.
- The controls should be run daily in parallel with the patient samples, after each calibration, or after the reagent lot has been changed.
- Please perform maintenance and standard operation as scheduled, including calibration and system analysis to assure the performance of measurement.
- The controls were tested with CE-approved kits for HIV-1/2 Antibody, HCV Antibody, Treponema pallidum Antibody and HBsAg.

The results are all negative. However, as no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazard.¹

- If there is any evidence of microbial contamination or excessive turbidity in the product, contact Wiener lab. Customer Service for help.
- Material safety data sheet is available upon request.
- Disposal of all waste material should be in accordance with local guidelines.
- The target values and ranges of Wiener lab. hs-cTnI Control are for reference only. Each laboratory can establish its own internal quality control procedure and set up its own target values and ranges.
- Please confirm the integrity of the package before use. Do not use the controls with damaged packages. Inappropriate storage will lead to incorrect results.
- If the controls are opened unintentionally before use, they shall be tightly capped and stored at 2-8°C, and then used within in-use stability time.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

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Prevention:

P261 Avoid breathing dust/fume/ gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response:

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
P302 + P352 IF ON SKIN: Wash with plenty of water.
P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/ container in accordance with local regulation.

Test Method

hs-cTnI Controls have similar components as test samples. Therefore, the controls can be tested by following the same procedure as tested samples, the target values and ranges are then calculated and compared to see if the controls meet specifications.

Preparation

- Refer to Wiener lab. System Operator's Manual and / or Customer Service to acquire information about control principle, control use, and control validity evaluation.
- Take the vial out of refrigerator, and let it stand for no less than 30 minutes at room temperature.
- Mix contents by gently inverting the vial at least 10 times before use. Avoid bubble formation.
- Because of possible evaporation effects, the control on the analyzers should be measured as soon as possible. After use, close the bottle immediately and store upright at 2-8°C.
- Run QC test you have requested or have it analyzed in the same way as patient samples.

Storage and Stability

hs-cTnI Control is stable up to the expiration date as indicated on the label when stored in unopened vial at 2-8°C. Once opened, it is stable for 30 days at 2-8°C.

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- Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer, reagent kits, and matched calibrators.
- General laboratory equipment.

Applicable Instrument

Wiener lab. CLIA series Chemiluminescence Immunoassay.

Assay Values

The control values (the target value and range) determined by Wiener lab. routine measurement procedure are listed in the target value sheet. The target value was obtained from Wiener lab. measurement system and the range was calculated as the target values $\pm$ 3 standard deviations. The control values are lot-specific. Please check the lot number of each control before use.

Quality Control

Use Wiener lab. controls as the samples to monitor the performance of system. The values of controls are lot-specific on Wiener lab. CLIA series Chemiluminescence Immunoassay Analyzer. The control results should be within the specified range as shown in the control target value sheet. If the controls are out of specified range, the measurement system and reagents should be checked, for instance, expiration date or storage condition of the calibrators, reagents and controls; possible contamination of reagents; position of reagents, samples placement on the analyzer; the parameter setting of software; or operation status of the analyzer. Each laboratory should establish its own internal quality control scheme and procedure for corrective and preventive actions.

Graphical Symbols

References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R. China

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WIENER Laboratorios S.A.I.C.
Referencia: ROTULOS E INSTRUCCIONES DE USO WIENER Laboratorios S.A.I.C.