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
BOY: 60,5mm

EN: 90mm

Handwritten signatures and stamps:

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

Valgothi R. Nicolas
Socio Gerente
NG MED S.R.L.
30-71687909-3

INTRODUCTION

The ABO system was the first human blood group system discovered by Landsteiner⁽¹⁾, and is still the most important system in transfusion practice.

The determination of Rh (D) is defined by the presence or absence of the D (RH1) antigen in the red blood cells. Anti-A, Anti-B and Anti-D reagents are used in ABO and Rh blood group typing.

INDICATION

The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping). Across Gel Double ABO / D^{VI} card is used for determination of A (ABO1), B (ABO2) and D (RH1) antigens.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in red blood cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples. Across Gel card contains polymerized dextran microspheres (gel) which function as a filter within the buffer solution in all microtubes on the card. Dextran and specific antibody containing reagent(s) or buffer solution are mixed. During the centrifugation of the cards, red blood cells agglutinates collected on the gel surface and throughout the length of the column according to their sizes. Non-agglutinated red blood cells go to the bottom of the microtube.

CONTENT

Microtubes of Across Gel Double ABO / D^{VI} card are identified by the front labels of the card:

- **A (ABO1) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- **B (ABO2) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- **D^{VI} (RH1) Microtube:** Monoclonal Anti-D (mixture of IgM and IgG antibodies, RUM-1, P3x61, ESD1M clone and MS-26)
- **Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

STORAGE CONDITIONS

Store cards upright at 2-25°C. Do not place cards in front of and/or next to any source of heat (places under direct sunlight, radiators, air conditioners, big devices which radiate heat when operating) or air ducts. Do not freeze cards. When the box is kept in intact package, it is stable until the expiry date specified on the outer label.

MATERIAL REQUIRED FOR TEST BUT NOT PROVIDED IN THE KIT

- Across Auto System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Reader
- Across System Dispenser
- Across System Pipette
- Across System Workstation
- Across LISS
- Disposable pipette tips
- Tube for preparation of suspension

COLLECTION AND PREPARATION OF SAMPLES

Samples

Blood samples with anticoagulant (containing plasma) can be used.
Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

Determination of the antigens of ABO/Rh System:

- Fresh blood samples collected in tubes containing anticoagulants like EDTA or citrate can be used.
- Blood samples stored at 2-8°C up to 48 hours after extraction can be used.
- Blood samples stored at 2-8°C, collected from bags containing ACD, CPD, CPDA and SAG-Mannitol which are not beyond the expiry date specified on the label must be used.

Sample Preparation

a) Preparing erythrocyte (red blood cell) suspension (To detect ABO/Rh system antigens)

Prepare 5% erythrocyte suspension (Add 25 µl erythrocyte to 0.5 ml Across LISS).

Valgini R. Nicolas
Socio Gerente
NG MED S.R.L.
30-7157908-3

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

TEST PROCEDURE

Across Gel Double ABO / D^{VI} card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18-25 °C). Check card condition before use (Do not use the card if you detect bubbles in microtube, lack of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Double ABO / D^{VI} test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments see the instrument manual.

MANUAL METHOD

ABO / D^{VI} Antigen Determination (A-B-D^{VI}-Ctl)

1. Each card is individually identified by a barcode. Unless a barcode reader is used identify them manually.
2. Carefully peel off the aluminium foil that covers the microtubes and keep the card upright to avoid reagent cross contamination.
3. Ensure the resuspension of the red blood cells before use.
4. Pipette 10 µl 5% red blood cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).
5. Centrifuge card or cards for 9 mins at 990 rpm.
6. Read the results.

RESULTS

Results Reading

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position.

Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8°C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

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+ = positive
0 = negative

Notes:

- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In ABO/Rh system: In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- D^{VI} microtube may not react with all D variants so negative reactions with the D^{VI} microtube must be verified using other reagents and techniques which may detect different variants of the D antigen.

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Double ABO/D^{VI} card have been analyzed with representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-D ^{VI}	1061	98.5%	100%

(I) Sensitivity: [(number of true positive results) / (number of true positive results + number of false negative results)] x 100

(II) Specificity: [(number of true negative results) / (number of true negative results + number of false positive results)] x 100

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative or a false positive test result.
- Adding a volume other than that stated in the method may change the reaction.
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens, and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high gold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent.

REFERENCES

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 10th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005

PACKAGING


Across Gel Double ABO / D^{VI}

(2x A / B / D^{VI} / Ctl)

Ref No: 810206 ----50 Cards

Document No:TK3414

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or

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INDICATION

The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) Across Gel Double ABO / D^{Vi} is card used for determination of A (ABO1), B (ABO2) and D (RH1) antigens.

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CONTENT

Microtubes of Across Gel Double ABO / D^{Vi} card are identified by the front labels of the card:

- **A(ABO1) Microtube:** Monoclonal Anti-A (Mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- **B(ABO2) Microtube:** Monoclonal Anti-B (Mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- **D^{Vi}(RH1) Microtube:** Monoclonal Anti-D (Mixture of IgM antibodies of human origin, RUM-1, P3x61)
- **Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

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- Across Auto System semi-automatic or automatic instruments.
- Across System Card Centrifuge
- Across System Card Reader
- Across System Dispenser
- Across System Pipette
- Across System Workstation
- Across LISS
- Disposable pipette tips
- Tube for preparation of suspension

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Samples

Blood samples with anticoagulant (containing plasma) can be used.

Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

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TEST PROCEDURE

Across Gel Double ABO / D^{Vi} is used with manual system and Across Auto System semi-automatic or automatic instruments. Allow the samples and reagents to reach room temperature (18- 25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, lack of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Double ABO / D^{Vi} test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments see the instrument manual.

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Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8 °C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

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0 = negative

Notes:

- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In ABO/Rh system: In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- D^{Vi} microtube may not react with all D variants so negative reactions with the D^{Vi} microtube must be verified using other reagents and techniques which may detect different variants of the D antigen.

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www.diapro.com.tr

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Double ABO / D^{Vi} card have been analyzed with representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-D ^{Vi}	1061	98,5%	100%

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WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative or a false positive test result.
- Adding a volume other than that stated in the method may change the reaction
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens, and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized

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
Ref No: 810207 ---- 50 Cards

Document No: TK3415

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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