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

Therefore, while classification of ABO system, both forward and reverse grouping must be done. In ABO forward grouping, red blood cells are tested with anti-A, anti-B, anti-AB antibodies, whereas in reverse grouping, ABO antibodies are determined in the serum/plasma using A1 and B red blood cells Anti D reagent is used for the determination of presence or absence of D (RH1) antigen on the surface of red blood cell.

INTENDED PURPOSE

Blood grouping is life sustaining before blood transfusion. The blood group compatibility between the patient and the donor before the blood transfusion is one of the things that should be done before the blood transfusion. This card is also used for this purpose. It is not used for diagnostic purposes and is not specific to the diagnosis of a specific disease. The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) and by the presence of antibodies in the serum corresponding to the antigen or antigens missing on the surface of red blood cells (reverse grouping). Across Gel Forward & Reverse ABO with D^{VI-}/D^{VI+} card is used for determination of A (ABO1), B (ABO2), D (RH1) antigens and of possible anti-A and/or anti-B antibodies.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in red blood cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples. Across Gel card contains polymerized dextran microspheres (gel) which function as a filter within the buffer solution in all microtubes on the card. Dextran and specific antibody-containing reagent(s) or buffer solution are mixed. During the centrifugation of the cards, red blood cells agglutinate collected on the gel surface and throughout the length of the column according to their sizes. Non-agglutinated red blood cells go to the bottom of the microtube. The reaction score is determined according to the location of the collapsed cells.

CONTENT

Microtubes of Across Gel Forward & Reverse ABO with D^{VI-}/D^{VI+} card are identified by the front labels of the card:

- **A (ABO1) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- **B (ABO2) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- **AB (ABO3) Microtube:** Monoclonal Anti-AB (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15, LB-2, 9621A8 clone)
- **D^{VI-}(RH1) Microtube:** Monoclonal Anti-D (mixture of IgM antibodies of human origin, RUM-1, P3x61)
- **D^{VI+}(RH1) Microtube:** Monoclonal Anti-D (mixture of IgM and IgG antibodies of human origin, clones RUM-1, P3x61, ESD1M, MS-26)
- **Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- **N Microtubes:** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

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TEST PROCEDURE

Across Gel Forward & Reverse ABO with D^{VI-}/D^{VI+} card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18-25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, absence of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Forward & Reverse ABO with D^{VI-}/D^{VI+} test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments see the instrument manual.

MANUAL METHOD

- Each card is individually identified by a barcode. Unless a barcode reader is used identify them manually.
- Carefully peel off the aluminium foil that covers the microtubes and keep the card upright to avoid reagent cross contamination.
- Reverse ABO Grouping (N_{A1} - N_B)**
- Homogenise Across Red Blood Cell reagents (A₁/B₁).
- Pipette 50 µl Across Red Blood Cell reagents (A₁/B₁) into the N_{A1} / N_B microtubes.
- Pipette 50 µl serum or plasma to microtube N_{A1} / N_B.
- Forward ABO Grouping + D^{VI-}/D^{VI+} Antigen Determination (A-B-AB-D^{VI-}-D^{VI+}-Ctl.)**
- Ensure the resuspension of the red blood cells before use.
- Pipette 10 µl 5% red blood cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).
- Centrifuge card or cards for 9 mins at 990 rpm.
- Read the results.

RESULTS

Results Reading

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scare small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

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0 = negative

Notes:

- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The meaning of Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In case of any discrepancy between forward-reverse groups, investigation must be made before giving the samples' test results (patient's age, clinical information, history, etc.).
- In ABO/Rh system : In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- D^{VI-}/D^{VI+} microtubes does not react with all D variants so negative reactions with the D^{VI-}/D^{VI+} microtubes must be verified using other reagents and techniques which may detect different variants of the D antigen.
- Any discrepancy between results obtained in D^{VI-} and D^{VI+} should be interpreted as weak or partial D antigen. Analyzing antigen expression is recommended.

- The O group cells should obtain negative results. Positive results may be caused by unexpected irregular antibodies reacting at room temperature with the corresponding antigens. If positive results are obtained, invalidate the test and proceed to investigate the discrepancy.
- Some individuals with A2 or A2B subgroups may present antibodies anti-A1, where a forward-reverse discrepancy can be seen.

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Forward & Reverse ABO with D^{VI-}/D^{VI+} card have been studied in a representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-AB	1061	100%	100%
anti-D ^{VI-}	1061	98,5%	100%
anti-D ^{VI+}	1061	98,5%	100%

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Reverse ABO group

There was compatibility between the results of forward and reverse grouping tests of Across Gel Forward & Reverse ABO with D^{VI-}/D^{VI+} card.

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Adding volumes other than those specified in the method can alter the reaction.
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent.

REFERENCES

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 10th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005

PACKAGING

Across Gel Forward & Reverse ABO With D^{VI-} / D^{VI+}

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INTENDED PURPOSE

Blood grouping is life sustaining before blood transfusion. The blood group compatibility between the patient and the donor before the blood transfusion is one of the things that should be done before the blood transfusion. This card is also used for this purpose. It is not used for diagnostic purposes and is not specific to the diagnosis of a specific disease. The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) and by the presence of antibodies in the serum corresponding to the antigen or antigens missing on the surface of red blood cells (reverse grouping). Across Gel Forward & Reverse ABO with D^{VI}- / CDE card is used for determination of (ABO1), B (ABO2), D (RH1), C (RH2) and E (RH3) antigens and of possible anti-A and/or anti-B antibodies.

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CONTENT

Microtubes of Across Gel Forward&Reverse ABO with D^{VI}- / CDE card are identified by the front labels of the card:

- A (ABO1) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- B (ABO2) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- AB (ABO3) Microtube:** Monoclonal Anti-AB (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15, LB-2, 9621A8 clone)
- D^{VI}-(RH1) Microtube:** Monoclonal Anti-D (mixture of IgM antibodies of human origin, RUM 1, P3x61)
- CDE (RH1,2,3) Microtubes:** Monoclonal Anti-CDE (mixture of IgM/IgG antibodies of human origin, clones MS-24, MS-201, MS-26, MS-80)
- Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- N Microtube:** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

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Prepare 5% erythrocyte suspension (Add 25 µl erythrocyte to 0.5 ml Across LISS)

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TEST PROCEDURE

Across Gel Forward & Reverse ABO with D^{VI}- / CDE card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18- 25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, absence of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Forward & Reverse ABO with D^{VI}- / CDE test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments see the instrument manual.

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Reverse ABO Grouping (N_{A1} - N_B)

- Homogenise Across Red Blood Cell reagents (A₁/B₁).
- Pipette 50 µl Across Red Blood Cell reagents (A₁/B₁) into the N_{A1} / N_B microtubes.
- Pipette 50 µl serum or plasma to microtube N_{A1} / N_B.

Forward ABO Grouping + D^{VI}- / CDE Antigen Determination (A-B-AB-D^{VI}-CDE-Ctl)

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It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position. Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8°C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

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- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The meaning of Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In case of any discrepancy between forward-reverse groups, investigation must be made before giving the samples' test results (patient's age, clinical information, history, etc.).
- In ABO/Rh system : In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- Negative results in D^{VI}- microtube must be verified using other reagents and techniques which may detect different variants of the D antigen.

- If a negative result is obtained in D^{VI}- microtube and a positive result in CDE microtube, determination of Rh system antigens is recommended.
- Some individuals with A2 or A2B subgroups may present antibodies anti-A1, where a forward-reverse discrepancy can be seen.

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Forward & Reverse ABO with D^{VI}- / CDE card have been studied in a representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-AB	1061	100%	100%
anti-D ^{VI} -	1061	98,5%	100%
anti-CDE	1047	100%	100%

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Reverse ABO group

There was compatibility between the results of forward and reverse grouping tests of Across Gel Forward & Reverse ABO with D^{VI}- / CDE card.

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Adding volumes other than those specified in the method can alter the reaction.
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent.

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PACKAGING

Across Gel Forward & Reverse ABO With D^{VI}- / CDE

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Handwritten signatures and stamps are visible on the right side of the page.

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Blv. No.25/1 Gebze / Kocaeli / Türkiye

www.diapro.com.tr

(A / B / AB / D^{VI+} / CDE / Ctl / N_{A1} / N_B) | For 'in vitro' diagnostic use only. Read carefully before using the product. |

INTRODUCTION

The ABO system was the first human blood group system discovered by Landsteiner⁽¹⁾ and is still the most important system in transfusion practice.

Therefore, while classification of ABO system, both forward and reverse grouping must be done. In ABO forward grouping, red blood cells are tested with anti-A, anti-B, anti-AB antibodies, whereas in reverse grouping, ABO antibodies are determined in the serum/plasma using A1 and B red blood cells. Anti D reagent is used for the determination of presence or absence of D (RH1) antigen on the surface of red blood cell. anti-CDE reagent is used for the determination of presence or absence of D (RH1), C (RH2) and E (RH3) antigens for the purpose of establishing Rh negative red blood cells.

INDICATION

Blood grouping is life sustaining before blood transfusion. The blood group compatibility between the patient and the donor before the blood transfusion is one of the things that should be done before the blood transfusion. This card is also used for this purpose. It is not used for diagnostic purposes and is not specific to the diagnosis of a specific disease. The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) and by the presence of antibodies in the serum corresponding to the antigen or antigens missing on the surface of red blood cells (reverse grouping). Across Gel Forward & Reverse ABO with D^{VI+} / CDE card is used for determination of (ABO1), B (ABO2), D (RH1), C (RH2) and E (RH3) antigens and of possible anti-A and/or anti-B antibodies.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in red blood cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples.

Across Gel card contains polymerized dextran microspheres (gel) which function as a filter within the buffer solution in all microtubes on the card. Dextran and specific antibody-containing reagent(s) or buffer solution are mixed.

During the centrifugation of the cards, red blood cells agglutinate collected on the gel surface and throughout the length of the column according to their sizes. Non-agglutinated red blood cells go to the bottom of the microtube.

CONTENT

Microtubes of Across Gel Forward&Reverse ABO with D^{VI+}/CDE card are identified by the front labels of the card:

- **A (ABO1) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- **B (ABO2) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- **AB (ABO3) Microtube:** Monoclonal Anti-AB (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15, LB-2, 9621A8 clone)
- **D^{VI+} (RH1) Microtube:** Monoclonal Anti-D (mixture of human IgM and IgG antibodies, clones RUM-1, P3x61, ESD1M ve MS-26)
- **CDE/ RH1,2,3 Microtubes:** Monoclonal Anti-CDE (mixture of IgM/IgG antibodies of human origin, clones MS-24, MS-201, MS-26, MS-80)
- **Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- **N Microtubes:** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

STORAGE CONDITIONS

Store cards upright at 2-25°C. Do not place cards in front of and/or next to any source of heat (places under direct sunlight, radiators, air conditioners, big devices which radiate heat when operating) or air ducts. Do not freeze cards. When the box is kept in intact package, it is stable until the expiry date specified on the outer label.

MATERIAL REQUIRED FOR TEST BUT NOT PROVIDED IN THE KIT

- Across Auto System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Reader
- Across System Dispenser
- Across System Pipette
- Across System Workstation
- Across LISS
- Disposable pipette tips
- Tube for preparation of suspension
- For reverse group test, Across RBC A/B red blood cell reagent

COLLECTION AND PREPARATION OF SAMPLES

Samples

Blood samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used. Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

Determination of the antigens of ABO/Rh System:

- Fresh blood samples collected in tubes containing anticoagulants like EDTA or citrate can be used.
- Blood samples stored at 2-8°C up to 48 hours after extraction can be used.
- Blood samples stored at 2-8°C, collected from bags containing ACD, CPD, CPDA and SAG-Mannitol which are not beyond the expiry date specified on the label must be used.

Determination of the reverse ABO group:

- Serum or plasma can be used for the determination of Reverse ABO group.
- Serum or plasma samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used.
- Serum or plasma samples stored at 2-8 °C and taken 48 hours ago can be used.
- Serum or plasma samples stored between -20 / -80 °C can be used.

Sample Preparation

a) Preparing erythrocyte (red blood cell) suspension (To detect ABO/Rh system antigens)

Prepare 5% erythrocyte suspension (Add 25 µl erythrocyte to 0.5 ml Across LISS).

b) Plasma or Serum (Reverse ABO grouping)

Perform necessary centrifuging procedure for clean plasma or serum.

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

TEST PROCEDURE

Across Gel Forward & Reverse ABO with D^{VI+}/CDE card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18- 25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, lack of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Forward & Reverse ABO with D^{VI+}/CDE test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments see the instrumental manual.

MANUAL METHOD

- Each card is individually identified by a barcode. Unless a barcode reader is used identify them manually.
- Carefully peel off the aluminium foil that covers the microtubes and keep the card upright to avoid reagent cross contamination.

Reverse ABO Grouping (N_{A1} - N_B)

- Homogenise Across Red Blood Cell reagents (A₁/B).
- Pipette 50 µl Across Red Blood Cell reagents (A₁/B) into the N_{A1} / N_B microtubes.
- Pipette 50 µl serum or plasma to microtube N_{A1} / N_B.

Forward ABO Grouping + D^{VI+}/CDE Antigen Determination (A-B-AB-D^{VI+}-CDE-Ctl)

- Ensure the resuspension of the red blood cells before use.
- Pipette 10 µl 5% red blood cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).
- Centrifuge card or cards for 9 mins at 990 rpm.
- Read the results.

RESULTS

Results Reading

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position. Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8 °C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results ABO System

Forward ABO Group				Reverse ABO Group		
A Microtube	B Microtube	AB Microtube	Ctl. Microtube	Microtube N+ A1 Red Blood Cells Reagent	Microtube N+ B Red Blood Cells Reagent	ABO Group
0	0	0	0	+	+	0
+	0	+	0	0	+	A
0	+	+	0	+	0	B
+	+	+	0	0	0	AB

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Notes:

Rh System (D,C and E antigens)

A positive result in CDE microtube represents presence of Rh system D, C or E antigens.

- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The meaning of Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In case of any discrepancy between forward-reverse groups, investigation must be made before giving the samples' test results (patient's age, clinical information, history, etc.).
- In ABO/Rh system : In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.

- Negative reactions with the D^{VI+} reagent must be verified using other reagents and techniques which may detect different variants of the D antigen.
- If a negative result is obtained in D^{VI+} microtube and a positive result in CDE microtube, determination of Rh system antigens is recommended.
- Some individuals with A2 or A2B subgroups may present antibodies anti-A1, where a forward-reverse discrepancy can be seen.

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Forward & Reverse ABO with D^{VI+}/CDE card have been studied in a representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-AB	1061	100%	100%
anti-D ^{VI+}	1061	98,5%	100%
anti-CDE	1047	100%	100%

(I) Sensitivity: [(number of true positive results) / (number of true positive results + number of false negative results)] x 100

(II) Specificity: [(number of true negative results) / (number of true negative results + number of false positive results)] x 100

Reverse ABO group

There was compatibility between the results of forward and reverse grouping tests of Across Gel Forward & Reverse ABO with D^{VI+}/CDE card.

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Adding a volume other than that stated in the method may change the reaction.
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr +90(262)751-53-34) and the competent.

REFERENCES

- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 10th edition, Blackwell Scientific Publications, Oxford, 1993.
- Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005

PACKAGING

Across Gel Forward & Reverse ABO With D^{VI+} / CDE

(A / B / AB / D^{VI+} / CDE / Ctl / N_{A1} / N_B)

Ref No: 810203 ---- 50 Cards

Document No: TK3411

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

IVD	"In Vitro" Diagnostic Medical Device
LOT	Batch Code
Use by date	
Temperature Limit	
Manufacturer	
Consult Instructions for Use	
REF	Catalogue Number
Do not re-use	
Do not use if package is damaged and consult	
Contains biological material of animal origin	
Contains biological material of human origin	

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www.diapro.com.tr

INTRODUCTION

The ABO system was the first human blood group system discovered by Landsteiner⁽¹⁾, and is still the most important system in transfusion practice.

Therefore, while classification of ABO system, both forward and reverse grouping must be done. In ABO forward grouping, red blood cells are tested with anti-A, anti-B, anti-AB antibodies, whereas in reverse grouping, ABO antibodies are determined in the serum/plasma using A1 and B red blood cells in the case of forward and reverse incompatibilities, further determination must be made with Lectin A 1 and anti-H reagents. In ABO blood group typing, anti-A, anti-B, anti-AB, anti-H and Lectin A1 reagents are used, and A1/B reverse group test are used as complementary.

INDICATION

The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) and by the presence of antibodies in the serum corresponding to the antigen or antigens missing on the surface of red blood cells (reverse grouping).

Across Gel ABO Typing card is used for determination of A (ABO1), B (ABO2), A1 (Lectin), H and of possible anti-A and/or anti-B antibodies.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in red blood cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples. Across Gel card contains polymerized dextran microspheres (gel) which function as a filter within the buffer solution in all microtubes on the card. Dextran and specific antibody-containing reagent(s) or buffer solution are mixed. During the centrifugation of the cards, red blood cells agglutinates collected on the gel surface and throughout the length of the column according to their sizes. Non-agglutinated red blood cells go to the bottom of the microtube.

CONTENT

The microtubes of Across Gel ABO Typing card are identified by the front labels of the card:

- A (ABO1) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- B (ABO2) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- AB (ABO3) Microtube:** Monoclonal Anti-AB (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15, LB-2, 9621A8 clone)
- A1 Lectin Microtube:** Dolichos Biflorus
H Microtube: Monoclonal IgM (IgM antibody of murine origin, LH1 clone)
- Ctl. Microtube:** Neutral gel (Buffered solution without antibodies)
- N microtubes:** Neutral gel (Buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

STORAGE CONDITIONS

Store cards upright at 2-25°C. Do not place cards in front of and/or next to any source of heat (places under direct sunlight, radiators, air conditioners, big devices which radiate heat when operating) or air ducts. Do not freeze cards. When the box is kept in intact package, it is stable until the expiry date specified on the outer label.

MATERIAL REQUIRED FOR TEST BUT NOT PROVIDED IN THE KIT

- Across Auto System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Reader
- Across System Dispenser
- Across System Pipette
- Across System Workstation
- Across LISS
- Disposable pipette tips
- Tube for preparation of suspension
- For reverse group test, Across RBC A/B red blood cell reagent

COLLECTION AND PREPARATION OF SAMPLES**Samples**

Blood samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used.

Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

Determination of the antigens of ABO/Rh System:

- Fresh blood samples collected in tubes containing anticoagulants like EDTA or citrate can be used.
- Blood samples stored at 2-8°C up to 48 hours after extraction can be used.
- Blood samples stored at 2-8°C, collected from bags containing ACD, CPD, CPDA and SAG-Mannitol which are not beyond the expiry date specified on the label must be used.

Determination of the reverse ABO group:

- Serum or plasma can be used for the determination of Reverse ABO group.
- Serum or plasma samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used.
- Serum or plasma samples stored at 2-8 °C and taken 48 hours ago can be used.
- Serum or plasma samples stored between -20 / -80 °C can be used.

Sample Preparation

a) Preparing erythrocyte (red blood cell) suspension (To detect ABO/Rh system antigens)
Prepare 5% erythrocyte suspension (Add 25 µl erythrocyte to 0.5 ml Across LISS).

b) Plasma or Serum (Reverse ABO grouping)
Perform necessary centrifuging procedure for clean plasma or serum.

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

TEST PROCEDURE

Across Gel ABO Typing card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18- 25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, lack of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel ABO Typing test method is the same principle for manual system and Across Auto System Auto System semi-automatic or automatic instruments see the instrumental manual.

MANUAL METHOD

- Each card is individually identified by a barcode. Unless a barcode reader is used identify them manually.
- Carefully peel off the aluminium foil that covers the microtubes and keep the card upright to avoid reagent cross contamination.

Reverse ABO Grouping (N_{A1} - N_B)

- Homogenise Across Red Blood Cell reagents (A₁/B).
- Pipette 50 µl Across Red Blood Cell reagents (A₁/B) into the N_{A1} / N_B microtubes.
- Pipette 50 µl serum or plasma to microtube N_{A1} / N_B.

Forward ABO Grouping + Antigen Determination (A-B-AB-A1 Lectin-H-Ctl)

- Ensure the resuspension of the red blood cells before use.
- Pipette 10 µl 5% red blood cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).
- Centrifuge card or cards for 9 mins at 990 rpm.
- Read the results.

RESULTS**Results Reading**

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position.

Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8 °C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results

Serological reactions observe in individuals⁽²⁾:

Forward ABO Group						Reverse ABO Group		Phenotype ABO Group
A	B	AB	A1 Lectin	H	Ctl.	A1	B	
+4	0	+4	+4	0	0	0	+4	A ₁
+4	0	+4	0	+2	0	†	+4	A ₂
+2/mf	0	+2/mf	0	+3	0	†	+4	A ₃
0/±	0	0/±	0	+4	0	0	+4	A _m
0/±	0	+1/+2	0	+4	0	0/+2	+4	A _x
0	0	0	0	+4	0	0/+2	+4	A _d
0	+4	+4	0	0	0	+4	0	B
0	+1/mf	+2/mf	0	+4	0	+4	0	B ₃
0	0	0/±	0	+4	0	+4	0	B _m
0	0/±	0/+2	0	+4	0	+4	0	B _x
+4	+4	+4	+3	0	0	0	0	AB
+4	+4	+4	0	0	0	†	0	A ₂ B
0	0	0	0	+4	0	+4	+4	O

† In these phenotypes, presence of Anti-A1 may vary; while most of the individuals with A2 group carry anti-A1 antibody, most of those with A3 group do not have this antibody.

+ = positive
mf = weak positive; a mixed field, similar to double population in the column
0 = negative

Notes:

- Results are not alone eligible for diagnostic use. They should be evaluated with other clinical findings of the patient and other data.
- The Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid.
- Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In ABO/Rh system : In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- In case of any disagreement between Forward-reverse groups, investigation must be made before giving the samples' test results (patient's age, clinical information, history, etc.).

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel ABO Typing card have been analyzed with representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-AB	1061	100%	100%

(I) Sensitivity: [(number of true positive results) / (number of true positive results + number of false negative results)] x 100

(II) Specificity: [(number of true negative results) / (number of true negative results + number of false positive results)] x 100

Reverse ABO group

There was agreement between the results of forward and reverse grouping tests of Across Gel ABO Typing card.

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Adding a volume other than that stated in the method may change the reaction.
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.

REFERENCES

- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 10th edition, Blackwell Scientific Publications, Oxford, 1993.
- Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005

PACKAGING**Across Gel ABO Typing**

(A / B / AB / A1 (Lectin) / H / Ctl / N_{A1} / N_B)

Ref No: 810208 ---- 50 Cards


810208 ---- 10 Cards

Document No: TK3416

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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