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
INDICATION

Across LISS is used for preparing erythrocyte suspensions, dilution of serum/plasma used in gel techniques.

PRINCIPLE

LISS (Low Ionic Strength Solution) significantly increases the antigen-antibody interaction by lowering the ion density in the environment. The use of this solution facilitates the binding of the antibody to red blood cells, reducing the time to form the antigen-antibody complex.

CONTENT

- Contains buffer solution with low ionic strength and preservatives.
- Reagents are ready for use.

STORAGE CONDITIONS

- Store Across LISS Solutions at 2-25°C.
- Do not freeze.
- Do not store at high temperature.
- If stored properly, it is stable until expiry date provided that contamination is prevented.

MATERIAL REQUIRED BUT NOT PROVIDED IN THE KIT

- Across Gel cards
- Follow the Instructions for use of the Across Gel cards to be used.

COLLECTION AND PREPARATION OF SAMPLES

Samples

Blood samples with anticoagulant (containing plasma) can be used. Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

To detect antigens on erythrocyte surface, work with:

- Fresh blood samples collected in tubes containing anticoagulants like EDTA, citrate etc. can be used.
- Blood samples stored at 2-8°C up to 48 hours after extraction can be used.
- Blood samples stored at 2-8°C, collected from bags containing ACD, CPD, CPDA and SAG-Mannitol which are not beyond the expiry date specified on the label can be used.

Sample Preparation:

a) Preparing erythrocyte (red blood cell) suspension

(To detect ABO/Rh- Kell system antigens)

To prepare 5% erythrocyte suspension, add 25 µl erythrocyte to 0,5 ml Across LISS.

b) Preparing erythrocyte (red blood cell) suspension (To screen unexpected antibodies such as Direct Coombs, Cross-Match)

To prepare 1% erythrocyte suspension, add 10 µl erythrocyte to 1 ml Across LISS.

TEST PROCEDURE

It can be worked with Across LISS manual equipment or Across System Auto semi-automatic and automatic instruments. For Across System Auto semi-automatic and automatic instruments, see the instrument's user manual. The working principles of Across LISS solution in manual equipment and Across Auto System semi-automatic and automatic instruments are equivalent.

WARNINGS AND PRECAUTIONS

- The product must be used by a qualified staff member.
- Preparing the red blood cell suspension in a concentration other than the recommended concentration may give a false negative or a false positive test result.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Adding a volume other than that stated in the methods in the Across Gel cards may change the reaction.
- Do not use a solution with past expiry date.
- Do not use the product in the presence of microbial or fungal contamination.
- Do not use the product if the package is damaged.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr/ +90(262)751-53-34) and the competent.

PACKAGING

Across LISS

Ref No : 830100 ---- 100 ml

830500 ---- 500 ml

Document No: TK3430

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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Valgoni R. Nicolas
Bacio Gruber
NG MED S.R.L.
30-71687909-3

INDICATION

Bromeline can be used in gel techniques to complement standard immunohematological testing procedures, as a reaction enhancer of certain antibodies with antigen, and in Crossmatch auto / allo compatibility testing.

PRINCIPLE

Bromeline is a proteolytic enzyme that reduces electronegative charges on the surfaces of red blood cells. Treatment of erythrocytes with proteolytic enzymes alters some cell surface components, allowing antibodies that do not normally bind in saline to bind to antigens on the cell surface. Some erythrocyte surface antigens become more accessible for serum antibodies after enzyme treatment, while others are destroyed. Such enzymes can increase the reactivity of certain antibodies, such as cold agglutinins and rhesus antibodies, while inhibiting the reactions of those specifically directed to MNS group antigens and Duffy blood group antigens.

CONTENT

Ready-to-use liquid bromeline solution with preservatives

STORAGE CONDITIONS

- Do not expose to extreme temperatures.
- Do not freeze.
- When not in use, store at 2 - 8 °C.
- If the product is stored properly at 2 - 8 °C, it is stable until the stated expiry date after opening.

TEST PROCEDURE

Adjust to room temperature before use. Follow the instructions for use on the relevant Across Gel cards.

WARNINGS AND PRECAUTIONS

- The product must be used by a qualified staff.
- Due to high sensitivity of enzymatic techniques, it is crucial that the recommended incubation times are not exceeded to avoid false positive results.
- The reaction strength of positive results depends on the age of the used
- The presence of auto-agglutinins in test specimens may cause false reactions.
- Erythrocytes coated with allo- or auto-antibodies (DAT positive) may produce false reactions blood sample.
- In normal serum/plasma, enzymatic treatment of red blood cells may cause the non-specific agglutination of the red blood cells.
- Do not use the product in the presence of microbial or fungal contamination.
- Use of decomposed samples and reagent or those contaminated by bacteria or otherwise can cause false positive or false negative results.
- As other proteolytic enzymes do, Across Bromelinee partially or completely breaks down Fya, Fyb, M, N, S, s, Xga antigens.
- Adding a volume other than that recommended in the method may change the reaction.
- Do not use a solution with past expiry date.
- Do not use the product if the package is damaged.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr/ +90(262)751-53-34) and the competent.

PACKAGING**Across Bromeline**

Ref No : 830400 ---- 25 ml


830200 ---- 100 ml

Document No: TK3431

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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VERONICA E. GIANPHERI
FARMACÉUTICA
M.N. 16215

Valgani R. Nicolas
Socio Gerente
RG MED S.R.L.
30-71697909-3

For 'in vitro' diagnostic use only. Read carefully before using the product.

INDICATION

It is used for screening and identification of unexpected antibodies, autocontrol and cross-compatibility (Crossmatch) tests by enzymatic method.

PRINCIPLE

Papaine is a proteolytic enzyme that reduces electronegative charges on the surfaces of red blood cells. Treatment of erythrocytes with proteolytic enzymes alters some cell surface components, allowing antibodies that do not normally bind in saline to bind to antigens on the cell surface. Some erythrocyte surface antigens become more accessible for serum antibodies after enzyme treatment, while others are destroyed. Such enzymes can increase the reactivity of certain antibodies, such as cold agglutinins and rhesus antibodies, while inhibiting the reactions of those specifically directed to MNS group antigens and Duffy blood group antigens.

CONTENT

Ready-to-use liquid papaine solution.

STORAGE CONDITIONS

- Do not expose to extreme temperatures.
- Do not freeze.
- When not in use, store at 2 - 8 °C.
- If the product is stored properly at 2 - 8 °C, it is stable until the stated expiry date after opening.

TEST PROCEDURE

Adjust to room temperature before use. Follow the instructions for use on the relevant Across Gel cards.

WARNINGS AND PRECAUTIONS

- The product must be used by a qualified staff member.
- Due to high sensitivity of enzymatic techniques, it is crucial that the recommended incubation times are not exceeded to avoid false positive results.
- The reaction strength of positive results depends on the age of the used.
- The presence of auto-agglutinins in test specimens may cause false reactions.
- Erythrocytes coated with allo- or auto-antibodies (DAT positive) may produce false reactions blood sample.
- In normal serum/plasma, enzymatic treatment of red blood cells may cause the non-specific agglutination of the red blood cells.
- Do not use the product in the presence of microbial or fungal contamination.
- Use of decomposed samples and reagent or those contaminated by bacteria or otherwise can cause false positive or false negative results.
- As other proteolytic enzymes do, Across Papaine partially or completely breaks down Fya, Fyb, M, N, S, s, Xga antigens.
- Adding a volume other than that recommended in the method may change the reaction.
- Do not use a solution with past expiry date.
- Do not use the product if the package is damaged.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr/ +90(262)751-53-34) and the competent.

PACKAGING

Across Papaine

Ref No : 830300 ---- 10 ml

Document No: TK3432

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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INDICATION

Internal washing of "In vitro" diagnostic devices

INTRODUCTION

- Across Wash O washing solution is used to prevent bacteria and other contaminants that can occur on Across Auto System semi-automatic or automatic instruments, analyzer's inside / outside surfaces of the fluidic system and probes and to eliminate these risks.
- Across Auto System semi-automatic or automatic instruments' analyzer's "Maintenance" and "System" fluid reservoirs are filled with 125 ml Across
- Wash O washing solution and complete to 10 liters with distilled water. After thorough mixing of the "Maintenance" and "System" liquid reservoirs, the caps are closed.

CONTENT

Saline-based concentrated washing solution used in fluidics system of automatic systems contains colourant and preservative. Wash O product to be used should be prepared with 8% distilled water. Only diluted solutions are used in devices.

STORAGE CONDITIONS

Store Across Wash O solutions at 18-25°C. Do not freeze. The expiry date of the product is on the label. The product is stable until the expiration date on the label if stored in suitable storage conditions and kept away from the agents causing the contamination.

Across Wash O solution should be used within maximum 7 days after being added to the "Maintenance" and "System" fluid reservoirs. If not used, it should be replaced with new solution.

WARNINGS AND PRECAUTIONS

- The product should only be used by trained and certified staff about Across System.
- Do not use if there is microbial or fungal contamination. The product can not be used for purposes other than the stated purpose. No guarantee can be given against problems that may arise if used for any other purpose.
- The "Maintenance" and "System" fluid reservoirs must be cleaned before the Across Wash O solution is ready for use.
- If washing solutions other than the specified product are used, the Across Auto System semi-automatic or automatic instruments' analyzer may cause problems on the inside / outside surfaces of the fluidic system and probes. (E.g. contamination due to bacteria and other agents, hemolysis in blood samples, fluid perception and sensitivity problems)
- Do not use the product in case of fractures, cracks, holes or dirt which may be on the packaging material.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr/ +90(262)751-53-34) and the competent.

PACKAGING

Across Wash O

Ref No: 830800 --- 500 ml

Document No: TK3433

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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IVD	"In Vitro" Diagnostic Medical Device
LOT	Batch Code
	Use by date
	Manufacturer
	Temperature Limit
	Consult Instructions for Use
REF	Catalogue Number

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INDICATION

Neutralizing agent of anti-CD38 antibody.

INTRODUCTION

Anti-CD38 antibody used in the treatment of multiple myeloma causes false positives in Cross-Match and antibody screening/identification tests applied before blood transfusion. Across Dara Remove is an agent that eliminates this false positive.

CONTENT

Across Dara Remove 1x500µl
Across S-LISS 1x 10ml

TEST PROCEDURE

1. A 1% erythrocyte suspension is prepared using the donor's erythrocytes and Liss solution. (10µl of donor erythrocytes are added to 1ml of Liss solution.)
2. Dilute Across Dara Remove solution to Across S-LISS (Special Liss) solution. (200µl Across S-LISS and 50µl Across Dara Remove for 250µl, 400µl Across S-LISS and 100µl Across Dara Remove for 500µl, 800µl Across S-LISS and 200µl Across Dara Remove for 1 ml)
3. Add 50µl of 1% erythrocytes onto the Across AHG card and add 25µl of Across Dara Remove solution. Gently pipette 5 times inside the card microtube to mix the erythrocytes and Across Dara Remove. The card is then incubated at 37°C for 30 minutes. (At this stage, Across Dara Remove washes CD38 antigens on the erythrocyte surface.)
4. 25µl of the patient's serum is added onto the erythrocytes on the card after the incubation is completed and gently pipetted 5 times again to mix. It is then incubated at 37°C for 15 minutes. (Cross-Match)
5. The card, whose incubation is finished, is centrifuged 2 times for 9 minutes at 990 rpm. It is centrifuged for 18 minutes in total. At the end of the centrifuge, the card is examined for the result.

WARNINGS AND PRECAUTIONS

- Across Dara Remove is an agent that denaturates the Kell antigen. Kell antigen should be checked before and after Dara Remove application.
- Across Dara Remove is stable up to the expiry date stated on the label.
- Across Dara Remove diluted with Across S-Liss is stable for 24 hours.
- Keep away from light.
- Store at 2-8°C.
- Wash skin thoroughly after handling.
- Wear protective gloves/ eye protection/ face protection.
- If in eyes, rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- If on skin, wash with plenty of water.
- If skin irritation occurs, get medical advice/ attention.
- Take off contaminated clothing and wash it before reuse.

PACKAGING

Across Dara Remove Kit

Ref No: 820205 --- 1x10ml , 1x500µl
Ref No 820206 Across Dara Remove
Ref No 820207 Across S-LISS

Document No: TK2358

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

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Therefore, while classification of ABO system, both forward and reverse grouping must be done. In ABO forward grouping, red blood cells are tested with Anti-A, Anti-B, Anti-AB antibodies, whereas in reverse grouping, ABO antibodies are determined in the serum/plasma using A1 and B red blood cells. Unlike ABO system antibodies, Rh and Kell antibodies develop in the antigens missing individuals due to the antigens present blood transfusion or pregnancy. In these systems, antigens with the highest immunogenicity are RH1 (D) and KEL1 (Kell).

INDICATION

The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) and by the presence of antibodies in the serum corresponding to the antigen or antigens missing on the surface of red blood cells (reverse grouping). It is used for determination of A (AB01), B (AB02), D (RH1) and Kell (KEL 1) antigens and of possible Anti-A and/or Anti-B antibodies.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in red blood cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples.

Across Gel card contains polymerized dextran microspheres (gel) which function as a filter within the buffer solution in all microtubes on the card. Dextran and specific antibody-containing reagent(s) or buffer solution are mixed.

During the centrifugation of the cards, red blood cells agglutinates collected on the gel surface and throughout the length of the column according to their sizes. Non-agglutinated red blood cells go to the bottom of the microtube.

CONTENT

Microtubes of Across Gel Forward&Reverse ABO with D^W/KELL card are identified by the front labels of the card:

- **A (AB01) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- **B (AB02) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- **AB(AB03) Microtube:** Monoclonal Anti-AB (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15, LB-2, 9621A8 clone)
- **D^W(RH1) Microtube:** Monoclonal Anti-D (mixture of human IgM antibodies, RUM-1, P3x61 clone)
- **Kell (KEL1) Microtube:** Monoclonal Anti-K (Human IgM antibodies, MS-56 clone)
- **Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- **N Microtubes :** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

STORAGE CONDITIONS

Store cards upright at 2-25°C. Do not place cards in front of and/or next to any source of heat (produced under direct sunlight, radiators, air conditioners, big devices which radiate heat when operating) or air ducts. Do not freeze cards. When the box is kept in intact package, it is stable until the expiry date specified on the outer label.

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Samples

Blood samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used.

Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

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- Serum or plasma can be used for the determination of Reverse ABO group.
- Serum or plasma samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used.
- Serum or plasma samples stored at 2-8 °C and taken 48 hours ago can be used.
- Serum or plasma samples stored between -20 / -80 °C can be used.

Sample Preparation

a) Preparing erythrocyte (red blood cell) suspension (To detect ABO/Rh/Kell system antigens)

Prepare 5% erythrocyte suspension (Add 25 µl erythrocyte to 0.5 ml Across LISS).

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TEST PROCEDURE

Across Gel Forward & Reverse ABO with D^W/Kell card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18-25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, lack of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Forward & Reverse ABO with D^W/Kell test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments see the instrumental manual.

MANUAL METHOD

- Each card is individually identified by a barcode. Unless a barcode reader is used identify them manually.
- Carefully peel off the aluminium foil that covers the microtubes and keep the card upright to avoid reagent cross contamination.
- Reverse ABO Grouping (N_{A1} - N_B)**
- Homogenise Across Red Blood Cell reagents (A₁/B).
- Pipette 50 µl Across Red Blood Cell reagents (A₁/B) into the N_{A1}/N_B microtubes.
- Pipette 50 µl serum or plasma to microtube N_{A1} /N_B.
- Forward ABO Grouping + D^W / KELL Antigen Determination (A-B-AB-D^W-Kell-Ctl)**
- Ensure the resuspension of the red blood cells before use.
- Pipette 10 µl 5% red blood cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).
- Centrifuge card or cards for 9 mins at 990 rpm.
- Read the results.

RESULTS

Results Reading

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

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- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The meaning of Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In case of any discrepancy between forward-reverse groups, investigation must be made before giving the samples' test results (patient's age, clinical information, history, etc.).
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- In ABO/Rh system : In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- Negative reactions with the D^W microtubes must be verified using other reagents and techniques which may detect different variants of the D antigen.
- Some individuals with A2 or A2B subgroups may present antibodies anti-A1, where a forward-reverse discrepancy can be seen.

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Forward & Reverse ABO with D^W/Kell card have been studied in a representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-AB	1061	100%	100%
anti-D ^W	1061	98,5%	100%
anti-Kell	220	100%	100%

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Reverse ABO group

There was compatibility between the results of forward and reverse grouping tests of Across Gel Forward & Reverse ABO with D^W/Kell card.

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Adding a volume other than that stated in the method may change the reaction.
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent.

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PACKAGING

Across Gel Forward & Reverse ABO With D^W/ KELL

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FORWARD & REVERSE ABO WITH D^{VI} / KELL

across GEL

(A / B / AB / D^{VI} / Kell / Ctl / N_{A1} / N_B) | For 'in vitro' diagnostic use only. Read carefully before using the product.

INTRODUCTION

The ABO system was the first human blood group system discovered by Landsteiner⁽¹⁾, and is still the most important system in transfusion practice.

Therefore, while classification of ABO system, both forward and reverse grouping must be done. In ABO forward grouping, red blood cells are tested with anti-A, anti-B, anti-AB antibodies, whereas in reverse grouping, ABO antibodies are determined in the serum/plasma using A1 and B red blood cells. Unlike ABO system antibodies, Rh and Kell antibodies develop in the antigens missing individuals due to the antigens present blood transfusion or pregnancy. In these systems, antigens with the highest immunogenicity are RH1 (D) and KEL1 (Kell).

INDICATION

The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) and by the presence of antibodies in the serum corresponding to the antigen or antigens missing on the surface of red blood cells (reverse grouping). Across Gel Forward & Reverse ABO with D^{VI}/KELL card is used for determination of A (AB01), B (AB02), D (RH1) and Kell (KEL 1) antigens and of possible anti-A and/or anti-B antibodies.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in red blood cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples. Across Gel card contains polymerized dextran microspheres (gel) which function as a filter within the buffer solution in all microtubes on the card. Dextran and specific antibody-containing reagent(s) or buffer solution are mixed. During the centrifugation of the cards, red blood cells agglutinate collected on the gel surface and throughout the length of the column according to their sizes. Non-agglutinated red blood cells go to the bottom of the microtube.

CONTENT

Microtubes of Across Gel Forward&Reverse ABO with D^{VI}/KELL card are identified by the front labels of the card:

- **A (ABO1) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- **B (ABO2) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- **AB (ABO3) Microtube:** Monoclonal Anti-AB (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15, LB-2, 9621A8 clone)
- **D^{VI} (RH1) Microtube:** Monoclonal Anti-D (mixture of IgG and IgM antibodies of human origin, clones RUM-1, P3x61, ESD1M and MS-26 clone)
- **Kell (KEL1) Microtube:** Monoclonal Anti-K (IgM antibodies of human origin, MS-56 clone)
- **Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- **N Microtubes:** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

STORAGE CONDITIONS

Store cards upright at 2-25°C. Do not place cards in front of and/or next to any source of heat (places under direct sunlight, radiators, air conditioners, big devices which radiate heat when operating) or air ducts. Do not freeze cards. When the box is kept in intact package, it is stable until the expiry date specified on the outer label.

MATERIAL REQUIRED FOR TEST BUT NOT PROVIDED IN THE KIT

- Across Auto System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Reader
- Across System Dispenser
- Across System Pipette
- Across System Workstation
- Across LISS
- Disposable pipette tips
- Tube for preparation of suspension
- For reverse group test, Across RBC A/B red blood cell reagent

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Samples

Blood samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used. Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

Determination of the antigens of ABO/Rh System:

- Fresh blood samples collected in tubes containing anticoagulants like EDTA or citrate can be used.
- Blood samples stored at 2-8°C up to 48 hours after extraction can be used.
- Blood samples stored at 2-8°C, collected from bags containing ACD, CPD, CPDA and SAG-Mannitol which are not beyond the expiry date specified on the label must be used.

Determination of the reverse ABO group:

- Serum or plasma can be used for the determination of Reverse ABO group.
- Serum or plasma samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used.
- Serum or plasma samples stored at 2-8°C and taken 48 hours ago can be used.
- Serum or plasma samples stored between -20 / -80 °C can be used.

Sample Preparation

a) Preparing erythrocyte (red blood cell) suspension (To detect ABO/Rh system antigens)
Prepare 5% erythrocyte suspension (Add 25 µl erythrocyte to 0.5 ml Across LISS).

b) Plasma or Serum (Reverse ABO grouping)

Perform necessary centrifuging procedure for clean plasma or serum.

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

TEST PROCEDURE

Across Gel Forward & Reverse ABO with D^{VI}/KELL card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18-25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, lack of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Forward & Reverse ABO with D^{VI}/KELL test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments see the instrumental manual.

MANUAL METHOD

- Each card is individually identified by a barcode. Unless a barcode reader is used identify them manually.
- Carefully peel off the aluminium foil that covers the microtubes and keep the card upright to avoid reagent cross contamination.
- Reverse ABO Grouping (N_{A1} - N_B)**
Homogenise Across Red Blood Cell reagents (A1/B).
- Pipette 50 µl Across Red Blood Cell reagents (A1/B) into the N_{A1} / N_B microtubes.
- Pipette 50 µl serum or plasma to microtube N_{A1} / N_B.
- Forward ABO Grouping + D^{VI} / KELL Antigen Determination (A-B-AB-D^{VI}-Kell-Ctl)**
- Ensure the resuspension of the red blood cells before use.
- Pipette 10 µl 5% red blood cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).
- Centrifuge card or cards for 9 mins at 990 rpm.
- Read the results.

RESULTS

Results Reading

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
Negative	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP	-	Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position. Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8°C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results ABO System

Forward ABO Group				Reverse ABO Group		ABO Group
A Microtube	B Microtube	AB Microtube	Ctl. Microtube	Microtube N+ A1 Red Blood Cells Reagent	Microtube N+ B Red Blood Cells Reagent	
0	0	0	0	+	+	0
+	0	+	0	0	+	A
0	+	+	0	+	0	B
+	+	+	0	0	0	AB

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Notes:

KELL System

A positive result in Kell microtube represents presence of Kell antigen. In the weaker than +3 reactions, presence of a weak antigen must be investigated.

- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In case of any discrepancy between Forward-reverse groups, investigation must be made before giving the samples' test results (patient's age, clinical information, history, etc.).
- In ABO/Rh system : In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Negative reactions with the D^{VI} reagent must be verified using other reagents and techniques which may detect different variants of the D antigen.
- Some individuals with A2 or A2B subgroups may present antibodies anti-A1, where a forward-reverse discrepancy can be seen.

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Forward & Reverse ABO with D^{VI}/KELL card have been studied in a representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-AB	1061	100%	100%
anti-D ^{VI}	1061	98,5%	100%
anti-Kell	220	100%	100%

(I) Sensitivity: [(number of true positive results) / (number of true positive results + number of false negative results)] x 100

(II) Specificity: [(number of true negative results) / (number of true negative results + number of false positive results)] x 100

Reverse ABO group

There was compatibility between the results of forward and reverse grouping tests of Across Gel Forward & Reverse ABO with D^{VI}/KELL card.

WARNINGS AND PRECAUTIONS

- It is used by the professional health team (health technician, biologist, health technician, etc.) working in health institutions such as hospital, blood centers where blood transfusion treatment is applied.
- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Adding a volume other than that stated in the method may change the reaction.
- Do not use a card with past expiry date.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens and may not be detected.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr/ +90(262)751-53-34) and the competent.

REFERENCES

- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 10th edition, Blackwell Scientific Publications, Oxford, 1993.
- Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005

PACKAGING

Across Gel Forward & Reverse ABO With D^{VI} / KELL

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(A / B / AB / D^{VI} / KELL / Ctl / N_{A1} / N_B)