	Convafiber	PM 2391-37
		Legajo N° 2391

Rótulo

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

 Convatec Limited First Avenue, Deeside Industrial Park Deeside,
Flintshire, CH5 2NU, UK

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INSTRUCTIONS FOR USE / INSTRUCCIONES DE USO / MODE D'EMPLOI / 사용 설명서

www.convatec.com

PRODUCT DESCRIPTION & CLINICAL BENEFITS

Aquacel ConvaFiber™ Dressings are soft, sterile, dressings made from sodium carboxymethylcellulose blended with strengthening cellulose fibres. This dressing absorbs high amounts of wound fluid and bacteria, creates a soft, cohesive gel that intimately conforms to the wound surface, manages excess exudate levels, provides a barrier to protect the wound bed from contamination, maintains a moist wound healing environment to support the body's healing process and aids in the removal of non-viable tissue from the wound (autolytic debridement)

INTENDED USE, INDICATIONS & INTENDED USERS

Aquacel ConvaFiber™ Dressings have been designed to be used as a primary dressing by Healthcare Professionals, carers and patients under the direction of a Healthcare Professional as per the indications below.

- Leg ulcers, including:
 - Venous ulcers
 - Arterial ulcers
- Diabetic ulcers including diabetic foot ulcers
- Pressure ulcers/injuries
- Surgical wounds
- Traumatic wounds
- Malignant wounds
- Partial thickness burns

CONTRAINDICATIONS

Aquacel ConvaFiber™ Dressings should not be used on individuals who are sensitive to or who had an allergic reaction to the dressing or any of its components listed in the Product Description.

PRECAUTIONS AND OBSERVATIONS

- Aquacel ConvaFiber™ Dressings
 - are for single-use only and should not be re-used. Reuse may lead to increased risk of infection, cross-contamination and delayed healing.
 - are not compatible with petroleum-based products.
 - are designed to be MRI safe
- Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open prior to use and dispose of the device according to local regulations.
- Due to the sterilisation process, there may be a slight odour on opening the primary packaging.
- During the body's normal healing process, non-viable tissue is removed from the wound (autolytic debridement), which could initially make the wound appear larger.
- Newly formed blood vessels may occasionally produce bloodstained wound fluid following removal of the dressing.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and

federal laws and regulations.

- If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the Competent Authority of the Member State in which the user and/or patient is established.

ASSOCIATED PRODUCT & THERAPIES

Aquacel ConvaFiber™ Dressings may be used alone or in combination with other wound care products such as; primary dressings, secondary dressings, securement bandages, compression therapy devices or medicated topical treatments. Use of additional products and therapies should be undertaken with the direction of a Healthcare Professional.

CONTACT AND DURATION OF USE

Aquacel ConvaFiber™ Dressings can be worn for up to 7 days, dressings should be changed earlier if clinically indicated.

DIRECTIONS FOR USE

- Before applying the dressing, cleanse the wound area with an appropriate wound cleanser.
- Dressings should overlap at least 1 cm onto the skin surrounding the wound. Dressings can be cut to size. Aquacel™ Ribbon should be used for cavity wounds
- Discard any unused portion of the dressing.
- A secondary dressing is required to hold the dressing in place, such as:
 - a moisture retentive cover dressing such as DuoDerm™ Extra Thin in lightly exuding wounds;
 - a secondary cover dressing such as ConvaFoam™ or Aquacel™ Foam or ConvaMax™, in moderate to heavily exuding wounds;
- See individual cover dressing package inserts for complete instruction for use.
- All wounds should be inspected regularly. Remove the dressing when clinically indicated (i.e., leakage, excessive bleeding, increased pain) or after a maximum of seven days.

FOR DRY WOUNDS

- Pre-moisten the dressings with sterile saline prior to dressing the wound.
- Fully cover the dressing with an occlusive dressing such as DuoDerm™ Extra Thin to avoid drying out of the dressing and subsequent dressing adherence to the wound.
- If on removal, the dressing is dry, hydrate with sterile saline prior to removal to reduce risk of trauma.

Store in cool dry place.

If further information or guidance is needed, please contact Convatec Professional Services.

Made in UK

© 2024 Convatec

™/® indicates a trademark of the Convatec Group

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
Effective

ARTWORK SPECIFICATION	
Artcode	1749828
Artwork Description	INS AQUACEL CONVAFIBER NAI
ICC*	423989, 423992, 423993, 423996
SAP*	1749830, 1749831, 1749832, 1749833
KOM/REF*	N/A
IS Code*	N/A
Component Type	Insert
LR#	LR-005151
Existing Art Code being Revised	N/A
REQUIREMENT	SPECIFICATION
Bucket	NAI
Local Code(s)	N/A
Size Code / Dieline	N/A
Dimensions	Flat: 270x450mm; Folded 90x90mm
Color(s) / Colour(s)	Black
Specification Reference(s)	MT56-001
Bar Code: UPC (A) or GTIN**	N/A
Bar Code: PZN – Code 39**	N/A
Bar Code: JAN-13**	N/A
* NOTE: Bar Code: See artwork for positioning, magnification and truncation. Scanability verified to ANSI grade “C” or better.	
HISTORY	
07-Nov-2024 /1749828: New artwork created.	

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