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
Beckman Coulter, Argentina

Confidencial

Confidential - Company Proprietary

Rolando F. Bornia
APODERADO

- 35 -

System Guide

AQUIOS STEM

Rx Only in the U.S.A.

B75978AA
March 2022

HUNOTECH S.A.S. Beckman Coulter
13 J, avenue de Lattre de Tassigny, BP 177
13276 Marseille cedex 9, France
(33) 941 172 727

Gabriela A. Cividino
Beckman Coulter Argentina S.A.
FARMACEUTICA
M.N. 15202/ M.P. 18093

**BECKMAN
COULTER**

Rolando F. Bornia
APODERADO

Beckman Coulter STEM Software System Guide
B75 78AA March 2022

© 2022 Beckman Coulter, Inc.
All Rights Reserved

Glossary of Symbols is available at
beckman.com/techdocs (PN C24689).

Original Instructions

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

Immunotech and the Immunotech product marks mentioned herein are trademarks or registered trademarks of Immunotech SAS, in the United States and other countries. Immunotech is a Beckman Coulter company.

All other trademarks are the property of their respective owners.

Contact Us

If you have any questions, contact our Customer Support Center.

Worldwide, find us via our website at
www.beckman.com/support/technical.

In the USA and Canada, call us at 1-800-369-0333.

In Austria, call us at +43 0810 300484.

In Germany, call us at +49 02151 333999.

In Sweden, call us at +46 (0)8 564 859 14.

In Netherlands, call us at +31 348 799 815.

In France, call us at +33 0825838306 6.

In the UK, call us at +44 845 600 1345.


In Ireland, call us at +353 (01) 4073082.

In Italy, call us at +39 0295392 456.

In other locales, contact your local Beckman Coulter Representative.

If you are a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

May be covered by one or more patents. - see
www.beckman.com/patents

Gabriela A. Cividino
Beckman Coulter Argentina S.A.
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Revision History

B75978AA, March 2022

Initial Issue

Software version 1.0

This document applies to the latest software listed and higher versions. When a subsequent software version affects the information in this document, a new issue will be released to the Beckman Coulter Web site. For labeling updates, go to www.beckman.com and download the latest version of the manual for your instrument.

Rolando F. Bornia
APODERADO


PN B75978AA


Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

iii

Rolando F. Bornia
APODERADO

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Safety Notice

Read all product manuals and consult with Beckman Coulter-trained personnel before attempting to operate instrument. Do not attempt to perform any procedure before carefully reading all instructions. Always follow product labeling and manufacturer's recommendations. If in doubt as to how to proceed in any situation, [Contact Us](#).

Beckman Coulter, Inc. urges its customers and employees to comply with all national health and safety standards such as the use of barrier protection. This may include, but it is not limited to, protective eyewear, gloves, and suitable laboratory attire when operating or maintaining this or any other automated laboratory analyzer.

Alerts for Warning and Caution

WARNINGS and CAUTIONS alert you as follows:

WARNING

WARNING indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury. May be used to indicate the possibility of erroneous data that could result in an incorrect diagnosis.

CAUTION

CAUTION indicates a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices. May be used to indicate the possibility of erroneous data that could result in an incorrect diagnosis.

CAUTION

System integrity could be compromised and operational failures could occur if:

- This equipment is used in a manner other than specified. Operate the instrument as instructed in the Product Manuals.
- You introduce software that is not authorized by Beckman Coulter into your computer. Only operate your system's computer with software authorized by Beckman Coulter.
- You install software that is not an original copyrighted version. Only use software that is an original copyrighted version to prevent virus contamination.

 CAUTION

If you purchased this product from anyone other than Beckman Coulter or an authorized Beckman Coulter distributor, and, it is not presently under a Beckman Coulter service maintenance agreement, Beckman Coulter cannot guarantee that the product is fitted with the most current mandatory engineering revisions or that you will receive the most current information bulletins concerning the product. If you purchased this product from a third party and would like further information concerning this topic, [Contact Us](#).

 WARNING

Risk of operator injury if:

- All doors, covers, and panels are not closed and secured in place prior to and during instrument operation.
- The integrity of safety interlocks and sensors is compromised.
- Instrument alarms and error messages are not acknowledged and acted upon.
- You contact moving parts.
- You mishandle broken parts.
- Doors, covers, and panels are not opened, closed, removed, and/or replaced with care.
- Improper tools are used for troubleshooting.

To avoid injury:

- Keep doors, covers, and panels closed and secured in place while the instrument is in use.
- Take full advantage of the safety features of the instrument. Do not defeat safety interlocks and sensors.
- Acknowledge and act upon instrument alarms and error messages.
- Keep away from moving parts.
- Report any broken parts to your Beckman Coulter Representative.
- Open/remove and close/replace doors, covers, and panels with care.
- Use the proper tools when troubleshooting.

Contents

Revision History, iii

Safety Notice, v

Alerts for Warning and Caution, v

Introduction, xv

Overview, xv

How to Use Your AQUIOS Manuals, xv

Before You Begin, xvi

About Your Cytometer Software, xvi

About this System Guide, xvi

Conventions, xvii

Symbols, xviii

Graphics Icons, xviii

CHAPTER 1: Use and Function, 1-1

Overview, 1-1

Intended Use, 1-1

Intended User, 1-2

System Components, 1-2

AQUIOS STEM System, 1-2

Unique Load & Go Feature, 1-3

AQUIOS STEM Test, 1-3

Consumables and Supplies, 1-3

Preparation, 1-4

AQUIOS STEM Lysing Solution, 1-4

Monoclonal Antibodies, 1-4

AQUIOS STEM-Count Fluorospheres, 1-5

Acquisition, 1-5

AQUIOS Sheath Fluid, 1-5


Quality Control, 1-5


AQUIOS STEM CD34 Control Cells, 1-5

Flow-Check Fluorospheres, 1-6

Limitations, 1-6

Material Safety Data Sheets (SDS/MSDS), 1-6

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

CHAPTER 2: Operating Principles, 2-1

Overview, 2-1

Principles of AQUIOS STEM Tests, 2-1

Testing with the AQUIOS STEM Kit Reagents, 2-2

Light Source Configuration, 2-2

CD34+HPC Enumeration Panels, 2-3

AQUIOS STEM Tests, 2-3

Analysis Algorithm Overview, 2-4

Introduction to Gating Strategy, 2-4

Gating Strategy: AQUIOS STEM Panel, AQUIOS STEM Duplicate, and

AQUIOS STEM Single, 2-5

Forward Scatter (FS) versus Side Scatter (SS) Plot, 2-6

CD34-PE versus CD45-FITC Plot, 2-6

BEADS versus TIME Plot, 2-7

Side Scatter (SS) versus 7-AAD Plot, 2-7

Side Scatter (SS) versus CD45-FITC Plot, 2-8

Side Scatter (SS) versus CD34-PE Plot, 2-8

Side Scatter (SS) versus CD45-FITC Plot, 2-9

CD45 DIM Side Scatter (SS) versus Forward Scatter (FS) Plot, 2-9

LYMPHS Side Scatter (SS) versus Forward Scatter (FS) Plot, 2-10

Side Scatter (SS) versus CD45-FITC Plot, 2-10

Side Scatter (SS) versus CD34-PE Plot, 2-11

Side Scatter (SS) versus CD45-FITC Plot, 2-11

Side Scatter (SS) versus Forward Scatter (FS) Plot, 2-12

Side Scatter (SS) versus Forward Scatter (FS) Plot, 2-12

Side Scatter (SS) versus 7-AAD Plot, 2-13

CHAPTER 3: Quality Control, 3-1

Overview, 3-1

Quality Control Checks, 3-1

AQUIOS STEM Quality Control Workflow, 3-1

QC Materials, 3-2

Daily QC Overview, 3-2

QC Setup, 3-3

QC Process, 3-3

Flow-Check Fluorospheres - Single Tube Loader, 3-4

STEM Single (AutoLoader) - AQUIOS STEM CD34 L1 and L2 Control Cells, 3-7

Adding Reagents, 3-7

Adding and Using AQUIOS STEM CD34 Control Lots, 3-12

Adding New Lots of AQUIOS STEM CD34 L1 and L2 Control Cells using QR Code, 3-12

Adding New Lots of AQUIOS STEM CD34 L1 and L2 Control Cells Manually, 3-12

Starting AQUIOS STEM CD34 L1 and L2 Control Cell QC With Pre-Existing Lot, 3-15

QC Screen, 3-16

QC Reports, 3-18

AQUIOS STEM (Single Tube Loader) Flow-Check Fluorospheres QC Result Examples, 3-18

AQUIOS STEM Single (Auto Loader) - AQUIOS STEM CD34 Control Cells Level 1 & 2 QC Result Examples, 3-19

AQUIOS STEM Single (Auto Loader) - AQUIOS STEM CD34 Control Cells Level 1 & 2 QC Result Examples, 3-20

Interlaboratory Quality Assurance (IQAP), 3-24

Before Setting Up IQAP Export, 3-24

Setting Up IQAP Export, 3-25

CHAPTER 4: Sample Analysis, 4-1

Overview, 4-1

Test Description, 4-2

AQUIOS STEM Application Workflow, 4-2

Sample preparation, 4-3

Storage Conditions and Stability, 4-4

Non-Frozen Specimen, 4-4

Frozen Specimen Type, 4-4

Reagents, 4-4

AQUIOS STEM Testing Workflow, 4-4

Create a Request, 4-6

Load Sample in the Instrument, 4-8

Options for Running Specimens, 4-10

Logging In Controls, 4-11

Running Samples Using the Autoloader, 4-11

Running Samples Using the Single-tube Loader, 4-11

Running Screen, 4-12

CHAPTER 5: Data Review, 5-1

Overview, 5-1

Review Screen, 5-1

Review Screen Display Areas (left to right), 5-2

Graphic Options, 5-6

AQUIOS STEM Panel, 5-7


Action Toolbox, 5-9


Data Review, 5-9

Review Data Displays, 5-9

Edit Regions for a Review Sample, 5-10

Population Statistics, 5-11

Gabriela A. Clividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N: 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Adjust Compensation Manually, 5-11

Result Screen, 5-15

Using the Result Screen to Select Data for Review, 5-15

Flags And Notifications, 5-16

Flags, 5-17

Notifications, 5-19

AQUIOS STEM Software Observations, 5-20

CHAPTER 6: Performance Characteristics , 6-1

Overview, 6-1

Linearity , 6-1

Limit of Blank/Limit of Detection/Limit of Quantitation, 6-2

Analytical Measuring Ranges, 6-3

Method Comparison, 6-4

Reference Interval, 6-6

Precision, 6-7

Using Control Cells, 6-7

Using Clinical Specimens, 6-7

Specificity, 6-8

Reagent Specificity, 6-8

Quality Control, 6-8

CHAPTER 7: Application Specifications, 7-1

Overview, 7-1

AQUIOS STEM Preparation, 7-1

AQUIOS STEM Sample Preparation, 7-1

AQUIOS STEM Throughput, 7-1

CHAPTER 8: Adjustment and Replacement Procedures, 8-1

Overview, 8-1

AQUIOS STEM Replacement Procedures, 8-1

Replacing the AQUIOS STEM Lysing Solution, 8-1

Replacing the AQUIOS STEM-Count Fluorospheres, 8-3

APPENDIX A: Reports, A-1

Overview, A-1


Reports, A-1


Report Description, A-8

Report Heading, A-8

Displays: Generating Reports, A-8

Result Table Heading, A-8

Gabriela A. Clividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Statistics: AQUIOS STEM Report, A-8

QC Report, A-8

QC Report Description, A-9

QC Result Report Heading, QC Instrument Report Heading, QC
Instrument Drift Report Heading, A-9

QC Result Displays & Statistics (Depend on Test Selected), A-9

QC Result First Table Heading, A-9

QC Result Second Table Heading, QC Instrument Second Table
Heading, QC Instrument Drift Second Table Heading, A-9

QC Reports, A-10

References

Index


Beckman Coulter, Inc.

Customer End User License Agreement

Related Documents

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO

Illustrations

- 2.1 Optical Filter Configuration, 2-2
- 2.2 FS versus SS Plot, 2-6
- 2.3 CD34-PE versus CD45-FITC Plot, 2-6
- 2.4 BEADS versus TIME Plot, 2-7
- 2.5 SS versus 7-AAD Plot, 2-7
- 2.6 SS versus CD45-FITC Plot, 2-8
- 2.7 SS versus CD34-PE Plot, 2-8
- 2.8 SS versus CD45-FITC Plot, 2-9
- 2.9 CD45 SS versus FS Plot, 2-9
- 2.10 LYMPHS SS versus FS Plot, 2-10
- 2.11 SS versus CD45-FITC Plot, 2-10
- 2.12 SS versus CD34-PE Plot, 2-11
- 2.13 SS versus CD45-FITC Plot, 2-11
- 2.14 SS versus FS Plot, 2-12
- 2.15 SS versus FS Plot, 2-12
- 2.16 SS versus 7-AAD Plot, 2-13
- 3.1 Results Tab Screen: Flow-Check Fluorospheres, 3-17
- 3.2 Results Tab: Flow-Check Fluorospheres Report Page 1, 3-18
- 3.3 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 1 (Run with AQUIOS STEM Single), 3-19
- 3.4 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 2 (Run with AQUIOS STEM Single), 3-19
- 3.5 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 1, Page 1 (Run with AQUIOS STEM Single), 3-20
- 3.6 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 1, Page 2 (Run with AQUIOS STEM Single), 3-21
- 3.7 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 2, Page 1 (Run with AQUIOS STEM Single), 3-22
- 3.8 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 2, Page 2 (Run with AQUIOS STEM Single), 3-23
- 3.9 QC Results Tab Screen: Export IQAP, 3-25
- 3.10 Export IQAP control data file, 3-25
- 4.1 Running Screen, 4-12
- 5.1 Review Screen: Display Areas, 5-2
- 5.2 Instrument Designation, 5-2

- 5.3 Auto Refresh, 5-2
- 5.4 Review Screen: Completed Runs List, 5-3
- 5.5 Review Screen: Sample Identifiers and Run Information Area, 5-3
- 5.6 Review Screen: Result Table Area (3 Wells), 5-4
- 5.7 Review Screen: Graphic Options Selection Menu, 5-4
- 5.8 Review Screen: Page Number, 5-4
- 5.9 Review Screen: Well#1 Page 1 (Mobilized Whole Blood), 5-7
- 5.10 Review Screen: Well#1 Page 2 (Mobilized Whole Blood), 5-8
- 5.11 Review Screen: Summary Page 1, 5-8
- 5.12 Action Toolbox, 5-9
- 5.13 Adjust Compensation, 5-13
- 5.14 Over Compensation, 5-14
- 5.15 Under Compensation, 5-14
- 6.1 Linearity: CD34+ cells/ μ L, 6-2
- 6.2 Linearity: CD45+ cells/ μ L, 6-2
- 6.3 Regression Analysis: CD34+ cells/ μ L, 6-5
- 6.4 Regression Analysis: CD34+ %, 6-5
- 6.5 Regression Analysis: CD45+ cells/ μ L, 6-5
- A.1 Report: STEM Panel - Apheresis - Well #1 Page 1, A-2
- A.2 Report: STEM Panel - Apheresis - Well #1 Page 2, A-3
- A.3 Report: STEM Panel - Apheresis - Well #1 Page 3, A-4
- A.4 Report: STEM Panel - Apheresis - Summary Page 1, A-5
- A.5 Report: STEM Panel - Apheresis - Summary Page 2, A-6
- A.6 Report: STEM Panel - Apheresis - Summary Page 3, A-7
- A.7 QC Results: Results View of AQUIOS STEM CD34 Control Cells Level 1, AQUIOS STEM Single, A-10
- A.8 QC Results: AQUIOS STEM CD34 Control Cells Level 1, AQUIOS STEM Single (Page 1), A-11
- A.9 QC Results: AQUIOS STEM CD34 Control Cells Level 1, AQUIOS STEM Single (Page 2), A-11
- A.10 QC Results: Results View of AQUIOS STEM CD34 Control Cells Level 2, AQUIOS STEM Single, A-12
- A.11 QC Results: AQUIOS STEM CD34 Control Cells Level 2, AQUIOS STEM Single (Page 1), A-13
- A.12 QC Results: AQUIOS STEM CD34 Control Cells Level 2, AQUIOS STEM Single (page 2), A-13

 Gabriela A. Cividino
 Beckman Coulter Argentina S.A.
 FARMACÉUTICA
 M.N. 15202/ M.P. 18093

 Rolando F. Bornia
 APODERADO

Tables

1.1	System Parameters, 1-1
2.1	Parameter Results Provided by STEM Panel , STEM Duplicate, or STEM Single Analysis, 2-3
2.2	Equivalence Between Primary Region Name and Replicated Region Name, 2-5
3.1	AQUIOS STEM QC Workflow, 3-1
3.2	Handling AQUIOS STEM Reagent Caps, 3-7
3.3	General Table on QC Screens: Defined, 3-17
4.1	Test Descriptions, 4-2
4.2	Processing Samples Using the STEM Application, 4-3
4.3	Minimum Specimen Volume Required per Specimen Tube, 4-3
5.1	Review Screen: Summary Area Display (AQUIOS STEM, Mobilized Whole Blood), 5-5
5.2	Review Screen: Summary Area Display, Page Number 1 (AQUIOS STEM Panel, Mobilized Whole Blood), 5-6
5.3	Population Statistics Recommended Use, 5-11
5.4	AQUIOS STEM Specific Flags, 5-17
5.5	AQUIOS STEM Specific Notifications, 5-19
6.1	LoB, LoD and LLoQ Results, 6-3
6.2	Analytical Measuring Ranges, 6-3
6.3	Regression statistics for all specimen types combined, 6-4
6.4	Method Comparison study results of AQUIOS STEM system compared with predicate method, 6-4
6.5	Regression statistics by specimen type, 6-6
6.6	Reference Intervals, 6-6
6.7	AQUIOS STEM CD34 Control Cells Level 1, 6-7
6.8	AQUIOS STEM CD34 Control Cells Level 2, 6-7
6.9	Repeatability Precision with clinical specimens, 6-8

Introduction

Overview

This introduction contains the following information:

- [How to Use Your AQUIOS Manuals](#)
- [Before You Begin](#)
- [About this System Guide](#)
- [Conventions](#)
- [Symbols](#)
- [Graphics Icons](#)

How to Use Your AQUIOS Manuals

Your AQUIOS CL Flow Cytometry system includes the manuals listed below:

- The **AQUIOS CL Flow Cytometer Instructions For Use** manual provides information for the day-to-day running of your AQUIOS CL Flow Cytometer and AQUIOS system software. It also provides procedures for cleaning the instrument and replacement procedures. You can find detailed step-by-step procedures for daily startup for running and reviewing quality control (QC) data, running samples, analyzing data, printing reports, and shutting the instrument down. The Instructions for Use manual also contains information on safety, troubleshooting, error messages, and operation principles. It contains a glossary of terms, a list of abbreviations and acronyms, and the warranty information for the instrument.
- This **AQUIOS STEM System Guide** provides reference information and instructions for using the AQUIOS System with the STEM Test. Refer to the glossary of terms in the AQUIOS CL Flow Cytometer Instructions for Use.
- The **AQUIOS STEM Host Transmission** manual provides the information needed to program the transmission interface between the AQUIOS CL Flow Cytometer and the laboratory's host computer.

To quickly determine which manual to read for the information needed, refer to [Related Documents](#).

Before You Begin

IMPORTANT The AQUIOS STEM System Guide is not an inclusive manual. It is imperative that an Operator is familiar with the information in the AQUIOS CL Flow Cytometer Instructions for Use manual, and the AQUIOS STEM Reagent and Consumables Instructions for Use and is capable of performing basic operations.

About Your Cytometer Software

Your local Beckman Coulter representative will install AQUIOS System Software and AQUIOS STEM Software.

About this System Guide

This AQUIOS STEM System Guide provides reference information and instructions for using the AQUIOS STEM System.

Use this guide in conjunction with the AQUIOS CL Flow Cytometer Instructions for Use manual.

Use this guide with the applicable reagent Instructions for Use listed under [Consumables and Supplies](#). These Instructions for Use contain reagent-specific information not contained in this guide.

The information in the AQUIOS STEM System Guide is organized as follows:

[CHAPTER 1, Use and Function](#)

Provides the intended use of the AQUIOS CL Flow Cytometer system as well as information concerning system components (Cytometer, STEM Test, reagents and consumables) and are referred to as AQUIOS STEM System application. Limitations are also provided.

[CHAPTER 2, Operating Principles](#)

Contains a description of the principles used in AQUIOS STEM testing, and an overview of the analysis algorithm.

[CHAPTER 3, Quality Control](#)

Provides information on running QC samples to verify the performance of the instrument before analyzing non-QC samples.

[CHAPTER 4, Sample Analysis](#)

Provides information for storage conditions and stability for specimens and reagents, information on the preparation required to successfully run samples using AQUIOS STEM System reagents, running samples using the autoloader, and running samples using the Single-tube Loader.

[CHAPTER 5, Data Review](#)

Provides information on reviewing sample results, including flagged results.

[CHAPTER 6, Performance Characteristics](#)

Provides information on clinical performance characteristics such as expected values, reference

ranges, linearity, accuracy, precision, and the specificity of AQUIOS STEM monoclonal antibodies and quality control materials.

CHAPTER 7, Application Specifications

Provides information for sample preparation, throughput, and sample rate.

CHAPTER 8, Adjustment and Replacement Procedures

Provides information for adjustment or replacement procedures to maintain AQUIOS STEM functionality.

This manual also includes an appendix, references, and an index.

Conventions

This guide uses the following conventions:

- Throughout this manual, the AQUIOS STEM System application (Cytometer, STEM Test, reagents and consumables) may also be referred to as system components.
- **Blue** text indicates that you can select on the text to access related information.
- **Bold font** indicates a screen icon, menu item, or software option on the Workstation screen.
- *Italics font* indicates screen text displayed on the Workstation.
- The term 'select' is used to indicate either one or both of the following actions:
 - to tap or touch your finger to the touch screen of the computer
 - to click with a mouse
- Tabs in the software are blue when not selected and inactive. Tabs in the software turn green when selected and active.

IMPORTANT IMPORTANT is used for comments that add value to the step or procedure being performed. Following the advice in the IMPORTANT adds benefit to the performance of a piece of equipment or to a process.

NOTE NOTE is used to call attention to notable information that should be followed during use or maintenance of this equipment.

Symbols

Safety symbols alert you to potentially dangerous conditions. These symbols, together with text, apply to specific procedures and appear as needed throughout this guide.

Symbol	Warning Condition	Action
	Consider all materials (specimens, reagents, controls, and monoclonal antibodies) and areas these materials come into contact with as being potentially infectious.	Wear appropriate barrier protection and follow safe laboratory procedures when handling any material in the laboratory.

Graphics Icons

All graphics, including screens and printouts, are for illustration purposes only and must not be used for any other purpose.

CHAPTER 1

Use and Function

Overview

This chapter contains information about:

- [Intended Use](#)
- [System Components](#)
- [Limitations](#)

Intended Use

The AQUIOS STEM System, comprised of AQUIOS STEM Software for the AQUIOS CL Flow Cytometry System, AQUIOS STEM-Kit Reagents, AQUIOS STEM CD34 Control Cells and Flow-Check Fluorospheres, is an in vitro diagnostic medical device intended to be used by laboratory professionals for the detection of the parameters in the specimen types shown in [Table 1.1, System Parameters](#) on the automated AQUIOS CL Flow Cytometry System. This is a quantitative assay used for the identification and enumeration of the cell populations specified in [Table 1.1](#) in patients who are preparing for autologous or allogeneic hematopoietic progenitor cell transplantation, and in donors who undergo mobilization or collection schemes for autologous or allogeneic hematopoietic progenitor cell transplantations. Refer to [Table 1.1](#) for details of the supported specimen types and anticoagulants.

Table 1.1 System Parameters

Specimen Type	Anticoagulant						Kits/Parameters
	K ₂ EDTA	K ₃ EDTA	Heparin	Heparin ACD-A Mix	ACD-A	CPD	AQUIOS STEM Kit
Whole Blood (WB)	X	X	X				CD34 Percent CD34 Absolute Count CD45 Absolute Count
Mobilized Whole Blood (MB)	X	X	X				
Fresh or Thawed Bone Marrow (BM)	X	X	X	X			
Fresh or Thawed Apheresis Product (AP)				X	X		
Fresh or Thawed Cord Blood (CB)					X	X	

AQUIOS STEM-Kit reagents consist of a CD45-FITC/CD34-PE murine monoclonal antibody reagent, a corresponding negative control (CD45-FITC/CD34-CTRL), an absolute count reagent (AQUIOS STEM-Count Fluorospheres), a cell viability reagent (7-AAD), and a ready-to-use lysing reagent (AQUIOS STEM Lysing Solution).

AQUIOS STEM CD34 Control Cells are liquid preparations of stabilized human leukocytes for the verification of the parameters CD34 and CD45 as part of the AQUIOS STEM System. Each kit contains 2 levels of CD34 with approximately 10 CD34+ cells/ μ L (level 1) and approximately 30 CD34+ cells/ μ L (level 2).

AQUIOS CL Flow Cytometry system software may be run on an AQUIOS workstation for off-line analysis of results generated by the AQUIOS CL Flow Cytometer with the monoclonal antibody reagents listed above. The off-line analysis must be performed in accordance with the product labeling. The independent computer workstation must be approved to be used with the AQUIOS software.

NOTE An off-line analysis workstation cannot be used for acquisition.

Intended User

This product is intended for laboratory professional use.

System Components

AQUIOS STEM System

IMPORTANT [Contact us](#) if a product is damaged upon receipt.

The AQUIOS STEM system applies the principles of flow cytometry to analyze a stained and lysed sample to identify various cellular populations as determined by the specific monoclonal antibodies and fluorochromes used.

Automatic operation is initiated by adding requests and loading a cassette containing specimen tubes in the autoloader or a specimen tube in the Single-tube Loader. Samples are automatically processed according to these requests. The sample is stained and incubated, the red blood cells are lysed using the AQUIOS STEM Lysing Solution. The white blood cells are analyzed on the AQUIOS CL Flow Cytometry system with the AQUIOS STEM Tests.

The AQUIOS Smart Scheduler tracks requests as they proceed from sample preparation to analysis. All AQUIOS reagents are continuously monitored by the system to track reagent consumption and location changes. The acquired data is stored and provides immediate data access for laboratories.

AQUIOS STEM reagents use a unique barcode identity for tracking reagent type, lot number, container number, expiration dates, and remaining reagent levels. The reagent consumption is monitored by the system as the samples are processed.

Unique Load & Go Feature

The AQUIOS CL Flow Cytometer is a quantitative automated analyzer that performs the STEM diagnostic applications in a “no-wash” sample preparation process. Since this system is intended to be an automated analyzer with hands-off processing of samples from specimen introduction to results reports, it is referred to as a Load & Go flow cytometer. The AQUIOS System Software and AQUIOS STEM Tests and Quality Control Reagents do not require user verification of standardization of light scatter, and fluorescence intensities. AQUIOS CL instrument settings are monitored specifically for the AQUIOS STEM System by running Flow-Check Fluorospheres.

This “closed” system will only run tests preconfigured by Beckman Coulter and does not allow for user-designed protocols or panels.

AQUIOS STEM Test

The AQUIOS STEM Tests are used:

- with quality control reagents on an AQUIOS CL Flow Cytometry system with AQUIOS System Software and AQUIOS STEM Software for system quality control, and
- with AQUIOS STEM Kit Reagent and AQUIOS System Software with AQUIOS STEM Software for automated sample analysis.

The reagents used in the AQUIOS STEM System allow for automated sample preparation and analysis at room temperature without a negative effect on cell viability within the time window required for automated sample preparation and analysis.

For more information, refer to the associated Reagent Instructions for Use.

Consumables and Supplies

The AQUIOS CL Flow Cytometer is optimized to operate with AQUIOS reagents. The stated analytical characteristics and specifications cited in this manual can only be guaranteed by using AQUIOS reagents and the AQUIOS 96-Deep Well Plate.

STEM sample preparation is optimized to operate using barcoded 96-deep conical-shaped well plates. Each well holds up to 600 µL. The plates are made of virgin polypropylene.

Additionally, the AQUIOS reagents (with the exception of the AQUIOS Sodium Hypochlorite Solution and Flow-Check Fluorospheres) and the AQUIOS Deep Well Plate are labeled with bar codes for smart-track monitoring to help minimize potential errors. The system will only run with AQUIOS reagents and the AQUIOS 96-Deep Well Plate.

NOTE AQUIOS reagents are for use on the AQUIOS CL Flow Cytometer only. Do not dilute, aliquot, or freeze the reagents. The product should be used according to labeled instructions.

NOTE Do not use reagents beyond the expiration date printed on the associated labeling and always use good laboratory practices when handling these reagents.

NOTE The AQUIOS CL Flow Cytometer system tracks the AQUIOS STEM-Count on-board time and expiration date. Ensure enough stability time is remaining to complete the day's runs by reviewing the information in the reagents window after startup or when new reagents are added.

NOTE The AQUIOS CL Flow Cytometer system tracks the AQUIOS Cleaning Agent for the percent remaining and expiration date.

NOTE The AQUIOS CL Flow Cytometer system does not track the AQUIOS Sodium Hypochlorite Solution. You must track the expiration date of the AQUIOS Sodium Hypochlorite Solution upon use.

NOTE The AQUIOS CL Flow Cytometer system does not track the Flow-Check Fluorospheres for the percent remaining.

For more information, refer to the associated reagent Instructions for Use on the Beckman Coulter Web site at www.beckman.com or [Contact us](#).

Preparation

AQUIOS STEM Lysing Solution

The AQUIOS STEM Lysing Solution is a highly specific, ready to use and gentle lysing agent. The AQUIOS STEM Lysing Solution lyses red blood cells in samples sufficiently to perform monoclonal antibody panel analysis without interference from red blood cells and allows analysis without damage to white blood cells.

NOTE The AQUIOS STEM Lysing Solution does not contain fixatives.

AQUIOS STEM Lysing Solution is housed inside the Cytometer behind the reagent door (the larger of the two doors on the Cytometer front cover).

When the AQUIOS STEM Lysing Solution bottle is replaced, the system looks for a specific bar code to ensure that the correct reagent is being loaded and tracked.

NOTE To replace the AQUIOS STEM Lysing Solution, refer to [Replacing the AQUIOS STEM Lysing Solution](#).

Risk of erroneous results. To prevent erroneous results, ensure the lot number on the AQUIOS STEM Lysing Solution is the same as the lot number of the AQUIOS STEM Kit.

Monoclonal Antibodies

The AQUIOS STEM Kit comprises two separate murine monoclonal antibody panel reagents as well as a viability dye. Each murine monoclonal antibody is conjugated to a specific fluorochrome and specific for a different cell surface antigen:

- CD45-FITC/CD34-PE reagent
- CD45-FITC Isoclonic Control-PE reagent
- 7-AAD Viability Dye

Specific cell surface staining is accomplished by incubating samples of a biological specimen with the two-color CD45-FITC / CD34-PE reagent.

An additional test of the same sample is stained with CD45-FITC / Isoclonic Control-PE reagent to check the non-specific binding of the CD34-PE monoclonal antibody. 7-AAD Viability Dye, a nucleic

acid dye that binds to accessible base pairs (cellular DNA), is used to distinguish between viable and nonviable cells.

The test depends on the ability of a monoclonal antibody to bind to the surface of cells expressing discrete antigenic determinants. Specific cell staining is accomplished by incubating the sample with the monoclonal antibody reagents.

AQUIOS STEM-Count Fluorospheres

AQUIOS STEM-Count Fluorospheres is a fluorescent microsphere reagent for direct determination of the absolute counts of white blood cell subsets in biological specimens. For more information, refer to [Replacing the AQUIOS STEM-Count Fluorospheres](#).

Acquisition

AQUIOS Sheath Fluid

AQUIOS Sheath Fluid is a cyanide-free, isotonic buffered saline solution. It is used to rinse components between sample analysis, and provide a sheath stream to transport the cells in single file through the flow cell.

For more information, refer to Chapter 1, System Overview of the AQUIOS CL Flow Cytometer Instructions for Use manual.

Quality Control

AQUIOS STEM CD34 Control Cells

AQUIOS STEM CD34 Control Cells are liquid preparations of stabilized human leukocytes (lymphocytes, monocytes and granulocytes) and erythrocytes that have lysing, light scatter, antigen expression and antibody staining properties representative of those found in human whole blood specimens.

AQUIOS STEM CD34 Control Cells Level 1 and Level 2 are in vitro diagnostic medical devices intended to be used by laboratory professionals as process controls for the verification of the parameters CD34 and CD45 as part of the AQUIOS STEM System on the automated AQUIOS CL Flow Cytometry System. They provide positive cell controls that are processed in the same manner as a biological sample. They also verify the methods used for staining targeted cells, the lysing of red blood cells, and the analysis of samples by the AQUIOS CL Flow Cytometer.

NOTE When these controls are used, scan the barcode to ensure the correct control is being loaded and tracked. Select the test as a QC request to ensure you are tracking results via the Levey-Jennings charts.

For more information, refer to the associated reagent Instructions for Use.

Flow-Check Fluorospheres

Flow-Check Fluorospheres are an assayed suspension of fluorospheres (fluorescent microspheres) used for daily verification of a flow cytometer's optical alignment and fluidics system.

Limitations

1. Accurate and reproducible results will be obtained, provided the procedures used are in accordance with the AQUIOS CL Flow Cytometer Instructions for Use manual, the AQUIOS STEM System Guide, associated reagent Instructions for Use, and good laboratory practices.
2. When modifying the algorithm-generated results through user intervention (regions and quadrants), it is the user's responsibility to follow the instructions appropriately when adjusting the compensation and the region's algorithm generated results. Results require review by a qualified laboratory professional.
3. Results obtained with flow cytometry may be erroneous if the laser is misaligned or the gates and regions are improperly set.
4. Certain patients may present special problems due to altered or very low numbers of certain cellular populations.
5. Specimens must be maintained at the appropriate temperature prior to placement on the system. See [CHAPTER 4, Storage Conditions and Stability](#).

Material Safety Data Sheets (SDS/MSDS)

To obtain an SDS or MSDS for Beckman Coulter reagents used on the AQUIOS system:

- On the Internet, go to www.beckman.com:
 1. Select **Safety Data Sheets (SDS/MSDS)** from the Support menu.
 2. Follow the instructions on the screen.
 3. [Contact us](#) if you have difficulty locating the information.
- If you do not have Internet access, [contact us](#).

CHAPTER 2

Operating Principles

Overview

This chapter contains information about:

- Principles of AQUIOS STEM Tests
- Light Source Configuration
- CD34+HPC Enumeration Panels
- Analysis Algorithm Overview

Principles of AQUIOS STEM Tests

AQUIOS STEM tests are based on the ability of monoclonal antibodies to bind to the surface of cells expressing discrete, or a combination of discrete, antigenic determinants.

The flow cytometer measures light scatter, and the fluorescence of cells. This allows for the separation of the population of interest within the electronic window defined on various histograms, which correlates with different parameters.

A series of dual histograms are used in the gating stage. The fluorescence of the delimited cells is analyzed in order to distinguish the positively-stained events from the unstained ones.

The AQUIOS STEM Software establishes the gating regions for the negative control/CD45 sample in an identical manner as the gating regions are established for the CD34/CD45 samples.

The gating strategy allows detection of leukocytes which are further categorized into various parameters. Refer to Parameter Results in [Table 2.1](#) and [Table 2.2](#), for a listing of parameters.

The parameters, triggers, discriminators, and regions are specific to the test.

AQUIOS STEM Software has stop criteria for 100 CD34+ cells AND 75,000 CD45+ cells OR 300 seconds. In the case where cell concentration prevents the collection of 100 CD34+ cells, a notification will inform the user that less than 100 CD34+ cells are acquired.

75000 is the maximum number of events and 30000 cells/ μ L is the maximum WBC concentration.

AQUIOS STEM Software uses a total of 4 discriminators in order to ensure all CD34+ cells are captured: Forward Scatter, Side Scatter, Fluorescence 1 (CD45-FITC), and Fluorescence 2 (CD34-PE). An event is acquired if the value of only one parameter is above the defined discriminator.

CD34+HPC Enumeration Panels

AQUIOS STEM Tests

AQUIOS STEM Tests are using the AQUIOS STEM Kit reagents containing monoclonal antibodies for the simultaneous identification and enumeration of viable absolute count of CD34+ HPC/ μL , viable absolute count of CD45+/ μL , and percentage of viable CD34+ HPC.

Three AQUIOS STEM tests are available as part of the AQUIOS STEM menu:

- Stem Panel: Run tests in duplicate plus negative control (3 wells)
- Stem Duplicate: Run tests in duplicate (no negative control; 2 wells)
- Stem Single: Single test (1 well)

Refer to [Table 4.1](#) for a description of each test.

[Table 2.1](#) contains the parameter results generated when a STEM Single, STEM Duplicate, or STEM Panel is used for analysis.

Table 2.1 Parameter Results Provided by STEM Panel , STEM Duplicate, or STEM Single Analysis

STEM Panel	STEM Duplicate	STEM Single
Viable CD34+ HPC Count Well 1 (cells/ μL)	Viable CD34+ HPC Count Well 1 (cells/ μL)	Viable CD34+ HPC Count (cells/ μL)
Viable CD34+ HPC Count Well 2 (cells/ μL)	Viable CD34+ HPC Count Well 2 (cells/ μL)	Viable Leukocytes Count (cells/ μL)
AVG Viable CD34+ HPC Count (cells/ μL)	AVG Viable CD34+ HPC Count (cells/ μL)	Viable CD34+ HPC % (of Viable Leukocytes)
Viable CD34+ HPC % Diff (Well 1 vs Well 2)	Viable CD34+ HPC % Diff (Well 1 vs Well 2)	DF: Dilution Factor
Viable Leukocytes Count Well 1 (cells/ μL)	Viable Leukocytes Count Well 1 (cells/ μL)	Viable CD34+ HPC Count (cells/ μL) (x DF)
Viable Leukocytes Count Well 2 (cells/ μL)	Viable Leukocytes Count Well 2 (cells/ μL)	Viable Leukocyte Count (cells/ μL) (x DF)
AVG Viable Leukocytes Count (cells/ μL)	AVG Viable Leukocytes Count (cells/ μL)	HV: Harvest Volume (mL)
Viable CD34+ HPC Well 1 % (of Viable Leukocytes)	Viable CD34+ HPC Well 1 % (of Viable Leukocytes)	Viable CD34+ HPC Count (x DF x HV)
Viable CD34+ HPC Well 2 % (of Viable Leukocytes)	Viable CD34+ HPC Well 2 % (of Viable Leukocytes)	BW: Body Weight (kg)
AVG Viable CD34+ HPC % (of viable Leukocytes)	AVG Viable CD34+ HPC % (of viable Leukocytes)	Viable CD34+ HPC Count per kg (x DF x HV / BW)
DF: Dilution Factor	DF: Dilution Factor	CAL
AVG Viable CD34+ HPC Count (cells/ μL) (x DF)	AVG Viable CD34+ HPC Count (cells/ μL) (x DF)	

Table 2.1 Parameter Results Provided by STEM Panel , STEM Duplicate, or STEM Single Analysis (Continued)

STEM Panel	STEM Duplicate	STEM Single
AVG Viable Leukocyte Count (cells/ μ L) (x DF)	AVG Viable Leukocyte Count (cells/ μ L) (x DF)	
HV: Harvest Volume (mL)	HV: Harvest Volume (mL)	
AVG Viable CD34+ HPC Count (x DF x HV)	AVG Viable CD34+ HPC Count (x DF x HV)	
BW: Body Weight (kg)	BW: Body Weight (kg)	
AVG Viable CD34+ HPC Count per kg (x DF x HV / BW)	AVG Viable CD34+ HPC Count per kg (x DF x HV / BW)	
CAL	CAL	
Viable CD34+ HPC Count Negative Ctrl (cells/ μ L)		

NOTE Refer to [Table 5.3](#) for Population Statistics.

Analysis Algorithm Overview

The AQUIOS STEM algorithm is the software analysis component that works in conjunction with AQUIOS STEM Kit reagents to automatically identify and enumerate viable Leukocytes (CD45+ populations) and viable CD34+ HPC cells according to ISHAGE guidelines. The algorithm combines information from AQUIOS parameters (FS, SS and Fluorescence FL1-FL4 to automatically generate gates and regions to accurately identify populations of interests.

Introduction to Gating Strategy

Software algorithms automatically detect the subset populations of interest to calculate and report the results for the predefined tests. The software generates flags and notifications after analyzing samples.

The hallmark of the ISHAGE Guidelines is a sequential gating strategy that derives the number of CD34+ cells from viable leukocytes. While the starting point of the sequential ISHAGE gating strategy, i.e. viable leukocytes, facilitates the correct quantification of rare CD34+ cell populations, it leads to the fact that as an end point all enumerated CD34+ cells are automatically viable, and it can be difficult to directly calculate the percentage of viable CD34+ cells out of all CD34+ cells. In addition, it may be challenging to retrieve the total number of CD45+ leukocytes from the applied gating strategy in order to assess overall specimen viability.

In order to overcome this challenge, AQUIOS STEM Software employs an additional background gating strategy starting with all leukocytes, including dead cells. This “replicated” gating strategy follows the same gating logic as the classical ISHAGE gating, but enables for the direct enumeration of all cells, as well as a direct comparison of total vs viable cells for selected cell populations. The additional parameters following this gating logic are marked with double asterisks (**) as Population Statistics.

To maintain the consistency of region placement, certain regions are replicated in more than one plot. Refer to [Table 2.2](#) for the matching plots and regions (primary vs. replicated) for three STEM tests.

Table 2.2 Equivalence Between Primary Region Name and Replicated Region Name

Primary Regions Name	Plot / Figure	Replicated Region Name	Plot / Figure
LEUKS	Figure 2.6	WBC	Figure 2.11
LYMPHS		ALL LYMPH	
CD34 POS	Figure 2.7	ALL CD34POS	Figure 2.12
CD45 DIM	Figure 2.8	ALL CD45 DIM	Figure 2.13
!FS LYMPH	Figure 2.10	ALL FS LYMPH	Figure 2.15
!CD34+ HPC	Figure 2.10	!CD34+ HPC	Figure 2.9
		ALL CD34+ HPC	Figure 2.14
!VIABLE	Figure 2.5	!CD34 VIABILITY	Figure 2.16

Gating Strategy: AQUIOS STEM Panel, AQUIOS STEM Duplicate, and AQUIOS STEM Single

The gating strategy is the same for the three AQUIOS STEM tests. The graphs below ([Figure 2.2](#) through [Figure 2.16](#)) show the gating of a fresh peripheral mobilized whole blood stained with CD45-FITC, CD34-PE, and 7AAD. The order of the graphs is aligned with the sequential gating strategy.


The first series of plots ([Figure 2.2](#) through [Figure 2.10](#)) show the identification and quantification of the viable cellular population.

The second series of plots ([Figure 2.11](#) through [Figure 2.16](#)) show the identification and quantification of the total (viable + non-viable) cellular population.

Forward Scatter (FS) versus Side Scatter (SS) Plot

The All Data Points gate is used to identify all cells of interest and eliminate debris. This FS/SS dot plot identifies the 'Cells and Beads' region consisting of leukocytes, and removes debris (lower left). Refer to [Figure 2.2](#).


Figure 2.2 FS versus SS Plot

CD34-PE versus CD45-FITC Plot

The All Data Points gate is used to identify all cells of interest, and eliminate debris. This CD34-PE versus CD45-FITC dot plot identifies the 'ALL BEADS' region consisting all AQUIOS STEM-Count Fluorospheres. Refer to [Figure 2.3](#).


Figure 2.3 CD34-PE versus CD45-FITC Plot

BEADS versus TIME Plot

From the *ALL BEADS* region, a BEADS/TIME dot plot is created to accurately identify AQUIOS STEM-Count Fluorospheres single events needed for the absolute count calculation of cells. The duplicates are removed. Refer to [Figure 2.4](#).

Figure 2.4 BEADS versus TIME Plot


Side Scatter (SS) versus 7-AAD Plot

The *Cells and No Beads* region is defined as *Cells and Beads* AND NOT *All Beads*. This dot plot identifies the viable events *VIABLE* gate. The calculation of the ****Sample Viability %** parameter result (population statistic) is derived from this region. Refer to [Figure 2.5](#).

NOTE AQUIOS STEM CD34 Control Cells are analyzed with the same gating scheme as the actual sample material. However, for AQUIOS STEM CD34 Control Cells, since they are stabilized cells, the rectangular region '*VIABLE*' on this plot is enlarged to gate total events stained by 7-AAD. Refer to the AQUIOS STEM CD34 Control Cells Instructions for Use.


Figure 2.5 SS versus 7-AAD Plot

Side Scatter (SS) versus CD45-FITC Plot

From the *!Viable* region, a SS/CD45-FITC dot plot is created to define the viable Leukocytes *LEUKS* gate, and the viable Lymphocytes *LYMPHS* gate. The calculation of Viable Leukocytes Count (cells/ μ L) parameter result is derived from this region. The red events are CD34+HPC+ back-scattered from the *!CD34+ HPC* region. Refer to [Figure 2.6](#).


Figure 2.6 SS versus CD45-FITC Plot

Side Scatter (SS) versus CD34-PE Plot

From the *LEUKS* region, a SS/CD34-PE dot plot is created to define the viable CD34 positive cells *CD34 POS* gate. Refer to [Figure 2.7](#).

Figure 2.7 SS versus CD34-PE Plot

G.A.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


R.F.B.

Rolando F. Bornia
APODERADO

Side Scatter (SS) versus CD45-FITC Plot

From the *CD34 POS* region, a SS/CD45-FITC dot plot is created to distinguish the *CD34+HPC CD45 DIM* from the lymphocytes populations by their CD45 expression level. The *CD34+HPC* shows a lower expression of CD45 compared to Lymphocytes cells. Refer to [Figure 2.8](#).


Figure 2.8 SS versus CD45-FITC Plot

CD45 DIM Side Scatter (SS) versus Forward Scatter (FS) Plot

From the *CD45 DIM* region, a SS/FS dot plot is created to identify a cluster of viable events *!CD34+ HPC* gate, meeting all the fluorescence and light scatter criteria of ISHAGE Guidelines for viable *CD34+ HPC*. The calculation of Viable *CD34+ HPC* Count (cells/ μ L) parameter result is derived from this region. The *!CD34+ HPC* region is linked to the *!FS LYMPH* region. Refer to [Figure 2.9](#).

Figure 2.9 CD45 SS versus FS Plot


NOTE The Viable *CD34+ HPC* % (of Viable Leukocytes) parameter result is calculated as follows:

$$\frac{\text{Viable CD34 + HPC Count (cells}/\mu\text{L)}}{\text{Viable Leukocytes Count (cells}/\mu\text{L)}} = (\text{Viable CD34 + HPC\% (of Viable Leukocytes)})$$

LYMPHS Side Scatter (SS) versus Forward Scatter (FS) Plot

From the *LYMPHS* region, a SS/FS dot plot is created to identify a cluster of viable events *!FS LYMPH* gate, meeting all the fluorescence and light scatter criteria for viable Lymphocyte population. The *!FS LYMPH* region is linked to the *!CD34+ HPC* region. Refer to [Figure 2.10](#).


Figure 2.10 LYMPHS SS versus FS Plot

Side Scatter (SS) versus CD45-FITC Plot

The *Cells and No Beads* region is defined as *Cells and Beads* AND NOT *All Beads*. This dot plot identifies the total (viable + non-viable) leukocytes *WBC* gate, and the total (viable + non-viable) Lymphocytes *ALL LYMPH* gate. The calculation of the ****WBC Count (cells/ μ L)** parameter result (population statistic) is derived from this region. Refer to [Figure 2.11](#).


Figure 2.11 SS versus CD45-FITC Plot

Side Scatter (SS) versus CD34-PE Plot

From the WBC region, a SS/CD34-PE dot plot is created to define the total (viable + non-viable) CD34 positive cells *ALL CD34 POS* gate. Refer to [Figure 2.12](#).


Figure 2.12 SS versus CD34-PE Plot

Side Scatter (SS) versus CD45-FITC Plot

From the *ALL CD34 POS* region, a SS/CD45-FITC dot plot is created to distinguish the total (live + non-viable) CD34+HPC *ALL CD45 DIM* gate from the total (live + non-viable) lymphocytes populations by their CD45 expression level. The CD34+HPC shows a lower expression of CD45 compared to Lymphocytes cells. Refer to [Figure 2.13](#).


Figure 2.13 SS versus CD45-FITC Plot

Side Scatter (SS) versus Forward Scatter (FS) Plot

From the *ALL CD45 DIM* region, a SS/FS dot plot is created to identify a cluster of total (viable + non-viable) events *ALL CD34+ HPC* gate, meeting all the fluorescence and light scatter criteria of ISHAGE Guidelines for total (viable + non-viable) CD34+ HPC. The calculation of ****Total CD34+ HPC Count** (cells/ μ L) parameter result (population statistic) is derived from this region. Refer to [Figure 2.14](#).


Figure 2.14 SS versus FS Plot

Side Scatter (SS) versus Forward Scatter (FS) Plot

From the *ALL LYMPH* region, a SS/FS dot plot is created to identify a cluster of total (viable + non-viable) events *ALL FS LYMPH* gate, meeting all the fluorescence and light scatter criteria for total (viable + non-viable) Lymphocyte population. Refer to [Figure 2.15](#).


Figure 2.15 SS versus FS Plot

Side Scatter (SS) versus 7-AAD Plot

From the *ALL CD34+ HPC* region, a SS/7AAD dot plot is created to distinguish viable CD34+ HPC *!CD34 Viability* gate, out of the total (viable + non-viable) CD34+ HPC. The calculation of the ***CD34+ Viability%* (of total CD34+) parameter result (population statistic) is derived from this region. Refer to [Figure 2.16](#).

Figure 2.16 SS versus 7-AAD Plot

Rolando F. Bornia
APODERADO

Gabriela A. Cividino
Beckman Coulter Argentina S.A.
FARMACÉUTICA
M.N. 15202/ M.P. 18093

CHAPTER 3

Quality Control

Overview

Perform daily QC samples to verify the performance of the instrument before analyzing non-QC samples. This includes both Flow-Check Fluorospheres and running AQUIOS STEM CD34 Control Cells. Daily Quality Control is a critical component to ensure the system's performance for the AQUIOS STEM Application.

This chapter contains information about:

- Quality Control Checks
- AQUIOS STEM Quality Control Workflow
- QC Materials
- Daily QC Overview
- QC Setup
- QC Screen

Quality Control Checks

AQUIOS STEM Quality Control Workflow

Table 3.1 AQUIOS STEM QC Workflow

Flow-Check Fluorospheres - Single Tube Loader	
Flow Check Fluorospheres - Alignment	Flow-Check Fluorospheres
Daily QC - AQUIOS STEM CD34 L1 and L2 Control Cells	
Add Reagents	AQUIOS STEM Lysing Solution
	AQUIOS STEM-Count Fluorospheres
	CD45-FITC/CD34-PE and 7-AAD Viability Dye
Load AQUIOS STEM CD34 L1 and L2 Control Cells	AQUIOS STEM CD34 Level 1 Control Cells
	AQUIOS STEM CD34 Level 2 Control Cells

The Quality Control workflow is as follows:

1. Run the Flow-Check Fluorospheres
2. Run the AQUIOS STEM CD34 Control Cells Level 1 and Level 2

The control data is automatically populated into the Levey-Jennings charts of the **QC Screen > Results**.

The AQUIOS system results are displayed and flagged based on the generated values. In general, the system displays fewer digits than that which is generated.

Results are rounded using standard rounding procedures. As a reminder, to round off decimals:

1. Find the place value you want (the 'rounding digit') and look at the digit to the immediate right of it.
2. If the digit to the right of the place value is less than 5, keep the original place value, and drop all digits to the right of the place value.
3. If the digit to the right of the place value is 5 or greater, add one to the place value, and drop all digits to the right of the place value.

For control results, if the assay range is, for example, 20-25, a result of 20 could be flagged as being out of range low, if the generated result is actually 19.9.

The same applies when comparing sample results to the normal and action ranges.

QC Materials

The following Beckman Coulter reagents are required to perform the daily QC for the AQUIOS STEM application:


- AQUIOS STEM Kit
 - AQUIOS STEM Lysing Solution
 - AQUIOS STEM-Count Fluorosphere
 - AQUIOS STEM CD45-FITC/CD34-PE
 - AQUIOS STEM CD45-FITC / IsoClonic Control-PE
 - AQUIOS STEM 7-AAD Viability Dye
- AQUIOS STEM CD34 Control Cells
 - Level 1 (L1) Control Cells
 - Level 2 (L2) Control Cells
- Flow-Check Fluorospheres

Daily QC Overview

Daily Quality Control is a critical component of ensuring the system's performance for the AQUIOS STEM application. Refer to below for instructions on how to run the Daily QC procedures for the

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Borna
APODERADO

AQUIOS STEM application. Refer to Daily QC Procedure in the AQUIOS CL Flow Cytometer Instructions for Use manual to QC the AQUIOS CL instrument itself.

IMPORTANT When running QC, if a failure occurs, the system will issue a flag. Check the individual runs for STEM to determine which QC run failed and repeat that QC run.

QC Setup

CAUTION

Risk of erroneous results. When using AQUIOS STEM, each control cell vial is limited to 15 piercings. Cap-piercing beyond the allotted number is not recommended.

Before running the AQUIOS STEM CD34 Control Cells, verify that:

1. The Daily Startup was performed. Refer to Startup in Chapter 3, Daily Startup in AQUIOS CL Flow Cytometer Instructions for Use manual for details.
2. The Beckman Coulter reagents and the reagent levels are suitable and sufficient for the QC procedure. To verify reagent availability, select  on the Status Bar. To replace a reagent, refer to Replacing Reagents and Consumables in Chapter 11, Replacement Procedures in the AQUIOS CL Flow Cytometer Instructions for Use manual.
3. The control assay values have been previously registered into the system. Refer to Reagent Info Screen in Chapter 8, Setup in the AQUIOS CL Flow Cytometer Instructions for Use manual for more information. If the control lot number information has not been registered in the system, scan the control assay sheet before adding a QC request. Refer to Starting a New Lot of Controls in Chapter 11, Replacement Procedures in the AQUIOS CL Flow Cytometer Instructions for Use manual.
4. If a scan error occurs, readjust the barcode and re-scan.

QC Process

A plastic polypropylene 13x75 mm test tube is required to run the Flow-Check Fluorospheres. Beckman Coulter does not recommend or guarantee the use of a specific tube to produce quality results. If you need information on a tube, [contact us](#).

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Flow-Check Fluorospheres - Single Tube Loader

Risk of spill. Do not use the autoloader position to run the tube containing the Flow-Check Fluorospheres.

Risk of needle crash. Do not use a non-pierceable cap to close the tube.

Follow the steps below to Add a new lot of Flow-Check fluorospheres, otherwise skip to step 5.

1 Select **Setup** on the workstation main screen and then select **Reagent Info**.


2 Select the following parameters on the screen:

- **Test:** FlowCheck_Alignment
- **Reagent/Sample:** Flow-Check

NOTE If the desired Lot already exists in the system, skip to step 5.

3 Select . On the **Add Reagent Lot** popup screen, input the lot number and expiration date.

IMPORTANT Verify input before selecting **OK**. Once saved in the system, The Lot number and Expiration Date cannot be edited. A table with the ranges for each field is displayed. You should accept the default ranges (but you have the option to modify if you want).

4 Select **OK**.

5 Select  to open the Single-tube door.


6 The message *Please press on the single tube door, and then wait for the system to release it is displayed.*

7 Open the door by pressing the top-middle part of the door in one single, solid motion.

NOTE If there is a tube present in the single tube loader, the message *Please remove any tubes from the Single-tube loader* is displayed. Remove the tube and select **OK**.

8 Gently invert the Flow-Check Fluorospheres bottle repeatedly to completely suspend (no visible pellet) fluorospheres and add >0.5 mL of undiluted Flow-Check Fluorospheres into a new plastic polypropylene 13x75 mm tube.

9 Place the 13x75 mm tube (without cap) in Position 1 of the single-tube loader.

10 On the **Sample Information Screen**, select **QC** radio button.

11 Enter **Flow Check** under **Sample ID** and select the following parameters,

- **Sample ID:** Flow Check
- **Test:** FlowCheck_Alignment
- **Type:** QC
- **Control:** Flow Check
- **Lot:** Select lot number

Sample Information + Request will be added

☐ Patient ☒ QC

Sample ID: Flow Check

Test: FlowCheck_Alignment

Type: QC

Control: Flow Check

Lot: 9435208

Lookup

12 Close the single-tube loader door. The request is automatically added.

13 Follow the sample preparation progression on the **Running** tab of the workstation main screen.

14 When the Flow-Check run is complete, view the results in the **Review** tab.

Select **stat** to open the Single-tube door to retrieve the control tubes.

15 Prior to running AQUIOS STEM CD34 Control Cells ensure Flow-Check HPCV and count meet specifications.

STEM Single (AutoLoader) - AQUIOS STEM CD34 L1 and L2 Control Cells

CAUTION

Risk of erroneous results. There are 15 tests per level of AQUIOS STEM CD34 Control Cells. Each cap should have no more than 15 piercings.

CAUTION

Make sure to correctly resuspend the cellular suspension by performing several manual mixing steps per inversion (6-7 times). Visually verify that there is no visible pellet of cells remaining in the tube.

NOTE The Control Cells qualify the system for use with all specimen types as outlined in [Table 1.1](#).

IMPORTANT Use the Type 1 cassette (13x75mm) for AQUIOS STEM CD34 Control Cells Level 1 and Level 2.

IMPORTANT QC test failure or failure to run QC tests (Flow-Check Fluorospheres and AQUIOS STEM CD34 L1 and L2 Control Cells) will result in a *Missing QC* flag being triggered for the patient test.

AQUIOS STEM CD34 Control Cells Level 1 and Level 2 are in vitro diagnostic medical reagents intended to be used by laboratory professionals as process controls for the verification of the parameters CD34 and CD45 as part of the AQUIOS STEM System on the automated AQUIOS CL Flow Cytometry System. Refer to [Table 1.1, System Parameters](#) for a list of supported specimen types.

Adding Reagents

Refer to [Table 3.1, AQUIOS STEM QC Workflow](#).

CAUTION

Risk of losing on-board samples. Replacing reagents and opening the reagent door while operating the Cytometer puts on-board samples at risk. Replacing reagents while the system is running should be performed as quickly as possible to avoid flagged results due to over lysing and over incubation.

1 Select Open Reagent Door. The open reagent door prompt will appear.

2 Press the reagent door to open it.

Table 3.2 Handling AQUIOS STEM Reagent Caps

Reagent	Cap Action
AQUIOS STEM Lysing Solution	Remove cap before loading into system.
AQUIOS STEM-Count Fluorospheres	Replace regular cap with screw cap with hole.
Monoclonal Antibodies & 7AAD Viability Dye Reagents	Replace shipping caps with pierceable caps.

 CAUTION

Risk of instrument damage. For the monoclonal antibodies and the 7AAD viability dye reagents, ensure the shipping caps have been replaced with the pierceable caps provided with the reagent packaging. For the AQUIOS STEM Lysing Solution, remove the cap before loading into the system. For the AQUIOS STEM-Count, the regular cap and seal foil have to be replaced with a screw cap with hole before loading into the system. (see further instructions below and in reagent IFU). Otherwise, the probe can be damaged.

IMPORTANT AQUIOS STEM-Count Fluorospheres should be vortexed for 10 to 12 seconds before placing in the carousel.

IMPORTANT When prompted to re-mix AQUIOS STEM-Count Fluorospheres, remove the screw cap with hole from the bottle and replace it with the standard cap before vortex mixing. After mixing, replace the standard cap with the screw cap with hole before returning to the instrument.

At the end of the day, remove the AQUIOS STEM-Count Fluorospheres bottle from the system and replace the screw cap with hole with the standard cap. Dry the inner surface of the screw cap with hole with low lint tissue paper before storing in the kit box. Store the kit box at 2-8°C until the next use.

IMPORTANT An additional label on the side of the bottle allows the association between AQUIOS STEM Lysing Solution and the kit lot number.

- 3** Add the required reagents to the reagent carousel, refer to [Table 3.1, AQUIOS STEM QC Workflow](#), and insert an AQUIOS 96-deep well plate in the correct orientation. Make sure reagent levels and number of unused wells are sufficient for the number of runs requested. If needed, add or replace with new reagents as appropriate before proceeding to analysis.

Refer to Heading , [AQUIOS STEM Replacement Procedures](#) in [CHAPTER 8, Adjustment and Replacement Procedures](#).

- a.** Place the vials of CD45-FITC /CD34-PE and 7AAD viability dye reagents in the carousel. Mount the carousel containing the vials of monoclonal antibody reagents on the carousel platform inside the reagent compartment.


NOTE When loading the carousel into the instrument, be sure to place it on the wheel and gently turn it until it locks into position on the wheel.

- b.** Remove the cap and place the AQUIOS STEM Lysing Solution bottle in any of the three available positions in the reagent compartment. Make sure the bar code label face outward.
- c.** Before testing or when a reminder pops up to remind mixing, with the original white cap on, mix the bottle of AQUIOS STEM-Count Fluorospheres thoroughly. Remove the original white cap and the aluminum heat seal (when applicable), and replace the original white cap with the screw cap with hole that is provided with the kit. Place the AQUIOS STEM-Count Fluorospheres in the bottle holder with the adapter. Make sure the barcode labels face outward.

NOTE Refer to the appropriate reagent IFU for details on how to properly mix reagents.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093



Rolando F. Bornia
APODERADO

- 4 Close the reagent door. The instrument scans the reagent barcodes.

A window appears to ask if Flow-Count Fluorospheres have been mixed. Select **Yes** if Flow-count Fluorospheres were mixed. Select **No** if Flow-count Fluorospheres were not mixed and open the reagent door to mix the beads.

NOTE For the AQUIOS STEM application you are using AQUIOS STEM-Count Fluorospheres, not Flow-Count Fluorospheres even though the message says Flow-Count Fluorospheres.



-
- 5 Use a pre-existing lot or add a new lot of AQUIOS STEM CD34 L1 and L2 Control cells and add a request.
- Refer to [Starting AQUIOS STEM CD34 L1 and L2 Control Cell QC With Pre-Existing Lot](#) to use a pre-existing lot.
 - Refer to [Adding New Lots of AQUIOS STEM CD34 L1 and L2 Control Cells using QR Code](#) to automatically add a lot with a barcode.
 - Refer to [Adding New Lots of AQUIOS STEM CD34 L1 and L2 Control Cells Manually](#) to manually add a lot.
-
- 6 If you select Pause for QC, once QC samples have completed analysis, the autoloader remains paused if QC fails so that the problem can be corrected before running subsequent samples. If QC passes, the autoloader continues processing additional cassettes in the autoloader following the list of samples to run in the Request folder.
-
- 7 Gently roll/invert the control reagent to mix the product in a manner similar to the handling of a whole blood specimen (mixing with vortex or a rotator is not recommended).
-

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

- 8 Insert control tubes into cassette #1 in the proper configuration with the barcode facing outward.

NOTE Refer to AQUIOS CL Flow Cytometer IFU, Chapter 5, for details on how to properly insert tubes into the cassette.

- 9 Place the cassette onto the cassette loading area (above autoloader).
- 10 Follow the sample preparation progression on the **Running** tab of the workstation main screen.

- 11 After the control tubes have ejected from the cytometer, the cassette will be moved to the unloading area below the autoloader. Remove and store the control tubes properly until next use.


Adding and Using AQUIOS STEM CD34 Control Lots


Adding New Lots of AQUIOS STEM CD34 L1 and L2 Control Cells using QR Code

- 1 Select  on the right side of the workstation main screen and the *Barcode Reagents* screen will appear.
- 2 After the system scanner light turns on, scan the QR code for the L1 and L2 control cells on the corresponding assay sheets.
- 3 Select **OK** on the *Barcode Reagents* screen.

Adding New Lots of AQUIOS STEM CD34 L1 and L2 Control Cells Manually

- 1 Select **Setup** on the workstation main screen.
- 2 Select **Reagent Info** tab on the left side of the main screen.
- 3 Select the following parameters:
 - **Test:** STEM Single.
 - **Reagent/Sample:** STEM CD34 – Level 1 or STEM CD34 – Level 2
- 4 Select . On **Add Reagent Lot** popup screen, input the lot number and expiration date.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

IMPORTANT Verify input before selecting **OK**. Once saved in the system, it cannot be edited.

Add Reagent Lot

Lot Number: Expiration Date:
(mm/dd/yyyy)

5 Select **OK**.

6 Select **Edit**. The *In Edit Mode* screen is displayed.

Test:

Reagent/Sample:

Lot Number: 

Expiration: CAL Factor:
(mm/dd/yyyy)

Note: Percents must be entered as decimal. Example 20% = .2

Result	Min	Max
CD34+ Leukocytes Percent		
CD34+ Count/uL		
CD45+ Count/uL		

7 Using the assay sheet, enter assay values for the following **MIN** and **MAX** parameters as shown in the example of assay sheet shown below:

- Percent Positive CD34+

NOTE For Percent Positive, ensure that the values within the parentheses from the assay sheet are used. For example, ("0.0020").

- Absolute Number

NOTE From CD34+ and CD45+.

IVD

AQUIOS STEM CD34 Control Cells

Table of Expected Results

REF C43667

LOT	11300063
	2021-08-07

Level 1

	Percent Positive Cell ^a	Min ^b	Max ^b	Absolute Number (cells / μ L)	Min	Max
CD34+	0.21 (.0021)	0.14 (.0014)	0.28 (.0028)	12.9	7.7	18.1
CD45+				6250	5250	7250

^a - based on leukocytes


^b - in case of manual entry, values must be entered as shown in parentheses

IMPORTANT: The Expiration date is displayed as YYYY-MM-DD. If entering the Lot and Expiration information manually, the expiration date must be entered in the software in the following format: MM/DD/YYYY.

8 Select **Save**.

Starting AQUIOS STEM CD34 L1 and L2 Control Cell QC With Pre-Existing Lot

- 1 Select  on the workstation main screen and then select  .
The *Add Request* screen appears.

- 2 Select **QC**.

- 3 Select the test **STEM Single** .

- 4 Select the  icon on the **Add Request** screen. After the system scanner light turns on scan the barcode on the control cell vial.

NOTE Screen will automatically switch from **Patient** to **QC** if the lot number already exists in the system.

- 5 **Sample ID** field will automatically populate with the control lot number on the vial when scanned.

NOTE Refer to AQUIOS CL Instructions for Use for more information on barcodes.

- 6 Verify and confirm the following parameters:

- **Test:** .
- **Type:** QC
- **Control:** STEM CD34 – Level 1 or STEM CD34 – Level 2
- **Lot:** Confirm the lot number

- 7 If you select Pause for QC, once QC samples have completed analysis, the autoloader remains paused if QC fails so that the problem can be corrected before running subsequent samples. If QC passes, the autoloader continues processing additional cassettes in the autoloader following the list of samples to run in the Request folder.

QC Screen

The Results tab on the QC screen tracks the changes over time in HPCV (Half Peak coefficient of variation) for FS, FL1, FL2, FL3, FL4 and FL5 parameters using Flow-Check Fluorospheres. This data reflects performance and stability of the instrument over the specified period.

NOTE Under **Test**, switch between **FlowCheck_Alignment** for Flow-Check results and **Stem Single** for Control Cells results.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Figure 3.1 Results Tab Screen: Flow-Check Fluorospheres


Table 3.3 General Table on QC Screens: Defined

1	Exclude Checkbox: Use the checkbox to eliminate this run from the QC range calculation and select  .	7	Min: The lowest value of all the non-excluded runs listed.
2	Lot: Control reagent lot number used for the specified QC run.	8	Max: The highest value of all the non-excluded runs listed.
3	Date run: The date on which the selected QC run was performed.	9	CV: The coefficient of variation of all the non-excluded runs listed (Standard Deviation/Average) x 100.
4	Results: The QC results quantitatively for the selected QC run. The results displayed depend on the selected test.	10	StDev: The standard deviation of all the non-excluded runs listed.
5	Lot: The lot numbers of the control reagent for all of the runs listed.	11	Period: Select the number of months to display in each graph.
6	Average: The average of all the non-excluded runs listed.	12	elQAP: Refer to Setting Up IQAP Export .


NOTE Summary statistics are listed at the bottom of the graphs. Use the vertical scroll bar to access all statistics. The customer can move the region.

QC Reports

AQUIOS STEM (Single Tube Loader) Flow-Check Fluorospheres QC Result Examples

Figure 3.2 Results Tab: Flow-Check Fluorospheres Report Page 1

Sample ID: 89352285090 Test: FlowCheck_Alignment v1.0
Run Date: 11 Jun 2021 08:10:01 Analysis Date: 11 Jun 2021 08:10:01
User: Admin Collect Date:
Specimen Type: Control Patient Name: Flow Check
Patient ID: Lot: 9435283 Gov. ID:
Status: Location:
Gender: Date Of Birth:
Physician: Physician Code:
Instrument Serial #: BA44086
Run Flags:
Run Notifications:
Comments:

Result	Value	Flag	Normal Range	Action Range
FS Count	5,000			4,000-6,000
FS Mean	399.25			
FS HPCV	0.78%			0.01%-4.00%
FL1 Mean	419.99			
FL1 HPCV	1.71%			0.01%-4.00%
FL2 Mean	412.61			
FL2 HPCV	1.56%			0.01%-4.00%
FL3 Mean	386.97			
FL3 HPCV	1.64%			0.01%-4.00%
FL4 Mean	395.90			
FL4 HPCV	1.81%			0.01%-4.00%
FL5 Mean	450.41			
FL5 HPCV	2.09%			0.01%-7.00%

Signature:

NOTE The region can be manually adjusted in order to capture all the fluorospheres.

AQUIOS STEM Single (Auto Loader) - AQUIOS STEM CD34 Control Cells Level 1 & 2 QC Result Examples

Figure 3.3 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 1 (Run with AQUIOS STEM Single)


Figure 3.4 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 2 (Run with AQUIOS STEM Single)


Figure 3.6 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 1, Page 2 (Run with AQUIOS STEM Single)

Sample ID: 10180063 Test: Stem Single v1.0
Run Date: 18 Mar 2021 09:16:45 Analysis Date: 18 Mar 2021 09:16:45


Result	Value	Flag	Normal Range	Action Range
Viable CD34+ HPC Count (cells/uL)	12.12			7.20-16.80
Viable Leukocytes Count (cells/uL)	6,526.36			5,451.00-7,527.00
Viable CD34+ HPC % (of Viable Leukocytes)	0.19%			0.11%-0.25%
DF: Dilution Factor	0.00			
Viable CD34+ HPC Count (cells/uL) (x DF)	0.00			
Viable Leukocyte Count (cells/uL) (x DF)	0.00			
HV: Harvest Volume (mL)	0.00			
Viable CD34+ HPC Count (x DF x HV)	0.00			
BW: Body Weight (kg)	0.00			
Viable CD34+ HPC Count per kg (x DF x HV / BW)				
CAL	989			

Statistic	Value	Flag	Normal Range	Action Range
**CD34+ Viability % (of total CD34+)	100.00%			
**Total CD34+ HPC Count (cells/uL)	12.12			
**WBC Count (cells/uL)	6,526.36			
**Sample Viability % (of WBC)	100.00%			
**DF: Dilution Factor	0.00			
**Total CD34+ HPC Count (cells/uL) (x DF)	0.00			
**WBC Count (cells/uL) (x DF)	0.00			
**Viable Events	76,508			
**Viable Leukocytes Events	76,495			
**Viable CD34+ HPC Events	142			
**Single Beads Events	11,592			

Signature:


Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Figure 3.7 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 2, Page 1 (Run with AQUIOS STEM Single)

Sample ID:	10180064	Test:	Stem Single v1.0
Run Date:	18 Mar 2021 09:23:47	Analysis Date:	18 Mar 2021 09:23:47
User:	Admin	Collect Date:	
Specimen Type:	Blood Control	Patient Name:	STEM CD34 - Level 2
Patient ID:	Lot: 10180064	Gov. ID:	
Status:		Location:	
Gender:		Date Of Birth:	
Physician:		Physician Code:	
Instrument Serial #:	BA44086	Harvest Volume (mL):	0
Report:	Well 1	Weight (kg):	0
		Dilution:	0
Run Flags:			
Run Notifications:			
Comments:			

** Population Statistics

Page #1 of 2

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO


Figure 3.8 Results Tab Screen: AQUIOS STEM CD34 Control Cells Level 2, Page 2 (Run with AQUIOS STEM Single)

Sample ID: 10180064 Test: Stem Single v1.0
Run Date: 18 Mar 2021 09:23:47 Analysis Date: 18 Mar 2021 09:23:47

Result	Value	Flag	Normal Range	Action Range
Viable CD34+ HPC Count (cells/uL)	37.63			23.90-43.50
Viable Leukocytes Count (cells/uL)	6,585.61			4,978.00-7,466.00
Viable CD34+ HPC % (of Viable Leukocytes)	0.57%			0.42%-0.66%
DF: Dilution Factor	0.00			
Viable CD34+ HPC Count (cells/uL) (x DF)	0.00			
Viable Leukocyte Count (cells/uL) (x DF)	0.00			
HV: Harvest Volume (mL)	0.00			
Viable CD34+ HPC Count (x DF x HV)	0.00			
BW: Body Weight (kg)	0.00			
Viable CD34+ HPC Count per kg (x DF x HV / BW)				
CAL	989			

Statistic	Value	Flag	Normal Range	Action Range
**CD34+ Viability % (of total CD34+)	100.00%			
**Total CD34+ HPC Count (cells/uL)	37.63			
**WBC Count (cells/uL)	6,585.61			
**Sample Viability % (of WBC)	100.00%			
**DF: Dilution Factor	0.00			
**Total CD34+ HPC Count (cells/uL) (x DF)	0.00			
**WBC Count (cells/uL) (x DF)	0.00			
**Viable Events	76,492			
**Viable Leukocytes Events	76,477			
**Viable CD34+ HPC Events	437			
**Single Beads Events	11,485			

Signature:

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Borna
APODERADO

Interlaboratory Quality Assurance (IQAP)

The Interlaboratory Quality Assurance Program (IQAP) is a Beckman Coulter program available to you through enrollment that complements and enhances your laboratory's in-house quality control. IQAP lets you submit your control recovery data to Beckman Coulter and in return, receive a personalized report that summarizes your results and compares them to those of your peer group (pool).

Before Setting Up IQAP Export


Before you export an IQAP file, you must register your AQUIOS instrument with IQAP. The IQAP group will assign an IQAP ID number for your AQUIOS.

- 1 Using an internet connected computer, navigate to website,
<https://qspweb.beckmancoulter.com>
and follow the instructions on the website to create an account and register your instrument.

- 2 On the AQUIOS main screen, select .


- 3 On the Setup screen, select .

- 4 On the System Setup screen, enter your IQAP ID number and select .

The screenshot shows a section titled "IQAP Settings" with a label "IQAP ID:" followed by a text input field. A red arrow points to the input field.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Setting Up IQAP Export

1 On the QC Results Tab Screen select .


NOTE Use the eIQAP button to export IQAP control data files.

Figure 3.9 QC Results Tab Screen: Export IQAP

2 Choose a location to export (save) the IQAP control data file.

Figure 3.10 Export IQAP control data file

Rolando F. Bornia
APODERADO


Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

- 3 Make a note of the location of the saved IQAP control data file, you will need to know the file location when you upload it to the IQAP website.

Refer to the Cytometry Interlaboratory Quality Assurance Program Procedure Manual for instructions for uploading your control data files. See [Related Documents](#).

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

CHAPTER 4

Sample Analysis

Overview

IMPORTANT Before performing an analysis using the AQUIOS STEM application, ensure you are thoroughly familiar with the operating procedures in the AQUIOS CL Flow Cytometer Instructions for Use manual.

This chapter contains the following sections:

- Test Description
- AQUIOS STEM Application Workflow
- Sample preparation
- Storage Conditions and Stability
- AQUIOS STEM Testing Workflow
- Create a Request
- Load Sample in the Instrument
- Options for Running Specimens
- Logging In Controls
- Running Samples Using the Autoloader
- Running Samples Using the Single-tube Loader

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Test Description

Three AQUIOS STEM tests are available in the AQUIOS STEM application:

- Stem Panel: Run test in duplicate plus negative control (3 wells)
- Stem Duplicate: Run test in duplicate, no negative control (2 wells)
- Stem Single: Run test in single, no negative control (1 well)

Refer to [Table 4.1, Test Descriptions](#) for the list of AQUIOS STEM kit reagents necessary for each test.

Table 4.1 Test Descriptions

Test Name	Reagent					Test Description
	CD45-FITC / CD34-PE	7-AAD Viability Dye	CD45-FITC / Isoclonic Control	Lysing Solution	STEM-Count Fluorospheres	
Stem Panel	X	X	X	X	X	Run tests in duplicate plus negative control (3 wells)
Stem Duplicate	X	X		X	X	Run tests in duplicate , no negative control (2 wells)
Stem Single	X	X		X	X	Single test (1 well)

AQUIOS STEM Application Workflow

When running samples using the AQUIOS STEM application with Single, Duplicate, and Panel tests, each test goes through four stages:

- Preparing** The specimen, monoclonal antibodies, and 7AAD viable dye are aspirated, dispensed and then mixed in the designated well of the 96-deep well plate. Refer to [Table 4.2](#).
- Incubating** An approximate 15-minute incubation allows proper staining of the specimen by the monoclonal antibodies. Refer to [Table 4.2](#).
- Lysing** The lysing system is dispensed into the well containing the stained specimen to lyse the RBCs and preserve the WBCs. Refer to [Table 4.2](#).
- Analyzing** The AQUIOS STEM-Count Fluorospheres are added to the prepared sample into the well and the mixture is passed through the flow cell to obtain the sample results. Refer to [Table 4.2](#).

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Table 4.2 Processing Samples Using the STEM Application

PREPARING	INCUBATING	LYSING	ANALYZING
1. User loads specimen^a. 2. Specimen is mixed and aspirated from the collection tube. 3. Specimen is dispensed into a microtiter plate well. 4. Antibodies are added to the well. 5. Suspension is mixed.	1. About 15 minutes.	1. AQUIOS STEM Lysing Solution is added to the well. 2. Solution is mixed. 3. Lysing, including incubation, for about 15 minutes.	1. AQUIOS STEM Kit STEM-Count Fluorospheres are added into the well. 2. Solution is mixed. 3. Sample is aspirated. 4. Sample is analyzed by flow cytometry up to 5 minutes.


a. When using the autoloader, specimens are automatically mixed. When using the Single-tube Loader, mix the collection tube as instructed before loading the specimen tube. (Refer to Chapter 5 in the AQUIOS CL Flow Cytometer Instructions for Use manual for details.)


Sample preparation

- Retain fresh specimens in blood collection tubes at recommended temperature depending on the specimen type, between 18°C and 25°C (64.4°F and 77°F) or between 2°C and 8°C (35.6°F and 46.4°F).
- Fresh samples are ideally analyzed as soon as possible after drawn, refer to [Storage Conditions and Stability](#) for storage conditions for each sample type.
- Make sure each patient specimen tube to be loaded complies with the minimum running volume. When using the autoloader, ensure that the specimen tube has a barcode on it before loading it to the system. Refer to [Table 4.3](#). Refer to AQUIOS CL Instructions for Use for complete list of specimen tubes for the Autoloader and Single-tube loader.

Table 4.3 Minimum Specimen Volume Required per Specimen Tube

Specimen Tube Size	Minimum Volume
13 mm x 75 mm	750 µL
13 mm x 100 mm	750 µL
16 mm x 100 mm	1000 µL
10.25 mm x 50 mm	650 µL
Sarstedt 13 mm x 75 mm	750 µL
Sarstedt 13 mm x 90 mm	750 µL
Sarstedt 11 mm x 66 mm	750 µL
Sarstedt 15 mm x 92 mm	1500 µL
Sarstedt 13 mm x 65 mm	750 µL

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Storage Conditions and Stability

Non-Frozen Specimen

Fresh samples are ideally analyzed as soon as possible after drawn, but no later than described below.

- Normal whole blood samples are stable for up to 20 hours if stored between 18-25°C (64-77°F).
- Mobilized whole blood samples are stable for up to 20 hours if stored between 18-25°C (64-77°F).
- Cord blood samples are stable for up to 24 hours if stored between 18-25°C (64-77°F).
- Apheresis samples are stable for up to 24 hours if stored between 2-8°C (35-46°F).
- Bone marrow samples are stable for up to 24 hours if stored between 2-8°C (35-46°F).

For best staining results, the following rules should be applied:

- Viable Leukocytes below 30000 Cells/ μ L
- Viable CD34+ HPC cells below 1000 Cells/ μ L

Frozen Specimen Type

Frozen samples should be processed within 15 minutes after thawing and kept on ice.

Beckman Coulter recommends using the Single Tube Loader function to analyze thawed samples.

Frozen samples should not be rerun.

In case of rerun, a new aliquot should be thawed.

Reagents

Refer to the appropriate package inserts for information on reagent stability.

- AQUIOS Cleaning Agent
- AQUIOS Sodium Hypochlorite Solution
- AQUIOS Sheath Solution
- AQUIOS STEM Kit Reagents
- AQUIOS STEM CD34 Control Cells
- Flow-Check Fluorospheres

AQUIOS STEM Testing Workflow

1 Perform AQUIOS Startup.

Refer to the AQUIOS CL Flow Cytometer Instructions for Use manual, as needed. Make sure AQUIOS CL Startup is completed and acceptable. Refer to Chapter 3, Daily Startup in the AQUIOS CL Flow Cytometer Instructions for Use manual.

- 2 Load AQUIOS STEM Kit reagents. The application specific reagents are automatically scanned by the system once they are loaded in their respective holders.

CAUTION

Risk of losing on-board samples. Replacing reagents and opening the reagent door while operating the Cytometer puts on-board samples at risk. Replacing reagents while the system is running should be performed as quickly as possible to avoid flagged results due to over lysing and over incubation.

CAUTION

Risk of instrument damage. For the monoclonal antibodies and the 7AAD viability dye reagents, ensure the shipping caps have been replaced with the pierceable caps provided with the reagent packaging. For the AQUIOS STEM Lysing Solution, remove the cap before loading into the system. For the AQUIOS STEM-Count, the regular cap and seal foil have to be replaced with a screw cap with hole before loading into the system. (see further instructions below and in reagent IFU). Otherwise, the probe can be damaged.

IMPORTANT AQUIOS STEM-Count Fluorospheres should be vortexed for 10 to 12 seconds before placing in the carousel.

IMPORTANT When prompted to re-mix AQUIOS STEM-Count Fluorospheres, remove the screw cap with hole from the bottle and replace it with the standard cap before vortex mixing. After mixing, replace the standard cap with the screw cap with hole before returning to the instrument.

Make sure AQUIOS STEM kit reagents were loaded and scanned correctly.

- Place the vials of CD45-FITC /CD34-PE and 7AAD viability dye reagents in the carousel.

NOTE When using STEM Panel test, you must also place a vial of CD45-FITC/ Isoclonic Control in the carousel. Refer to [Table 4.1, Test Descriptions](#).

- Mount the carousel containing the vials of monoclonal antibody reagents on the carousel platform inside the reagent compartment

NOTE When loading the carousel into the instrument, be sure to place it on the wheel and gently turn it until it locks into position on the wheel.

- Remove the cap and place the AQUIOS STEM Lysing Solution bottle in any of the three available positions in the reagent compartment. Ensure the bar code label faces outward.
- Before testing or when a reminder pops up to remind mixing, with the original white cap on, mix the bottle of AQUIOS STEM-Count Fluorospheres thoroughly. Remove the original white cap and the aluminum heat seal (when applicable), and replace the original white cap with the screw cap with hole that is provided with the kit. Place the AQUIOS STEM-Count Fluorospheres in the bottle holder with the adapter. Make sure the bar code label faces outward.

NOTE Refer to the appropriate reagent IFU for details on how to properly mix reagents.


- Insert an AQUIOS 96-deep well plate in the correct orientation. Make sure reagent levels and number of unused wells are sufficient for the number of runs requested. If needed, add or replace with new reagents as appropriate before proceeding to analysis.

At the end of the test day, remove the STEM-Count bottle from the system and replace the screw cap with hole with the standard cap. Dry the inner surface of the screw cap with hole with low lint tissue paper before storing in the kit box. Store the kit box at 2-8°C until the next use. Store the AQUIOS STEM Lysing Solution bottle at room temperature.

- 3 Ensure Quality Control checks have been completed and are acceptable within 24 hours using AQUIOS STEM CD34 Control Cells Level 1 and level 2. Refer to Chapter 4, Quality Control in the AQUIOS CL Flow Cytometer Instructions for Use manual as needed

Create a Request

- 1 Select  on the workstation main screen and then select . The *Add Request* screen appears.

- 2 Select **Patient**.

- 3 Select the  icon on the **Add Request** screen. After the system scanner light turns on scan the barcode on the control cell vial.

NOTE It is possible to manually populate the Sample ID field.

NOTE Refer to AQUIOS CL Instructions for Use for more information on barcodes.

4 Select an AQUIOS STEM test.

Three AQUIOS STEM tests are available: Stem Single, Stem Duplicate or Stem Panel.

Add Request

☒ Patient ☐ QC

Sample ID:

Test: **Select Test**

Specimen Type: **Select Test**

Patient Name:

Patient ID:

Status:

Location:

Gender: **Select** Date of Birth:

Physician Name:

Physician Code:

OK Cancel

5 Select the Specimen Type.

Five specimen types are available:

- Apheresis
- Bone Marrow
- Cord Blood
- Mobilized Blood
- Whole Blood

Add Request

☒ Patient ☐ QC

Sample ID:

Test: **Stem Panel**

Specimen Type: **Select Type**

Patient Name:

Patient ID:

Status:

Location:

Gender: **Select** Date of Birth:

Physician Name:

Physician Code:

Harvest Vol. (mL): Weight (kg): Dilution:

OK Cancel

-
- 6 Enter Harvest Volume, Patient Weight and Dilution Factor information.

 CAUTION

Risk of erroneous results. If no values are entered for these parameters, the instrument will default them to 0. This will cause some calculations in the results page to equal 0. Input "1" for Harvest Volume, Weight, and Dilution if these fields would have otherwise been left blank.

The harvest volume (**Harvest Vol. (mL)**) of the sample, the patient weight (**Weight (kg)**) and the dilution factor of the tested sample (**Dilution**) are available as shown above in the **Add Request** screen. The values can be integer or with decimals, the maximum number of digits is 5.

- **Harvest Vol. (mL)** should be expressed in milliliter unit.
- **Weight (Kg)** should be expressed in Kilogram unit.
- If the **Dilution** of the sample is 1:10, the entered value will be 10 (Note: 1:10 equals 1 part of sample plus 9 parts of diluent).
- Enter additional information if necessary, in the Add Request screen: Collect Date, Patient Name, Status, etc.


-
- 7 When finished entering information on the **Add Request** screen, select **OK**.
-

Load Sample in the Instrument

-
- 1 Gently roll/invert the sample tube to mix the product (mixing with vortex or a rotator is not recommended).



Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Borna
APODERADO

- 2 Insert sample tubes into a cassette that matches the tube size in the proper orientation with the barcode facing outward.

NOTE Refer to AQUIOS CL Flow Cytometer IFU, Chapter 5, for details on how to select the appropriate cassette and properly insert tubes into the cassette.

- 3 Place the cassette onto the cassette loading area (above autoloader).
- 4 Follow the sample preparation progression on the **Running** tab of the workstation main screen.

- 5 After the sample tubes have ejected from the cytometer, the cassette will be moved to the unloading area below the autoloader. Remove the sample tube and store properly if it might be used again. Refer to [Storage Conditions and Stability](#)

Options for Running Specimens

WARNING

Risk of personal injury or contamination:

- If not properly shielded while using the instrument, it is possible to become injured or contaminated. To prevent possible injury or biological contamination, you must wear proper laboratory attire, including gloves, a laboratory coat, and eye protection.
- Operating the instrument while the waste level sensor is disconnected can cause biological contamination. To prevent biological contamination, do not operate the instrument while the waste level sensor is disconnected.

CAUTION


Risk of clogging or erroneous results if samples contain clots, fibrin strands, or cellular debris, or if samples are lipemic. Use good laboratory practices and follow your laboratory procedures for specimen handling and/or inspection prior to placing on the system.

Specimen tubes can be loaded in the autoloader or the Single-tube Loader. The Single-tube Loader functions as a stat loader. A higher priority is given to a specimen tube loaded into the Single-tube Loader; therefore, it will be processed ahead of the specimen tubes loaded in an autoloader cassette.

For the AQUIOS STEM application, the autoloader has continuous loading capability with up to 15 runs (15 specimens in Single, 7 in Duplicate, and 5 in Panel) in order to avoid overincubation with the lysing reagent. The specimens are automatically processed according to the unprocessed stack cassette position (bottom to top) at the top of the autoloader holder and according to the sample position within a cassette (position 1 to 5). Place the cassette that you want to run first at the bottom of the unprocessed stack and place the specimen you want to run first in position one (1) of the cassette. Refer to Running Samples Using the Autoloader in Chapter 5, Sample Processing in the AQUIOS CL Flow Cytometer Instructions for Use manual.

The AQUIOS Tetra tests can be simultaneously loaded on the AQUIOS CL Flow Cytometry system with AQUIOS STEM Tests. The AQUIOS Smart Scheduler tracks requests for both tests as they proceed from sample preparation to analysis. Ensure that the appropriate QC has been run for the

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

AQUIOS Tetra tests prior to loading patient samples. Ensure that the appropriate and sufficient AQUIOS Tetra reagents are loaded. Refer to AQUIOS Tetra System Guide for more information.

NOTE Wait until the AQUIOS instrument displays green light on autoloader and status bar to process additional samples.

Logging In Controls


For detailed information on Logging In Controls, refer to Chapter 11, Replacement Procedures in the AQUIOS CL Flow Cytometer Instructions for Use manual.

Running Samples Using the Autoloader


For detailed information on Running Samples Using the Autoloader, refer to Chapter 5, Sample Processing, in the AQUIOS CL Flow Cytometer Instructions for Use manual.

Running Samples Using the Single-tube Loader

For detailed information on Running Samples Using the Single-tube Loader, refer to Chapter 5, Sample Processing in the AQUIOS CL Flow Cytometer Instructions for Use manual.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Running Screen

The Running screen (Figure 4.1) tracks any samples currently being processed by the Cytometer. The sample workflow is displayed on the screen (**Running**) as the sample progresses from one stage to the next (from Preparing to Analyzing). The **Details** portion of the screen provides visibility and the specifics of each stage, including timing.

Figure 4.1 Running Screen

For more information about the Running Screen, refer to Appendix A in the AQUIOS CL Flow Cytometer Instructions for Use manual.

Gabriela A. Cividino
Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
Rolando F. Bornia
APODERADO

Overview

This chapter contains information about:

- Review Screen
- Data Review
- Adjust Compensation Manually
- Result Screen
- Flags And Notifications

The AQUIOS Workstation provides the user with the ability to review both QC data and completed patient tests.

The review process provides the user with the ability to review data based on the Review/Editor and Administrator rights function to edit completed runs. The following are the items which can be edited:

- Sample information
- Gate and Region Placement
- Compensation

Review Screen

IMPORTANT There is a risk of reporting incorrect results. Data displays used to arrive at the test result for light scatter patterns, antibody staining profiles, and all gates and boundaries should be reviewed by a laboratory professional when interpreting the data. If results are suspect, follow your laboratory's procedures to resolve.

For a more detailed explanation of this screen and its components, refer to Chapter 6, Data Review in the AQUIOS CL Flow Cytometer Instructions for Use manual.


The Review screen displays analyzed samples that were previously run and that are currently pending review and/or transmission. The Review screen can be customized using a drop down menu with results from each well (Well#1, Well#2 or Control) or with a pre-defined graph combinations of the 3 wells (Panel test) or 2 wells (Duplicate test) called Summary. The access to the different well is dependent of the Test option selected. These screens allow you to review, edit, and save changes to regions, gates, compensation and/or patient information prior to transmitting results to the LIS. For easy editing, you can zoom any of these displays as needed. Also, the results can be printed or exported as a PDF, as a Flow Cytometry Standard listmode (.lmd) or native format (.crd) file. Both sample history changes and reagent information are available for review. All the displays can be zoomed in for easier editing.

The following Review screen (Figure 5.1) illustrates the different display areas associated with a typical Review screen. This is an example; individual results vary.

NOTE The Review Screen can be enabled to automatically refresh to display the last completed run.

NOTE If results for a run are auto-transmitted, they will not appear under the Review Tab. In this case, you will need to reference these results under the Results Tab. Refer to Results Screen.

Figure 5.1 Review Screen: Display Areas

1. Instrument Designation
2. Auto Refresh
3. Completed Run List
4. Sample Identifiers and Run Information Area
5. Result Table Area
6. Graphic Options Selection (drop-down box)
7. Action Toolbar
8. Selected Graphic Display
9. Pages

Review Screen Display Areas (left to right)

Figure 5.2 Instrument Designation

Instrument: All

Allows for selection of data based on the instrument that was used to collect the data. This is important when multiple instruments are networked. This feature allows the Operator to view data generated by other instruments from any networked Workstation.

Figure 5.3 Auto Refresh

☒ Auto Refresh

Automatically refreshes to display the last completed run.

Figure 5.4 Review Screen: Completed Runs List

Sample ID	Test Name	Run Date
11300064	Stem Single	7/21/2021 1
89352285182	FlowCheck_Alignment	7/21/2021 1
89352285179	Stem Panel	7/12/2021 3
89352285178	Stem Panel	7/12/2021 3
11300063	Stem Single	7/12/2021 1
11300064	Stem Single	7/12/2021 1
89352285171	FlowCheck_Alignment	7/12/2021 9
89352285171	FlowCheck_Alignment	7/7/2021 9:1
89352285171	FlowCheck_Alignment	7/6/2021 9:1
11300064	Stem Single	6/29/2021 1
11300063	Stem Single	6/29/2021 9
89352285171	FlowCheck_Alignment	6/29/2021 9
89352285172	Stem Panel	6/14/2021 2
11300064	Stem Single	6/14/2021 1
11300063	Stem Single	6/14/2021 1
89352285171	FlowCheck_Alignment	6/14/2021 9
89352285147	FlowCheck_Alignment	6/14/2021 9
11300063	Stem Single	6/11/2021 8
89352285147	FlowCheck_Alignment	6/11/2021 8
89352285147	FlowCheck_Alignment	6/10/2021 8

This area presents completed runs that are awaiting review and release. The last run completed is located at the top of the list. The runs are identified by:

- Sample ID
- Test Name
- Run Date
- Patient ID
- Patient Name
- ID

Any test that generated one or more Notification is highlighted in orange.

Any test that resulted in the generation of one or more Flags is highlighted in red.

Any test having both notification and flags is highlighted in red.

Samples not highlighted have no flags or notifications.

NOTE Use the horizontal scroll bar to see all fields.

Figure 5.5 Review Screen: Sample Identifiers and Run Information Area

Sample ID:	89352285172	Test:	Stem Panel
Run Date:	6/14/2021 2:38 PM	Collect Date:	06/14/2021
Specimen Type:	Mobilized Blood	Inst. SN:	AZ23062
Patient Name:	FreshCD34verifMWB		
Patient ID:	FreshCD34verifMWB df12	Gov. ID:	
Status:		Location:	
Gender:	Male	Date of Birth:	05/16/2002
Physician Name:			
Physician Code:		H.V.: 5	Weight: 70
		Dil.: 2	
Run Flags:			
Run Notifications:	Sample ID reused		
Comments:			

This area provides the sample information that was entered when the test request was generated. It also provides the information regarding when the test was run and displays any Flag or Notification messages generated during the run. Refer to [Flags And Notifications](#).

Next to the comments field is . Select this to see details about your sample, such as population statistics (non-IVD results). If a result is modified, the comment field is available for documentation regarding why the modification was performed.

 Gabriela A. Cividino
 Beckman Coulter Argentina S.A
 FARMACÉUTICA
 M.N. 15202/ M.P. 18093


 Rolando F. Bornia
 APODERADO

Figure 5.6 Review Screen: Result Table Area (3 Wells)


Result	Value	Flag	Normal Range	Action Range
Viable CD34+ HPC Count Well 1 (cells/uL)	25.40			
Viable CD34+ HPC Count Well 2 (cells/uL)	26.57			
AVG Viable CD34+ HPC Count (cells/uL)	25.99			
Viable CD34+ HPC % Diff (Well 1 vs Well 2)	-2.25%			
Viable Leukocytes Count Well 1 (cells/uL)	28,166.21			
Viable Leukocytes Count Well 2 (cells/uL)	29,594.76			
AVG Viable Leukocytes Count (cells/uL)	28,880.48			
Viable CD34+ HPC Well 1 % (of Viable Leukocytes)	0.09%			
Viable CD34+ HPC Well 2 % (of Viable Leukocytes)	0.09%			
AVG Viable CD34+ HPC % (of Viable Leukocytes)	0.09%			
DF: Dilution Factor	2.00			
AVG Viable CD34+ HPC Count (cells/uL) (x DF)	51.97			
AVG Viable Leukocyte Count (cells/uL) (x DF)	57,760.97			
HV: Harvest Volume (mL)	5.00			
AVG Viable CD34+ HPC Count (x DF x HV)	259,851.21			
BW: Body Weight (kg)	70.00			
AVG Viable CD34+ HPC Count per kg (x DF x HV / BW)	3,712.16			
CAL	989			
Viable CD34+ HPC Count Negative Ctrl (cells/uL)	0.28			

This area provides the results generated from the test run, in tabular format. Refer to [Quality Control Checks](#) for information on how results are rounded.

NOTE Parameters listed depend on the Test option requested: 1 well, 2 wells or 3 wells (as shown in [Figure 5.6](#)).

NOTE Results can be shown or hidden. Refer to Test Setup in Chapter 8 Setup in the AQUIOS CL Flow Cytometer Instructions For Use manual.

Figure 5.7 Review Screen: Graphic Options Selection Menu

Graphic options are accessed by selecting the desired view from a drop-down menu.

Results from each well (Well#1, Well#2 or Control) or one pre-defined graph combinations (Summary) are accessed by selecting the desired view from a drop-down menu:

- Well #1
- Well #2
- Control
- Summary

On [Figure 5.7](#) and [Table 5.1](#), the graphic options drop-down box specifies **Summary**.

Figure 5.8 Review Screen: Page Number

Allows you to switch between page 1 and page 2 to review all the graphics.


[Table 5.2](#) shows examples of graphics displayed on Page 1 and on Page 2.


NOTE Graphics shown are dependent upon the AQUIOS STEM test selected.

GABRIELA A. CIVIDINO
Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

ROLANDO F. BORNIÁ
Rolando F. Bornia
APODERADO

Table 5.1 Review Screen: Summary Area Display (AQUIOS STEM, Mobilized Whole Blood)

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093



Rolando F. Borna
APODERADO

Table 5.2 Review Screen: Summary Area Display, Page Number 1 (AQUIOS STEM Panel, Mobilized Whole Blood)

Graphics from Page1

Graphic Options

The following Review screens illustrate the different graphic options available for each individual patient result. There are 2 graphic options, one is for each tested well separately (such as Well #1, Well #2, and Control if the STEM Panel analysis option is used), and the other is combined wells Summary. These are examples only. Results may vary.

NOTE When a sample is selected for review, the Summary Review screen appears for that sample. Select the Well #1, Well #2 or Control option from the drop-down box.

Figure 5.10 Review Screen: Well#1 Page 2 (Mobilized Whole Blood)



Figure 5.11 Review Screen: Summary Page 1

Gabriela A. Cividino
Beckman Coulter Argentina S.A.
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Action Toolbox


Use the Action Toolbox (Figure 5.12) on the  screen to review, edit, and save changes to both regions and/or patient information prior to transmitting results to the LIS. Also, the results can be printed or exported in a PDF or CRD or LMD format (CRD and LMD files can be exported from the Review screen only). Both the sample history changes and the reagent information are available for review. All the displays can be zoomed in for easier editing. Refer to Chapter 6, Data Review in the AQUIOS CL Flow Cytometer Instructions for Use manual.

Figure 5.12 Action Toolbox


Data Review


Review Data Displays

IMPORTANT There is a risk of reporting incorrect results. Data displays used to arrive at the test result for light scatter patterns, antibody staining profiles, and all gates and boundaries should be reviewed by a laboratory professional when interpreting the data. If results are suspect, follow your laboratory's procedures to resolve.

Data review should include the following:

1. Review the light scatter.
 - a. Verify any population of interest is resolved as much as possible from other populations that are present (within biological constraints).
2. Review the antibody staining patterns.
 - a. Verify overall expected staining patterns are observed (the fluorescence staining correlates to the antibody used in sample preparation).
 - b. Verify abnormally dim or negative sample staining patterns are reviewed to ensure agreement with clinical presentation and other diagnostic data.
3. Review Gate and Region boundaries.
 - a. Verify gates and regions encompass any populations of interest and exclude any undesired events.
 - b. All red blood cells may not be lysed under the following conditions: presence of nucleated red blood cells, abnormal protein concentration, or hemoglobinopathies. Use of CD45 as a gating condition for leukocytes ensures that RBCs, as CD45 negative events, are not included in the leukocytes gate.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Edit Regions for a Review Sample


NOTE The viable region is automatically extended to capture all the 7AAD- cells.

IMPORTANT Do NOT move the regions for control runs.


IMPORTANT Sample ID, Test, Run Date, and Serial # cannot be edited.

IMPORTANT Administrator, Reviewer, or Editor role required. Editing or adjusting regions is not a routine adjustment.

- 1 Select **Review**

- 2 Select the sample of interest from the Review Screen Completed Run List. Select the run to edit from the menu below.


- 3 Select . The Edit action toolbox appears, and the fields for Sample ID, Test and Run Date and Instrument Serial Number turn grey to indicate that this is information that cannot be changed.


- 4 Select any other field to add or change patient information.

- 5 Select the plot that you want to edit.
- Double-click the graph to enlarge.
 - Select on the region you want to edit.

- 6 Once adjustment is done, select  to refresh the region statistics.

- 7 If you need to correct your modification, select  to reject the change. Repeat the steps above.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Borna
APODERADO

- 8 Select  to save the changes to the file.

Population Statistics

The statistics presented in the [Table 5.3](#) are used as indicators of pre-analytical sample handling issues and to assist in troubleshooting instrument performance issues.

Table 5.3 Population Statistics Recommended Use

Population Statistics	Recommended Use
CD34+ Viability % (of total CD34+)	Can be used as numeric aid for the direct assessment of CD34+ viability out of all CD34+ cells.
Sample Viability % (of WBC)	Provides a numeric aid for the direct assessment of overall sample viability out of all CD45+ white blood cells.
Total CD34+ HPC Count (cells/ μ L)	Can be used as numeric aid for the direct assessment of all CD34+ cells (both viable and non-viable).
WBC Count (cells/ μ L)	Can be used as numeric aid for the direct assessment of all CD45+ white blood cells (both viable and non-viable).


Adjust Compensation Manually

IMPORTANT Adjustment to compensation may be required in rare cases. Circumstances may include:

- At installation when compensation is performed the first time.
- After service is performed.
- If after reviewing a run, the operator determines that compensation needs to be adjusted.


IMPORTANT Compensation may only be adjusted by an Administrator, Reviewer, or Editor

If it is determined that compensation needs to be adjusted, follow the steps described below:

- 1 Select .
- 2 Select the sample of interest from the Review Details List.
- 3 Select . The Edit action toolbox appears.


4 Select .


5 Select the run to edit from the menu below. 1 is Well#1, 2 is Well#2 and #3 is Control:

6 Select  to edit compensation or select  to cancel the adjust compensation action.

7 Adjust from Display tab as show below:

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO

Figure 5.14 Over Compensation



Figure 5.15 Under Compensation



d. Select .

- 8 Select , select to save the changes to Future Runs or to cancel the changes.

NOTE Selecting Future Runs only applies to all tests that are run in the future that use this compensation. After selecting Future Runs, re-run of controls is recommended.

NOTE Compensation changes are recorded in the Maintenance Log. Refer to Maintenance Screen in Chapter 8, Setup, in the AQUIOS CL Flow Cytometer Instructions for Use manual.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

9 Select:

NOTE To revert to the previous compensation, go to the previous run and resave the compensation.

Result Screen

NOTE For a more detailed explanation of this screen and its components, refer to Chapter 6, Data Review in the AQUIOS CL Flow Cytometer Instructions for Use manual.

The Result screen provides a search menu to navigate through the results database. The Search criteria is entered under  and the results for the query are displayed under .

NOTE Results can also be accessed and viewed through  or .

To review the status of a sample, refer to the **Review Option** icons below (under the action toolbox).

NOTE Once the Review Option is selected, the sample is removed from the Review screen but it continues to be accessible through the Result screen.

indicates that the sample was reviewed.


indicates that the sample was reviewed and rejected.

indicates that the sample was not reviewed.


Using the Result Screen to Select Data for Review

Once you have selected , you will have access to the following options:

- Search results by the following criteria:
 - Sample ID
 - Patient ID
 - Gov. ID
 - Patient Name
 - Test

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

- Instrument Start Date
- Instrument End Date
- QC Results
- Reviewed Only
- Clear search criteria.
- Export results to .crd, .lmd, .pdf or .csv file.
- Do advanced searches using .
- Fields to be displayed as columns in View Results screen.
- Exclude run from View Results screen.
- Load a Query request.
- Append to a current search.
- Save a query result.

For detailed instructions for using any of these options, refer to Chapter 6 in the AQUIOS CL Flow Cytometer Instructions manual.

Flags And Notifications

CAUTION

Beckman Coulter Inc. does not claim to identify every abnormality in all samples. Beckman Coulter suggests using all available options to optimize the sensitivity of instrument results. Options include:

- **Flags**
- **Notifications**

The AQUIOS System Software generates a set of algorithm flags and notifications. The AQUIOS STEM Tests also generate flags and notifications that are specific to the AQUIOS STEM application.

The flag and notification messages appear on the Review and Result screen. Samples with flags are highlighted in red, while samples with notifications are highlighted in orange.

All results with flags or notifications are automatically sent to the review tab for further review. Only samples with no flags or notifications can be setup for automatic transmission. A laboratory professional should review all results, regardless of flags and notifications, prior to data release. To setup automatic transmission, refer to the **Auto Transmit** heading in Chapter 8 in the AQUIOS CL Flow Cytometer Instructions for Use manual.

- For the list of AQUIOS System Software flags, notifications, and specific result notifications, refer to Chapter 6 in the AQUIOS CL Flow Cytometer Instructions for Use manual.
- For the list of AQUIOS STEM specific flags and criteria, refer to [Table 5.4](#). For AQUIOS STEM specific notifications and criteria, refer to [Table 5.5](#).

Flags

A Flag indicates an issue that affects the results, follow the Recommended User Action in the table below.

NOTE Bleach is defined as high-quality, fragrance-free, gel-free bleach (5 to 6% solution of sodium hypochlorite - available chlorine).

Risk of chemical injury from bleach. To avoid contact with the bleach, use barrier protection, including protective eyewear, gloves, and suitable laboratory attire. Refer to the Safety Data Sheet for details about chemical exposure before using the chemical.

Table 5.4 AQUIOS STEM Specific Flags

Flag	Definition	Recommended User Action
Missing QC	Occurs if the Setup option QC Tracking was enabled and QC was not run successfully.	1. Confirm the Setup Option QC Tracking was set as expected. 2. Perform Daily QC procedure. 3. Reject, Re-Request from the Review Options menu. 4. If the problem persists, contact us.
General Error	An error occurred that did not fall into one of the categories below.	1. Run a Bleach Cycle followed by a Flush Cycle. 2. Reject, Re-Request from the Review Options menu. 3. If the problem persists, contact us.
High Count Rate	Based on events/second.	1. Confirm the sample WBC concentration that is within Beckman Coulter recommendation. Refer to Principles of AQUIOS STEM Tests. 2. Run a Bleach Cycle followed by a Flush Cycle. 3. Reject, Re-Request from the Review Options menu. 4. If the problem persists on the same sample, the sample is limited to run the application. Reject the sample. 5. If the problem persists, contact us.
Over Incubation	Sample staining incubation time exceeded.	1. Run a Bleach Cycle followed by a Flush Cycle. 2. Reject, Re-Request from the Review Options menu. 3. If the problem persists, contact us.
Over Lyse Incubation	Lyse incubation time exceeded.	1. Run a Bleach Cycle followed by a Flush Cycle. 2. Reject, Re-Request from the Review Options menu. 3. If the problem persists, contact us.
Clog or Bubble	The system detected unstable sample flow while analyzing the sample.	1. Run a Bleach Cycle followed by a Flush Cycle. 2. Reject, Re-Request from the Review Options menu. 3. If the problem persists, contact us.

Table 5.4 AQUIOS STEM Specific Flags

Flag	Definition	Recommended User Action
Control Range Failure	Only applies to daily QC control samples. Control results are outside of the acceptable ranges.	1. Run a Bleach Cycle followed by a Flush Cycle. 2. Check the controls expiration date. 3. Check if the controls blood level is sufficient. 4. Repeat Quality Control. 5. If the problem persists, contact us.
AQUIOS STEM-Count Mix Time Expired	Time since last mixing of flow counts exceeded (6 hours).	1. Mix the AQUIOS STEM-Count Fluorospheres, refer to the AQUIOS STEM Kit Instructions for Use. 2. Reject, Re-Request from the Review Options menu.
General algorithm Error	Algorithm did not complete the analysis of current sample.	1. Reject, Re-Request from the Review Options menu. 2. If the problem persists, contact us.
Not Enough Events	Not enough data points in the sample.	1. Verify that there is sufficient specimen in the tube. 2. Reject, Re-Request from the Review Options menu. 3. If the problem persists, contact us.
Bead Events Error	Number of events in the "fluorospheres" region is too low or too high.	1. Check the volume of AQUIOS STEM-Count Fluorospheres. 2. Run a Bleach Cycle followed by a Flush Cycle. 3. Reject, Re-Request from the Review Options menu. Mix the AQUIOS STEM-Count Fluorospheres before the repeated run.
Laser Power Error	Insufficient laser power.	1. Select Result Rejected from the Review Options Menu. 2. Perform a Shutdown per CHAPTER 7 in the AQUIOS CL Flow Cytometer Instructions For Use. 3. Perform a startup of the Startup per CHAPTER 3 in the AQUIOS CL Flow Cytometer Instructions For Use. 4. Review maintenance log and if the problem persists, contact us.
Vacuum Low	Insufficient vacuum.	1. Select Result Rejected from the Review Options Menu. 2. Perform a Shutdown per CHAPTER 7 in the AQUIOS CL Flow Cytometer Instructions For Use. 3. Perform a startup of the Startup per CHAPTER 3 in the AQUIOS CL Flow Cytometer Instructions For Use. 4. Review maintenance log and if the problem persists, contact us.
Power Supply Error	Detected by measuring the power supply voltages.	1. Select Result Rejected from the Review Options Menu. 2. Perform a Shutdown per CHAPTER 7 in the AQUIOS CL Flow Cytometer Instructions For Use. 3. Perform a startup of the Startup per CHAPTER 3 in the AQUIOS CL Flow Cytometer Instructions For Use. 4. Review maintenance log and if the problem persists, contact us.

Table 5.4 AQUIOS STEM Specific Flags


Flag	Definition	Recommended User Action
Onboard Expiration	A reagent has exceeded the maximum allowed time on board the system.	1. Replace onboard expired reagent with a non-expired reagent. 2. Reject, Re-Request from the Review Options menu. 3. If the problem persists, contact us.
Expired reagent	Based on barcode expiration dates (open, on-board, and closed vial), The reagent is run after the expiration date has passed.	1. Replace expired reagent with a non-expired reagent. 2. Reject, Re-Request from the Review Options menu. 3. If the problem persists, contact us.

Notifications

Notifications are generated to convey an abnormal cell distribution or population which may be specimen specific. Follow the Recommended User Action in the table below.

Table 5.5 AQUIOS STEM Specific Notifications

Notification	Definition	Recommended User Action
High CD34 count in control run	CD34 counts in the negative control run is above the requirements: 1. Average of CD34 cells is less than 50 cells/μL , IsoClonic Control of CD34 shall be less than or equal to 5 cells/μL 2. Average of CD34 cells is greater or equal to 50 cells/μL , IsoClonic Control of CD34 shall be less than or equal to 10% of the average of CD34 cells count.	1. Review the gating of CD34+HPC, CD45+ and fluorospheres from the Review screen. 2. Adjust the gating if determined to be necessary. 3. If the problem persists after gate adjustment, lab professionals shall evaluate the results based on their assessment and standard practice in lab. NOTE Notification is not removed if the count goes back to normal range after adjusting the gating.
Action Range	Result was outside the user defined action range.	1. Follow your laboratory's policies and procedures for out-of-range results for the specific Action Range.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO


Table 5.5 AQUIOS STEM Specific Notifications

Notification	Definition	Recommended User Action
Sample ID Reused	Only patient runs are identified with this notification. This occurs when: 1. The Sample ID for this specimen run has been used in a prior run (patient or control) within the last 24 hours. 2. There is a duplicate Sample ID in the Running tab for this run.	1. Ensure this is a re-run of a previous run. 2. Check Sample ID of the specimen to ensure it is correct for the specimen.
User Modified	The demographics or the gating for this run were modified by a user.	1. Review the changes and determine if run is to be rejected and/or re-run (Rejected, Re-Request option) or accepted (Reviewed).
Less than 100 CD34+ events acquired	Less than 100 CD34+ events acquired	1. This notification is for informational use only, no further action is required. Lab professionals shall evaluate the results based on their assessment and standard practice in laboratory.

AQUIOS STEM Software Observations

While running QC for AQUIOS STEM CD34 L1 and L2 Control Cells, the software asks if Flow-Count Fluorospheres have been mixed. For the AQUIOS STEM application you are using AQUIOS STEM-Count Fluorospheres, not Flow-Count Fluorospheres even though the message says Flow-Count Fluorospheres.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Performance Characteristics

Overview


The performance characteristics presented below were collected on the AQUIOS CL Flow Cytometry System with AQUIOS System software and AQUIOS STEM application with no modification of the algorithm generated results.

This chapter contains the following information:


- Linearity
- Method Comparison
- Reference Interval
- Precision
- Specificity
- Quality Control

Linearity

Three replicates were made at each of 13 dilutions for CD34+ and CD45+. Ranges for CD34+ and CD45+ were established by staining an Apheresis sample with AQUIOS STEM Kit reagents, and analyzing by flow cytometry using AQUIOS system software. Statistical determination of linearity was performed using the guideline CLSI EP06-Ed2 Nov 2020, Evaluation of Linearity of Quantitative Measurement Procedures. See [Figure 6.1, Linearity: CD34+ cells/μL](#) and [Figure 6.2, Linearity: CD45+ cells/μL](#). Values are expressed in terms of absolute count (cells/μL).

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Figure 6.1 Linearity: CD34+ cells/ μ L



Figure 6.2 Linearity: CD45+ cells/ μ L

Limit of Blank/Limit of Detection/Limit of Quantitation

Limit of Blank (LoB) is the highest measurement result that is likely to be observed (with a stated probability or type I error) for a blank sample per the CLSI EP17-A3 guideline, Protocols for Determination of Limits of Detection and Limits of Quantitation.

Limit of Detection (LoD) is the lowest number of cells in a sample that can be detected with (stated) probability, although perhaps not quantified as an exact value, per the CLSI-EP17-A3 guideline, Protocols for Determination of Limits of Detection and Limits of Quantitation.

Lower Limit of Quantitation (LLoQ) is the lowest number of cells in a sample that can be quantitatively determined with stated acceptable precision and trueness, under stated

experimental conditions per the CLSI EP17-A3 guideline, Protocols for Determination of Limits of Detection and Limits of Quantitation.

Table 6.1 LoB, LoD and LLoQ Results

Analyte	LoB	LoD	LLoQ	Upper 95% Confidence Limit of LLoQ
CD34c	0.49	1.49	2.00	4.00
CD45c	31.33	35.33	36.00	36.00

Analytical Measuring Ranges


The Analytical Measuring Ranges provided indicate the AQUIOS STEM System's analytical capabilities to enumerate CD34+ cells and CD45+ cells. The AQUIOS STEM Analytical Measuring Ranges are derived from Linearity and Lower Limit of Quantitation values:

- **Linearity:** Ranges over which results are directly proportional to the amount of analyte in a sample, per the CLSI EP06-Ed2 Nov 2020, Evaluation of Linearity of Quantitative Measurement Procedures.
- **Lower Limit of Quantitation:** The lowest amount of cell number in a sample that can be quantitatively determined with stated acceptable precision and trueness, under stated experimental conditions per the CLSI EP17-A3 guideline, Protocols for Determination of Limits of Detection and Limits of Quantitation.

Table 6.2 Analytical Measuring Ranges

Parameters	Units	AQUIOS STEM Analytical Measuring Ranges
CD 34+	cells/ μ L	4-1000
CD 45+	cells/ μ L	500-30000

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Method Comparison

The comparability of AQUIOS STEM System was assessed by comparing the results with the predicate method using the estimation of difference as described in CLSI EP9-A3 Method Comparison and Bias Estimation Using Patient Samples. For viable CD34+ stem cells (count and %) and CD45+ cells (count). The comparator method is the stemCXP System. A total of 442 specimens were collected from a multi-center clinical study, including 161 Whole Blood, 53 Mobilized Whole Blood, 101 Bone Marrow, 97 Apheresis Product, and 30 Cord Blood.

Table 6.3 Regression statistics for all specimen types combined

Specimen Type	Parameter	N*	Correlation	Slope	Intercept
Overall	Viable CD34+ (cells/uL)	289	0.990	1.00	0.04
	% viable CD34+ of viable CD45+	442	0.966	0.92	0.01
	Viable CD45+ (cells/uL)	440	0.991	1.00	161.57

* CD34+ and CD45+ counts below LLOQ have been excluded from the respective analyses

Table 6.4 Method Comparison study results of AQUIOS STEM system compared with predicate method

Parameter	N*	Percentile	Level	Bias	95% Confidence Limit of Bias
Viable CD34+ (cells/uL)	289	25th	3.00	0.06	-0.87 – 1.00
		50th	25.75	0.06	-0.54 – 0.66
		75th	98.50	0.06	-2.25 – 2.37
% viable CD34+ of viable CD45+	442	25th	0.05	0.01	-0.10 – 0.12
		50th	0.22	-0.00	-0.08 – 0.07
		75th	0.86	-0.05	-0.13 – 0.02
Viable CD45+ (cells/uL)	440	25th	5701.00	169.78	94.31 – 245.26
		50th	9104.50	173.50	92-73 – 254.27
		75th	16486.00	181.55	-39.43 – 402.54

* CD34+ and CD45+ counts below LLOQ have been excluded from the respective analyses

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO

Figure 6.3 Regression Analysis: CD34+ cells/ μ L


Figure 6.4 Regression Analysis: CD34+ %



Figure 6.5 Regression Analysis: CD45+ cells/ μ L

PN B75978AA

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

6-5

Table 6.5 Regression statistics by specimen type

Specimen Type	Parameter	N*	Correlation	Slope	Intercept
Peripheral blood	Viable CD34+ (cells/uL)	64 (normal: 18, mobilized: 46)	0.996	0.94	0.47
	% viable CD34+ of viable CD45+	214 (normal: 161, mobilized: 53)	0.999	0.99	80.95
	Viable CD45+ (cells/uL)	213 (normal: 161, Mobilized: 52)	0.994	0.97	-0.01
Bone Marrow	Viable CD34+ (cells/uL)	98	0.982	1.00	1.66
	% viable CD34+ of viable CD45+	101	0.942	0.90	241.61
	Viable CD45+ (cells/uL)	100	0.988	1.03	0.05
Apheresis	Viable CD34+ (cells/uL)	97	0.991	1.04	-2.07
	% viable CD34+ of viable CD45+	97	0.995	1.06	8.05
	Viable CD45+ (cells/uL)	97	0.972	0.96	-0.00
Cord Blood	Viable CD34+ (cells/uL)	30	0.990	0.89	0.28
	% viable CD34+ of viable CD45+	30	0.980	0.95	319.15
	Viable CD45+ (cells/uL)	30	0.978	0.96	-0.01

Reference Interval

A total of 120 normal Whole Blood specimens from apparently healthy adult donors were collected from 3 geographically diverse sites in order to assess the reference interval using AQUIOS STEM System. The selection of donors is consistent with guidelines stated in CLSI EP28-A3C Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory.

Data analysis was performed in accordance to CLSI EP28-A3c. Non-parametric calculation of the 95% upper and lower reference limits for CD34 and CD45, with calculation of 90% confidence intervals were summarized in the table below. These values are intended to be representative only. Each laboratory must establish its own.

Table 6.6 Reference Intervals

Parameter	Units	N	Mean	Lower Limit (90% confidence bounds)	Upper Limit (90% confidence bounds)
CD34+	Absolute viable count (cells/uL)	120	2.75	0.88 (0.54 – 1.02)	6.16 (5.33 – 10.80)
% CD34+	% viable CD34+ of viable CD45+	120	0.05	0.01 (0.01 – 0.02)	0.09 (0.08 – 0.14)
CD45+	Absolute viable count (cells/uL)	120	6203.35	3746.91 (2660.36 – 4363.74)	8901.82 (8262.02 – 9883.23)

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Precision

Using Control Cells

The imprecision performance of the AQUIOS STEM System was demonstrated using AQUIOS STEM Control Cells (Level 1 and Level 2). One of three sites ran Control Cells in duplicate, twice a day for 20 days. The other two sites ran Control Cells in triplicate, twice a day for 5 days. Absolute count and percentage of CD34+ cells and absolute count of CD45+ cells were measured.

Data analysis was performed in accordance with EP05-A3 and summarized below in Table 6.7 and Table 6.8.

Table 6.7 AQUIOS STEM CD34 Control Cells Level 1

Name	N	Mean	Repeatability		Between Runs		Between Days		Between Sites		Reproducibility	
			SD	CV%	SD	CV%	SD	CV%	SD	CV%	SD	CV%
CD34+	140	8.71	1.206	13.85	0.000	0.00	0.305	3.50	0.286	3.28	1.276	14.66
% CD34+	140	0.14	0.019	13.71	0.000	0.00	0.004	2.94	0.01	6.90	0.022	15.62
CD45+	140	6173.88	203.567	3.30	0.000	0.00	69.142	1.12	266.13	4.31	342.119	5.54

Table 6.8 AQUIOS STEM CD34 Control Cells Level 2

Name	N	Mean	Repeatability		Between Runs		Between Days		Between Sites		Reproducibility	
			SD	CV%	SD	CV%	SD	CV%	SD	CV%	SD	CV%
Name	N	Mean	SD	CV%	SD	CV%	SD	CV%	SD	CV%	SD	CV%
CD34+	140	28.85	2.203	7.64	0.375	1.30	0.091	0.32	0.000	0.00	2.236	7.75
% CD34+	140	0.46	0.030	6.56	0.000	0.00	0.000	0.00	0.016	3.40	0.034	7.39
CD45+	140	6311.97	216.117	3.42	105.999	1.68	80.976	1.28	304.309	4.82	396.363	6.28

Using Clinical Specimens

In order to demonstrate the repeatability precision performance of the AQUIOS STEM System for clinical specimens, a total of 63 samples per intended use with a mix of specimen and anticoagulant types were collected from multiple clinical sites. Each clinical sample was analyzed with minimum three, maximum five replicates, and the absolute count and percentage of CD34+ cells and absolute count of CD45+ cells were measured.

Data analysis was performed in accordance with CLSI H26-A2. Mean, standard deviation and coefficient of variations were calculated for each specimen. Profiles for CD34 and CD45 were created based on the relationship between means and %CV. Percent CV was estimated at percentiles and relevant medical decision levels (5, 10, 20 CD34+ cells/ μ L).

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO

Table 6.9 Repeatability Precision with clinical specimens

Parameter	Percentile	Mean	%CV
Viable CD34+ (cells/uL)	25th	14.57	7.55
	50th	24.12	6.67
	75th	84.77	4.89
% viable CD34+ of viable CD45+	25th	0.07	10.88
	50th	0.28	7.03
	75th	1.09	4.59
Viable CD45+ (cells/uL)	25th	5873.96	2.28
	50th	9576.57	2.30
	75th	21054.44	2.35

Specificity

Reagent Specificity

The CD45 antigen is expressed on every type of hematopoietic cell except mature erythrocytes and their immediate progenitors.^{8, 9} It has not been detected in differentiated nonhematopoietic tissue.^{8, 9, 10, 11}


The CD34 antigen is expressed on virtually all hematopoietic precursor cells including multipotent stem cells.


CD34 epitopes have been separated into three classes based upon sensitivity to neuraminidase and to *Pasteurella haemolytica* glycoprotease.⁴² The CD34 antibody recognizes a class III epitope of the CD34 molecule.⁴⁴

Quality Control

Ensure that controls and Flow-Check Fluorospheres are run according to manufacturer's recommendations as listed in the instrument manuals. Target value ranges are provided for each lot of CD34 Control Cells, Level 1 and Level 2.

Before samples are analyzed, assayed controls need to be run to verify antibody reactivity.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Application Specifications

Overview

This chapter contains the following information:

- AQUIOS STEM Preparation
- AQUIOS STEM Throughput

AQUIOS STEM Preparation

All samples are prepared and analyzed onboard using a 96-deep well plate.

AQUIOS STEM Sample Preparation

Refer to [Table 4.1](#) for test descriptions and reagents used accordingly.

The following description applies to the AQUIOS STEM Single test (= 1 well from the 96-Deep Well Plate).

The test uses 43 µL of samples stained with 13 µL of each reagent (monoclonal antibodies and viability dye). After 15 minutes of incubation, the sample is lysed using 430 µL Lysing Reagent, then AQUIOS STEM-Count is added after the lyse incubation of approximately 15 minutes. The sample is then mixed and aspirated for analysis.

For AQUIOS STEM Duplicate or AQUIOS STEM Panel, 2 or 3 wells will be used respectively

AQUIOS STEM Throughput


Time from loading a specimen to result for the first sample is approximately 45 minutes for single (1 well) option, 52 minutes for duplicate (2 well) option, and 59 minutes for panel (3-well) option.

Each additional well results take 7 minutes, so results for single (1 well) option will be 7 minutes, whereas duplicate (2 well) option will be 14 minutes, and panel (3 well) option will be 21 minutes for final result.

NOTE AQUIOS STEM System uses a flow rate of 70 µL/minute.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Borna
APODERADO

Rolando F. Bornia
APODERADO

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Adjustment and Replacement Procedures

Overview

This chapter contains the following information:


- AQUIOS STEM Replacement Procedures

AQUIOS STEM Replacement Procedures

Replacing the AQUIOS STEM Lysing Solution

NOTE You should discard a reagent bottle (on board) that is expiring and use a new reagent bottle to avoid errors.


- 1 Select  on the status bar. When prompted, open the reagent door by pressing the top middle part of the door in one single, solid motion.

- 2 Remove the empty AQUIOS STEM Lysing Solution bottle.
- 3 Remove the caps on the new AQUIOS STEM Lysing Solution bottle before loading.


NOTE These caps are not pierceable.

- 4 Place the new bottles in the appropriate holders with the barcodes facing outwards. The AQUIOS STEM Lysing Solution bottles are interchangeable and can be placed in any of the three available positions.

NOTE Remove all reagents having invalid barcode reading before proceeding.

- 5 Close the reagent door. The reagents and plate are scanned, then the Reagent Levels screen appears.


- 6 Close the Reagent Levels screen and the Main screen appears.


Replacing the AQUIOS STEM-Count Fluorospheres

IMPORTANT Do not use the AQUIOS Flow-Count Fluorospheres reagent in replacement of the AQUIOS STEM-Count Fluorospheres. The system will not start.


- 1 Select  on the status bar. When prompted, open the reagent door by pressing the top middle part of the door in one single, solid motion.

- 2 Remove the empty AQUIOS STEM-Count Fluorospheres bottle, pull out the bottle and adapter as a single unit. Hold the adapter in one hand and pull out the AQUIOS STEM-Count Fluorospheres bottle with the other hand.


- 3 Properly mix the new AQUIOS STEM-Count Fluorospheres bottle (vortex 10 to 12 seconds is required).
 - 4 Remove the cap on the new AQUIOS STEM-Count Fluorospheres bottle, and completely remove the aluminum seal before loading.
NOTE This cap and aluminum seal are not pierceable.
- IMPORTANT** Place the specific holder for AQUIOS STEM-Count Fluorospheres in the instrument prior to placing the AQUIOS STEM-Count Fluorospheres bottle in the holder.
- 5 Place the screw cap with hole on the AQUIOS STEM-Count Fluorospheres bottle.
 - 6 Place the new AQUIOS STEM-Count Fluorospheres bottle in the adapter with the barcode facing outward and ensure that there is no space between the bottom of the AQUIOS STEM-Count Fluorospheres bottle and the adapter. The fit will be snug.


NOTE If you need to rotate the bottle in the holder, ensure that the barcode or label is not scratched by the securing metal spring.


Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

- 7 Place the adapter into the AQUIOS STEM-Count Fluorospheres holder. The adapter can only be inserted one way, the key in the adapter must fit in the groove of the holder as shown.

- 8 Close the reagent door. The reagents and plate are scanned, and the Reagent Levels screen appears.


NOTE On the Reagent Levels screen, if a sample is running, the AQUIOS STEM-Count Fluorospheres column displays the number of tests scheduled to have AQUIOS STEM-Count Fluorospheres added in parentheses. AQUIOS STEM-Count Fluorospheres information will also be available on the Reagent Levels screen including name, date, and number of tests remaining.

9 Close the Reagent Levels screen and the Main screen appears.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Overview

This appendix contains the following information:

- [Reports](#)
- [QC Report](#)

Reports

Reports can be printed or saved (PDF) from the Review or the Result screens. The report will be generated for each well individually, nevertheless all the results (IVD parameters and population statistics) will be included. For instance, in the case of a sample run with AQUIOS STEM Panel, the report for the well#1 will include only the graphics from the well#1, but all the results from the 3 wells (Well#1, Well#2 and Control).

Refer to below examples of Reports for various specimens.

NOTE The report can be customized under System Setup to display the name, address, and facility logo.


- Example report formats for STEM
 - [Figure A.1, Report: STEM Panel - Apheresis - Well #1 Page 1](#)
 - [Figure A.2, Report: STEM Panel - Apheresis - Well #1 Page 2](#)
 - [Figure A.3, Report: STEM Panel - Apheresis - Well #1 Page 3](#)
 - [Figure A.4, Report: STEM Panel - Apheresis - Summary Page 1](#)
 - [Figure A.5, Report: STEM Panel - Apheresis - Summary Page 2](#)
 - [Figure A.6, Report: STEM Panel - Apheresis - Summary Page 3](#)

NOTE The parameters with ** are population statistics.

Also refer to [Report Description](#) and [QC Report Description](#).

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO

Figure A.1 Report: STEM Panel - Apheresis - Well #1 Page 1

Sample ID: 89352211984
Run Date: 15 Mar 2022 15:56:15
User: Admin
Specimen Type: Apheresis
Patient ID:
Status:
Gender: F
Physician:
Instrument Serial #: AZ23062
Report: Well 1
Run Flags:
Run Notifications:
Comments:

Test: Stem Panel v1.0
Analysis Date: 15 Mar 2022 16:14:49
Collect Date: 14 Mar 2022
Patient Name:
Gov. ID:
Location:
Date Of Birth: 12 Oct 1977
Physician Code:
Harvest Volume (mL): 50
Weight (kg): 75 Dilution: 20

** Population Statistics

Page #1 of 3

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Figure A.2 Report: STEM Panel - Apheresis - Well #1 Page 2


Sample ID: 89352211984

Test: Stem Panel v1.0


Run Date: 15 Mar 2022 15:56:15

Analysis Date: 15 Mar 2022 16:14:49

Result	Value	Flag	Normal Range	Action Range
Viable CD34+ HPC Count Well 1 (cells/uL)	31.61			
Viable CD34+ HPC Count Well 2 (cells/uL)	29.72			
AVG Viable CD34+ HPC Count (cells/uL)	30.66			
Viable CD34+ HPC % Diff (Well 1 vs Well 2)	3.08%			
Viable Leukocytes Count Well 1 (cells/uL)	9,308.80			
Viable Leukocytes Count Well 2 (cells/uL)	9,146.02			
AVG Viable Leukocytes Count (cells/uL)	9,227.41			
Viable CD34+ HPC Well 1 % (of Viable Leukocytes)	0.34%			
Viable CD34+ HPC Well 2 % (of Viable Leukocytes)	0.32%			
AVG Viable CD34+ HPC % (of Viable Leukocytes)	0.33%			
DF: Dilution Factor	20.00			
AVG Viable CD34+ HPC Count (cells/uL) (x DF)	613.25			
AVG Viable Leukocyte Count (cells/uL) (x DF)	184,548.21			
HV: Harvest Volume (mL)	50.00			
AVG Viable CD34+ HPC Count (x DF x HV)	30,662,346.20			
BW: Body Weight (kg)	75.00			
AVG Viable CD34+ HPC Count per kg (x DF x HV / BW)	408,831.28			
CAL	962			
Viable CD34+ HPC Count Negative Ctrl (cells/uL)	0.61			

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Figure A.3 Report: STEM Panel - Apheresis - Well #1 Page 3

Sample ID: 89352211984 Test: Stem Panel v1.0
Run Date: 15 Mar 2022 15:56:15 Analysis Date: 15 Mar 2022 16:14:49

Statistic	Value	Flag	Normal Range	Action Range
**CD34+ Viability Well 1 % (of total CD34+)	98.41%			
**CD34+ Viability Well 2 % (of total CD34+)	96.31%			
**AVG CD34+ Viability % (of total CD34+)	97.36%			
**Total CD34+ HPC Count Well 1 (cells/uL)	32.12			
**Total CD34+ HPC Count Well 2 (cells/uL)	30.86			
**AVG Total CD34+ HPC Count (cells/uL)	31.49			
**WBC Count Well 1 (cells/uL)	9,601.19			
**WBC Count Well 2 (cells/uL)	9,493.66			
**AVG WBC Count (cells/uL)	9,547.43			
**Sample Viability Well 1 % (of WBC)	96.96%			
**Sample Viability Well 2 % (of WBC)	96.34%			
**AVG Sample Viability % (of WBC)	96.65%			
**Total CD34+ HPC Count Negative Ctrl (cells/uL)	0.70			
**DF: Dilution Factor	20.00			
**AVG Total CD34+ HPC Count (cells/uL) (x DF)	629.75			
**AVG WBC Count (cells/uL) (x DF)	190,948.54			
**Viable Events Well 1	72,803			
**Viable Events Well 2	72,364			
**Viable Leukocytes Events Well 1	72,748			
**Viable Leukocytes Events Well 2	72,322			
**Viable CD34+ HPC Events Well 1	247			
**Viable CD34+ HPC Events Well 2	235			
**Viable CD34+ HPC Events Negative Ctrl	13			
**Single Beads Events Well 1	7,518			
**Single Beads Events Well 2	7,607			
**Single Beads Events Negative Ctrl	20,535			

Signature:

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093



Rolando F. Bornia
APODERADO

Figure A.4 Report: STEM Panel - Apheresis - Summary Page 1

Sample ID:	89352211984	Test:	Stem Panel v1.0
Run Date:	15 Mar 2022 15:56:15	Analysis Date:	15 Mar 2022 16:14:49
User:	Admin	Collect Date:	14 Mar 2022
Specimen Type:	Apheresis	Patient Name:	
Patient ID:		Gov. ID:	
Status:		Location:	
Gender:	F	Date Of Birth:	12 Oct 1977
Physician:		Physician Code:	
Instrument Serial #:	AZ23062	Harvest Volume (mL):	50
Report:	Summary	Weight (kg):	75
Run Flags:		Dilution:	20
Run Notifications:			
Comments:			

** Population Statistics

Page #1 of 3

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Figure A.5 Report: STEM Panel - Apheresis - Summary Page 2

Sample ID: 89352211984 Test: Stem Panel v1.0
Run Date: 15 Mar 2022 15:56:15 Analysis Date: 15 Mar 2022 16:14:49

Result	Value	Flag	Normal Range	Action Range
Viable CD34+ HPC Count Well 1 (cells/uL)	31.61			
Viable CD34+ HPC Count Well 2 (cells/uL)	29.72			
AVG Viable CD34+ HPC Count (cells/uL)	30.66			
Viable CD34+ HPC % Diff (Well 1 vs Well 2)	3.08%			
Viable Leukocytes Count Well 1 (cells/uL)	9,308.80			
Viable Leukocytes Count Well 2 (cells/uL)	9,146.02			
AVG Viable Leukocytes Count (cells/uL)	9,227.41			
Viable CD34+ HPC Well 1 % (of Viable Leukocytes)	0.34%			
Viable CD34+ HPC Well 2 % (of Viable Leukocytes)	0.32%			
AVG Viable CD34+ HPC % (of Viable Leukocytes)	0.33%			
DF: Dilution Factor	20.00			
AVG Viable CD34+ HPC Count (cells/uL) (x DF)	613.25			
AVG Viable Leukocyte Count (cells/uL) (x DF)	184,548.21			
HV: Harvest Volume (mL)	50.00			
AVG Viable CD34+ HPC Count (x DF x HV)	30,662,346.20			
BW: Body Weight (kg)	75.00			
AVG Viable CD34+ HPC Count per kg (x DF x HV / BW)	408,831.28			
CAL	962			
Viable CD34+ HPC Count Negative Ctrl (cells/uL)	0.61			

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO

Figure A.6 Report: STEM Panel - Apheresis - Summary Page 3

Sample ID: 89352211984

Test:

Stem Panel v1.0

Run Date: 15 Mar 2022 15:56:15


Analysis Date: 15 Mar 2022 16:14:49

Statistic	Value	Flag	Normal Range	Action Range
**CD34+ Viability Well 1 % (of total CD34+)	98.41%			
**CD34+ Viability Well 2 % (of total CD34+)	96.31%			
**AVG CD34+ Viability % (of total CD34+)	97.36%			
**Total CD34+ HPC Count Well 1 (cells/uL)	32.12			
**Total CD34+ HPC Count Well 2 (cells/uL)	30.86			
**AVG Total CD34+ HPC Count (cells/uL)	31.49			
**WBC Count Well 1 (cells/uL)	9,601.19			
**WBC Count Well 2 (cells/uL)	9,493.66			
**AVG WBC Count (cells/uL)	9,547.43			
**Sample Viability Well 1 % (of WBC)	96.96%			
**Sample Viability Well 2 % (of WBC)	96.34%			
**AVG Sample Viability % (of WBC)	96.65%			
**Total CD34+ HPC Count Negative Ctrl (cells/uL)	0.70			
**DF: Dilution Factor	20.00			
**AVG Total CD34+ HPC Count (cells/uL) (x DF)	629.75			
**AVG WBC Count (cells/uL) (x DF)	190,948.54			
**Viable Events Well 1	72,803			
**Viable Events Well 2	72,364			
**Viable Leukocytes Events Well 1	72,748			
**Viable Leukocytes Events Well 2	72,322			
**Viable CD34+ HPC Events Well 1	247			
**Viable CD34+ HPC Events Well 2	235			
**Viable CD34+ HPC Events Negative Ctrl	13			
**Single Beads Events Well 1	7,518			
**Single Beads Events Well 2	7,607			
**Single Beads Events Negative Ctrl	20,535			

Signature:

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Report Description

Report Heading

For a description of how to set up report headings, refer to Appendix A, Software Screens in the AQUIOS CL Flow Cytometer Instructions for Use manual.

Displays: Generating Reports

Reports can be generated using the available pre-defined formats for data display. Use the Graph Format drop-down menu under the Review or Result screens to select the display desired for the report.

Result Table Heading

- **Result** - use to view the assay statistics
- **Value** - use to view the patient's result obtained with the AQUIOS system
- **Flag** - Results are flagged as H or L if they are outside the normal ranges, and HH or LL if outside of the action ranges
- **Normal Range** - Reference values defined by the user in the Setup window, Test Settings tab
- **Action Range** - Reference values defined by the user in the Setup window, Test Settings tab

Statistics: AQUIOS STEM Report

- Viable CD34+ HPC Count (cells/ μ L)
- Viable Leukocytes Count (cells/ μ L)
- Viable CD34+ HPC % (of Viable Leukocytes)


QC Report

The QC reports can be printed, saved (PDF) or created as a spreadsheet (CSV) file from any of the three QC screens. All three screens generate the same all-inclusive report containing the QC data in the following order: Results, Instrument, and Instrument Drift. Note that the QC report statistics depend on the assay selected under the Test drop-down menu. Also, the report can be customized under System Setup to display the name, address and facility logo. For generating a QC report, refer to the AQUIOS CL Flow Cytometer Instructions for Use manual.

The QC report can show all lot numbers or a single lot number. QC can be displayed for up to four months using the "Period" option. If lot (All) is selected, available lot numbers under this time frame will be displayed in the same graph identified with a different color.

NOTE Users with Administrator privileges or Reviewer/Editor privileges can exclude/restore data points.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

QC Report Description

QC Result Report Heading, QC Instrument Report Heading, QC Instrument Drift Report Heading

- **Facility Name** - location where test was performed.
- **Test** - assay used for the control.
- **Control** - AQUIOS control selected for the report.
- **Lot number** - single lot selected or (All) if plotting all available lots.
- **Date** - month, year, and period selected for the report.

NOTE QC can be displayed for up to four months using the "Period" option.

QC Result Displays & Statistics (Depend on Test Selected)

- Viable CD34+ HPC Count (cells/ μ L)
- Viable Leukocytes Count (cells/ μ L)
- Viable CD34+ HPC % (of Viable Leukocytes)


QC Result First Table Heading

- **Lot** - lot number used for the specified QC run.
- **Average** - average of all listed runs.
- **Min** - lowest value of all listed runs.
- **Max** - highest value of all listed runs.
- **CV** - coefficient of variation of all listed runs (Standard Deviation/Average) x 100
- **StDev** - standard deviation of all listed runs.

QC Result Second Table Heading, QC Instrument Second Table Heading, QC Instrument Drift Second Table Heading


- **Exclude** - data points included or excluded in the report.
- **Lot** - single lot selected or (All) if plotting all available lots.
- **Date run** - date in which control was analyzed.
- **Pass/Fail** - If Pass, data point is within expected range. If Fail, data point is not within the expected range.
- **Assay statistics** - percent or count/ μ L of selected assay.

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

QC Reports

Figure A.7 QC Results: Results View of AQUIOS STEM CD34 Control Cells Level 1, AQUIOS STEM Single

Gabriela A. Cividino
Gabriela A. Cividino
Beckman Coulter Argentina S.A.
FARMACEUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
Rolando F. Bornia
APODERADO

Figure A.8 QC Results: AQUIOS STEM CD34 Control Cells Level 1, AQUIOS STEM Single (Page 1)


Figure A.9 QC Results: AQUIOS STEM CD34 Control Cells Level 1, AQUIOS STEM Single (Page 2)


Result	Lot	Avg	Min	Max	CV	StDev
Viable CD34+ HPC Count (cells/uL)	10180063	8.83	7.87	9.91	8.38 %	0.74
Viable Leukocytes Count (cells/uL)	10180063	6,534.89	6,121.83	7,223.83	5.78 %	377.72
Viable CD34+ HPC % (of Viable Leukocytes)	10180063	0.14%	0.12%	0.14%	6.24 %	0.01%

Exclude (Yes/No)	Lot	DateRun	Instrument	Pass/Fail	Viable CD34+ HPC Count (cells/uL)	Viable Leukocytes Count (cells/uL)	Viable CD34+ HPC % (of Viable Leukocytes)
No	10180063	4/7/2021	5	Pass	8.82	6,121.83	0.14%
No	10180063	4/8/2021	5	Pass	9.40	6,603.64	0.14%
No	10180063	4/9/2021	5	Pass	7.87	6,511.14	0.12%
No	10180063	4/12/2021	5	Pass	8.73	6,449.65	0.14%
No	10180063	4/13/2021	5	Pass	9.91	7,223.83	0.14%
No	10180063	4/15/2021	5	Pass	8.26	6,299.23	0.13%

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Borna
APODERADO

Figure A.10 QC Results: Results View of AQUIOS STEM CD34 Control Cells Level 2, AQUIOS STEM Single

Gabriela A. Cividino
Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
Rolando F. Bornia
APODERADO

Figure A.11 QC Results: AQUIOS STEM CD34 Control Cells Level 2, AQUIOS STEM Single (Page 1)

Results:



Figure A.12 QC Results: AQUIOS STEM CD34 Control Cells Level 2, AQUIOS STEM Single (page 2)

Result	Lot	Avg	Min	Max	CV	StDev
Viable CD34+ HPC Count (cells/uL)	10180064	31.57	28.90	35.13	7.19 %	2.27
Viable Leukocytes Count (cells/uL)	10180064	6,697.32	6,301.14	7,166.41	5.20 %	348.27
Viable CD34+ HPC % (of Viable Leukocytes)	10180064	0.47%	0.43%	0.50%	5.56 %	0.03%


Exclude (Yes/No)	Lot	DateRun	Instrument	Pass/Fail	Viable CD34+ HPC Count (cells/uL)	Viable Leukocytes Count (cells/uL)	Viable CD34+ HPC % (of Viable Leukocytes)
No	10180064	4/7/2021	5	Pass	30.89	6,301.14	0.49%
No	10180064	4/8/2021	5	Pass	33.08	6,652.57	0.50%
No	10180064	4/9/2021	5	Pass	35.13	7,166.41	0.49%
No	10180064	4/12/2021	5	Pass	28.90	6,307.89	0.46%
No	10180064	4/13/2021	5	Pass	31.60	6,805.36	0.46%
No	10180064	4/15/2021	5	Pass	29.79	6,950.55	0.43%

Gabriela A. Clivdino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Borna
APODERADO

Rolando F. Bornia
APODERADO

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093


References



1. Centers for Disease Control and Prevention. 2003 Guidelines for performing single-platform absolute CD4+ T-cell determinations with CD45 gating for persons infected with Human Immunodeficiency Virus (HIV). MMWR 2003; 57 (No. RR. 2): 1-13.
2. Nicholson JKA, Hubbard M and Jones BM. Use of CD45 fluorescence and side scatter characteristics for gating lymphocytes when using the whole blood lysis procedure and flow cytometry. Cytometry (Communications in Clinical Cytometry), 1996; 26: 16-21.
3. Centers for Disease Control and Prevention. Guidelines for performing single-platform absolute CD4+ T-cell determinations with CD45 gating for persons infected with human immunodeficiency virus. MMWR, 2003;52 (No. RR-2):1-13.
4. CDC: 1997 Revised Guidelines for Performing CD4+ T-Cell Determinations in Persons Infected with Human Immunodeficiency Virus (HIV). MMWR 1997; 46 (RR-2); 1-29.
5. H42-A2 CLSI: Enumeration of Immunologically Defined Cell Populations by Flow Cytometry: Approved Guideline - 2nd Edition 2007.
6. Zloza, et al. Multiple populations of T lymphocytes are distinguished by the level of CD4 and CD8 coexpression and require individual consideration. Journal of Leukocyte Biology 2006; 79 (4-6).
7. Alter G, Malenfant JM, Delabre RM, Burgett NC, Yu XG, Lichterfeld M, Zaunders J, Altfeld M. Increased Natural Killer Cell Activity in Viremic HIV-1 Infection. The Journal of Immunology 2004(173):5305-5311.
8. Coffman RL and Weissman IL. A B-cell-specific member of the T200 glycoprotein family. Nature, 1981, B220; 289: 681-683.
9. Dalchau R and Fabre JW. Identification with a monoclonal antibody of a predominantly B lymphocyte-specific determinant of the human leukocyte common antigen. J Exp Med, 1980; 153: 753-765.
10. Omary MB, Trowbridge IS and Battifora HA. Human homologue of murine T200 glycoprotein. J Exp Med, 1980; 152: 842-852.
11. Dalchau R, Kirkley J and Fabre JW. Monoclonal antibody to a human leukocyte-specific membrane glycoprotein probably homologous to the leukocyte-common (L-C) antigen of the rat. Eur J Immunol, 1981; 10: 737-744.
12. Reinherz EL and Schlossman SF. The differentiation and function of human T lymphocytes. Cell, 1980; 19:821-827.
13. Reinherz EL, Meuer S, Fitzgerald KA, Hussey RE, Levine H, and Schlossman SF. Antigen recognition by human T lymphocytes is linked to surface expression of the T3 molecular complex. Cell, 1982; 30: 735-743.
14. Meuer SC, Acuto O, Hussey RE, Hodgdon JC, Fitzgerald KA, Schlossman SF, and Reinherz EL. Evidence for the T3-associated 90K heterodimer as the T-cell antigen receptor. Nature, 1983; 303: 808-810.
15. De Martini RM and Parker JW. Immunologic alterations in human immunodeficiency virus infection: A review. J Clin Lab Anal, 1989; 3:56-70.
16. Reinherz EL, Meuer SC, and Schlossman SF. The delineation of antigen receptors on human T lymphocytes. Immunol Today, 1983; 4: 5-8.

17. Caligiuri M, Murray C, Buchwald D, Levine H, Cheney P, Peterson D, Komaroff AL, and Ritz J. Phenotypic and functional deficiency of natural killer cells in patients with chronic fatigue syndrome. *J Immunol*, 1987; 139: 3306-3313.
18. Ravetch, J.V., Perussia, B., "Alternative membrane forms of FcγRIII (CD16) on human natural killer cells and neutrophils. Cell type-specific expression of two genes that differ in single nucleotide substitutions." 1989, *J. Exp. Med.*, 170(2), 481-497.
19. Huizinga, T.W., Roos, D., von dem Borne, A.E., "Neutrophil Fc-gamma receptors: A two-way bridge in the immune system", 1990, Mar 15, *Blood*, 75(6), 1211-1214.
20. Lanier, L.L., Yu, G., Phillips, J.H., "Coassociation of CD3 zeta with a receptor (CD16) for IgG Fc on human natural killer cells", 1989 Dec 14, *Nature*, 342(6251), 803- 805.
21. Hibbs, M., L., Selvaraj, P., Carpen, O., Springer, T.A., Kuster, H., Jouvin, M. H., Kinet, J.P., "Mechanisms for regulating expression of membrane isoforms of Fc gamma RIII (CD16)", 1989 Dec 22, *Science*, 246(4937), 1608-1611.
22. Tamm, A., Schmidt, R.E., "The binding epitopes of human CD16 (Fc gamma RIII) monoclonal antibodies: Implications for ligand binding", 1996 Aug 15, *J. Immunol.*, 157(4), 1576-1581.
23. Tamm, A., Bassmann, Schmidt, R.E., "Natural killer cell structural studies: Localization of the epitopes of human CD16 (FcRIII) monoclonal antibodies on the molecular model of CD16", 1997, *Leucocyte Typing VI, White Cell Differentiation Antigens*, Kishimoto, T., et al, Eds., Garland Publishing, Inc., 324-326.
24. Ritz, J., Trinchieri, G., Lanier, L.L., "NKcell antigens: section report", 1995, *Leucocyte Typing V, White Cell Differentiation Antigens*. Schlossman, S.F., et al., Eds., Oxford University Press, 1367-1372.
25. McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, Hogg N, Horton M, Ling N, MacLennan ICM, Mason DY, Milstein C, Spiegelhalter D and Waldman H, eds. 1987. *Leukocyte Typing III*. Oxford University Press, Oxford, UK. pp. 38, 40, 42, 43, 116, 167, 170-172, 176, 199, 202, 206, 302-308, 315, 475, 796-801.
26. Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, Ritz J, Shaw S, Silverstein R, Springer R, Tedder TF, and Todd RF, eds. 1995. *Leukocyte Typing V*. Oxford University Press, Oxford, UK. Volume 1 pp. 262, 263, 268, 270, 507-511, Volume 2 pp. 1398-1400.
27. Barclay AN, Birkeland ML, Brown MH, Beyers AD, Davis SJ, Somoza C, Williams AF, eds. *The Leukocyte Antigen Facts Book*. Academic Press, London, 1993; 106-109, 142-143, 228-229.
28. Tedder TF and Isaacs CM. 1989. Isolation of a cDNAs encoding the CD19 antigen of human and mouse B lymphocytes: a new member of the immunoglobulin superfamily. *J Immunol* 143:712-717.
29. Griffin JD, Hercend T, Beveridge R and Schlossman SF. 1983. Characterization of an antigen expressed by human natural killer cells. *J Immunol* 130:2947-2951.
30. Hercend T, Griffin JD, Bensussan A, Schmidt RE, Edson MA, Brennan A, Murray C, Daley JF, Schlossman SF and Ritz J. 1985. Generation of monoclonal antibodies to a human natural killer clone. Characterization of two natural killer-associated antigens, NKH1A and NKH2, expressed on subsets of large granular lymphocytes. *J Clin Invest* 75:932-943.
31. Schmidt RE, Murray C, Daley JF, Schlossman SF and Ritz J. 1986. A subset of natural killer cells in peripheral blood displays a mature T-cell phenotype. *J Exp Med* 164:351-356.
32. Schmidt RE, Michon JM, Woronicz J, Schlossman SF, Reinhertz EL, and Ritz J. 1987. Enhancement of natural killer function through activation of the T11 E-rosette receptor. *J Clin Invest* 79:305-308.

33. Caligiuri M, Murray C, Buchwald D, Levine H, Cheney P, Peterson D, Komaroff AL and Ritz J. 1987. Phenotypic and functional deficiency of natural killer cells in patients with chronic fatigue syndrome. *J Immunol* 139:3306-3313.
34. Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman SF, eds: 1984. *Leukocyte Typing*. Springer-Verlag, New York, NY, pp. 28, 41-42, 44, 196.
35. Knapp W, Dorken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, and Kr. von dem Borne AEG, eds. 1989. *Leukocyte Typing IV, White Cell Differentiation Antigens*. Oxford University Press, Oxford, UK. pp. 22, 36-38, 295, 296, 536, 541, 699-702.
36. Garfall AL, Fraietta JA, Maus MV. Immunotherapy with chimeric antigen receptors for multiple myeloma. *Discov Med*. 2014 Jan;17(91):37-46.
37. Sutherland DR, Anderson L, Keeney M, Nayar R, Chin-Yee I: The ISHAGE guidelines for CD34+ cell determination by flow cytometry. *International Society of Hematotherapy and Graft Engineering. J Hematother*. 1996 Jun;5(3):213-26.
38. Keeney M, Chin-Yee I, Weir K, Popma J, Nayar R, Sutherland DR: Single platform flow cytometric absolute CD34+ cell counts based on the ISHAGE guidelines. *International Society of Hematotherapy and Graft Engineering. Cytometry*. 1998 Apr 15;34(2):61-70.
39. EDQM Council of Europe: *European Pharmacopoeia (Ph. Eur.)*. 10th Edition 2019-2022. Website: <https://pheur.edqm.eu/home>
40. Foundation for the Accreditation of Cellular Therapy (FACT), Joint Accreditation Committee – ISCT and EBMT (JACIE): *International standards for hematopoietic cellular therapy product collection, processing, and administration*. 8th Edition, Version 8.0, 2021. Website: <http://www.factwebsite.org/>
41. Foundation for the Accreditation of Cellular Therapy (FACT): *NetCord-FACT international standards for cord blood collection, banking, and release for administration*. 7th Edition, 2019. Website: <http://www.factwebsite.org/>
42. Whitby A, Whitby L, Fletcher M, Reilly JT, Sutherland DR, Keeney M, Barnett D: ISHAGE protocol: are we doing it correctly? *Cytometry B Clin Cytom*. 2012 Jan;82B:9-17.
43. Sutherland DR, Marsh JCW, Davidson J, Baker MA, Keating A, and Mellors A. Differential sensitivity of CD34 epitopes to cleavage by *Pasteurella haemolytica* glycoprotease: Implications for purification of CD34-positive progenitor cells. *Exp Hematol*, 1992; 20:590-599.
44. Greaves MF, Titley I, Colman SM, Buhning H-J, Campos L, Castoldi GL, Garrido F,
45. Gaudernack G, Girard J-P, Ingles-Esteve J, Invernizzi R, Knapp W, Lansdorp PM, Lanza F, Merle-Beral H, Parravicini C, Razak K, Ruiz-Cabello F, Springer TA, van der Schoot CE, and Sutherland DR. CD34 cluster workshop report, published in *Leukocyte Typing V*, Schlossman SF et al., eds. Oxford University Press, Oxford, UK. 1995; 840-846.

 Gabriela A. Cividino
 Beckman Coulter Argentina S.A
 FARMACÉUTICA
 M.N. 15202/ M.P. 18093

 Rolando F. Bornia
 APODERADO

Rolando F. Bornia
APODERADO

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

A

- Action Toolbox
 - Review screen, 5-9
- adjust compensation
 - manually, 5-11
- Analysis Algorithm
 - 3-dimensional lymphocyte gating, 2-4
 - overview, 2-4
- Analysis reagent
 - Sheath Fluid, 1-5
- AQUIOS Manuals
 - Use, xv
- AQUIOS STEM Panel monoclonal antibody
 - reagents, 1-4
- AQUIOS Tetra System software
 - used with, 1-3

B

- bar codes
 - smart-track monitoring, 1-3
- Beckman Coulter Customer Support Center,
 - contacting, ii
- blood volume
 - minimum required per specimen tube, 4-3

C

- call center, contact information, ii
- compensation
 - to manually adjust, 5-11
- contact information, Beckman Coulter
 - Customer Support Center, ii
- control lot number, registering. *See* Starting a New Lot of Beckman Coulter Control in Chapter 9 of AQUIOS CL Flow Cytometer Instructions for Use

D

- Daily Startup. *See* Startup in Chapter 3 of AQUIOS CL Flow Cytometer Instructions for Use

 Gabriela A. Cividino
 Beckman Coulter Argentina S.A
 FARMACÉUTICA
 M.N. 15202/ M.P. 18093

- data sheets, material safety
 - how to order, 1-6
- disclaimer
 - graphics, xviii

F

- Flags and Notifications, 5-16

G

- graphic
 - options, 5-6
- graphics disclaimer, xviii

H

- help, Beckman Coulter Customer Support Center, ii

I


- intended use, 1-1

L

- limitations, 1-6
- linearity, 6-1
- Lymphocyte Immunophenotyping Panels
 - Tetra-1 Panel, 2-3

M

- manual overview, xv
- manual, updates (primary entry), -iii
- manually adjust
 - compensation, 5-11
- material safety data sheets (MSDS/SDS)
 - how to order, 1-6
- method comparison, 6-4
- MSDS (material safety data sheets)
 - how to order, 1-6

 Rolando F. Bornia
 APODERADO

O

- options
 - graphic, 5-6
 - running specimens, 4-10
- overview
 - 3-dimensional lymphocyte gating, 2-4
 - Analysis Algorithm, 2-4

P

- Parameter results of STEM Single, Duplicate, and Panel
 - parameter results, 2-3
- population statistics, 5-11
 - recommended use, 3-1, 3-7, 5-11
- precision
 - Precision Performance of Quantitative Measurement Methods, 6-7
- preparation
 - monoclonal antibodies, 1-4
 - Tetra-1 sample, 7-1
- Principles of Tetra Tests, 2-1

Q

- Q C workflow, 3-1
- QC
 - report, A-8
 - report description, A-9
 - Results screen, 3-16
- QC materials
 - required, 3-2
- QC overview, 3-2
- QC process, 3-3
- Quality Control, 6-8
- Quality Control reagents
 - AQUIOS STEM CD34 Control Cells, 1-5

R

- reagent
 - specificity, 6-8
- reference ranges, 6-6
- registering control lot number. *See* Starting a New Lot of Beckman Coulter Control in Chapter 9 of AQUIOS CL Flow Cytometer Instructions for Use, 3-3
- replacing reagents and consumables. *See* AQUIOS CL Flow Cytometer Instructions for Use

- report
 - QC, A-8
- report description
 - displays
 - generating reports, A-8
 - QC, A-9
 - report heading, A-8
 - result table heading, A-8
 - statistics, Tetra-1 Report, A-8
- reports, A-1
- Results screen
 - search menu, 5-15
 - selecting data for review, 5-15
- Review screen
 - Action Toolbox, 5-9
 - display areas, 5-2
 - samples pending transmission, 5-1
- Running screen, 4-12
- running specimens
 - options, 4-10

S

- safety symbols, xviii
 - potentially infectious, xviii
- sample populations
 - to automatically identify and enumerate, 2-4
- SDS (safety data sheets)
 - how to order, 1-6
- service, contact information, ii
- specificity
 - reagent, 6-8
- STEM
 - throughput, 7-1
- STEM Application
 - Workflow, 4-2
- STEM Single, Duplicate, Panel Results, 2-3
- storage conditions and stability
 - reagents, 4-4
 - Whole Blood Specimen, 4-4
- support, Beckman Coulter customer, ii
- System Components
 - "closed" system, 1-3
 - autoloader, 1-2
 - fully automatic operation, 1-2
 - load & go feature, 1-3
 - Single-tube Loader, 1-2
- System Guide
 - about, xvi

system limitations, 1-6

T


throughput
STEM, 7-1

W

Workflow
STEM Application, 4-2


Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACEUTICA
M.N. 15202/ M.P. 18093


Rolando F. Bornia
APODERADO

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18993

Rolando F. Bornia
APODERADO

Beckman Coulter, Inc.

Customer End User License Agreement

This Product contains software that is owned by Beckman Coulter, Inc. or its suppliers and is protected by United States and international copyright laws and international trade provisions. You must treat the software contained in this Product like any other copyrighted material. This license and your right to use the Product terminate automatically if you violate any part of this agreement.

This is a license agreement and not an agreement for sale. Beckman Coulter hereby licenses this Software to you under the following terms and conditions:

You May:

1. Use this software in the computer supplied to you by Beckman Coulter;
2. Maintain one copy of this software for backup purposes (the backup copy shall be supplied by Beckman Coulter);
3. After written notification to Beckman Coulter, transfer the entire Product to another person or entity, provided you retain no copies of the Product software and the transferee agrees to the terms of this license agreement.

You May Not:

1. Use, copy or transfer copies of this Software except as provided in this license agreement;
2. Alter, merge, modify or adapt this Software in any way including disassembling or decompiling;
3. Loan, rent, lease, or sublicense this Software or any copy.

Limited Warranty

Beckman Coulter warrants that the software will substantially conform to the published specifications for the Product in which it is contained, provided that it is used on the computer hardware and in the operating system environment for which it was designed. Should the media on which your software arrives prove defective, Beckman Coulter will replace said media free of charge within 90 days of delivery of the Product. This is your sole remedy for any breach of warranty for this software.

Except as specifically noted above, Beckman Coulter makes no warranty or representation, either expressed or implied, with respect to this software or its documentation including quality, performance, merchantability, or fitness for a particular purpose.

No Liability for Consequential Damages

In no event shall Beckman Coulter or its suppliers be liable for any damages whatsoever (including, without limitation, damages for loss of profits, business interruption, loss of information, or other pecuniary loss) arising out of the use of or inability to use the Beckman Coulter Product software. Because some states do not allow the exclusion or limitation of liability for consequential damages, the above limitation might not apply to you.


General

This agreement constitutes the entire agreement between you and Beckman Coulter and supersedes any prior agreement concerning this Product software. It shall not be modified except by written agreement dated subsequent to the date of this agreement signed by an authorized Beckman Coulter representative. Beckman Coulter is not bound by any provision of any purchase order, receipt, acceptance, confirmation, correspondence, or otherwise, unless Beckman Coulter specifically agrees to the provision in writing. This agreement is governed by the laws of the State of California.


Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

PN B/5978AA

Rolando F. Bornia
APODERADO


Warranty-1

Gabriela A. Cividino
Beckman Coulter Argentina S.A
FARMACÉUTICA
M.N. 15202/ M.P. 18093

Rolando F. Bornia
APODERADO

Related Documents

Your AQUIOS documentation can be found on our Web site at www.beckman.com

AQUIOS STEM System Guide

PN B75978

AQUIOS CL Flow Cytometer Instructions for Use

PN B21896

AQUIOS STEM LIS Host Transmission

PN C62645

Cytometry Interlaboratory Quality Assurance Program Procedure Manual

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INFORMATION FOR USA ONLY
 Reagent for the identification and enumeration of human hematopoietic progenitor cells.
 Reactive ingredients: Formaldehyde <1%, Pyrrolidine <5%, EDTA <1%, Disodium Salt <1%, Potassium Bicarbonate <1%,
 7-Aminoclonomycin D 0.005%, Carbonic Anhydrase <1%, Bovine Serum Albumine <1%, Polystyrene Latex Beads <5%,
 Monoclonal Antibodies (Murine Source) <1%.

May be covered by one or more pat. - see www.beckmancoulter.com/patents.

BECKMAN COULTER **IMMUNOTECH SAS**
 a Beckman Coulter Company
 13276 Marseille Cedex 9 France, + (33) 4 91 17 27 27

DANGER
 H315 H317 H319 H341 H350
 Formaldehyde <3%, Dimethyl
 Sulfoxide 1 - 2%, Pyrrolidina
 1.0 - 1.5%.

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
RÓTULO ORIGINAL DEL PRODUCTO

3) C43667 AQUIOS STEM CD34 Control Cells Kit

AQUIOS STEM CD34 Control Cells IVD Info. for USA Only Rx Only Intended as an assayed, lysable whole blood quality control product for immunophenotyping analysis using monoclonal antibody reagents for flow cytometry. Reactive Ingredients: Stabilized Human Erythrocyte and Leukocyte Preparation 2 x 3 mL May be covered by one or more pat. – see www.beckman.com/patents. Manufactured for IMMUNOTECH S.A.S. a Beckman Coulter Company  beckman.com/techdocs  IMMUNOTECH S.A.S. a Beckman Coulter Company 130, avenue de Lattre de Tassigny, BP 177 13276 Marseille cedex 9, France Tél. (33) 491 172 727 C43555-AA Made in U.S.A.	AQUIOS STEM CD34 Control Cells REF C43667 CONTENTS 1x 3 mL Level 1 CONTENTS 1x 3 mL Level 2 Package Insert 2°C 8°C  BECKMAN COULTER C43550-AA
---	---

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1. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
2. Sutherland DR, Anderson L, Keeney M, Nayar R, Chin-Yee I: The ISHAGE guidelines for CD34+ cell determination by flow cytometry. International Society of Hematotherapy and Graft Engineering. *J Hematother.* 1996 Jun;5(3):213-26.
3. Keeney M, Chin-Yee I, Weir K, Popma J, Nayar R, Sutherland DR: Single platform flow cytometric absolute CD34+ cell counts based on the ISHAGE guidelines. International Society of Hematotherapy and Graft Engineering. *Cytometry.* 1998 Apr 15;34(2):61-70.
4. Whitby A, Whitby L, Fletcher M, Reilly JT, Sutherland DR, Keeney M, Barnett D: ISHAGE protocol: are we

doing it correctly? Cytometry B Clin Cytom. 2012
Jan;82B:9-17.

5. Centers for Disease Control and Prevention, National
Institutes of Health: Biosafety in Microbiological and
Biomedical Laboratories. 5th Edition. HHS Publication
No. (CDC) 21-1112, Dec. 2009.

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Fabricado para IMMUNOTECH S.A.S.
IMMUNOTECH S.A.S. es una compañía de Beckman Coulter
Referencia: BECKMAN COULTER ARGENTINA S.A. rótulos e instrucciones de uso