ROTULOS

1.

5010004

AK-Calibrant 4x1 ml

LABELS

p. 1 | 1


confidential

Forma 1 - 10/2016 - 04

Vial Labels

Box Label

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6.

5005100

TECHNOCLOT PT Owren Capillary Calibration Set

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

7.

5005102

TECHNOCLOT® PT Owren Capillary Control Set

TECHNOCLOT® PT Owren
Capillary CaCl₂ 17 mM 1mL

TECHNOCLOT® PT Owren
Capillary Control AK 1mL

TECHNOCLOT® PT Owren
Capillary Control N 1mL

TECHNOCLOT® PT Owren
Capillary Water 1mL

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Technoclone
TECHNOCLOT® PT Owren
Capillary Control Set

Kit Contents:

- 1x Capillary Control N 1mL 1P37800
- 1x Capillary Control AK 1mL 2B35800
- 2x CaCl₂ 17mM 1mL 9U428B0
- 2x Water 1mL W1418B0

TECHNOCLOT® PT Owren
Capillary Control Set

Capillary Control Set for accuracy and precision control of TECHNOCLOT® PT Owren manual test.

REF 5005102

TECHNOCLOT® PT Owren
Capillary Control Set

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



9.

5006011
Ceveron TGA RB Kit

TC technoclone
LABELS,
TVB300.01
Date
p. 1 | 1 confidential

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


IF-2026-81838177-APN-INPM#ANMAT

5006013
Ceveron TGA RC Low Kit

technoclone
LABELS,
TVB400.01
Date
p. 1 | 1 confidential

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11.

5006015
Ceveron TGA RC High Kit

technoclone
LABELS,
TVB500.01
Date
p. 1 | 1 confidential

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12.

5006347
Ceveron TGA CAL Set
1Set

technoclone
LABELS,
TVB0700.01
Date
p. 1 | 1 confidential

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PRODUCT DESCRIPTION

INTENDEND USE

The AK calibration set can be used for the following purposes:

- For direct INR determination of a patient sample (INR reference curve)
- For establishing laboratory, instrument and reagent specific ISI and PT normal plasma values
- For constructing a reference curve in percentage of the normal.

COMPOSITION

The calibration set contains 4 calibration plasmas for standardising the PT test. The AK-Calibrant A is a lyophilised normal plasma pool which is produced from the selected citrated plasmas (Austrian source plasma) of at least 100 blood samples of healthy donors. AK-Calibrant B, C, and D are produced exclusively from the plasma of donors on long-term oral anticoagulant treatment. As a result, the prothrombin complex factors (II, VII, IX, X, Protein C and Protein S) are present at the level normally expected in plasma from patients under oral anticoagulant treatment, together with the Proteins induced by vitamin K absence (PIVKA inhibitors).

No preservatives are added to the lyophilised plasmas.

mL	Reagent	Description
1 x 1	AK-Calibrant A	lyophilised normal plasma, INR approx.1.0
1 x 1	AK-Calibrant B	lyophilised AK-plasma, INR about 2.0
1 x 1	AK-Calibrant C	lyophilised AK-plasma, INR about 3.0
1 x 1	AK-Calibrant D	lyophilised AK-plasma, INR about 4.0

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water
- Control Plasma and reagents*

REF	5020040	Coagulation Control N	5 x 1 mL
REF	5020050	Coagulation Control N	50 x 1 mL
REF	5011050	Coagulation Control AK	5 x 1 mL
REF	5011060	Coagulation Control AK	50 x 1 mL
REF	5003128	TECHNOCLLOT® PT	10 x 4 mL
REF	5003129	TECHNOCLLOT® PT	10 x 10 mL
REF	5003220	TECHNOCLLOT® PT Plus	5 x 4 mL
REF	5003222	TECHNOCLLOT® PT Plus	10 x 10 mL
REF	5005032	TECHNOCLLOT® PT Owren	10 x 4 mL
REF	5005037	TECHNOCLLOT® PT Owren	10 x 10 mL

* or any other package sizes, special Ceveron® alpha or TECHNOCLLOT® Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of AK-Calibrant are tested and found negative for Hb_sAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution:

RT**	+15°C	+4°C
4 hours	6 hours	8 hours

The vials can be only frozen once.

Upon storage, caps should be screwed tightly.

**=room temperature

TEST PROCEDURE

PREPARATION AND PERFORMANCE OF THE TEST

CEVERON

Technoclone provides application sheets Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

- Reconstitute the lyophilised plasma in 1 mL distilled water (RT) before performing the test and allow to stand for 10 minutes***. For standardising tests a reconstitution time of 30 minutes is recommended.
- The reconstituted plasmas should be treated in the same way as patient samples, using the instructions given for the test.

ANALYSES RESULTS

EVALUATION

Evaluations in item 1 and 2 are based on the INR definition.

1. Establishing an INR Reference Curve

Determine the coagulation times of AK-Calibrant A, B, C and D in duplicate and plot the mean values on double log paper (x-axis: INR values; y-axis: Coagulation time in seconds see Fig. 1). From this reference curve the INR's of patient plasmas can be directly read off.

2. Establishing Laboratory Specific ISI's

A linear relationship exists between the log INR value and log PT (sec) as expressed by the equation $[\log(\text{PT patient plasma}) = 1/\text{ISI} \times \log(\text{PT normal plasma})]$.

A local ISI value can either be calculated by linear regression analysis or read from the graphical plot of the results using the enclosed graph sheet. This sheet contains a pre-set ISI scale printed on the Y2 axis which is derived mathematically from the equation given above.

3. Graphical Evaluation

- The INR calibration curve is established as described above in Section 1 using the enclosed graph evaluation sheet and the curve extended to intersect the Y1 and the Y2 axis.
- A second line is then drawn parallel to the calibration curve starting from the co-ordinates INR 1.0/10 sec on the x-axis and continuing through to the Y2 axis. This line is referred to as the standardised straight-line (see Figure 2). The local ISI value for the thromboplastin is read on the printed ISI scale on the Y2 axis where the standardised straight-line intersects this axis (see Figure 2).
- An alternative method to the standardised straight-line is to measure the distance (a) between the calibration curve and x-axis at the co-ordinates INR 1.0/10 sec and then measure the same distance down on the Y2 axis starting from the intersect of the calibration curve on this axis (see Figure 2)

4. Establishing a Laboratory Specific PT Normal Plasma

Using the following equation:

$$\text{PT normal plasma (sec)} = \frac{\text{PT AK-Calibrant A (sec)}}{\text{INR AK-Calibrant A}}$$

5. Setting up the Reference (Calibration) Curve as a Percentage of the Norm

For thromboplastins, e. g. TECHNOCLLOT® PT Plus, the percentages of norm values may be read off from the enclosed table. These values are specific to each thromboplastin and cannot be transferred to other reagents. In order to establish the reference curve, the percentages of the norm are plotted against the mean time values in seconds and the points joined by a straight line (Fig. 3).

STANDARDISATION

In order to monitor oral anticoagulant treatment the Prothrombin time (PT) should, according to the WHO, be reported as an International Normalised Ratio (INR). The INR is calculated using the ISI (International Sensitivity Index) by the following formula:

$$\text{INR} = \left[\frac{\text{TPZ patient plasma}}{\text{TPZ normal plasma}} \right]^{\text{ISI}}$$

The ISI of a thromboplastin is assigned by the manufacturer, and is batch specific. The ISI is also influenced by the type of instrument used and variations can be clearly seen when different instruments are used.

LIMITATION OF THE TEST

A new calibration is required for each batch of thromboplastin and for each instrument used. Also, if software changes are introduced or following a major service of either instruments or equipment, recalibration is recommended.

The figures in the table are only applicable to the indicated lot number of AK-Calibrant and the method has to adhere strictly as prescribed for the respective reagent.

PERFORMANCE CHARACTERISTICS

Within-run precision was assessed over multiple runs using a specific lot of reagent and AK-Calibrant plasmas. The coefficients of variation obtained in this study were < 3% for PT sec for Level A, and < 6% for Levels B-D.

ISTH and CLSI documents state that the difference between duplicate PT (seconds) as well as the difference between the mean of the run (day- to-day) results must be within 10%.

TRACEABILITY

The INR's for AK-Calibrant are determined according to WHO (WHO_TRS_889_A3_1999) guidelines by Technoclone and/or by external quality assessment.

LITERATURE

Please contact Technoclone or your local distributor.

*** When modifying methods, the reconstitution volumes might change.

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RU

SRB

REF	5003009	TECHNOPLASTIN® HIS	12 x 2 mL
REF	5003010	TECHNOPLASTIN® HIS f. Ceveron	12 x 4 mL
REF	5003030	TECHNOPLASTIN® HIS	2 x 10 mL
REF	5003026	TECHNOPLASTIN® HIS	6 x 10 mL
REF	5003028	TECHNOPLASTIN® HIS f. Ceveron	6 x 20 mL
REF	5003021	TECHNOPLASTIN® HIS	20 x 10 mL
REF	5003022	TECHNOPLASTIN® HIS f. Ceveron	20 x 20 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Kľúčová slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / дистрибуирана вода / destilovaná voda / Destilisaná Voda		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg bruksvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fidiť se instrukciami / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / počet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дистрибуирана вода / destilovaná voda / Destilisaná Voda	LOT	lot / Charge / lotto / lote / lote / lot / sats / série / Parti / партиа / партида номер / šarže / лот / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionsbuffer / Reaktionspuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	MTP	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροπιλοδότησης / Микротительна пласка / Микропланшет / Mikrotitracijska destička / Mikrotitracione ploče
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / färdig att användas / færdig til brug / klar til bruk/ έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / разтворете или разрежете с / zředíť anebo rozpustit v / разбавить или растворить в / naředite nebo rozpustite v / razrediti ili rastvoriti u	SUB	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer / Inkubationsbuffer / Inkubationspuffer / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
IVD	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico in vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рпш диагностичкѣм εντόβ σωλήνα / за ин витро диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnostický prostředek in vitro / in vitro dijagnostika	IF-2026-81838177-APN-INPM#ANMAT	

TECHNOPLASTIN® HIS

Ca-thromboplastin for the determination of the prothrombin time (PT, Quick)

GB

PRODUCT DESCRIPTION INTENDEND USE

TECHNOPLASTIN® HIS is used for the determination of the prothrombin time (Quick test). This screening test is used

- for the control of blood coagulation disorders of the extrinsic system,
- for the surveillance of oral anticoagulant therapy,
- for separate determination of individual coagulation factors of the extrinsic system,
- to check the liver's function of synthesizing coagulation factors in hepatopathies.

A heparin neutralizing agent in the reagent permits the assessment of plasma obtained from patients on conventional heparin therapy (0.2 - 0.8 U. hep./mL).

COMPOSITION

TECHNOPLASTIN® HIS (HIS = HeparinInSensitive) is a standardised Ca-thromboplastin reagent obtained from rabbit brain which is characterized by sensitivity to coagulation factors II, V, VII and X. In addition, the reagent contains a heparin neutralizing agent.

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- Pipettes
- Distilled water
- Control Plasmas and Calibrators

REF	5020040	Coagulation Control N **	5 x 1 mL
REF	5011050	Coagulation Control AK **	5 x 1 mL
REF	5010004	AK-Calibrant	4 x 1 mL

** or any other package sizes, special Ceveron® alpha or Technoclott Control and Calibration reagents of Technoclone.

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STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...+8 °C.

Stability after reconstitution:

37 °C 6 hours	RT* 24 hours	12 °C (Ceveron) 5 days	+2...+8 °C 6 days	-20 °C 2 months
------------------	-----------------	---------------------------	----------------------	--------------------

Upon storage, caps should be screwed tightly. * = room temperature

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Plasma separation: Mix 9 parts of venous blood and 1 part of sodium citrate solution (0.11 mol/L) and centrifuge for 15 min. at an RCF of at least 2500 (corresponding to DIN 58905). Keep the plasma at room temperature (up to 4 hours).

PREPARATION OF REAGENT

The lyophilized reagent is to be reconstituted with the indicated volume of distilled water at room temperature and may be used immediately. For standardisation a reconstitution time of 30 min is recommendend.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides Application sheets for Ceveron® alpha. The Application sheets contain analyser/assay specific handling and performance information which may differ from that provided in these instruction for use. In this case the information contained in the Application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Pipette the following substances into a test tube preheated to 37°C:

	Citrated plasma	0.1 mL
	Incubate for 1 min at + 37°C	
+	TECHNOPLASTIN® HIS (37°C)	0.2 mL

Upon addition of TECHNOPLASTIN® HIS, start the stopwatch or the coagulometer resp. and determine the clotting time.

ANALYSES RESULTS

CALCULATION OF THE RESULTS

The prothrombin time is indicated in seconds, in % of normal (Quick), or in INR. Values can be converted by means of the enclosed reference curve in form of a table, considering the indicated method. Reference curves have to be checked for accuracy by means of control plasma prior to use. If the values are beyond the confidence range given for control plasma, a separate reference curve has to be prepared.

For the preparation of the reference curve Technoclone's calibration plasma or the following dilutions of a normal plasma pool (Coagulation Control N) with isotonic saline solution may be used.

- % of normal dilution	100 % undil.	50 % 1+1	25 % 1+3	12.5 % 1+7	10 % 1+9
---------------------------	-----------------	-------------	-------------	---------------	-------------

These dilutions are to be determined like patients samples. The mean values of triple determinations are plotted on a reciprocal curve sheet and joined linearly.

The values indicated in % of normal can be transformed to INR according to the enclosed reference curve table. The indication of the prothrombin time in INR, as recommended by the WHO, is based on the following formula:

$$INR = \left(\frac{PT \text{ patient plasma (sec)}}{PT \text{ normal plasma (sec)}} \right)^{ISI}$$

REFERENCE RANGE

Normal time – approximate value: 13 – 17 sec.

Normal range in % of normal: 70 - > 100 %

Therapeutic range of oral anticoagulants: INR of 2.0 – 4.5 of 20 – 42% of normal, resp.

STANDARDISATION

The reference batch is calibrated against the Secondary Int. Reference Preparation for rabbit brain thromboplastin reagent (plain), RBT/79.

The standardisation is performed by use of the "Simplified calibration procedure" (Loeliger et al. Thrombosis Diathesis Haemorrh 33-172, 1975 and Loeliger et al. Thrombosis and Haemostasis 40,272, 1978).

LIMITATION OF THE TEST

Inhibitors of the Lupus type anticoagulant can influence prothrombin time and lead to INRs that do not accurately reflect the true level of anticoagulation.

The choice of anticoagulant (i.e. oxalate instead of citrate) and the condition of the specimen (e.g. hemolyzed, lipemic, parenteral feeding, etc.) may affect results. The latter is particularly true of optical instrumentation measurement of the PT.

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Following correlation (INR) was obtained in comparing TECHNOPLASTIN HIS with :

Recombiplastin 2G	y = 1.073x	R² = 0.9837
Thromborel S	y = 1.0082x	R² = 0.9927

LITERATURE

Moll, St and Ortel, Th. L. Monitoring Warfarin Therapy in Patients with Lupus Anticoagulants. Ann Int Med 127, 177-185 (1997).

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Coagulation Control AK - English

INTENDED USE

Coagulation Control AK can be used for

- Accuracy control of prothrombin time determination in the therapeutic range of oral anticoagulant therapy.
 - in INR
 - in % of normal
- Precision control of prothrombin time in the therapeutic range of oral anticoagulant therapy.
 - in INR
 - in % of normal

SUMMARY

Coagulation Control AK is prepared from plasma of patients stabilized on long-term oral anticoagulant therapy. Therefore the factors of the prothrombin complex (II, VII, IX and X) have the same coagulation activity as plasma from patients on oral anticoagulant therapy. The PIVKA inhibitors (Proteins induced by Vitamin K Absence) are present in Coagulation Control AK as they are in patient's plasma. Coagulation Control AK contains stabilizers but no bactericide additives.

REAGENTS

The [REF](#) 5011050 Coagulation Control AK contains:

	Reagent / Content	Description
5 x 1 mL	Coagulation Control AK	Human plasma, lyophilized

The [REF](#) 5011060 Coagulation Control AK contains:

	Reagent / Content	Description
50 x 1 mL	Coagulation Control AK	Human plasma, lyophilized

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Warning and precautions

- IVD for in vitro diagnostics use.
- This kit is intended for use by personnel trained in laboratory procedures and universal precautions for the use of chemicals and potentially biohazardous substances.
- All human blood or plasma products as well as samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Calibrators and control plasmas are made from human blood and any individual plasma involved in the procedure are tested HbsAg, HIV 1/2 Ab and HCV-Ab-negative. However, all human blood products should be handled as potentially infectious material.
- Get a Material Safety Data Sheet for this product from [www.technoclone.com](#).

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TEST PROCEDURE

Preparation of Coagulation Control AK

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- Coagulation Control AK*: Reconstitute each control vial with 1.0 mL of distilled water. Allow the reconstituted material to stand at room temperature for 10 minutes followed by swirling of the vial before use.
 - For standardization a reconstitution time of 30 minutes is recommended.

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LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of Coagulation Control AK and the method has to adhere strictly as prescribed for the respective reagent.

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INTERPRETATION OF RESULTS

The figures in the table are only applicable to the indicated lot number of Coagulation Control AK and the method has to adhere strictly as prescribed for the respective reagent.

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LITERATURE

Please contact Technoclone ([www.technoclone.com](#)) or your local distributor.

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Coagulation Control N - English

INTENDED USE

For use as accuracy control of screening tests (PT/aPTT/TT/Fibrinogen) and determination of various coagulation factors within normal coagulation activity.

Precision control, i.e. the recording of the day to day variation, is an important procedure in coagulation testing; checking the accuracy of the test system is particularly important.

Coagulation Control N can be used for:

- Accuracy control of all screening tests, i.e. prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT), fibrinogen determination (Clauss), and of various coagulation factors in the normal coagulation activity.

- Precision control of all parameters listed in the table.

SUMMARY

Coagulation Control N is prepared from selected citrated plasma donations of healthy donors. The clotting activity is normally distributed. Thereby the "average" presence of all coagulation factors and inhibitors is secured. Coagulation Control N contains stabilizers but no bactericide additives.

REAGENTS

The REF 5020040 Coagulation Control N contains:

	Reagent / Content	Description
5 x 1 mL	Coagulation Control N	Human plasma, lyophilized

The REF 5020050 Coagulation Control N contains:

	Reagent / Content	Description
50 x 1 mL	Coagulation Control N	Human plasma, lyophilized

Material required (not supplied with the kit)

- Distilled water
- Precision pipettes
- Laboratory timer

Warning and precautions

- IVD for in vitro diagnostics use.
- This kit is intended for use by personnel trained in laboratory procedures and universal precautions for the use of chemicals and potentially biohazardous substances.
- All human blood or plasma products as well as samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Calibrators and control plasmas are made from human blood and any individual plasma involved in the procedure are tested HbsAg, HIV 1/2 Ab and HCV-Ab-negative. However, all human blood products should be handled as potentially infectious material.
- Get a Material Safety Data Sheet for this product from www.technoclone.com.

Stability and storage

The expiry date printed on the labels is only applicable to storage of the unopened containers at +2...8 °C.

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Use the reconstituted plasma in single factor determinations within 2 hours, stored at room temperature.

TEST PROCEDURE

Preparation of Coagulation Control N

Before starting the test, all the required components are to be brought to room temperature.

When reconstituting plasmas, mixing reagents or buffers avoid foaming.

- Coagulation Control N*: Reconstitute each control vial with 1.0 mL of distilled water. Allow the reconstituted material to stand at room temperature for 10 minutes followed by swirling of the vial before use.

For standardization a reconstitution time of 30 minutes is recommended.

Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of Coagulation Control N and the method has to adhere strictly as prescribed for the respective reagent.

The guidelines for evaluation given in this leaflet have been established assuming that a constant sensitivity of different batches is guaranteed by the manufactures of each reagent.

INTERPRETATION OF RESULTS

The figures in the table are only applicable to the indicated lot number of Coagulation Control N and the method has to adhere strictly as prescribed for the respective reagent.

STANDARDISATION

The reported values were determined over multiple runs using a specific lot of reagent and against an Internal Reference Standard which is traceable to the current International Standards, identified in the acceptance range table.

For the tests where International Standards are not available, these parameters have been assigned against an Internal Reference Standard which is traceable to a frozen normal plasma pool of 100 donors.

LITERATURE

Please contact Technoclone (www.technoclone.com) or your local distributor.

EDITORIAL NOTE

This document is available in several languages. The translations have been done using the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.

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

IF-2026-81838177-APN-INPM#ANMAT

PRODUCT DESCRIPTION**INTENDED USE**

For use as accuracy control of screening tests (PT/aPTT/TT/Fibrinogen) within normal coagulation activity.

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TECHNOCLOT® Control N can be used for:

- Accuracy control for screening tests, prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT), fibrinogen determination (Clauss), in the normal coagulation activity.
- Precision control of all parameters listed in the table.

COMPOSITION

TECHNOCLOT® Control N is prepared from selected citrated plasma donations of healthy donors. The clotting activity is normally distributed. Thereby the "average" presence of all coagulation factors and inhibitors is secured. TECHNOCLOT® Control N contains stabilizers but no bactericide additives.

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- Pipettes
- Distilled water

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of TECHNOCLOT® Control N are tested and found negative for HbsAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

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*=room temperature ¹

TEST PROCEDURE**PREPARATION AND PERFORMANCE OF THE TEST**

- Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial until the product is completely reconstituted (avoid frothing).
- Allow the reconstituted TECHNOCLOT® Control N to stand for 10 minutes at room temperature before use ². Invert to mix before use (avoid foaming).
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

ANALYSES RESULTS**EVALUATION OF THE RESULTS**

The figures in the table are only applicable to the indicated lot number of TECHNOCLOT® Control N and the method has to adhere strictly as prescribed for the respective reagent.

STANDARDIZATION AND TRACEABILITY

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LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of TECHNOCLOT® Control N and the method has to adhere strictly as prescribed for the respective reagent.

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¹ Use the reconstituted plasma in single factor determinations within 2 hours stored at room temperature.

² For standardization a reconstitution time of 30 minutes is recommended.

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Coagulation Control A - English

INTENDED USE

For use as accuracy control of screening tests (PT/aPTT/TT/Fibrinogen) and determination of various coagulation factors in the interface between normal and abnormal coagulation activity.

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Coagulation Control A can be used for:

- Accuracy control of all screening tests, i.e. prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT), fibrinogen determination (Claus), and of various coagulation factors in the interface between normal and abnormal coagulation activity.

- Precision control of all parameters listed in the table.

SUMMARY

Coagulation Control A is a lyophilized abnormal human citrated plasma, where the clotting factors are reduced. Coagulation Control A contains stabilizers but no bactericide additives.

REAGENTS

The REF 5021055 Coagulation Control A contains:

	Reagent / Content	Description
5 x 1 mL	Coagulation Control A	Human plasma, lyophilized

The REF 5021060 Coagulation Control A contains:

	Reagent / Content	Description
50 x 1 mL	Coagulation Control A	Human plasma, lyophilized

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Stability opened/ in use:

Reagent	18...25 °C	2...8 °C	< -20 °C
Coagulation Control A	4 hours	8 hours	1 month

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TEST PROCEDURE

Preparation of Coagulation Control A

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- Coagulation Control A*: Reconstitute each control vial with 1.0 mL of distilled water. Allow the reconstituted material to stand at room temperature for 10 minutes followed by swirling of the vial before use.

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LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of Coagulation Control A and the method has to adhere strictly as prescribed for the respective reagent.

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INTERPRETATION OF RESULTS

The figures in the table are only applicable to the indicated lot number of Coagulation Control A and the method has to adhere strictly as prescribed for the respective reagent.

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
IF-2026-81838177-APN-INPM#ANMAT

PRODUCT DESCRIPTION**INTENDED USE**

For use as accuracy control of screening tests (PT/aPTT) in the interface between normal and abnormal coagulation activity.

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TECHNOCLOT® Control A can be used for:

- Accuracy control for prothrombin time (PT), partial thromboplastin time (aPTT), in the interface between normal and abnormal coagulation activity.
- Precision control of all parameters listed in the table.

COMPOSITION

TECHNOCLOT® Control A is an lyophilized abnormal human citrated plasma, where the clotting factors are reduced. TECHNOCLOT® Control A contains stabilizers but no bactericide additives.

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- Pipettes
- Distilled water

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of TECHNOCLOT® Control A are tested and found negative for HbsAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

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- Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial until the product is completely reconstituted (avoid frothing).
- Allow the reconstituted TECHNOCLOT® Control A to stand for 10 minutes at room temperature before use¹. Invert to mix before use (avoid foaming).
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

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The figures in the table are only applicable to the indicated lot number of TECHNOCLOT® Control A and the method has to adhere strictly as prescribed for the respective reagent.

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LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of TECHNOCLOT® Control A and the method has to adhere strictly as prescribed for the respective reagent.

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¹ For standardization a reconstitution time of 30 minutes is recommended.

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
IF-2026-81838177-APN-INPM#ANMAT

PRODUCT DESCRIPTION**INTENDED USE**

Siron LS is an activated partial thromboplastin time (aPTT) reagent to be used:

- as a screening test for the intrinsic coagulation system
 - in the diagnosis and treatment of haemophilia A and B
 - for specific determinations of factors VIII, IX, X, XI and XII when used in combination with the corresponding deficiency plasma, for control of heparin treatment
 - for detection of Lupus anticoagulant
- Siron LS measures the overall activity of the intrinsic coagulation system. Siron LS is sensitive to the factors XII, XI, X, IX, VIII, V and II, but insensitive to factors VII and XIII.

COMPOSITION

Siron LS is an activated partial thromboplastin time (aPTT) test. This liquid preparation is a stable aqueous suspension of phospholipids. FXII activation of Siron LS is achieved by addition of ellagic acid.

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- Buffer**

[REF] 5277015 CaCl₂ 25 mmol/L 100 mL

- Control Plasma Normal and Abnormal **

[REF] 5020040 Coagulation Control N 5 x 1 mL

[REF] 5021055 Coagulation Control A 5 x 1 mL

[REF] 5220130 Coagulation Reference 5 x 1 mL

** or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- RUO for research use only
- Siron LS contains no human material. FXII activation of Siron LS is achieved by ellagic acid and phospholipids from synthetic origins.
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8°C.

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*=room temperature (+18...25°C)

TEST PROCEDURE**PREPARATION OF PLASMA SAMPLES**

Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L or 0.13 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905) to obtain platelet poor plasma.

STORAGE OF PLASMA

Samples of plasma are stored in plastic tubes or siliconized vials at room temperature before testing:

- For patients that are not on heparin therapy, within **3 hours**
- For patients on heparin therapy, within **1 hour**
- For longer period - 20°C, 2 weeks
- 70°C, 6 months

PREPARATION OF REAGENT

Each vial contains 4 mL / 10 mL ready to use aPTT reagent. Before each use, let the reagent stabilize for 30 min at room temperature (+18...25°C) while shaking the vial from time to time.

PERFORMANCE OF THE TEST**CEVERON**

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

The CaCl₂ 25 mmol/L should be preheated to +37°C.

Pipetting Scheme:

100 µL plasma sample
+ 100 µL Siron LS in a pre-warmed tube
incubate for 3 minutes
+ 100 µL CaCl ₂ solution (25 mmol/L; +37°C) Start stop-watch simultaneously
Leave the tube for 20 sec. before tilting. Determine the clotting end point when clot is formed. The results are given in seconds.

ANALYSES RESULTS**REFERENCE RANGE**

Normal range: 25 - 38 seconds

The normal range refers to analyzers with optical detection (Ceveron® alpha). It is recommended that each laboratory establish and controls its normal range.

Siron LS in control of heparin therapy:

Siron LS seconds of 2 - 4 times the normal has been found to correspond to the accepted therapeutic level of heparin. However, the actual therapeutic level should be decided from the clinical situation.

STANDARDISATION

No international calibrator is available for the standardization of aPTT tests. Siron LS documentation is based on studies testing plasma samples from normal subjects and heparinised plasma samples. An in-house reference batch is established in order to avoid batch-to batch variation.

LIMITATION OF THE TEST

- The blood samples should be taken in non-activating samples tubes in order to avoid activation of the clotting factors in the plasma sample.
- Plasma samples should be stored at room temperature to avoid cold activation.
- Platelet poor plasma samples should be used.
- It is important to use the correct incubation time, which is 3 minutes for Siron LS when testing both heparinised and normal plasma samples.
- The calcium concentration used for re-calcification should be 25 mmol/L when testing normal, abnormal or heparin plasmas.
- The REFERENCE RANGE represents approximate values only. Each laboratory should determine its own normal range to take account of laboratory specific factors, such as the type of coagulometer used, reagent temperature, patient population, etc...
- The values found when testing control plasma normal and control plasma abnormal should be compared to the value given on the data-key for the corresponding lot of Control plasma.
- If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of Siron LS and additional reagents (e.g. deficient plasmas) where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.
- Clotting time is not affected in this assay by samples containing ≤ 1000 mg/dL Haemoglobin, ≤ 46 mg/dL unconjugated bilirubin and ≤ 1000 mg/dL triglycerides. However elevated levels of conjugated bilirubin causes a prolonged clotting time.

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	21	21	20	20
Mean (sec)	28.2	73.3	28.2	65.7
CV (%)	0.99	0.92	1.96	2.3

COMPARISON OF METHODS OR CORRELATION

Following correlation was obtained in comparing SIRON LS:

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
PRODUCT DESCRIPTION**INTENDEND USE**

Siron LIS (lupus insensitive) is a liquid activated partial thromboplastin time (aPTT) reagent with **low sensitivity to lupus anticoagulants** to be used:

- as a screening test for the intrinsic coagulation system
- for specific determinations of factors VIII, IX, XI and XII when used in combination with the corresponding deficiency plasma.
- for control of heparin treatment.

COMPOSITION

Siron LIS is a liquid preparation of a stable aqueous suspension of phospholipids. FXII activation of Siron LIS is achieved by addition ellagic acid.

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- Buffer*			
REF	5277015	CaCl ₂ 25 mmol/L	100 mL
- Control Plasma Normal and Abnormal *			
REF	5020040	Coagulation Control N	5 x 1 mL
REF	5021055	Coagulation Control A	5 x 1 mL

*or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- Siron LIS contains no human material. FXII activation of Siron is achieved by ellagic acid and phospholipids from natural origins.
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

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Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L or 0.13 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905).

STORAGE OF PLASMA SAMPLES

Plasma samples are stored in plastic tubes or siliconized vials at room temperature before testing:

- For patients that are not on heparin therapy, testing required within 3 hours
- For patients on heparin therapy, testing required within 1 hour
- For long term storage
 - 20°C, 2 weeks
 - 70°C, 6 months

PREPARATION OF REAGENT

Each vial contains 4mL / 10mL ready to use aPTT reagent. Before each use, let the reagent stabilize for 30 min at room temperature (18...25°C).

Some settling of the activator may occur during storage; thus, repeatedly invert Siron LIS gently just before use to ensure homogeneity. Continuous **stirring is not required** during use.

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MANUAL

The CaCl₂ 25 mmol/L should be preheated to 37°C.

Pipetting Scheme:

100 µL plasma sample
+ 100 µL Siron LIS in a prewarmed tube
incubate for 3 minutes
+ 100 µL CaCl ₂ solution (25 mmol/L; 37°C) Start stop-watch simultaneously
Leave the tube for 20 sec. before tilting. Determine the clotting end point when clot is formed. The results are given in seconds.

ANALYSES RESULTS**REFERENCE RANGE**

Normal range: 30 - 44 seconds

The normal range refers to manual method. It is recommended that each laboratory establish and controls its normal range.

Siron LIS in control of heparin therapy:

A Siron LIS time of 2-4 times of the normal has been found to correspond to the accepted therapeutic level of heparin. However, the actual therapeutic level should be decided from the clinical situation.

STANDARDISATION

No international calibrator is available for the standardization of aPTT tests. Siron LIS documentation is based on studies testing plasma samples from normal subjects and heparinised plasma samples. An in-house reference batch is established in order to avoid batch-to batch variation.

LIMITATION OF THE TEST

- The blood samples should be taken in non-activating samples tubes in order to avoid activation of the clotting factors in the plasma sample.
- Plasma samples should be stored at room temperature to avoid cold activation.
- Platelet poor plasma samples should be used.
- It is important to use the correct incubation time, which is 4 minutes for Siron LIS when testing both heparinised and normal plasma samples.
- The calcium concentration used for recalcification should be 25 mmol/L when testing normal, abnormal or heparin plasmas.
- The reference range represents approximate values only. Each laboratory should determine its own normal range to take account of laboratory specific factors, such as the type of coagulometer used, reagent temperature, patient population, etc...
- The values found when testing control plasma should be compared to the values given on the control batch table for the corresponding lot of Control plasma. If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of Siron LIS and the additional reagents (e.g. deficient plasmas) where a calibration curve is necessary. In addition a new calibration is recommended, if software changes are introduced or after a major service of either instruments or equipment.
- Clotting time is not affected in this assay by samples containing ≤ 1000 mg/dL Haemoglobin, ≤ 46 mg/dL unconjugated bilirubin and ≤ 1000 mg/ml triglycerides. However elevated levels of conjugated bilirubin causes a prolonged clotting time.

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	Intra assay		Inter assay	
Sample	Sample 1	Sample 2	Sample 3	Sample 4
n	21	21	15	15
MV (sec)	33.3	80.9	32.3	79.9
CV (%)	1.54	0.44	1.85	0.45

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Following correlation was obtained in comparing SIRON LIS with:

Actin® FS $y = 1.04x - 2.77$ $R^2 = 0.986$

APPLICATIONS FOR INSTRUMENTS

Application Sheets are available from Technoclone or your local distributor upon request.

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PRODUCT DESCRIPTION

INTENDEND USE

The Activated Partial Thromboplastin Time (aPTT) assay is a universally accepted screening procedure used to detect abnormalities in the intrinsic coagulation System. Dapttin® TC can be used to detect deficiencies of Factors II, V, VIII, IX, X, XI, and XII. In addition, Dapttin® TC can be used to detect the Lupus anticoagulant and is recommended when monitoring heparin therapy since it is sensitive to the presence of heparin.

COMPOSITION

Dapttin® TC (double-activated aPTT) is a standardized aPTT-(activated partial thromboplastin time) reagent consisting of silica and sulphatide as surface activators and a blend of highly purified phospholipids. Response to all factors and inhibitors is optimal.

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- Pipette: 100 µL
- Distilled water
- Buffers

REF	5277015	CaCl ₂ 25 mmol/L	100 mL
-	-	Control Plasma Normal and Abnormal **	
REF	5020040	Coagulation Control N	5 x 1 mL
REF	5021055	Coagulation Control A	5 x 1 mL

** or any other package sizes, special Ceveron® alpha or Technoclot Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- Dapttin® TC contains no human material
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

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Upon storage, caps should be screwed tightly. * = room temperature (20...25 °C)
Repeated freezing is not recommended

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905). Store the plasma at room temperature and use within 4 hours. Stability at -20°C: 1 month.

PREPARATION OF REAGENT

The lyophilised reagent should be dissolved in the volume of distilled water (room temperature) indicated. Allow the dissolved reagent to stand for 10 minutes before use.

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When using equipment with a mechanical measuring system, it is recommended to preheat Dapttin® TC to 37°C.

Pipetting Scheme:

0.10 mL plasma sample
+ 0.10 mL DAPTIN® TC
Shake briefly and incubate for two minutes at 37°C
+ 0.10 mL CaCl ₂ solution (25 mmol/L; 37°C)
Start timing after the addition of the CaCl ₂ solution and determine the clotting end point

ANALYSES RESULTS

REFERENCE RANGE

29 – 42 seconds

For quality control a normal control plasma (i.e. Coagulation Control N) and abnormal control plasma (i.e. Coagulation Control A) should be tested in each test series.

Dapttin® TC has a low opacity and is suitable for use in all commercially available automated coagulometers. The reagent can be used in conjunction with standard machine aPTT programs and the incubation time varied to suit the machine parameters (range 2 – 5 minutes).

STANDARDISATION

No international calibrator is available for the standardization of aPTT tests. Dapttin® TC documentations are based on studies testing plasma samples from normal subjects and heparinised plasma samples. A house reference batch is established in order to avoid batch-to-batch variation.

LIMITATION OF THE TEST

- Only citrate should be used as anticoagulant (neither EDTA nor heparin)
- The blood samples should be taken in non-activating samples tubes in order to avoid activation of the clotting factors in the plasma sample.
- Plasma samples should be stored at room temperature to avoid cold activation.
- Platelet Poor Plasma samples should be used.
- It is important to use the correct incubation time, when testing both heparinised and normal plasma samples.
- The calcium concentration used for recalcification should be 25 mmol/L when testing normal, abnormal or heparin plasmas.
- The REFERENCE RANGE represents approximate values only. It is recommended that each laboratory should determine its own normal range to take account of laboratory specific factors, such as the type of coagulometer used, reagent temperature, patient population, etc...
- The values found when testing Control Plasma Normal and Control plasma Abnormal should be compared to the value given on the data-key for the corresponding lot of Control plasma.
- If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of Dapttin® TC and additional reagents (e.g. deficient plasmas) where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.

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PRECISION

Reproducibility was determined with different samples (in series and day to day).

The following results were obtained:

	Intra assay		Inter assay	
Sample	Sample 1	Sample 2	Sample 1	Sample 2
n	12	12	10	10
MV sec	34.0	85.7	33.1	80.1
SD (%)	0.394	0.478	0.408	1.44
CV (%)	1.16	0.56	1.23	1.80

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Following correlation (%) was obtained in comparing DAPTIN TC (Technoclone) with Actin FS (Siemens) in aPTT method:

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
PRODUCT DESCRIPTION

INTENDENT USE

Thrombin Reagent is used for the determination of the Thrombin Time (TT) in order to check:

- the fibrinogen level (fibrin degradation products)
- heparin therapy
- fibrinolytic therapy

COMPOSITION

The standardised Thrombin Reagent contains bovine thrombin.

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- Pipette: 200 µL - Distilled water
- Control plasma normal and abnormal**

REF 5020040 Coagulation Control N 5 x 1 mL

REF 5021055 Coagulation Control A 5 x 1 mL

** or any other package sizes, special Ceveron® alpha or TECHNOCLLOT® Control and Calibration reagents of Technoclone.

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PREPARATION OF PLASMA SAMPLES

Plasma separation:

Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L) and centrifuge for 15 min at a RCF of at least 2500 g (corresponding to DIN 58905). Store the plasma at room temperature.

PREPARATION OF REAGENT

Carefully open the vial and reconstitute the lyophilised reagent in 6 mL (2 mL in case of samples containing heparin) distilled water by slowly rotating the vial. Allow the reconstituted reagent to stand for 10 minutes at room temperature. For standardisation a reconstitution time of 30 min is recommended.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheet supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Pipetting scheme:

0.20 mL citrated plasma (room temperature)
Incubate for 1 minute at 37°C ¹
+ 0.20 mL Thrombin Reagent (room temperature)
determine the point of coagulation

¹ When using plastic tubes, incubate for 2 minutes at 37°C.

LIMITATION OF THE TEST

The values found when testing normal control plasma (i.e. Coagulation Control N) and abnormal control plasma (i.e. Coagulation Control A) should be compared to the value given on the data-key for the corresponding lot of Control plasma. If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.

If coagulation analyzers with other measurement principles are used, the obtained coagulation times can deviate from the stated confidence range. In this case it is recommended that each laboratory establishes its own reference range.

ANALYSES RESULTS

CALCULATION OF THE RESULTS

The results of the thrombin time are assayed in seconds.

See "REFERENCE RANGE"

- In order to control the normal value, normal control plasma (i.e. Coagulation Control N) should be included in each batch of assays.
- For quality control in the abnormal range, the use of abnormal control plasma (i.e. Coagulation Control A) or a control plasma containing heparin is advisable.

REFERENCE RANGE

NORMAL AND THERAPEUTIC RANGES

	3.3 IU Thrombin/mL (reconstituted in 6 mL distilled water)	10 IU Thrombin/mL (reconstituted in 2 mL distilled water)
Normal range	< 22 sec.	6 – 11 sec.
Heparin therapy*	45 – 120 sec.	20 – 60 sec.
Fibrinolytic therapy (urokinase*)	> 120 sec.	20 – 60 sec.

STANDARDIZATION

No international calibrator is available for the standardization of thrombin time. An in-house reference batch is established in order to avoid batch-to batch variation.

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	Intra assay		Inter assay	
Sample	Sample 1	Sample 2	Sample 1	Sample 2
n	11	12	11	11
MV sec	11.7	15.1	12.1	14.9
SD (%)	0.234	0.273	0.519	0.379
CV (%)	2.0	1.81	4.29	2.55

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Following correlation (%) was obtained in comparing THROMBIN REAGENT (Technoclone) with THROMBIN REAGENT (Roche) in THROMBIN TIME methode:

Thrombin time method: n=94 y = 1.0083x + 3.2242 R² = 0.8123

DETECTION LIMIT

Heparin: LMW: ≤ 1E/ml (3,3IU/ml), ≤ 2E/ml (10IU/ml)

Triglyceride: no detection limit to 500mg/dl

Bilirubin: no detection limit to 0.4mg/dl

LITERATURE

* R. Zimmermann et al.: Laboratory Control of Heparin and Thrombolytic Therapy by Thrombin Clotting Time and Activated Partial Thromboplastin Time. Ärztl. Lab. 33: 121-126 (1987)

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LITERATURA

* R. Zimmermann et al.: Laboratory Control of Heparin and Thrombolytic Therapy by Thrombin Clotting Time and Activated Partial Thromboplastin Time. Arztl. Lab. 33: 121-126 (1987)

Fibrinogen Reagent

GB

IT

PT

SE

NO

RU

DE

ES

FR

DK

BG

GR

REF 5138080 Fibrinogen Reagent 5 x 5 mL

REF 5138085 Fibrinogen Reagent 5 x 2 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / výrobc / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / Θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg bruksvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj uputstvo pre upotrebu
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / počet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Destilisana Voda	LOT	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτία / партида номер / šarže / lot / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционный буфер / Рабочий буферный раствор / Reakční puf / Reakcioni puffer	MTP	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / плъка микротитлоδότησης / Микротитърна плака / Микропланшет / Mikrotitracijska destička / Mikrotitrationsplättchen
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конъюгат / Конъюгат / Konjugát / Konjugat	RTU	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контроль / Контрольный образец / Kontrola / Kontrola	STOP	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / разворете или разредете с / zředí anebo rozpustí v / разбавить или растворить в / nafedte nebo rozpustte v / razrediti ili rastvoriti u	SUB	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační puf / Inkubacioni puffer	WASH	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrat promývacieho roztoku / Koncentrat solucije za ispiranje
IVD	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico in vitro / diagnóstico in vitro / diagnostic in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рηση διαγνωστικηтм εντόб σωλήνα / за ин vitro диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnostický prostředek in vitro / in vitro dijagnostika		

IF-2026-81838177-APN-INPM#ANMAT

Fibrinogen Reagent

GB

PRODUCT DESCRIPTION

INTENDED USE

For the determination of Fibrinogen based on a modified Clauss method. The method for Fibrinogen determination is based on a modification of the method as described by Clauss. Citrated plasma is diluted, mixed with an excess of Thrombin (approx. 80 I.U./mL) and the coagulation time is determined. There is a linear correlation between the logarithm of the coagulation time and the logarithm of the fibrinogen concentration.

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- Solutions / Buffer: Distilled water
- REF 5410012 Imidazole Buffer 90 mL
- Control Plasmas and Calibrators**
- REF 5020040 Coagulation Control N 5 x 1 mL
- REF 5021055 Coagulation Control A 5 x 1 mL
- REF 5220110 Coagulation Reference 5 x 1 mL

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WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

RT*	+2...8 °C / +12°C (Ceveron)
5 days	12 days

Upon storage, caps should be screwed tightly. Avoid contamination with microorganisms.

*room temperature

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Mix 9 parts of venous blood and 1 part of sodium citrate solution (0.11 mol/L) and centrifuge for 15 min at the RCF of at least 2500g (corresponding to DIN 58905). Store the plasma at room temperature (up to 1 day). Keep the plasma in a plastic tube to avoid contact activation. For quality control a normal control plasma (e.g. Coagulation Control N) and abnormal control plasma (e.g. Coagulation Control A) should be tested in each test series.

PREPARATION OF REAGENT

The lyophilized reagents should be reconstituted in the volume of distilled water (room temperature) indicated. Allow the dissolved reagents to stand for 10 minutes at room temperature. For standardized tests a reconstitution time of 30 minutes is recommended.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides Application sheets for Ceveron® alpha. The Application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the Application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Patient plasma and Control Plasmas are diluted 1:10 with Imidazole Buffer (1 volume plasma and 9 volumes of buffer). The optimum range of determination is between 9 and 40 seconds. If the time value is below 9 seconds the determination should be repeated using a higher dilution, e.g. 1:20, and if the time value is above 40 seconds with a lower dilution, e.g. 1:5. Pipetting scheme:

0.20 mL	diluted citrated plasma
	incubate for 1 minute at 37°C
+ 0.20 mL	Fibrinogen reagent (room temperature)
	determine the point of coagulation


ANALYSES RESULTS

TABLE FOR DETERMINATION OF FIBRINOGEN

Each package contains a lot specific table (calibration curve) for determination of fibrinogen. The values for the above mentioned controls are read off this table and have to be in range for the test to be valid. Alternatively, an own calibration curve can be established.

CALIBRATION CURVE

An own calibration curve for a 1:10 dilution of patient's plasma is established using the Coagulation Reference, not included in the package: Reconstitute the Coagulation Reference with 1mL distilled water, prepare a geometric dilution of 1:5 (0.5 mL Coagulation Reference and 2.0 mL buffer) and from this dilution a 1:10, 1:20, and 1:40 dilution (see below) using imidazole Buffer.


- The coagulation times of these dilutions are determined.
- Calculate the fibrinogen concentration of the dilutions from the given value (mg Fibrinogen/100 mL) of the Coagulation Reference by multiplying by 10 and dividing by the appropriate dilution factor.
- Plot the values on double log paper x-axis: concentration in mg / 100 mL, y-axis: coagulation time in seconds

Example:

Dilutions	sec	mg Fibrinogen / 100 mL*
1:5	7.2	490
1:10	14.4	245
1:20	28.3	123
1:40	57.0	61.3

* mg Fibrinogen / 100 mL of the Coagulation Reference (e.g. 245) multiplied by 10 and divided by the dilution factors 5, 10, 20 and 40.

Example of a calibration curve for a 1:10 dilution of patient's plasma:

EVALUATION

Time values of the 1:10 diluted plasma samples can be read off the calibration curve directly. Using other dilutions the Fibrinogen values found have to be calculated as follows:

$$\frac{\text{fibrinogen value (table or calibration curve)}}{10} \times \text{dilution factor} = \frac{\text{mg fibrinogen}}{100 \text{ mL plasma}}$$

REFERENCE RANGE

180-450 mg Fibrinogen / 100mL

STANDARDISATION

The standardization of the lot specific table for determination was done by Technoclone, by using Coagulation Reference.

LIMITATION OF THE TEST

The tilting method is not suited for determination of Fibrinogen based on the Clauss method. We, therefore, recommend the hooking method or a coagulometer to determine the coagulation times. Some auto analyzers, based on optical clot detection, need special applications.

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	Intra assay		Inter assay	
Sample	Sample 1	Sample 2	Sample 1	Sample 2
n	12	12	12	12
MV %	284.2	101.7	267.4	122.7
SD (%)	13.79	7.177	12.27	4.437
CV (%)	4.85	7.06	4.59	3.62

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Following correlation (%) was obtained in comparing Fibrinogen determination based on a modified Clauss method with Fibrinogen (Technoclone) and Fibrinogen (Diagnostica STAGO):
Fibrinogen method: n=192 y = 0.986821x + 4.61536 R² = 0.911212

LINEARITY

Fibrinogen method: 60 - 700 mg/dl

INTERFERENCES

No interference for: Heparin: UFH: ≤ 2 IU/ml
CBn Fibrinogen Fragment: ≤ 10mg/ml
Triglyceride: ≤ 500mg/dl
Bilirubin: ≤ 0.4mg/dl

LITERATURE

Please contact Technoclone or your local distributor

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	Intra assay	Sample 1	Sample 2	Inter assay	Sample 1	Sample 2
Sample	Sample 1	Sample 2	Sample 1	Sample 2	Sample 1	Sample 2
n	12	12	12	12	12	12
MV %	284.2	101.7	267.4	122.7		
SD (%)	13.79	7.177	12.27	4.437		
CV (%)	4.85	7.06	4.59	3.62		

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Coagulation Reference - English

INTENDED USE

Calibration plasma for all screening tests and coagulation factors.

Coagulation Reference can be used for:

- Establishing calibration curves for screening tests.
- Establishing calibration curves for all coagulation factors given in the table.
- Accuracy control for all screening tests, i.e. prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT) and determination of fibrinogen (Claus) within the normal range.
- Precision control of all parameters listed in the table.

SUMMARY

Coagulation Reference is prepared from selected citrated plasma donations of healthy donors. The clotting activity is normally distributed. Thereby the "average" presence of all coagulation factors and inhibitors is secured. Coagulation Reference contains stabilizers but no bactericide additives.

REAGENTS

The REF 5220110 Coagulation Reference contains:

	Reagent / Content	Description
5 x 1 mL	Coagulation Reference	Human plasma, lyophilized

The REF 5220120 Coagulation Reference contains:

	Reagent / Content	Description
50 x 1 mL	Coagulation Reference	Human plasma, lyophilized

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Stability and storage

The expiry date printed on the labels is only applicable to storage of the unopened containers at +2...+8 °C.

Stability opened/ in use:

Reagent	18...25 °C
Coagulation Reference	2 hours

TEST PROCEDURE

Preparation of Coagulation Reference

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- Coagulation Reference*: Reconstitute each calibrator vial with 1.0 mL of distilled water. Allow the reconstituted material to stand at room temperature for 10 minutes followed by swirling of the vial before use.

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Freezing of the reconstituted reagent is not recommended due to possible loss of activity.

LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of Coagulation Reference and the method has to adhere strictly as prescribed for the respective reagent.

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INTERPRETATION OF RESULTS

The figures in the table are only applicable to the indicated lot number of Coagulation Reference and the method has to adhere strictly as prescribed for the respective reagent.

STANDARDISATION

The expected values in the batch table are determined according to the manufacturers instructions. Calibration according to WHO by Technoclone.

The Coagulation Reference is calibrated against the respective WHO International Standard or an internal house standard, where no international standard is available. Consult the batch table.

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