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2024-12-31

2°C 8°C

Ortho-Clinical Diagnostics
Penc oed Bridgend
CF35 6PZ, UK
Made in United Kingdom


2023-12-31

CE 0459

SN 0003

Vitros is a trademark
GEM. C295

T56/JEC295/100

(01)10758750032228
(17)241231(11)231231
(10)5555

Reactive Ingredients

1 VITROS Anti-HTLV I/II Calibrator (goat polyclonal anti-HTLV in anti-HTLV negative human plasma with antimicrobial agent, 2 mL).

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CH REP Ortho-Clinical Diagnostics Switzerland GmbH, Grabenstrasse 25, 6340 Baar, Switzerland

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6900448 LOTE DD.MM.AAAA

3) VITROS Immunodiagnostic Products Anti- HTLV I/II Controls

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Reactive Ingredients

3 sets of VITROS Anti-HTLV I/II Controls 1 and 2 (human plasma with antimicrobial agent, 1 mL).

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2024-12-31

2°C 8°C

Ortho-Clinical Diagnostics
Pencoe Bridgend
CF35 5PZ, UK
Made in United Kingdom

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CH REP

Ortho-Clinical Diagnostics Switzerland GmbH, Grabenstrasse 25, 6340 Baar, Switzerland

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IF-2026-90344634-APN-INPM#ANMAT

BIOARS S.A.
BIOQ. CLAUDIA ETCHÉVES
DIRECTOR TÉCNICO

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1) VITROS Immunodiagnostic Products Anti-HTLV I/II Reagent Pack

Embalaje de aluminio del reactivo VITROS^{IP} Products Immunodiagnostic IVD Anti-HTLV I/II Reagent Pack 2°C 8°C LOT 0010 2020-10-08 SN 0000 T50/LGEM1295FP/100 Ortho-Clinical Diagnostics UK	100 Coated Wells VITROS^{IP} Products Immunodiagnostic IVD Anti-HTLV I/II Reagent Pack LOT 0010 2020-10-08 SN 0000 Ortho-Clinical Diagnostics UK T50/LGEM1295RP/100
--	--

2) VITROS Immunodiagnostic Products Anti-HTLV I/II Calibrator

1 Vial de Calibrador

VITROS^{IP} Products
Immunodiagnostic IVD

Anti-HTLV I/II Calibrator

2°C 8°C

LOT 5555

2023-12-31

SN 0003

T52/G1295CL1/100

Ortho-Clinical Diagnostics UK

3) VITROS Immunodiagnostic Products Anti-HTLV I/II Controls

Control 1 (3 viales) VITROS^{IP} Products Immunodiagnostic IVD HTLV I/II Negative Control LOT 5555 2023-12-31 SN 0003 T54/G15295C1/100 Ortho-Clinical Diagnostics UK	Control 2 (3 viales) VITROS^{IP} Products Immunodiagnostic IVD HTLV I/II Positive Control LOT 5555 2023-12-31 SN 0003 T54/G15295C2/100 Ortho-Clinical Diagnostics UK
---	---

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aHTLV

VITROS Immunodiagnostic Products Anti-HTLV I/II Reagent Pack

REF 690 0446

VITROS Immunodiagnostic Products Anti-HTLV I/II Calibrator

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
- Robert-Guroff M, Weiss SH, Giron JA, et al. Prevalence of antibodies to HTLV-I, -II, and -III in intravenous drug abusers from an AIDS endemic region. *JAMA* 1986;255(2):3133-3137.
- Lee HH, Swanson P, Rosenblatt JD, et al. Relative prevalence and risk factors of HTLV-I and HTLV-II infection in US blood donors. *Lancet* 1991;337:1435-1439.
- Okochi K, Sato H, Hinuma Y. A retrospective study on transmission of adult T-cell leukemia virus by blood transfusion: seroconversion in recipients. *Vox Sang* 1984;46:245-253.
- Sato H, Okochi K. Transmission of human T-cell leukemia virus (HTLV-I) by blood transfusion: demonstration of proviral DNA in recipients' blood lymphocytes. *Int J Cancer* 1986;37:395-400.
- Kaplan JE, Litchfield B, Rouault C, et al. HTLV-I associated myelopathy associated with blood transfusion in the United States: epidemiologic and molecular evidence linking donor and recipient. *Neurology* 1991;41:192-197.
- Saxton EH, Lee H, Swanson P, et al. Detection of human T-cell leukemia/lymphoma virus type I in a transfusion recipient with chronic myelopathy. *Neurology* 1989;39:841-844.
- Donegan E, Pell P, Lee H, Shaw GM, Mosley JW. Transmission of human T-lymphotropic virus type I by blood components from a donor lacking anti-p24: a case report. The Transfusion Safety Study Group. *Transfusion* 1992;32:68-71.
- Hino S, Yamaguchi K, Katimine S, et al. Mother-to-child transmission of human T-cell leukemia virus type I. *Jpn J Cancer Res* 1985;76:474-480.
- Kajiyama W, Kashiwagi S, Ikematsu H, Hayashi J, Nomura H, Okochi K. Intrafamilial transmission of adult T cell leukemia virus. *J Infect Dis* 1986;154:851-857.
- Sullivan MT, Williams AE, Fang CT, Notari EP, Poiesz BJ, Ehrlich GD. Human T-lymphotropic virus (HTLV) types I and II infection in sexual contacts and family members of blood donors who are seropositive for HTLV type I or II. *Transfusion* 1993;33:585-590.
- Centers for Disease Control (CDC). Human T-lymphotropic virus type I screening in volunteer blood donors — United States, 1989. *MMWR Morb Mortal Wkly Rep* 1990;39:915-924.
- Blattner WA, Takatsuki K, Gallo RC. Human T-cell leukemia-lymphoma virus and adult T-cell leukemia. *JAMA* 1983;250:1074-1080.
- Hinuma Y, Komoda H, Chosa T, et al. Antibodies to adult T-cell leukemia-virus-associated antigen (ATLA) in sera from patients with ATL and controls in Japan: a nationwide sero-epidemiologic study. *Int J Cancer* 1982;29:631-635.
- Gessain A, Barin F, Vermant JC, et al. Antibodies to human T-lymphotropic virus type-I in patients with tropical spastic paraparesis. *Lancet* 1985;2:407-410.
- Rodgers-Johnson P, Gajdusek DC, Morgan OSC, Zaninovic V, Sarin PS, Graham DS. HTLV-I and HTLV-II antibodies and tropical spastic paraparesis. *Lancet* 1985;2:1247-1248.
- Osame M, Usuku K, Izumo S, et al. HTLV-I associated myelopathy, a new clinical entity. *Lancet* 1986;1:1031-1032.
- Gessain A, Gout O. Chronic myelopathy associated with human T-lymphotropic virus type I (HTLV-I). *Ann Intern Med* 1992;117:933-946.
- Wiley CA, Nerenberg M, Cros D, Soto-Aguilar MC. HTLV-I polymyositis in a patient also infected with the human immunodeficiency virus. *N Engl J Med* 1989;320:992-995.
- Blattner WA. Epidemiology of HTLV-I and associated diseases. In: Blattner WA, ed. Human retrovirology: HTLV. New York: Raven Press, 1990:251-66.
- Gibbs WN, Lofters WS, Campbell M, et al. Non-Hodgkin lymphoma in Jamaica and its relation to adult T-cell leukemia-lymphoma. *Ann Intern Med* 1987;106:361-368.
- Newton M, Cruickshank K, Miller D, et al. Antibody to human T-lymphotropic virus type-I in West-Indian-born UK residents with spastic paraparesis. *Lancet* 1987;1:415-416.
- The T- and B-cell malignancy study group. Statistical analyses of clinico-pathological, virological and epidemiological data on lymphoid malignancies with special references to adult T-cell leukemia/lymphoma: a report of the second nationwide study of Japan. *Jpn J Clin Oncol* 1985;15:517-535.
- Saxinger W, Blattner WA, Levine PH, et al. Human T-cell leukemia virus (HTLV-I) antibodies in Africa. *Science* 1984;225:1473-1476.
- Maloney EM, Biggar RJ, Neel JV, Taylor ME, et al. Endemic human T-cell lymphotropic virus type II infection among isolated Brazilian Amerindians. *J Infect Dis* 1992;166:100-107.
- Babona DV, Nurse GT. HTLV-I antibodies in Papua New Guinea. *Lancet* 1988;2:1148-1149.
- Yanagihara R, Jenkins CL, Alexander SS, Mora CA, Garruto RM. Human T lymphotropic virus type I infection in Papua New Guinea: high prevalence among the Hagahai confirmed by western analysis. *J Infect Dis* 1990;162:649-654.
- Bastian I, Gardner J, Webb D, Gardner I. Isolation of a human T-lymphotropic virus type I strain from Australian aborigines. *J Virol* 1993;67:843-851.
- Kirkland MA, Frasca J, Bastian I. Adult T-cell leukaemia; lymphoma in an aborigine. *Aust NZ J Med* 1991;21:739-741.
- Blayney DW, Blattner WA, Robert-Guroff M, et al. The human T-cell leukemia-lymphoma virus in the southeastern United States. *JAMA* 1983;250:1048-1052.

IF-2026-90344634-APN-INPM#ANMAT

30. Blattner WA, Nomura A, Clark JW, et al. Modes of transmission and evidence for viral latency from studies of human T-cell lymphotropic virus type I in Japanese migrant populations in Hawaii. *Proc Natl Acad Sci USA* 1986;83:4895-4898.
31. Williams AE, Fang CT, Slamon DJ, et al. Seroprevalence and epidemiological correlates of HTLV-I infection in US blood donors. *Science* 1988;240:643-666.
32. Rosenblatt JD, Golde DW, Wachsman W, et al. A second isolate of HTLV-II associated with atypical hairy-cell leukemia. *New Engl J Med* 1986;315:372-377.
33. Murphy EL, Engstrom JW, Miller K, Sacher RA, Busch MP, Hollingsworth CG. HTLV-II associated myelopathy in a 43-year-old woman. REDS Investigators. *Lancet* 1993;341:757-758.
34. Jacobson S, Lehky T, Nishimura M, Robinson S, McFarlin DE, Dhib-Jalbut S. Isolation of HTLV-II from a patient with chronic progressive neurological disease clinically indistinguishable from HTLV-I-associated myelopathy/tropical spastic paraparesis. *Ann Neurol* 1993;33:392-396.
35. Hjelle B, Appenzeller O, Mills R, Alexander S et al. Chronic neurodegenerative disease associated with HTLV-II infection. *Lancet* 1992;339:645-646.
36. Hjelle B, Mills R, Swenson S, Mertz G, Key C, Allen S. Incidence of hairy cell leukemia, mycosis fungoides, and chronic lymphocytic leukemia in first known HTLV-II-endemic population. *J Infect Dis* 1991;163:435-440.
37. Maloney EM, Biggar RJ, Neel JV, et al. Endemic human T-cell lymphotropic virus II infection among isolated Brazilian Amerindians. *J Infect Dis* 1992;166:100-107.
38. Wiktor SZ, Cannon RD, Atkinson WL, et al. Infection with human T-lymphotropic virus types I and II in sexually transmitted disease clinics in Baltimore and New Orleans. *J Infect Dis* 1992;165:920-924.
39. Zella D, Mori L, Sala M, et al. HTLV-II infection in Italian drug abusers. *Lancet* 1990;336:575-576.
40. Hjelle B, Wilson C, Cyrus S, et al. Human T-cell leukemia virus type II infection frequently goes undetected in contemporary US blood donors. *Blood* 1993;81:1641-1644.
41. Summers M et al. Luminogenic Reagent Using 3-Chloro 4-Hydroxy Acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. *Clin Chem*. 41:S73; 1995.
42. CLSI. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
43. Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.
44. Calam RR. Specimen Processing Separator Gels: An Update. *J Clin Immunoassay*. 11:86-90; 1988.
45. CLSI. *Collection of Diagnostic Venous Blood Specimens. 7th ed*. CLSI standard GP41. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
46. CLSI. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. CLSI document GP39-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
47. CLSI. *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline - Fourth Edition*. CLSI Guideline C24, Wayne, PA: Clinical and Laboratory Standards Institute; 2016.
48. CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline – Third Edition*. CLSI document EP07-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
49. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. ed. 4. Washington, D.C.: AACC Press, 1995.
50. Friedman RB, Young DS. *Effects of Disease on Clinical Laboratory Tests*. ed. 3. Washington, D.C.: AACC Press; 1997.
51. Tryding N, Tufvesson C, Sonntag O (eds). *Drug Effects in Clinical Chemistry*. Ed. 7. Stockholm: The National Corporation of Swedish Pharmacies, Pharmasoft AB, Swedish Society for Clinical Chemistry; 1996.
52. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

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
BIOARS S.A.
BIOQ CLAUDIA ETCHÉVES
DIRECTOR TÉCNICO

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aHTLV

VITROS Immunodiagnostic Products
Controles Anti-HTLV I/II

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VITROS Immunodiagnostic Products Controles Anti-HTLV I/II

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1. CLSI. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Cuarta edición*. CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
2. Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

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© Ortho Clinical Diagnostics, 2020–2025