Item# 3001-3906-70
rev. 03\25
Draft 10

402628-01

XXXXXXXXXX
CÓDIGO DE LOTE

YYYY-MM-DD
FECHA DE FABRICACIÓN

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5x4-2001U.XXX
REFERENCIA

BARCODE
FPO -
magenta does not print

ORAQUICK®
AUTOPRUEBA DE VIH

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ENGLISH

ORAQUICK® HIV SELF-TEST

For Outside USA Use Only • *In Vitro* Diagnostic Use • Do Not Reuse

INSTRUCTIONS FOR USE

The OraQuick® HIV Self-Test is an in-vitro diagnostic medical device (IVD) that is used for self-testing of antibodies for HIV-1 and HIV-2 in oral fluid. This test is intended as an aid to detect antibodies to HIV-1 and HIV-2 from infected individuals. **You must follow the test directions carefully to get an accurate result.** Guidance test for self-detection of HIV-1 and HIV-2 without diagnostic value.

WARNING: If you are on HIV treatment you may get a false result. Clinical data has not been collected to demonstrate the performance of OraQuick® HIV Self-Test in individuals that are undergoing PrEP.

Do not eat or drink for at least 15 minutes before you start the test or use mouth cleaning products 30 minutes before you start the test.

VIEW INSTRUCTIONS

oraquick.com

HOW TO USE THE ORAQUICK® HIV SELF-TEST KIT

1

YOU WILL NEED A WAY TO TIME THE TEST

2

Kit contains: **Test Stick**, **Test Tube** (containing 1 mL of fluid), **Test Stand**, **instructions for use**, and **disposal bag.** Remove these items to begin testing.

3

Your test kit contains two pouches.

4

Tear open the pouch containing the **tube.**

5

Remove the cap.

6

DO NOT pour out the liquid.
DO NOT drink.

7


Slide the tube into the **stand.**

8

Tear open pouch containing the **test device** and remove. **DO NOT** touch the flat pad with your fingers. **DO NOT** eat or swallow the preservative.

9


Press the **Flat Pad** firmly against your gum and swab it along your **upper gum once** (fig. 1) and **your lower gum once** (fig. 2).

10

Put the **flat pad** all the way into the tube until it touches the bottom.

11

LEAVE IT THERE for 20 MINUTES before reading the results.
DO NOT read the result after 40 minutes.

INTERPRETING RESULTS Read test results in a well-lit area

HIV POSITIVE RESULT

Two complete lines, even if the line is faint, means you may be HIV POSITIVE and you need to seek additional testing by a trained professional to confirm an HIV diagnosis.

99.4% of people (152 out of 153) correctly reported their result as positive. This means that 1 out of 153 people infected with HIV reported a negative test result. This is called a false negative.

As soon as possible . . .

Visit your nearest HIV Testing Center or Health Facility

HIV NEGATIVE RESULT

IF READ BEFORE 20 MINUTES, RESULT MAY NOT BE CORRECT

ONE LINE next to the "C" and NO line next to the "T", your result is HIV NEGATIVE.

99.0% of people (717/724) correctly reported their result as negative. This means that 7 out of 724 people not infected with HIV reported a positive test result. This is called a false positive.

Seek regular testing. If you may have been exposed to HIV, test again in 3 months.

INVALID RESULT

If there is no line next to the "C" (even when there is a line next to the "T"), the test line or control line are not complete (all the way across the window), or a red background makes it impossible to read the test, the test is not working and should be repeated. **You will need to obtain another test.**

1.8% of study subjects (16 out of 900) failed to obtain a test result.

The test did not work properly.

Visit your nearest HIV Testing Centre or Health Facility to test again.

NOT SURE OF RESULT

You do not know your result or you are unsure of your result. Visit your nearest HIV Testing Centre or Health Facility to test again.

DISPOSE

Remove the test stick, put the cap on the test tube, place in the disposal bag provided and throw away all contents in the normal trash.

PRODUCT INFORMATION 5X4-1000, 5X4-1001, 5X4-2001

Warnings and Precautions

- Most people feel a little bit anxious when taking an HIV test. But, if you feel very anxious about taking the test, you may want to wait until you are calmer to take it, or get tested by your doctor or local clinic.
- Make sure you are taking the test in a place with good lighting.
- If you are on antiretroviral drugs (ARVs), including for HIV treatment, pre-exposure prophylaxis (PrEP), or post-exposure prophylaxis (PEP), you may get a false result
- If you have participated in a HIV vaccine clinical trial, you may get a positive result using this test, but it may not mean that you are infected with HIV. You should seek follow-up with the research group.
- Do not use if the packaging is opened or any of the contents are opened or missing.
- If today is after the expiration date on the outside of the packaging, do not use this test.
- Do not open any of the packets until you are ready to begin your test.
- Do not eat or drink for at least 15 minutes before you start the test or use mouth cleaning products 30 minutes before you start the test..
- Remove dental products such as dentures or any other products that cover your gums prior to the oral collection.
- Do not use the test if it has been exposed to household cleaning products.
- Do not use this test if it has been stored outside the acceptable temperature of 36°-86° F (2°-30°C).

Biological principles of the test:

The OraQuick HIV Self-Test is a manually performed, visually read, 20 minute immunoassay for the qualitative detection of antibodies to HIV-1 and HIV-2. The assay test strip contains synthetic peptides representing the HIV envelope region and a goat-anti-human IgG procedural control immobilized onto a nitrocellulose membrane in the Test (T) zone and the Control (C) zone, respectively.

Limitations of the test

- Oral bleeding may result in an invalid result. If the test result is invalid, visit your nearest testing centre or healthcare facility.
- The OraQuick® HIV Self-Test may not detect HIV infections that have occurred within the last 3 months.
- For a positive result, the intensity of the test line does not necessarily equal the amount of antibody in the specimen.
- This product has not been evaluated for use in self-testing for individuals younger than 12 years of age. For children ages 2-11, testing must be performed by a trained health care worker.

Interfering substances and unrelated medical conditions

If you are HBV, HCV or HTLV (I/II) positive, you may get a false result. It is recommended that users observe a 15 minute wait period after food and drink and a 30 minute wait period after using oral care products.

Kit contents (one of each):

Test stick, test tube (containing 1 mL of fluid), test stand, Instructions for use and disposal bag.

Materials needed but not provided:

- A timer, watch or something that can time 20 to 40 minutes.
- Eyeglasses – If you use eyeglasses to read, you should wear them while reading your test results.

Questions & Answers

1. **What does the test do?**
 The OraQuick® HIV Self-Test is an in-vitro diagnostic home-use test for HIV (HIV-1 and HIV-2) in oral fluid. This test works by looking for your body's response (antibodies) to fighting the HIV virus. A positive result is preliminary and follow-up confirmatory testing is needed.
2. **What kinds of events put me at risk for HIV?**
 - Sex (vaginal, oral or anal) with multiple sex partners
 - Sex with someone who is HIV positive or whose HIV status you don't know
 - Sex between a man and another man
 - Using illegal injected drugs or steroids
 - Shared needles or syringes
 - Exchanged sex for money
 - Having been diagnosed or treated for hepatitis, tuberculosis or a sexually transmitted disease like syphilis
3. **How soon after a risk event can I test myself?**
 This test detects HIV infection if used **3 months after a risk event**. If you want to be tested before 3 months, you should go to your doctor or local healthcare provider.
4. **Why shouldn't I use this test right after a risk event?**
 When you have been infected with the HIV virus, your body tries to defend against the HIV virus by producing antibodies to it. These antibodies can be found in your blood or oral fluid. It takes your body up to 3 months to produce these antibodies at levels that can be detected by this test.
5. **How often should someone test for HIV?**
 How often you should get tested for HIV depends on your circumstances. If you have never been tested for HIV, you should be tested at least once.
 The Centers for Disease Control and Prevention (CDC) recommends being tested at least once a year if you do things that can result in HIV infection. These include:
 - Sex (vaginal, oral or anal) with multiple sex partners
 - Sex with someone who is HIV positive or whose HIV status you don't know
 - Sex between a man and another man
 - Using illegal injected drugs or steroids
 - Shared needles or syringes
 - Exchanged sex for money
 - Having been diagnosed or treated for hepatitis, tuberculosis or a sexually transmitted disease like syphilis

 If you have been exposed to events that could result in HIV infection, you should test 3 months after the event. If you feel you are at risk of being infected with HIV, it is a good idea to test regularly.
6. **What does a negative result mean?**
 A negative result means that the test has not detected HIV antibodies (your body's way of fighting the virus). If it has been at least 3 months since you have had a risk event and you followed the directions carefully, you likely do not have HIV. However, as with any test, there is a chance for false results. This is called a "false negative". To reduce the chance of false results, be sure to use the test correctly. If it has been less than 3 months since you have had a risk event, wait until the full 3 months have passed to take the test. If you did not follow the directions carefully, you should take the test again to be sure your result is accurate.
7. **What should I do if I get a negative result?**
 If you feel you are at risk of being infected with HIV, it is a good idea to test regularly.
8. **Will the results of this test tell me if it is safe to have unprotected sex?**
 No, you should not use this test to make decisions on behavior that may put you at increased risk for HIV, such as having unprotected sex.
9. **What is a false negative?**
 When a person who is actually infected with HIV receives a negative test result even though they are HIV positive.
10. **What can cause a false negative result?**
 A false negative result can occur for the following reasons:
 - **If you have had a risk event less than 3 months prior to taking the test**
 - Incorrectly reading test result as negative
 - Not following the test directions carefully
 - If you wore dental products such as dentures or any other products that cover your gums while swiping your gums
11. **What does a positive result mean?**
 A positive result means that you may have HIV. **A doctor, clinic or healthcare professional must confirm your OraQuick® HIV Self-Test result**.

LABORATORIOS RICHMOND S.A.C.I.F.
Farm. Mariana C. Rodriguez
Co - Directora Técnica
M.N. 15.581

12. **What should I do if I get a positive result?**
 You need to follow up with a doctor or healthcare provider to get a confirmatory test. At that time your doctor or healthcare provider will discuss the next steps that need to be taken.
13. **What is a false positive result?**
 When a person who is not infected with HIV receives a positive test result even though they are HIV negative.
14. **What can cause a false positive result?**
 A false positive result can occur for any of the following reasons:
 - Incorrectly reading test result as positive
 - Not following the test directions carefully
 - Not waiting 15 minutes after eating or drinking, or 30 minutes after using mouth cleaning products before taking the test.
 - If you have participated in a HIV vaccine clinical trial
 - Swiping each gum more than once during oral collection
15. **How can I tell that my test is working correctly?**
 If your test is working correctly you will see a complete line next to the "C" on your Test Stick. If there is no line or an incomplete line next to the "C" your test did not work. In this case, the result is considered 'invalid' and will need to be repeated with a new test device.
16. **Do any drugs or medications affect the test?**
 To date there is no evidence that the use of antibiotics or medications (not HIV related) affect the test results; however, individuals on treatment for HIV should not use this test.

LABORATORIOS RICHMOND S.A.C.I.F.
Lic. Elvira Zini
Apoderada

EXPLANATION OF SYMBOLS

LOT Batch Code	REF Catalog Number	Caution, Consult Accompanying Documents	Consult Instructions for Use
Do Not Reuse	IVD *In Vitro* Diagnostic Medical Device	Manufacturer	EXP Date of Expiration
Temperature Limitation	Use By	Age Restriction	Date of Manufacturing

Name of the manufacturing establishment:
OraSure Technologies, Inc.
2685 Opus Way, Bethlehem, PA 18020 - Estados Unidos
(+1) 610-882-1820 • www.OraSure.com

Imported and secondary packaged by:
Laboratorios Richmond S.A.C.I.F.
Legal Address:
Av. Elcano 4938 - CABA - República Argentina (TE: +54 11 5555 1600)

Technical Director:
Pablo D. Da Pos, MN 13296
Self-assessment product
Authorized by ANMAT N° PM-2997-1

IF-2025-108444612-APN-INPM#ANMAT

402625-01

© 2025 OraSure Technologies, Inc. • OraQuick®, logo design, and configuration are trademarks of OraSure Technologies, Inc.

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Importado y envasado secundariamente por: Laboratorios Richmond S.A.C.I.F.
Referencia: LABORATORIOS RICHMOND S.A.C.I.F. (2997-1) ROTULOS E INSTRUCCIONES DE USO