



Ministerio de Salud
Secretaría de Calidad en Salud
A.N.M.A.T.

CERTIFICADO DE AUTORIZACIÓN DE MODIFICACIÓN

N° rev: 1882-32#0004

El Instituto Nacional de Productos Médicos de la Administración Nacional de Medicamentos, Alimentos y Tecnología Médica (ANMAT) certifica que de acuerdo con lo solicitado por Bio Lap S.A , se autoriza la modificación en el Registro Nacional de Productores y Productos de Tecnología Médica (RPPTM), del PM 1882-32 aprobado según:

Disposición autorizante N° 7229/16 de fecha 05 julio 2016

Disposiciones modificatorias y reválidas N°: 1882-32#0001 de año 2022

Del siguiente/s dato/s característico/s:

DATO IDENTIFICATORIO A MODIFICAR	DATO AUTORIZADO HASTA LA FECHA	MODIFICACION/RECTIFICACION AUTORIZADA
Modelos	Modelos: CT2.8L30 Standard cortical screw - Ø2.8 mm - L 30 mm CT2.8L30-ST Standard cortical screw - Ø2.8 mm - L 30 mm - STERILE CT2.8L32 Standard cortical screw - Ø2.8 mm - L 32 mm CT2.8L32-ST Standard cortical screw - Ø2.8 mm - L 32 mm - STERILE CT2.8L34 Standard cortical screw - Ø2.8 mm - L 34 mm CT2.8L34-ST Standard cortical screw - Ø2.8 mm - L 34 mm - STERILE CT2.8L36 Standard cortical screw - Ø2.8 mm - L 36 mm	Modelos CT2.8L30 Ø2.8 mm non-locking screw - L30 mm CT2.8L30-ST Ø2.8 mm non-locking screw - L30 mm - STERILE CT2.8L32 Ø2.8 mm non-locking screw - L32 mm CT2.8L32-ST Ø2.8 mm non-locking screw - L32 mm - STERILE CT2.8L34 Ø2.8 mm non-locking screw - L34 mm CT2.8L34-ST Ø2.8 mm non-locking screw - L34 mm - STERILE CT2.8L36 Ø2.8 mm non-locking screw - L36 mm CT2.8L36-ST Ø2.8 mm non-locking screw - L36 mm - STERILE CT2.8L38 Ø2.8 mm non-locking screw - L38

CT2.8L36-ST Standard cortical screw - Ø2.8 mm - L 36 mm - STERILE	mm
CT2.8L38 Standard cortical screw - Ø2.8 mm - L 38 mm	CT2.8L38-ST Ø2.8 mm non-locking screw - L38 mm - STERILE
CT2.8L40-ST Standard cortical screw - Ø2.8 mm - L 40 mm - STERILE	CT2.8L40 Ø2.8 mm non-locking screw - L40 mm
CT2.8L42 Standard cortical screw - Ø2.8 mm - L 42 mm	CT2.8L40-ST Ø2.8 mm non-locking screw - L40 mm - STERILE
CT2.8L42-ST Standard cortical screw - Ø2.8 mm - L 42 mm - STERILE	CT2.8L42 Ø2.8 mm non-locking screw - L42 mm
CT2.8L44 Standard cortical screw - Ø2.8 mm - L 44 mm	CT2.8L42-ST Ø2.8 mm non-locking screw - L42 mm - STERILE
CT2.8L44-ST Standard cortical screw - Ø2.8 mm - L 44 mm - STERILE	CT2.8L44 Ø2.8 mm non-locking screw - L44 mm
CT2.8L46 Standard cortical screw - Ø2.8 mm - L 46 mm	CT2.8L44-ST Ø2.8 mm non-locking screw - L44 mm - STERILE
CT2.8L46-ST Standard cortical screw - Ø2.8 mm - L 46 mm - STERILE	CT2.8L46 Ø2.8 mm non-locking screw - L46 mm
CT2.8L48 Standard cortical screw - Ø2.8 mm - L 48 mm	CT2.8L46-ST Ø2.8 mm non-locking screw - L46 mm - STERILE
CT2.8L48-ST Standard cortical screw - Ø2.8 mm - L 48 mm - STERILE	CT2.8L48 Ø2.8 mm non-locking screw - L48 mm
CT2.8L50 Standard cortical screw - Ø2.8 mm - L 50 mm	CT2.8L48-ST Ø2.8 mm non-locking screw - L48 mm - STERILE
CT2.8L50-ST Standard cortical screw - Ø2.8 mm - L 50 mm - STERILE	CT2.8L50 Ø2.8 mm non-locking screw - L50 mm
CT2.8L55 Standard cortical screw - Ø2.8 mm - L 55 mm	CT2.8L50-ST Ø2.8 mm non-locking screw - L50 mm - STERILE
CT2.8L55-ST Standard cortical screw - Ø2.8 mm - L 55 mm - STERILE	CT2.8L55 Ø2.8 mm non-locking screw - L55 mm
CT2.8L60 Standard cortical screw - Ø2.8 mm - L 60 mm	CT2.8L55-ST Ø2.8 mm non-locking screw - L55 mm - STERILE
CT2.8L60-ST Standard cortical screw - Ø2.8 mm - L 60 mm - STERILE	CT2.8L60 Ø2.8 mm non-locking screw - L60 mm
CT3.5L10 Standard cortical screw - Ø3.5 mm - L 10 mm	CT2.8L60-ST Ø2.8 mm non-locking screw - L60 mm - STERILE
CT3.5L10-ST Standard cortical screw - Ø3.5 mm - L 10 mm - STERILE	CT3.5L10 Ø3.5 mm non-locking screw - L10 mm
	CT3.5L10-ST Ø3.5 mm non-locking screw - L10 mm - STERILE
	CT3.5L12 Ø3.5 mm non-locking screw - L12 mm
	CT3.5L12-ST Ø3.5 mm non-locking screw - L12 mm - STERILE
	CT3.5L14 Ø3.5 mm non-locking screw - L14 mm
	CT3.5L14-ST Ø3.5 mm non-locking screw - L14 mm - STERILE
	CT3.5L16 Ø3.5 mm non-locking screw - L16 mm
	CT3.5L16-ST Ø3.5 mm non-locking screw -

CT3.5L12 Standard cortical screw - Ø3.5 mm - L 12 mm	L16 mm - STERILE
CT3.5L12-ST Standard cortical screw - Ø3.5 mm - L 12 mm - STERILE	CT3.5L18 Ø3.5 mm non-locking screw - L18 mm
CT3.5L14 Standard cortical screw - Ø3.5 mm - L 14 mm	CT3.5L18-ST Ø3.5 mm non-locking screw - L18 mm - STERILE
CT3.5L14-ST Standard cortical screw - Ø3.5 mm - L 14 mm - STERILE	CT3.5L20 Ø3.5 mm non-locking screw - L20 mm
CT3.5L16 Standard cortical screw - Ø3.5 mm - L 16 mm	CT3.5L20-ST Ø3.5 mm non-locking screw - L20 mm - STERILE
CT3.5L16-ST Standard cortical screw - Ø3.5 mm - L 16 mm - STERILE	CT3.5L22 Ø3.5 mm non-locking screw - L22 mm
CT3.5L18 Standard cortical screw - Ø3.5 mm - L 18 mm	CT3.5L22-ST Ø3.5 mm non-locking screw - L22 mm - STERILE
CT3.5L18-ST Standard cortical screw - Ø3.5 mm - L 18 mm - STERILE	CT3.5L24 Ø3.5 mm non-locking screw - L24 mm
CT3.5L20 Standard cortical screw - Ø3.5 mm - L 20 mm	CT3.5L24-ST Ø3.5 mm non-locking screw - L24 mm - STERILE
CT3.5L20-ST Standard cortical screw - Ø3.5 mm - L 20 mm - STERILE	CT3.5L26 Ø3.5 mm non-locking screw - L26 mm
CT3.5L22 Standard cortical screw - Ø3.5 mm - L 22 mm	CT3.5L26-ST Ø3.5 mm non-locking screw - L26 mm - STERILE
CT3.5L22-ST Standard cortical screw - Ø3.5 mm - L 22 mm - STERILE	CT3.5L28 Ø3.5 mm non-locking screw - L28 mm
CT3.5L24 Standard cortical screw - Ø3.5 mm - L 24 mm	CT3.5L28-ST Ø3.5 mm non-locking screw - L28 mm - STERILE
CT3.5L24-ST Standard cortical screw - Ø3.5 mm - L 24 mm - STERILE	CT3.5L30 Ø3.5 mm non-locking screw - L30 mm
CT3.5L26 Standard cortical screw - Ø3.5 mm - L 26 mm	CT3.5L30-ST Ø3.5 mm non-locking screw - L30 mm - STERILE
CT3.5L26-ST Standard cortical screw - Ø3.5 mm - L 26 mm - STERILE	CT3.5L32 Ø3.5 mm non-locking screw - L32 mm
CT3.5L28 Standard cortical screw - Ø3.5 mm - L 28 mm	CT3.5L32-ST Ø3.5 mm non-locking screw - L32 mm - STERILE
CT3.5L28-ST Standard cortical screw - Ø3.5 mm - L 28 mm - STERILE	CT3.5L40 Ø3.5 mm non-locking screw - L40 mm
CT3.5L30 Standard cortical screw - Ø3.5 mm - L 30 mm	CT3.5L40-ST Ø3.5 mm non-locking screw - L40 mm - STERILE
CT3.5L30-ST Standard cortical screw - Ø3.5 mm - L 30 mm - STERILE	CT3.5L45 Ø3.5 mm non-locking screw - L45 mm
	CT3.5L45-ST Ø3.5 mm non-locking screw - L45 mm - STERILE
	CT3.5L50 Ø3.5 mm non-locking screw - L50 mm
	CT3.5L50-ST Ø3.5 mm non-locking screw - L50 mm - STERILE
	CT3.5L55 Ø3.5 mm non-locking screw - L55 mm
	CT3.5L55-ST Ø3.5 mm non-locking screw - L55 mm - STERILE
	CT3.5L60 Ø3.5 mm non-locking screw - L60 mm

CT3.5L32 Standard cortical screw - Ø3.5 mm - L 32 mm	mm
CT3.5L32-ST Standard cortical screw - Ø3.5 mm - L 32 mm - STERILE	CT3.5L60-ST Ø3.5 mm non-locking screw - L60 mm - STERILE
CT3.5L40 Standard cortical screw - Ø3.5 mm - L 40 mm	CT3.5L65 Ø3.5 mm non-locking screw - L65 mm
CT3.5L40-ST Standard cortical screw - Ø3.5 mm - L 40 mm - STERILE	CT3.5L65-ST Ø3.5 mm non-locking screw - L65 mm - STERILE
CT3.5L45 Standard cortical screw - Ø3.5 mm - L 45 mm	CT3.5L70 Ø3.5 mm non-locking screw - L70 mm
CT3.5L45-ST Standard cortical screw - Ø3.5 mm - L 45 mm - STERILE	CT3.5L70-ST Ø3.5 mm non-locking screw - L70 mm - STERILE
CT3.5L50 Standard cortical screw - Ø3.5 mm - L 50 mm	FTS1 Straight plate - Diaphyseal fibula - Symmetrical - Size 1
CT3.5L50-ST Standard cortical screw - Ø3.5 mm - L 50 mm - STERILE	FTS1-ST Straight plate - Diaphyseal fibula - Symmetrical - Size 1 - STERILE
CT3.5L55 Standard cortical screw - Ø3.5 mm - L 55 mm	
CT3.5L55-ST Standard cortical screw - Ø3.5 mm - L 55 mm - STERILE	H1.4QT4.0L40 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L40 mm
CT3.5L60 Standard cortical screw - Ø3.5 mm - L 60 mm	H1.4QT4.0L45 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L45 mm
CT3.5L60-ST Standard cortical screw - Ø3.5 mm - L 60 mm - STERILE	H1.4QT4.0L50 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L50 mm
CT3.5L65-ST Standard cortical screw - Ø3.5 mm - L 65 mm - STERILE	H1.4QT4.0L55 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L55 mm
	H1.4QT4.0L60 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L60 mm
CT3.5L70 Standard cortical screw - Ø3.5 mm - L 70 mm	H1.4QT4.0L65 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L65 mm
CT3.5L70-ST Standard cortical screw - Ø3.5 mm - L 70 mm - STERILE	H1.4QT4.0L70 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L70 mm
FTS1 Plate for diaphyseal fibula fracture - Size 1	QDT2.8L10 Ø2.8 mm non-locking screw - L10 mm
FTS1-ST Plate for diaphyseal fibula fracture - Size 1 - STERILE	QDT2.8L10-ST Ø2.8 mm non-locking screw - L10 mm - STERILE
H1.4QT4.0L30 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 30 mm	QDT2.8L12 Ø2.8 mm non-locking screw - L12 mm
H1.4QT4.0L30-ST Self-drilling compressive screw - Ø4.0 mm -	QDT2.8L12-ST Ø2.8 mm non-locking screw - L12 mm - STERILE
	QDT2.8L14 Ø2.8 mm non-locking screw - L14 mm

cannulated Ø1.4 mm - short thread - L 30 mm - STERILE	QDT2.8L14-ST Ø2.8 mm non-locking screw - L14 mm - STERILE
H1.4QT4.0L35 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 35 mm	QDT2.8L16 Ø2.8 mm non-locking screw - L16 mm
H1.4QT4.0L35-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 35 mm - STERILE	QDT2.8L16-ST Ø2.8 mm non-locking screw - L16 mm - STERILE
H1.4QT4.0L40 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 40 mm	QDT2.8L18 Ø2.8 mm non-locking screw - L18 mm
H1.4QT4.0L40-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 40 mm - STERILE	QDT2.8L18-ST Ø2.8 mm non-locking screw - L18 mm - STERILE
H1.4QT4.0L45 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 45 mm	QDT2.8L20 Ø2.8 mm non-locking screw - L20 mm
H1.4QT4.0L45-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 45 mm - STERILE	QDT2.8L20-ST Ø2.8 mm non-locking screw - L20 mm - STERILE
H1.4QT4.0L50 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 50 mm	QDT2.8L22 Ø2.8 mm non-locking screw - L22 mm
H1.4QT4.0L50-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 50 mm - STERILE	QDT2.8L22-ST Ø2.8 mm non-locking screw - L22 mm - STERILE
H1.4QT4.0L55 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 55 mm	QDT2.8L24 Ø2.8 mm non-locking screw - L24 mm
H1.4QT4.0L55-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 55 mm - STERILE	QDT2.8L24-ST Ø2.8 mm non-locking screw - L24 mm - STERILE
H1.4QT4.0L60 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 60 mm	QOT3.5L10 Ø3.5 mm non-locking screw - L10 mm
H1.4QT4.0L60-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 60 mm - STERILE	QOT3.5L10-ST Ø3.5 mm non-locking screw - L10 mm - STERILE
	QOT3.5L12 Ø3.5 mm non-locking screw - L12 mm
	QOT3.5L12-ST Ø3.5 mm non-locking screw - L12 mm - STERILE
	QOT3.5L14 Ø3.5 mm non-locking screw - L14 mm
	QOT3.5L14-ST Ø3.5 mm non-locking screw - L14 mm - STERILE
	QOT3.5L16 Ø3.5 mm non-locking screw - L16 mm
	QOT3.5L16-ST Ø3.5 mm non-locking screw - L16 mm - STERILE
	QOT3.5L18 Ø3.5 mm non-locking screw - L18 mm
	QOT3.5L18-ST Ø3.5 mm non-locking screw - L18 mm - STERILE
	QOT3.5L20 Ø3.5 mm non-locking screw - L20 mm
	QOT3.5L20-ST Ø3.5 mm non-locking screw - L20 mm - STERILE
	QOT3.5L22 Ø3.5 mm non-locking screw - L22 mm
	QOT3.5L22-ST Ø3.5 mm non-locking screw - L22 mm - STERILE

H1.4QT4.0L65 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 65 mm	QOT3.5L24 Ø3.5 mm non-locking screw - L24mm
H1.4QT4.0L65-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 65 mm - STERILE	QOT3.5L24-ST Ø3.5 mm non-locking screw - L24 mm - STERILE
H1.4QT4.0L70 Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 70 mm	QT4.0L40 Ø4.0 mm lag screw - L40 mm
H1.4QT4.0L70-ST Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L 70 mm - STERILE	QT4.0L40-ST Ø4.0 mm lag screw - L40 mm - STERILE
QDT2.8L10 Non locking screw - Ø2.8 mm - L 10 mm	QT4.0L45 Ø4.0 mm lag screw - L45 mm
QDT2.8L10-ST Non locking screw - Ø2.8 mm - L 10 mm - STERILE	QT4.0L45-ST Ø4.0 mm lag screw - L45 mm - STERILE
QDT2.8L12 Non locking screw - Ø2.8 mm - L 12 mm	QT4.0L50 Ø4.0 mm lag screw - L50 mm
QDT2.8L12-ST Non locking screw - Ø2.8 mm - L 12 mm - STERILE	QT4.0L50-ST Ø4.0 mm lag screw - L50 mm - STERILE
QDT2.8L14 Non locking screw - Ø2.8 mm - L 14 mm	QT4.0L55 Ø4.0 mm lag screw - L55 mm
QDT2.8L14-ST Non locking screw - Ø2.8 mm - L 14 mm - STERILE	QT4.0L55-ST Ø4.0 mm lag screw - L55 mm - STERILE
QDT2.8L16 Non locking screw - Ø2.8 mm - L 16 mm	QT4.0L60 Ø4.0 mm lag screw - L60 mm
QDT2.8L16-ST Non locking screw - Ø2.8 mm - L 16 mm - STERILE	QT4.0L60-ST Ø4.0 mm lag screw - L60 mm - STERILE
QDT2.8L18 Non locking screw - Ø2.8 mm - L 18 mm	QT4.0L65 Ø4.0 mm lag screw - L65 mm
	QT4.0L65-ST Ø4.0 mm lag screw - L65 mm - STERILE
QDT2.8L20 Non locking screw - Ø2.8 mm - L 20 mm	QT4.0L70 Ø4.0 mm lag screw - L70 mm
QDT2.8L20-ST Non locking screw - Ø2.8 mm - L 20 mm - STERILE	QT4.0L70-ST Ø4.0 mm lag screw - L70 mm - STERILE
QDT2.8L22 Non locking screw - Ø2.8 mm - L 22 mm	RITSM0 Standard plate - Distal tibia - Symmetrical - Size 0
QDT2.8L22-ST Non locking screw - Ø2.8 mm - L 22 mm - STERILE	RITSM0-ST Standard plate - Distal tibia - Symmetrical - Size 0 - STERILE
QDT2.8L24 Non locking screw -	RITSM1 Standard plate - Distal tibia - Symmetrical - Size 1
	RITSM1-ST Standard plate - Distal tibia - Symmetrical - Size 1 - STERILE
	RITSMH1 Hook plate - Distal tibia - Symmetrical - Size 1
	RITSMH1-ST Hook plate - Distal tibia - Symmetrical - Size 1 - STERILE
	RTDLS1 Lateral plate - Distal fibula - Standard - Right - Size 1
	RTDLS1-ST Lateral plate - Distal fibula - Standard - Right - Size 1 - STERILE
	RTDLS2 Lateral plate - Distal fibula - Standard - Right - Size 2
	RTDLS2-ST Lateral plate - Distal fibula - Standard - Right - Size 2 - STERILE
	RTDLS3 Lateral plate - Distal fibula - Standard - Right - Size 3
	RTDLS3-ST Lateral plate - Distal fibula - Standard - Right - Size 3 - STERILE
	RTDLS4 Lateral plate - Distal fibula -

Ø2.8 mm - L 24 mm QDT2.8L24-ST Non locking screw - Ø2.8 mm - L 24 mm - STERILE QOT3.5L10 Non locking screw - Ø3.5 mm - L10 mm QOT3.5L10-ST Non locking screw - Ø3.5 mm - L10 mm - STERILE QOT3.5L12 Non locking screw - Ø3.5 mm - L12 mm QOT3.5L12-ST Non locking screw - Ø3.5 mm - L12 mm - STERILE QOT3.5L14 Non locking screw - Ø3.5 mm - L14 mm QOT3.5L14-ST Non locking screw - Ø3.5 mm - L14 mm - STERILE QOT3.5L16 Non locking screw - Ø3.5 mm - L16 mm QOT3.5L16-ST Non locking screw - Ø3.5 mm - L16 mm - STERILE QOT3.5L18 Non locking screw - Ø3.5 mm - L18 mm QOT3.5L18-ST Non locking screw - Ø3.5 mm - L18 mm - STERILE QOT3.5L20 Non locking screw - Ø3.5 mm - L20 mm QOT3.5L20-ST Non locking screw - Ø3.5 mm - L20 mm - STERILE QOT3.5L22 Non locking screw - Ø3.5 mm - L22 mm QOT3.5L22-ST Non locking screw - Ø3.5 mm - L22 mm - STERILE QOT3.5L24 Non locking screw - Ø3.5 mm - L24 mm QOT3.5L24-ST Non locking screw - Ø3.5 mm - L24 mm - STERILE QT4.0L40 Lag screw - Ø4.0 mm - L 40 mm QT4.0L40-ST Lag screw - Ø4.0 mm - L 40 mm - STERILE QT4.0L45 Lag screw - Ø4.0 mm - L 45 mm	Standard - Right - Size 4 RTDLS4-ST Lateral plate - Distal fibula - Standard - Right - Size 4 - STERILE RTDLS5 Lateral plate - Distal fibula - Standard - Right - Size 5 RTDLS5-ST Lateral plate - Distal fibula - Standard - Right - Size 5 - STERILE RTDQ1 Posterolateral plate - Distal fibula - Right - Size 1 RTDQ1-ST Posterolateral plate - Distal fibula - Right - Size 1 - STERILE RTDQ2 Posterolateral plate - Distal fibula - Right - Size 2 RTDQ2-ST Posterolateral plate - Distal fibula - Right - Size 2 - STERILE RTDQ3 Posterolateral plate - Distal fibula - Right - Size 3 RTDQ3-ST Posterolateral plate - Distal fibula - Right - Size 3 - STERILE RTGLS1 Lateral plate - Distal fibula - Standard - Left - Size 1 RTGLS1-ST Lateral plate - Distal fibula - Standard - Left - Size 1 - STERILE RTGLS2 Lateral plate - Distal fibula - Standard - Left - Size 2 RTGLS2-ST Lateral plate - Distal fibula - Standard - Left - Size 2 - STERILE RTGLS3 Lateral plate - Distal fibula - Standard - Left - Size 3 RTGLS3-ST Lateral plate - Distal fibula - Standard - Left - Size 3 - STERILE RTGLS4 Lateral plate - Distal fibula - Standard - Left - Size 4 RTGLS4-ST Lateral plate - Distal fibula - Standard - Left - Size 4 - STERILE RTGLS5 Lateral plate - Distal fibula - Standard - Left - Size 5 RTGLS5-ST Lateral plate - Distal fibula - Standard - Left - Size 5 - STERILE RTGQ1 Posterolateral plate - Distal fibula - Left - Size 1 RTGQ1-ST Posterolateral plate - Distal fibula - Left - Size 1 - STERILE RTGQ2 Posterolateral plate - Distal fibula - Left - Size 2 RTGQ2-ST Posterolateral plate - Distal fibula - Left - Size 2 - STERILE RTGQ3 Posterolateral plate - Distal fibula - Left - Size 3 RTGQ3-ST Posterolateral plate - Distal fibula
---	--

QT4.0L45-ST Lag screw - Ø4.0 mm - L 45 mm - STERILE	- Left - Size 3 - STERILE
QT4.0L50 Lag screw - Ø4.0 mm - L 50 mm	RTSLN1 Lateral plate - Distal fibula - Narrow - Symmetrical - Size 1
QT4.0L50-ST Lag screw - Ø4.0 mm - L 50 mm - STERILE	RTSLN1-ST Lateral plate - Distal fibula - Narrow - Symmetrical - Size 1 - STERILE
QT4.0L55 Lag screw - Ø4.0 mm - L 55 mm	RTSLN2 Lateral plate - Distal fibula - Narrow - Symmetrical - Size 2
QT4.0L55-ST Lag screw - Ø4.0 mm - L 55 mm - STERILE	RTSLN2-ST Lateral plate - Distal fibula - Narrow - Symmetrical - Size 2 - STERILE
QT4.0L60 Lag screw - Ø4.0 mm - L 60 mm	SDT2.8L10 Ø2.8 mm locking screw - L10 mm
QT4.0L60-ST Lag screw - Ø4.0 mm - L 60 mm - STERILE	SDT2.8L10-STI Ø2.8 mm locking screw - L10 mm - STERILE
QT4.0L65 Lag screw - Ø4.0 mm - L 65 mm	SDT2.8L12 Ø2.8 mm locking screw - L12 mm
QT4.0L65-ST Lag screw - Ø4.0 mm - L 65 mm - STERILE	SDT2.8L12-STI Ø2.8 mm locking screw - L12 mm - STERILE
QT4.0L70 Lag screw - Ø4.0 mm - L 70 mm	SDT2.8L14 Ø2.8 mm locking screw - L14 mm
QT4.0L70-ST Lag screw - Ø4.0 mm - L 70 mm - STERILE	SDT2.8L14-STI Ø2.8 mm locking screw - L14 mm - STERILE
RITSM0 Standard plate for distal tibia - Size 0	SDT2.8L16 Ø2.8 mm locking screw - L16 mm
RITSM0-ST Standard plate for distal tibia - Size 0 - STERILE	SDT2.8L16-STI Ø2.8 mm locking screw - L16 mm - STERILE
RITSM1 Standard plate for distal tibia - Size 1	SDT2.8L18 Ø2.8 mm locking screw - L18 mm
RITSM1-ST Standard plate for distal tibia - Size 1 - STERILE	SDT2.8L18-STI Ø2.8 mm locking screw - L18 mm - STERILE
RITSMH1 Hook plate for distal tibia - Size 1	SDT2.8L20 Ø2.8 mm locking screw - L20 mm
RITSMH1-ST Hook plate for distal tibia - Size 1 - STERILE	SDT2.8L20-STI Ø2.8 mm locking screw - L20 mm - STERILE
RTDLS1 Lateral plate for distal fibula - Standard Right - Size 1	SDT2.8L22 Ø2.8 mm locking screw - L22 mm
RTDLS1-ST Lateral plate for distal fibula - Standard Right - Size 1 - STERILE	SDT2.8L22-STI Ø2.8 mm locking screw - L22 mm - STERILE
RTDLS2 Lateral plate for distal fibula - Standard Right - Size 2	SDT2.8L24 Ø2.8 mm locking screw - L24 mm
RTDLS2-ST Lateral plate for distal fibula - Standard Right - Size 2 - STERILE	SDT2.8L24-STI Ø2.8 mm locking screw - L24 mm - STERILE
RTDLS3 Lateral plate for distal fibula - Standard Right - Size 3	SDT2.8L30 Ø2.8 mm locking screw - L30 mm
RTDLS3-ST Lateral plate for distal fibula - Standard Right - Size 3 - STERILE	SDT2.8L30-STI Ø2.8 mm locking screw - L30 mm - STERILE
RTDLS4 Lateral plate for distal	SDT2.8L32 Ø2.8 mm locking screw - L32 mm
	SDT2.8L32-STI Ø2.8 mm locking screw - L32 mm - STERILE
	SDT2.8L34 Ø2.8 mm locking screw - L34 mm
	SDT2.8L34-STI Ø2.8 mm locking screw - L34 mm - STERILE
	SDT2.8L36 Ø2.8 mm locking screw - L36 mm
	SDT2.8L36-STI Ø2.8 mm locking screw - L36 mm - STERILE
	SDT2.8L38 Ø2.8 mm locking screw - L38 mm
	SDT2.8L38-STI Ø2.8 mm locking screw - L38 mm - STERILE
	SDT2.8L40 Ø2.8 mm locking screw - L40 mm
	SDT2.8L40-STI Ø2.8 mm locking screw - L40 mm - STERILE

fibula - Standard Right - Size 4 RTDLS4-ST Lateral plate for distal fibula - Standard Right - Size 4 - STERILE RTDLS5-ST Lateral plate for distal fibula - Standard Right - Size 5 - STERILE RTDQ1 Posterolateral plate for distal fibula - Right - Size 1 RTDQ1-ST Posterolateral plate for distal fibula - Right - Size 1 - STERILE RTDQ2 Posterolateral plate for distal fibula - Right - Size 2 RTDQ2-ST Posterolateral plate for distal fibula - Right - Size 2 - STERILE RTDQ3 Posterolateral plate for distal fibula - Right - Size 3 RTDQ3-ST Posterolateral plate for distal fibula - Right - Size 3 - STERILE RTGLS1 Lateral plate for distal fibula - Standard Left - Size 1 RTGLS1-ST Lateral plate for distal fibula - Standard Left - Size 1 - STERILE RTGLS2 Lateral plate for distal fibula - Standard Left - Size 2 RTGLS2-ST Lateral plate for distal fibula - Standard Left - Size 2 - STERILE RTGLS3 Lateral plate for distal fibula - Standard Left - Size 3 RTGLS3-ST Lateral plate for distal fibula - Standard Left - Size 3 - STERILE RTGLS4 Lateral plate for distal fibula - Standard Left - Size 4 RTGLS4-ST Lateral plate for distal fibula - Standard Left - Size 4 - STERILE RTGLS5 Lateral plate for distal fibula - Standard Left - Size 5 RTGLS5-ST Lateral plate for distal fibula - Standard Left - Size 5 - STERILE RTGQ1 Posterolateral plate for distal fibula - Left - Size 1	mm - STERILE SDT2.8L42 Ø2.8 mm locking screw - L42 mm SDT2.8L42-STI Ø2.8 mm locking screw - L42 mm - STERILE SDT2.8L44 Ø2.8 mm locking screw - L44 mm SDT2.8L44-STI Ø2.8 mm locking screw - L44 mm - STERILE SDT2.8L46 Ø2.8 mm locking screw - L46 mm SDT2.8L46-STI Ø2.8 mm locking screw - L46 mm - STERILE SDT2.8L48 Ø2.8 mm locking screw - L48 mm SDT2.8L48-STI Ø2.8 mm locking screw - L48 mm - STERILE SDT2.8L50 Ø2.8 mm locking screw - L50 mm SDT2.8L50-STI Ø2.8 mm locking screw - L50 mm - STERILE SDT2.8L55 Ø2.8 mm locking screw - L55 mm SDT2.8L55-STI Ø2.8 mm locking screw - L55 mm - STERILE SDT2.8L60 Ø2.8 mm locking screw - L60 mm SDT2.8L60-STI Ø2.8 mm locking screw - L60 mm - STERILE SOT3.5L10 Ø3.5 mm locking screw - L10 mm SOT3.5L10-ST Ø3.5 mm locking screw - L10 mm - STERILE SOT3.5L12 Ø3.5 mm locking screw - L12 mm SOT3.5L12-ST Ø3.5 mm locking screw - L12 mm - STERILE SOT3.5L14 Ø3.5 mm locking screw - L14 mm SOT3.5L14-ST Ø3.5 mm locking screw - L14 mm - STERILE SOT3.5L16 Ø3.5 mm locking screw - L16 mm SOT3.5L16-ST Ø3.5 mm locking screw - L16 mm - STERILE SOT3.5L18 Ø3.5 mm locking screw - L18 mm SOT3.5L18-ST Ø3.5 mm locking screw - L18 mm - STERILE SOT3.5L20 Ø3.5 mm locking screw - L20 mm SOT3.5L20-ST Ø3.5 mm locking screw - L20 mm - STERILE SOT3.5L22 Ø3.5 mm locking screw - L22 mm SOT3.5L22-ST Ø3.5 mm locking screw - L22 mm - STERILE SOT3.5L24 Ø3.5 mm locking screw - L24 mm SOT3.5L24-ST Ø3.5 mm locking screw - L24 mm - STERILE WASH-T4 Compression washer for Ø4.0 mm compressive screws WASH-T4-ST Compression washer for Ø4.0 mm compressive screws - STERILE
---	--

RTGQ1-ST Posterolateral plate for distal fibula - Left - Size 1 - STERILE	Instrumental
RTGQ2 Posterolateral plate for distal fibula - Left - Size 2	ANC016 2.5 mm hexagonal screwdriver
RTGQ2-ST Posterolateral plate for distal fibula - Left - Size 2 - STERILE	ANC042 Mini set - Base
RTGQ3 Posterolateral plate for distal fibula - Left - Size 3	ANC042/CB Mini set - Cambered lid
RTGQ3-ST Posterolateral plate for distal fibula - Left - Size 3 - STERILE	ANC046C Ø2.0 mm non threaded guide gauge for Ø2.8 mm screws
RTSLN1 Lateral plate for distal fibula - Narrow symmetrical - Size 1	ANC082E 2.0 mm quick coupling hexagonal prehensor screwdriver
RTSLN1-ST Lateral plate for distal fibula - Narrow symmetrical - Size 1 - STERILE	ANC083C 2 in 1 : 2.5 mm hexagonal prehensor screwdriver - Ø3.5 mm countersink
RTSLN2 Lateral plate for distal fibula - Narrow symmetrical - Size 2	ANC084 Ø2.8 mm quick coupling countersink
RTSLN2-ST Lateral plate for distal fibula - Narrow symmetrical - Size 2 - STERILE	ANC088C Ø2.0 mm quick coupling drill bit - L125 mm
SDT2.8L10 Locking screw - Ø2.8 mm - L 10 mm	ANC089C Ø2.7 mm quick coupling drill bit - L125 mm
SDT2.8L10-ST Locking screw - Ø2.8 mm - L 10 mm - STERILE	ANC102 Length gauge for Ø2.8 mm screws
SDT2.8L12 Locking screw - Ø2.8 mm - L 12 mm	ANC103 2.0 mm hexagonal non prehensor screwdriver
SDT2.8L12-ST Locking screw - Ø2.8 mm - L 12 mm - STERILE	ANC107 2.5 mm quick coupling hexagonal non prehensor screwdriver
SDT2.8L14 Locking screw - Ø2.8 mm - L 14 mm	ANC124 Length gauge for Ø3.5 mm screws
SDT2.8L14-ST Locking screw - Ø2.8 mm - L 14 mm - STERILE	ANC186 Ø2.7 mm threaded guide gauge for Ø3.5 mm screws
SDT2.8L16 Locking screw - Ø2.8 mm - L 16 mm	ANC191 Ø2.7 mm non threaded bent guide gauge for Ø3.5 mm screws
SDT2.8L16-ST Locking screw - Ø2.8 mm - L 16 mm - STERILE	ANC252 Fast drilling guide for lateral plate - Distal fibula - Standard - Left (RTGLSx)
SDT2.8L18 Locking screw - Ø2.8 mm - L 18 mm	ANC253 Fast drilling guide for lateral plate - Distal fibula - Standard - Right (RTDLSx)
SDT2.8L18-ST Locking screw - Ø2.8 mm - L 18 mm - STERILE	ANC254/B Activ Ankle set - Base
SDT2.8L20 Locking screw - Ø2.8 mm - L 20 mm	ANC254/C Activ Ankle set - Lid
SDT2.8L20-ST Locking screw - Ø2.8 mm - L 20 mm - STERILE	ANC254/I1 Activ Ankle set - Insert 1 - Implants
SDT2.8L22 Locking screw -	ANC254/I2 Activ Ankle set - Insert 2 - Instruments
	ANC254/I3 Activ Ankle set - Insert 3 - Syndesmosis / Medial malleolus
	ANC254/R Activ Ankle set - Screw rack
	ANC256M Ø2.7 mm quick coupling drill bit - L180 mm
	ANC261M Ø2.7 mm non threaded bent long guide gauge for Ø3.5 and Ø4.0 mm screws
	ANC268C Ø2.0 mm threaded guide gauge for Ø2.8 mm screws
	ANC349 15 cm verbrugge forceps
	ANC350 Ø4.5 mm AO quick coupling handle - Size 1

Ø2.8 mm - L 22 mm SDT2.8L22-ST Locking screw - Ø2.8 mm - L 22 mm - STERILE SDT2.8L24 Locking screw - Ø2.8 mm - L 24 mm SDT2.8L24-ST Locking screw - Ø2.8 mm - L 24 mm - STERILE SDT2.8L30 Locking screw - Ø2.8 mm - L 30 mm SDT2.8L30-ST Locking screw - Ø2.8 mm - L 30 mm - STERILE SDT2.8L32 Locking screw - Ø2.8 mm - L 32 mm SDT2.8L32-ST Locking screw - Ø2.8 mm - L 32 mm - STERILE SDT2.8L34 Locking screw - Ø2.8 mm - L 34 mm SDT2.8L34-ST Locking screw - Ø2.8 mm - L 34 mm - STERILE SDT2.8L36 Locking screw - Ø2.8 mm - L 36 mm SDT2.8L38 Locking screw - Ø2.8 mm - L 38 mm SDT2.8L38-ST Locking screw - Ø2.8 mm - L 38 mm - STERILE SDT2.8L40 Locking screw - Ø2.8 mm - L 40 mm SDT2.8L40-ST Locking screw - Ø2.8 mm - L 40 mm - STERILE SDT2.8L42 Locking screw - Ø2.8 mm - L 42 mm SDT2.8L42-ST Locking screw - Ø2.8 mm - L 42 mm - STERILE SDT2.8L44 Locking screw - Ø2.8 mm - L 44 mm SDT2.8L44-ST Locking screw - Ø2.8 mm - L 44 mm - STERILE SDT2.8L46 Locking screw - Ø2.8 mm - L 46 mm SDT2.8L46-ST Locking screw - Ø2.8 mm - L 46 mm - STERILE SDT2.8L48 Locking screw - Ø2.8 mm - L 48 mm SDT2.8L48-ST Locking screw - Ø2.8 mm - L 48 mm - STERILE SDT2.8L50 Locking screw - Ø2.8 mm - L 50 mm SDT2.8L50-ST Locking screw - Ø2.8 mm - L 50 mm - STERILE	ANC388 2.5 mm quick coupling hexagonal non prehensor screwdriver - cannula Ø1.4 mm ANC414M Quick coupling drill bit Ø2.9 mm - cannula 1.4 mm - L125 mm ANC452 Bending iron ANC454 Fast drilling guide for posterolateral plate - Distal fibula - Left (RTGQx) ANC455 Fast drilling guide for posterolateral plate - Distal fibula - Right (RTDQx) ANC456 Fast drilling guide for lateral plate - Distal fibula - Narrow - Symmetrical (RTSLNx) ANC463 Ø3.5 mm quick coupling countersink ANC503 150 mm reduction forceps ANC504 150 mm pointed reduction forceps ANC542 Ø3.5 mm quick coupling drill bit - L125 mm ANC689 Length gauge for Ø2.8 mm screws - Measures 10 - 60 mm ANC690 Ø2.0 mm quick coupling drill bit - L180 mm ANC691 Ø2.0 mm threaded long guide gauge for Ø2.8 mm screws ANC692 Ø2.0 mm non threaded bent guide gauge for Ø2.8 mm screws ANC1137 Ø4.5 mm cannulated AO quick coupling straight handle with ratchet - Size 2 330.213.120 Pin Ø1.3 - L120 mm 33.0214.120 Pin Ø1.4 - L120 mm
--	--

	SDT2.8L55 Locking screw - Ø2.8 mm - L 55 mm SDT2.8L55-ST Locking screw - Ø2.8 mm - L 55 mm - STERILE SDT2.8L60 Locking screw - Ø2.8 mm - L 60 mm SDT2.8L60-ST Locking screw - Ø2.8 mm - L 60 mm - STERILE SOT3.5L10 Locking screw - Ø3.5 mm - L 10 mm SOT3.5L10-ST Locking screw - Ø3.5 mm - L 10 mm - STERILE SOT3.5L12 Locking screw - Ø3.5 mm - L 12 mm SOT3.5L12-ST Locking screw - Ø3.5 mm - L 12 mm - STERILE SOT3.5L14 Locking screw - Ø3.5 mm - L 14 mm SOT3.5L14-ST Locking screw - Ø3.5 mm - L 14 mm - STERILE SOT3.5L16 Locking screw - Ø3.5 mm - L 16 mm SOT3.5L16-ST Locking screw - Ø3.5 mm - L 16 mm - STERILE SOT3.5L18 Locking screw - Ø3.5 mm - L 18 mm SOT3.5L18-ST Locking screw - Ø3.5 mm - L 18 mm - STERILE SOT3.5L20 Locking screw - Ø3.5 mm - L 20 mm SOT3.5L20-ST Locking screw - Ø3.5 mm - L 20 mm - STERILE SOT3.5L22 Locking screw - Ø3.5 mm - L 22 mm SOT3.5L22-ST Locking screw - Ø3.5 mm - L 22 mm - STERILE SOT3.5L24 Locking screw - Ø3.5 mm - L 24 mm SOT3.5L24-ST Locking screw - Ø3.5 mm - L 24 mm - STERILE WASH-T4 Washer WASH-T4-ST Washer - STERILE Instrumental ANC042 Alians instrumentation set (Mini box) ANC042/CB Mini box - Cambered top ANC046C Ø2.0 mm non	
--	---	--

	threaded guide gauge for Ø2.8 mm screws ANC082E 2.0 mm quick coupling hexagonal prehensor screwdriver ANC083C 2 in 1 : 2.5 mm hexagonal prehensor screwdriver - Ø3.5 mm countersink ANC088C Ø2.0 mm quick coupling drill bit - L 125 mm ANC089C Ø2.7 mm quick coupling drill bit - L 125 mm ANC102 Length gauge for Ø2.8 mm screws ANC103 2.0 mm hexagonal non prehensor screwdriver ANC107 2.5 mm quick coupling hexagonal non prehensor screwdriver ANC124 Length gauge for Ø3.5 mm cortical screws ANC186 Ø2.7 mm threaded guide gauge for Ø3.5 mm screws ANC191 Ø2.7 mm non threaded bent guide gauge for Ø3.5 mm screws ANC252 Fast drilling guide for RTGLSx plates ANC253 Fast drilling guide for RTDLSx plates ANC254/B Activ Ankle instrumentation set - Base ANC254/C Activ Ankle instrumentation set - Top ANC254/I1 Activ Ankle instrumentation set - Implants insert ANC254/I2 Activ Ankle instrumentation set - Instruments insert ANC254/I3 Activ Ankle instrumentation set - Syndesmosis / Medial malleollus insert ANC254/I3 Activ Ankle instrumentation set - Syndesmosis / Medial malleollus	
--	--	--

	insert ANC254/R Activ Ankle instrumentation set - Screw tray ANC256M Ø2.7 mm quick coupling drill bit - L 180 mm ANC261M Ø2.7 mm non threaded bent long guide gauge for Ø3.5 and Ø4.0 mm screws ANC268C Ø2.0 mm threaded guide gauge for Ø2.8 mm screws ANC279/C Activ instrumentation set - Top ANC279/I Activ instrumentation set - Insert ANC349 15 cm verbrugge forceps ANC350 Ø4.5 mm AO quick coupling handle - Size 1 ANC388 2.5 mm quick coupling hexagonal non prehensor screwdriver - cannula Ø1.4 mm ANC414M Quick coupling drill bit Ø2.9 mm - cannula 1.4 mm - L 125 mm ANC452 Bending iron ANC454 Fast drilling guide for RTGQx plates ANC455 Fast drilling guide for RTDQx plates ANC456 Fast drilling guide for RTSLNx plates ANC463 Ø3.5 mm quick coupling countersink ANC503 150 mm reduction forceps ANC504 150 mm pointed reduction forceps ANC542 Ø3.5 mm quick coupling drill bit - L 125 mm ANC654 Holder for fibula plate ANC689 Length gauge for Ø2.8 mm screw - Measures 10 - 60 mm ANC690 Ø2.0 mm quick coupling drill bit - L 180 mm ANC691 Ø2.0 mm threaded long guide gauge for Ø2.8 mm screws ANC692 Ø2.0 mm non threaded	
--	--	--

	bent guide gauge for Ø2.8 mm screws ANC1137 Ø11 mm cannulated quick coupling straight handle with ratchet – Taille 2 33.0213.120 Pin Ø1.3 L120 mm	
Rótulos y/o instrucciones de uso	Instrucciones de uso y rótulos aprobados por Disposición N°7229/16	Se incorpora la siguiente declaración de exclusividad: "Instrumental de uso exclusivo del PM 1882-32 "Sistema de Tornillos y Placas para osteosíntesis de Peroné"

Quedando el producto inscripto con los siguientes datos característicos:

Nombre descriptivo: SISTEMA DE TORNILLOS Y PLACAS PARA OSTEOSÍNTESIS DE PERONÉ E INSTRUMENTAL ASOCIADO

Código de identificación y nombre técnico del producto médico, (ECRI-UMDNS):
13-050-Placas, para Huesos

Marca(s) de (los) producto(s) médico(s): NEWCLIP TECHNICIS

Clase de Riesgo: III

Indicación/es autorizada/s: El sistema de placas y tornillos para osteosíntesis de peroné (distal y diáfisis), está indicado para la fijación de fracturas simples y complejas, osteotomías y pseudoartrosis de la porción distal y diáfisis del peroné

Modelos: Modelos

CT2.8L30 Ø2.8 mm non-locking screw - L30 mm
CT2.8L30-ST Ø2.8 mm non-locking screw - L30 mm - STERILE
CT2.8L32 Ø2.8 mm non-locking screw - L32 mm
CT2.8L32-ST Ø2.8 mm non-locking screw - L32 mm - STERILE
CT2.8L34 Ø2.8 mm non-locking screw - L34 mm
CT2.8L34-ST Ø2.8 mm non-locking screw - L34 mm - STERILE
CT2.8L36 Ø2.8 mm non-locking screw - L36 mm
CT2.8L36-ST Ø2.8 mm non-locking screw - L36 mm - STERILE
CT2.8L38 Ø2.8 mm non-locking screw - L38 mm
CT2.8L38-ST Ø2.8 mm non-locking screw - L38 mm - STERILE
CT2.8L40 Ø2.8 mm non-locking screw - L40 mm
CT2.8L40-ST Ø2.8 mm non-locking screw - L40 mm - STERILE
CT2.8L42 Ø2.8 mm non-locking screw - L42 mm
CT2.8L42-ST Ø2.8 mm non-locking screw - L42 mm - STERILE
CT2.8L44 Ø2.8 mm non-locking screw - L44 mm
CT2.8L44-ST Ø2.8 mm non-locking screw - L44 mm - STERILE
CT2.8L46 Ø2.8 mm non-locking screw - L46 mm
CT2.8L46-ST Ø2.8 mm non-locking screw - L46 mm - STERILE
CT2.8L48 Ø2.8 mm non-locking screw - L48 mm
CT2.8L48-ST Ø2.8 mm non-locking screw - L48 mm - STERILE
CT2.8L50 Ø2.8 mm non-locking screw - L50 mm
CT2.8L50-ST Ø2.8 mm non-locking screw - L50 mm - STERILE
CT2.8L55 Ø2.8 mm non-locking screw - L55 mm
CT2.8L55-ST Ø2.8 mm non-locking screw - L55 mm - STERILE
CT2.8L60 Ø2.8 mm non-locking screw - L60 mm

CT2.8L60-ST Ø2.8 mm non-locking screw - L60 mm - STERILE
 CT3.5L10 Ø3.5 mm non-locking screw - L10 mm
 CT3.5L10-ST Ø3.5 mm non-locking screw - L10 mm - STERILE
 CT3.5L12 Ø3.5 mm non-locking screw - L12 mm
 CT3.5L12-ST Ø3.5 mm non-locking screw - L12 mm - STERILE
 CT3.5L14 Ø3.5 mm non-locking screw - L14 mm
 CT3.5L14-ST Ø3.5 mm non-locking screw - L14 mm - STERILE
 CT3.5L16 Ø3.5 mm non-locking screw - L16 mm
 CT3.5L16-ST Ø3.5 mm non-locking screw - L16 mm - STERILE
 CT3.5L18 Ø3.5 mm non-locking screw - L18 mm
 CT3.5L18-ST Ø3.5 mm non-locking screw - L18 mm - STERILE
 CT3.5L20 Ø3.5 mm non-locking screw - L20 mm
 CT3.5L20-ST Ø3.5 mm non-locking screw - L20 mm - STERILE
 CT3.5L22 Ø3.5 mm non-locking screw - L22 mm
 CT3.5L22-ST Ø3.5 mm non-locking screw - L22 mm - STERILE
 CT3.5L24 Ø3.5 mm non-locking screw - L24 mm
 CT3.5L24-ST Ø3.5 mm non-locking screw - L24 mm - STERILE
 CT3.5L26 Ø3.5 mm non-locking screw - L26 mm
 CT3.5L26-ST Ø3.5 mm non-locking screw - L26 mm - STERILE
 CT3.5L28 Ø3.5 mm non-locking screw - L28 mm
 CT3.5L28-ST Ø3.5 mm non-locking screw - L28 mm - STERILE
 CT3.5L30 Ø3.5 mm non-locking screw - L30 mm
 CT3.5L30-ST Ø3.5 mm non-locking screw - L30 mm - STERILE
 CT3.5L32 Ø3.5 mm non-locking screw - L32 mm
 CT3.5L32-ST Ø3.5 mm non-locking screw - L32 mm - STERILE
 CT3.5L40 Ø3.5 mm non-locking screw - L40 mm
 CT3.5L40-ST Ø3.5 mm non-locking screw - L40 mm - STERILE
 CT3.5L45 Ø3.5 mm non-locking screw - L45 mm
 CT3.5L45-ST Ø3.5 mm non-locking screw - L45 mm - STERILE
 CT3.5L50 Ø3.5 mm non-locking screw - L50 mm
 CT3.5L50-ST Ø3.5 mm non-locking screw - L50 mm - STERILE
 CT3.5L55 Ø3.5 mm non-locking screw - L55 mm
 CT3.5L55-ST Ø3.5 mm non-locking screw - L55 mm - STERILE
 CT3.5L60 Ø3.5 mm non-locking screw - L60 mm
 CT3.5L60-ST Ø3.5 mm non-locking screw - L60 mm - STERILE
 CT3.5L65 Ø3.5 mm non-locking screw - L65 mm
 CT3.5L65-ST Ø3.5 mm non-locking screw - L65 mm - STERILE
 CT3.5L70 Ø3.5 mm non-locking screw - L70 mm
 CT3.5L70-ST Ø3.5 mm non-locking screw - L70 mm - STERILE
 FTS1 Straight plate - Diaphyseal fibula - Symmetrical - Size 1
 FTS1-ST Straight plate - Diaphyseal fibula - Symmetrical - Size 1 - STERILE

H1.4QT4.0L40 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L40 mm

H1.4QT4.0L45 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L45 mm

H1.4QT4.0L50 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L50 mm

H1.4QT4.0L55 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L55 mm

H1.4QT4.0L60 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L60 mm

H1.4QT4.0L65 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L65 mm

H1.4QT4.0L70 Ø4.0 mm compressive screw - cannula Ø1.4 - short thread - L70 mm

QDT2.8L10 Ø2.8 mm non-locking screw - L10 mm

QDT2.8L10-ST Ø2.8 mm non-locking screw - L10 mm - STERILE

QDT2.8L12 Ø2.8 mm non-locking screw - L12 mm

QDT2.8L12-ST Ø2.8 mm non-locking screw - L12 mm - STERILE

QDT2.8L14 Ø2.8 mm non-locking screw - L14 mm

QDT2.8L14-ST Ø2.8 mm non-locking screw - L14 mm - STERILE

QDT2.8L16 Ø2.8 mm non-locking screw - L16 mm

QDT2.8L16-ST Ø2.8 mm non-locking screw - L16 mm - STERILE

QDT2.8L18 Ø2.8 mm non-locking screw - L18 mm

QDT2.8L18-ST Ø2.8 mm non-locking screw - L18 mm - STERILE

QDT2.8L20 Ø2.8 mm non-locking screw - L20 mm

QDT2.8L20-ST Ø2.8 mm non-locking screw - L20 mm - STERILE

QDT2.8L22 Ø2.8 mm non-locking screw - L22 mm

QDT2.8L22-ST Ø2.8 mm non-locking screw - L22 mm - STERILE

QDT2.8L24 Ø2.8 mm non-locking screw - L24 mm

QDT2.8L24-ST Ø2.8 mm non-locking screw - L24 mm - STERILE

QOT3.5L10 Ø3.5 mm non-locking screw - L10 mm

QOT3.5L10-ST Ø3.5 mm non-locking screw - L10 mm - STERILE

QOT3.5L12 Ø3.5 mm non-locking screw - L12 mm

QOT3.5L12-ST Ø3.5 mm non-locking screw - L12 mm - STERILE

QOT3.5L14 Ø3.5 mm non-locking screw - L14 mm

QOT3.5L14-ST Ø3.5 mm non-locking screw - L14 mm - STERILE

QOT3.5L16 Ø3.5 mm non-locking screw - L16 mm

QOT3.5L16-ST Ø3.5 mm non-locking screw - L16 mm - STERILE

QOT3.5L18 Ø3.5 mm non-locking screw - L18 mm

QOT3.5L18-ST Ø3.5 mm non-locking screw - L18 mm - STERILE

QOT3.5L20 Ø3.5 mm non-locking screw - L20 mm

QOT3.5L20-ST Ø3.5 mm non-locking screw - L20 mm - STERILE

QOT3.5L22 Ø3.5 mm non-locking screw - L22 mm

QOT3.5L22-ST Ø3.5 mm non-locking screw - L22 mm - STERILE

QOT3.5L24 Ø3.5 mm non-locking screw - L24mm

QOT3.5L24-ST Ø3.5 mm non-locking screw - L24 mm - STERILE

QT4.0L40 Ø4.0 mm lag screw - L40 mm

QT4.0L40-ST Ø4.0 mm lag screw - L40 mm - STERILE

QT4.0L45 Ø4.0 mm lag screw - L45 mm

QT4.0L45-ST Ø4.0 mm lag screw - L45 mm - STERILE

QT4.0L50 Ø4.0 mm lag screw - L50 mm

QT4.0L50-ST Ø4.0 mm lag screw - L50 mm - STERILE

QT4.0L55 Ø4.0 mm lag screw - L55 mm

QT4.0L55-ST Ø4.0 mm lag screw - L55 mm - STERILE

QT4.0L60 Ø4.0 mm lag screw - L60 mm

QT4.0L60-ST Ø4.0 mm lag screw - L60 mm - STERILE

QT4.0L65 Ø4.0 mm lag screw - L65 mm

QT4.0L65-ST Ø4.0 mm lag screw - L65 mm - STERILE
 QT4.0L70 Ø4.0 mm lag screw - L70 mm
 QT4.0L70-ST Ø4.0 mm lag screw - L70 mm - STERILE
 RITSM0 Standard plate - Distal tibia - Symmetrical - Size 0
 RITSM0-ST Standard plate - Distal tibia - Symmetrical - Size 0 - STERILE
 RITSM1 Standard plate - Distal tibia - Symmetrical - Size 1
 RITSM1-ST Standard plate - Distal tibia - Symmetrical - Size 1 - STERILE
 RITSMH1 Hook plate - Distal tibia - Symmetrical - Size 1
 RITSMH1-ST Hook plate - Distal tibia - Symmetrical - Size 1 - STERILE
 RTDLS1 Lateral plate - Distal fibula - Standard - Right - Size 1
 RTDLS1-ST Lateral plate - Distal fibula - Standard - Right - Size 1 - STERILE
 RTDLS2 Lateral plate - Distal fibula - Standard - Right - Size 2
 RTDLS2-ST Lateral plate - Distal fibula - Standard - Right - Size 2 - STERILE
 RTDLS3 Lateral plate - Distal fibula - Standard - Right - Size 3
 RTDLS3-ST Lateral plate - Distal fibula - Standard - Right - Size 3 - STERILE
 RTDLS4 Lateral plate - Distal fibula - Standard - Right - Size 4
 RTDLS4-ST Lateral plate - Distal fibula - Standard - Right - Size 4 - STERILE
 RTDLS5 Lateral plate - Distal fibula - Standard - Right - Size 5
 RTDLS5-ST Lateral plate - Distal fibula - Standard - Right - Size 5 - STERILE
 RTDQ1 Posterolateral plate - Distal fibula - Right - Size 1
 RTDQ1-ST Posterolateral plate - Distal fibula - Right - Size 1 - STERILE
 RTDQ2 Posterolateral plate - Distal fibula - Right - Size 2
 RTDQ2-ST Posterolateral plate - Distal fibula - Right - Size 2 - STERILE
 RTDQ3 Posterolateral plate - Distal fibula - Right - Size 3
 RTDQ3-ST Posterolateral plate - Distal fibula - Right - Size 3 - STERILE
 RTGLS1 Lateral plate - Distal fibula - Standard - Left - Size 1
 RTGLS1-ST Lateral plate - Distal fibula - Standard - Left - Size 1 - STERILE
 RTGLS2 Lateral plate - Distal fibula - Standard - Left - Size 2
 RTGLS2-ST Lateral plate - Distal fibula - Standard - Left - Size 2 - STERILE
 RTGLS3 Lateral plate - Distal fibula - Standard - Left - Size 3
 RTGLS3-ST Lateral plate - Distal fibula - Standard - Left - Size 3 - STERILE
 RTGLS4 Lateral plate - Distal fibula - Standard - Left - Size 4
 RTGLS4-ST Lateral plate - Distal fibula - Standard - Left - Size 4 - STERILE
 RTGLS5 Lateral plate - Distal fibula - Standard - Left - Size 5
 RTGLS5-ST Lateral plate - Distal fibula - Standard - Left - Size 5 - STERILE
 RTGQ1 Posterolateral plate - Distal fibula - Left - Size 1
 RTGQ1-ST Posterolateral plate - Distal fibula - Left - Size 1 - STERILE
 RTGQ2 Posterolateral plate - Distal fibula - Left - Size 2
 RTGQ2-ST Posterolateral plate - Distal fibula - Left - Size 2 - STERILE
 RTGQ3 Posterolateral plate - Distal fibula - Left - Size 3
 RTGQ3-ST Posterolateral plate - Distal fibula - Left - Size 3 - STERILE
 RTSLN1 Lateral plate - Distal fibula - Narrow - Symmetrical - Size 1
 RTSLN1-ST Lateral plate - Distal fibula - Narrow - Symmetrical - Size 1 - STERILE
 RTSLN2 Lateral plate - Distal fibula - Narrow - Symmetrical - Size 2
 RTSLN2-ST Lateral plate - Distal fibula - Narrow - Symmetrical - Size 2 - STERILE
 SDT2.8L10 Ø2.8 mm locking screw - L10 mm
 SDT2.8L10-STI Ø2.8 mm locking screw - L10 mm - STERILE
 SDT2.8L12 Ø2.8 mm locking screw - L12 mm
 SDT2.8L12-STI Ø2.8 mm locking screw - L12 mm - STERILE
 SDT2.8L14 Ø2.8 mm locking screw - L14 mm
 SDT2.8L14-STI Ø2.8 mm locking screw - L14 mm - STERILE

SDT2.8L16 Ø2.8 mm locking screw - L16 mm
 SDT2.8L16-STI Ø2.8 mm locking screw - L16 mm - STERILE
 SDT2.8L18 Ø2.8 mm locking screw - L18 mm
 SDT2.8L18-STI Ø2.8 mm locking screw - L18 mm - STERILE
 SDT2.8L20 Ø2.8 mm locking screw - L20 mm
 SDT2.8L20-STI Ø2.8 mm locking screw - L20 mm - STERILE
 SDT2.8L22 Ø2.8 mm locking screw - L22 mm
 SDT2.8L22-STI Ø2.8 mm locking screw - L22 mm - STERILE
 SDT2.8L24 Ø2.8 mm locking screw - L24 mm
 SDT2.8L24-STI Ø2.8 mm locking screw - L24 mm - STERILE
 SDT2.8L30 Ø2.8 mm locking screw - L30 mm
 SDT2.8L30-STI Ø2.8 mm locking screw - L30 mm - STERILE
 SDT2.8L32 Ø2.8 mm locking screw - L32 mm
 SDT2.8L32-STI Ø2.8 mm locking screw - L32 mm - STERILE
 SDT2.8L34 Ø2.8 mm locking screw - L34 mm
 SDT2.8L34-STI Ø2.8 mm locking screw - L34 mm - STERILE
 SDT2.8L36 Ø2.8 mm locking screw - L36 mm
 SDT2.8L36-STI Ø2.8 mm locking screw - L36 mm - STERILE
 SDT2.8L38 Ø2.8 mm locking screw - L38 mm
 SDT2.8L38-STI Ø2.8 mm locking screw - L38 mm - STERILE
 SDT2.8L40 Ø2.8 mm locking screw - L40 mm
 SDT2.8L40-STI Ø2.8 mm locking screw - L40 mm - STERILE
 SDT2.8L42 Ø2.8 mm locking screw - L42 mm
 SDT2.8L42-STI Ø2.8 mm locking screw - L42 mm - STERILE
 SDT2.8L44 Ø2.8 mm locking screw - L44 mm
 SDT2.8L44-STI Ø2.8 mm locking screw - L44 mm - STERILE
 SDT2.8L46 Ø2.8 mm locking screw - L46 mm
 SDT2.8L46-STI Ø2.8 mm locking screw - L46 mm - STERILE
 SDT2.8L48 Ø2.8 mm locking screw - L48 mm
 SDT2.8L48-STI Ø2.8 mm locking screw - L48 mm - STERILE
 SDT2.8L50 Ø2.8 mm locking screw - L50 mm
 SDT2.8L50-STI Ø2.8 mm locking screw - L50 mm - STERILE
 SDT2.8L55 Ø2.8 mm locking screw - L55 mm
 SDT2.8L55-STI Ø2.8 mm locking screw - L55 mm - STERILE
 SDT2.8L60 Ø2.8 mm locking screw - L60 mm
 SDT2.8L60-STI Ø2.8 mm locking screw - L60 mm - STERILE
 SOT3.5L10 Ø3.5 mm locking screw - L10 mm
 SOT3.5L10-ST Ø3.5 mm locking screw - L10 mm - STERILE
 SOT3.5L12 Ø3.5 mm locking screw - L12 mm
 SOT3.5L12-ST Ø3.5 mm locking screw - L12 mm - STERILE
 SOT3.5L14 Ø3.5 mm locking screw - L14 mm
 SOT3.5L14-ST Ø3.5 mm locking screw - L14 mm - STERILE
 SOT3.5L16 Ø3.5 mm locking screw - L16 mm
 SOT3.5L16-ST Ø3.5 mm locking screw - L16 mm - STERILE
 SOT3.5L18 Ø3.5 mm locking screw - L18 mm
 SOT3.5L18-ST Ø3.5 mm locking screw - L18 mm - STERILE
 SOT3.5L20 Ø3.5 mm locking screw - L20 mm
 SOT3.5L20-ST Ø3.5 mm locking screw - L20 mm - STERILE
 SOT3.5L22 Ø3.5 mm locking screw - L22 mm
 SOT3.5L22-ST Ø3.5 mm locking screw - L22 mm - STERILE
 SOT3.5L24 Ø3.5 mm locking screw - L24 mm

SOT3.5L24-ST Ø3.5 mm locking screw - L24 mm - STERILE
 WASH-T4 Compression washer for Ø4.0 mm compressive screws
 WASH-T4-ST Compression washer for Ø4.0 mm compressive screws - STERILE
 Instrumental
 ANC016 2.5 mm hexagonal screwdriver
 ANC042 Mini set - Base
 ANC042/CB Mini set - Cambered lid
 ANC046C Ø2.0 mm non threaded guide gauge for Ø2.8 mm screws
 ANC082E 2.0 mm quick coupling hexagonal prehensor screwdriver
 ANC083C 2 in 1 : 2.5 mm hexagonal prehensor screwdriver - Ø3.5 mm countersink
 ANC084 Ø2.8 mm quick coupling countersink
 ANC088C Ø2.0 mm quick coupling drill bit - L125 mm
 ANC089C Ø2.7 mm quick coupling drill bit - L125 mm
 ANC102 Length gauge for Ø2.8 mm screws
 ANC103 2.0 mm hexagonal non prehensor screwdriver
 ANC107 2.5 mm quick coupling hexagonal non prehensor screwdriver
 ANC124 Length gauge for Ø3.5 mm screws
 ANC186 Ø2.7 mm threaded guide gauge for Ø3.5 mm screws
 ANC191 Ø2.7 mm non threaded bent guide gauge for Ø3.5 mm screws
 ANC252 Fast drilling guide for lateral plate - Distal fibula - Standard - Left (RTGLSx)
 ANC253 Fast drilling guide for lateral plate - Distal fibula - Standard - Right (RTDLSx)
 ANC254/B Activ Ankle set - Base
 ANC254/C Activ Ankle set - Lid
 ANC254/I1 Activ Ankle set - Insert 1 - Implants
 ANC254/I2 Activ Ankle set - Insert 2 - Instruments
 ANC254/I3 Activ Ankle set - Insert 3 - Syndesmosis / Medial malleolus

 ANC254/R Activ Ankle set - Screw rack
 ANC256M Ø2.7 mm quick coupling drill bit - L180 mm
 ANC261M Ø2.7 mm non threaded bent long guide gauge for Ø3.5 and Ø4.0 mm screws
 ANC268C Ø2.0 mm threaded guide gauge for Ø2.8 mm screws

 ANC349 15 cm verbrugge forceps
 ANC350 Ø4.5 mm AO quick coupling handle - Size 1
 ANC388 2.5 mm quick coupling hexagonal non prehensor screwdriver - cannula Ø1.4 mm
 ANC414M Quick coupling drill bit Ø2.9 mm - cannula 1.4 mm - L125 mm
 ANC452 Bending iron
 ANC454 Fast drilling guide for posterolateral plate - Distal fibula - Left (RTGQx)
 ANC455 Fast drilling guide for posterolateral plate - Distal fibula - Right (RTDQx)
 ANC456 Fast drilling guide for lateral plate - Distal fibula - Narrow - Symmetrical (RTSLNx)
 ANC463 Ø3.5 mm quick coupling countersink
 ANC503 150 mm reduction forceps
 ANC504 150 mm pointed reduction forceps
 ANC542 Ø3.5 mm quick coupling drill bit - L125 mm

 ANC689 Length gauge for Ø2.8 mm screws - Measures 10 - 60 mm
 ANC690 Ø2.0 mm quick coupling drill bit - L180 mm
 ANC691 Ø2.0 mm threaded long guide gauge for Ø2.8 mm screws
 ANC692 Ø2.0 mm non threaded bent guide gauge for Ø2.8 mm screws

ANC1137 Ø4.5 mm cannulated AO quick coupling straight handle with ratchet - Size 2
330.213.120 Pin Ø1.3 - L120 mm
33.0214.120 Pin Ø1.4 - L120 mm

Período de vida útil: Producto estéril: 5 años
Producto no estéril: No corresponde

Condición de uso: Uso exclusivo a profesionales e instituciones sanitarias

Fuente de obtención de la materia prima de origen biotecnológico: N/A

Forma de presentación: Envase individual estéril
Envase individual No estéril

Método de esterilización: Radiación gamma

Nombre del fabricante: NEWCLIP TECHNICS

Lugar de elaboración: 45 rue des Garottières- PA de la Lande Saint Martin-44115 Haute Goulaine. Francia.

AUTORIZADO

Habiéndose cumplimentado con lo previsto en la Disposición ANMAT N° 64/25 y 9688/19, esta Administración autoriza las modificaciones solicitadas en Argentina, Ciudad de Buenos Aires, 11 septiembre 2026.

**BODNAR
Andrea
Veronica**

Firmado digitalmente
por BODNAR Andrea
Veronica
Fecha: 2026.10.09
14:11:07 -03'00'

Dirección Evaluación y Registro de Productos
Médicos
Firma y Sello

Instituto Nacional de Productos Médicos
ANMAT
Firma y Sello

El presente certificado será válido únicamente cuando se presente junto con las Disposiciones previas del PM enunciadas anteriormente y sea verificado con su código QR a través de la página de ANMAT. Este certificado no modifica el período de vigencia.

Fecha de emisión: 11 septiembre 2026



ABRIOLA Leticia Adriana
CUIL 27323186494



La validez del presente documento deberá verificarse mediante el código QR.

Firmado digitalmente por: TERRIZZANO Maria
Lorena
Fecha y hora: 13.10.2026 09:32:47